JustRight™ 5mm Stapler

REF JR-ST25 -2.0

5 mm shaft diameter, 20 cm shaft length, 25 mm jaw length, 2.0 mm staple

JustRight™ 5mm Reload

REF JR-REL25-2.0

5 mm shaft diameter, 20 cm shaft length, 25 mm jaw length, 2.0 mm staple

- en - Instructions for Use
- da – Brugsanvisning
- de – Gebrauchsanweisung
- el – Οδηγίες χρήσης
- es – Instrucciones de uso
- fi – Käyttöohjeet
- fr – Notice d'utilisation
- is – Notkunarleiðbeiningar
- it – Istruzioni per l'uso
- nl – Gebruiksaanwijzing
- no – Bruksanvisning
- pt – Instruções de utilização
- ro – Instrucțiuni de utilizare
- sv – Bruksanvisning
- tr – Kullanım Kılavuzu

Caution

Read all warnings, cautions, and instructions provided with this instrument before use.

If any incidents occur while using the JustRight™ 5mm Stapler and Reload, users should report these issues directly to Bolder Surgical by phone (866.683.1743) or by email (complaints@boldsurg.com), and to your local competent authority.

Symbols

	Sterilized Using Ethylene Oxide		Do Not Use if Package is Opened or Damaged and Consult Instructions for Use
	Catalog, Reorder or Reference Number		Caution
	Consult Instructions for Use		Store in Temperatures Between 18°C – 28°C
	Do Not Re-use		Transport Conditions Must Not Exceed: -18°C - 60°C
	Manufacturer		Not Made with Natural Rubber Latex
	Do Not Resterilize		Medical Device
	Use-by Date		Keep Dry
	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician		Lot Number
	Date of Manufacture		European Community Authorized Representative
	This package forms the sterile barrier		Plastics - Other

Technical Characteristics for the JustRight™ 5.0mm Stapler and Reload

Size of Staple: 2.0 mm		Length of Shaft: 20 cm	
Compatible Trocar Sizes: ≥5 mm		Distal Jaw Opening: ≤15°	
Staple Pattern: Staggered Rows per Cartridge: 4		Compressed Tissue Range: 0.75 - 1.0 mm	
Indicated Firing Zone Length of Staple Line: 25 mm			
Maximum Number of Reloads per Stapler Handle: 15			
Pre-fire Compression Time: ≥10 seconds			

Description

The JustRight™ 5mm Stapler is a sterile, single patient use device that staples and simultaneously divides tissue. The JustRight™ 5mm Stapler fitted with a JustRight™ 5mm Reload places four 25 mm staggered rows of titanium staples while simultaneously dividing the tissue between the two inner staple rows. The unformed staple leg length is 2 mm and the staple back span dimension is 3 mm. The JustRight™ 5mm Stapler is provided with a stapler handle and one JustRight 5mm Reload staple cartridge. The JustRight™ 5mm Reload cartridges are also provided separately to allow for multiple staple line applications within one surgical procedure. One JustRight™ 5mm Stapler handle can be fired up to 15 times in a single procedure using 15 JustRight™ 5mm Reload cartridges. The JustRight™ 5mm Stapler is compatible for introduction and use through endoscopic cannulas as small as 5 mm. The stapler handle has a rotation knob to allow for 360° rotation of the JustRight™ 5mm Reload cartridge. The JustRight™ 5mm Stapler is used on tissue that compresses to a thickness of 0.75 to 1.0 mm.

Titanium, as found in JustRight Surgical surgical staples, is acceptable for use with MRI and NMR procedures as it is nonferromagnetic. The implanted staples are not subject to displacement or induced torque, and do not result in localized heating. The implantable staples are categorized as MR Safe.

Indications for Use

The JustRight™ 5mm Stapler is intended for use in abdominal, gynecologic, pediatric, and thoracic surgery for resection, transection, and creation of anastomosis.

Contraindications

- Do not use the JustRight™ 5mm Stapler 2.0 mm staples on any tissue that compresses to less than 0.75 mm in thickness, or on any tissue that cannot comfortably compress to 1.0 mm.
- Do not use the JustRight™ 5mm Stapler on tissue such as liver, hepatic vasculature, spleen, pancreas, or biliary structures where compressibility is such that closure of the instrument would be destructive.
- Do not use the JustRight™ 5mm Stapler where adequacy of hemostasis cannot be verified visually after application.
- Do not use the JustRight™ 5mm Stapler with a staple line reinforcement material. Adequate staple formation cannot be verified when using staple line buttressing materials such as PERI-STRIPS or SEAMGUARD products.
- Do not use the JustRight™ 5mm Stapler on the aorta or on tissues that are necrotic, friable, or have altered integrity.

Warnings and Cautions

Warning
This product is intended for single use only; it cannot be adequately cleaned or resterilized for safe reuse. Attempts to clean or sterilize may result in biocompatibility infection or product failure risks to the patient.
Place tissue in the center of the jaws. To avoid incomplete ligation, do not grasp structure beyond the staple line surface; do not place tissue in the jaw hinge. Any tissue extending beyond the cut mark will not be transected.
When using the stapler more than once during a SINGLE surgical procedure, be sure to remove the empty JustRight™ 5mm Reload single use unit and reload a new one. When the JustRight™ 5mm Reload is empty, the blue staple deployment and cut lever will not engage. This prevents inadvertent actuation of the integrated knife.
After staples have been deployed, do not reclamp an empty staple cartridge on tissue or use to manipulate structures. Inadvertent tissue damage may occur.
Preoperative radiotherapy may result in changes to tissue. These changes may, for example, cause the tissue thickness to exceed the indicated range for the staple size. Careful consideration with regard to staple size should be given to any pre-surgical treatment the patient may have undergone.
Device performance testing was conducted preclinically in animals. No human clinical testing was performed.
Always inspect the tissue thickness prior to application of the JustRight™ 5mm Stapler.
Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. Prior to performing any endoscopic procedures, consult the medical literature relative to technique, complications, and hazards.
When using endoscopic instruments and accessories from different manufacturers in a procedure, be sure to verify compatibility and ensure that electrical isolation or grounding is not compromised.
Do not use the JustRight™ 5mm Stapler unless properly trained.

Caution


- After firing, always inspect the staple line for hemostasis. Minor bleeding or oozing may be controlled by electrocautery or manual sutures.
- The anvil of the stapler must be completely visible (past the cannula sleeve) prior to opening the jaws within the body cavity.

The JustRight™ 5mm Stapler is intended to be used with a laparoscopic cannula with a minimum diameter of 5.6 mm.
Always close the blue Top Lever on the stapler handle prior to introducing and removing the stapler from the cannula sleeve.
If packaging is exposed to environmental conditions outside of those specified, do not use.
Inspect packaging for damage. If damaged, do not use.
Avoid placing fingers between the lever and handle, or in the jaws when clamping. Injury to the user may result.
This device is intended for use by medical professionals only.


Preparing the Stapler

Handle - JustRight™ 5mm Stapler

- To insert the staple reload, the top blue lever must be in the closed position.

Reload - JustRight™ 5mm Reload

- Remove the green ship lock from the JR-REL25-2.0 staple reload by depressing the ribbed grip area on the reload while rotating the green ship lock in the direction of the arrow and pulling on it.

NOTE: Do not remove the green ship lock on the Reload until you are ready to load into the handle.

- Align the blue dot on the reload connector with the blue dot on the handle's rotation knob and insert until fully seated.


NOTE: The reload can only be inserted one way; the tabs on the reload connector are different widths and match the slots in the rotation knob.

- Rotate the reload shaft clockwise until the blue dot on the outside of the reload shaft aligns with the blue dot on the rotation knob.


Using the Stapler

1. Insert the instrument into the cannula sleeve.

NOTE: The top blue lever (1) must be in the closed position before attempting to insert the stapler into cannula sleeve.

2. Lift the top blue lever (1) to open the jaws

3. Turn the rotation knob until the jaws are in the required position.


Warning

Establish proximal control of blood vessels prior to stapling where practical and implement methods of blood vessel control in the event of stapler failure

When positioning the stapler on the application site, ensure that no obstruction, such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.

Take care when positioning the stapler on the application site to avoid the unintended stapling of other anatomic structures.

4. Position the tissue to be transected between the black lines on the outside of the instrument jaws.


5. Pull back on the Jaw Clamp Lever (2) until the jaws lock on tissue to be transected.

Warning

If the jaws of the device do not lock closed, the target tissue is outside the compression specifications. Reposition the jaws on to tissue within compression specification range.

If the tissue does not compress to 0.75 to 1.0 mm in thickness it is outside the compression specifications. Do not staple or cut on this tissue. Reposition the jaws on to tissue within the compression specification range.

6. To release and reposition jaws, depress the blue Jaw Release button (5) on the back of the handle.

7. Repeat Step 5.

NOTE: After closing, pausing for 10 seconds prior to firing may result in better tissue compression and better staple formation.

8. Flip the green Safety Lock (3) down. Do not hold the safety lock down.

Safety Lock (green) - Up

Safety Lock (green) - Down

9. Pull back completely on the blue Staple Deployment and Cut Lever (4) and release to reset the tissue structure.

NOTE: Only pull back on the blue Staple Deployment and Cut Lever (4) once to completely fire a staple line.

10. Press the blue Jaw Release button (5) to open the jaws.

11. Push down on the blue Top Lever (1) (closed position) and remove the device from the cannula.

12. Remove and set aside a spent staple Reload after firing.

13. Reload the device if necessary:

- Remove the spent reload by depressing the ribbed grip area on the reload and rotating it until the alignment pins are disengaged, then pull the spent reload from the handle and discard it.
- Remove the green ship lock from a new JR-REL25-2.0 staple reload by depressing the ribbed grip area on the reload while rotating the green ship lock in the direction of the arrow and pulling on it.
- Insert a new reload by aligning the blue dot on the reload connector with the blue dot on the handle's rotation knob and inserting until fully seated.
- Rotate the reload shaft clockwise until the blue dot on the outside of the reload shaft aligns with the blue dot on the rotation knob.

14. The JustRight™ 5mm Stapler may be reloaded and fired up to 15 times in a single procedure.

After Surgery

Warning

Do not reuse or resterilize instruments and accessories labeled "disposable" or "single use only" and dispose according to the facility's policy for biohazards and sharps.

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All rights reserved.
Emergo Europe B.V.