

ThinPrep® Stain Protocol using Rinse II and Bluing II Solutions on the Medite TST 44 Stainer

Use this protocol ONLY when ThinPrep Rinse II and Bluing II Solutions (P/Ns ASY-04875, ASY-04876) are components of the ThinPrep Stain.

Use this protocol for staining slides fixed with 95% alcohol.

Maintain bath heights to completely cover the slides at full immersion. (Note: Bath fill volume = 450 ml)

Caution: Change the Distilled / Instrument Feed Water bath at Station 19 after every **four** racks of slides.

Required: Agitation "ON"				
No.	Bath	Reagent	Time	Cor*
1	Start	70% Reagent or Ethyl alcohol	----	----
2	16	50% Reagent or Ethyl alcohol	01:00	no
3	17	Distilled or Instrument Feed Water (IFW) ¹	01:00	no
4	18	ThinPrep Nuclear Stain	07:00	yes
5	19	Distilled or Instrument Feed Water ¹ (Replace after every 4 racks.)	00:10	yes
6	20	ThinPrep Rinse II Solution	01:00	yes
7	21	Distilled or Instrument Feed Water ¹	00:30	no
8	22	ThinPrep Bluing II Solution	00:30	yes
9	23	Distilled or Instrument Feed Water ¹	00:30	no
10	24	Distilled or Instrument Feed Water ¹	00:30	no
11	25	95% Reagent or Ethyl alcohol	00:30	no
12	26	ThinPrep Orange G Solution	02:00	yes
13	40	95% Reagent or Ethyl alcohol	00:15	yes
14	39	95% Reagent or Ethyl alcohol	00:15	yes
15	38	ThinPrep EA Solution	04:00	yes
16	37	95% Reagent or Ethyl alcohol	01:00	yes
17	36	95% Reagent or Ethyl alcohol	01:00	yes
18	35	100% Reagent or Ethyl alcohol	00:30	no
19	34	100% Reagent or Ethyl alcohol	00:30	no
20	33	100% Reagent or Ethyl alcohol	00:30	no
21	32	Xylene or other Hologic approved clearing agent ²	01:00	no
22	31	Xylene or other Hologic approved clearing agent ²	03:00	no
23	Unload	Xylene or other Hologic approved clearing agent ²	----	----

Remove slides to a separate clearing bath. Coverslip with an appropriate Hologic approved mountant. ²

¹ Hologic specification for IFW: ≥ 1.0 megohm-cm Resistivity **OR** ≤ 1.0 μ Siemens/cm Conductivity (Refer to Clinical and Laboratory Standards Institute (CLSI) document C3-A4, 2006)

² See ThinPrep Stain User's Manual, Section 1, "STAINING" or contact Hologic for current list of Hologic approved clearing agents and mounting media.



Hologic, Inc. • 250 Campus Drive • Marlborough, MA 01752 USA • 1-800-442-9892 • www.hologic.com

Hologic BV • Da Vincilaan 5 • 1930 Zaventem • Belgium

Solution Stations

Medite TST 44 Stainer

Using ThinPrep Rinse II and Bluing II Solutions for staining slides fixed with 95% alcohol

Step Number →

2

Solution →

50%
Alcohol

Bath Station Number →

9

22	21	20	19	18	17	16	15	14	13					
Xylene	Xylene	100% Alcohol	100% Alcohol	100% Alcohol	95% Alcohol	95% Alcohol	ThinPrep EA	95% Alcohol	95% Alcohol					
31	32	33	34	35	36	37	38	39	40	41	42	43	44	45

2	3	4	5	6	7	8	9	10	11	12				
50% Alcohol	Distilled Water	ThinPrep Nuclear	Distilled Water	ThinPrep Rinse II	Distilled Water	ThinPrep Bluing II	Distilled Water	Distilled Water	95% Alcohol	ThinPrep Orange G				
16	17	18	19	20	21	22	23	24	25	26	27	28	29	30

1	1	1	1	23	23	23	23							
Start 70% Alcohol	Start 70% Alcohol	Start 70% Alcohol	Start 70% Alcohol	Unload Xylene	Unload Xylene	Unload Xylene	Unload Xylene							
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15