celero[®]

Breast Biopsy Device
Breast Biopsy Introducer

Instructions For Use
Gebruiksaanwijzing
Mode d'emploi
Bedienungsanleitung
Istruzioni per l'uso
Instruções de Utilização
Instrucciones de uso

HOLOGIC[®]

Celero[®] Biopsy Device

Instructions for Use

Please read all information carefully. Failure to properly follow the instructions may lead to unintended consequences.

Important: This package insert is designed to provide instructions for use (IFU) for the Celero breast biopsy device. It is not a reference to surgical techniques.

Indications (Product codes: Celero-12 and Celero INTRO-12 (if applicable))

The Celero breast biopsy device is indicated for use to obtain percutaneous core biopsy samples from soft tissue and tumors of the breast. This product is intended for diagnostic use only.

Contraindications

Sound professional judgment should be used for patients where increased risk or complications may be associated with biopsy based on the physician's judgment, such as patients receiving anticoagulant therapy or those who have known bleeding disorders.

Device Description

The Celero breast biopsy device is a vacuum-assisted spring loaded core biopsy device which is used to obtain percutaneous core biopsy samples from soft tissue. The Celero breast biopsy device is designed for biopsy and diagnosis only. The device is for single-patient use and the entire handpiece is fully disposable. The device primarily consists of a hollow inner cannula with a side aperture and a sharpened outer cannula that extends at high speeds over the aperture to acquire targeted tissue. During this process, vacuum is generated inside the inner cannula to assist in pulling tissue into the aperture. An introducer system, including a stylet, specifically designed for use with the Celero breast biopsy device is provided separately. This introducer system provides single insertion access to the targeted area of interest and assists in deploying biopsy site identification systems. The Celero breast biopsy device is designed to acquire a maximum of 10 tissue cores.

Device preparation (all steps to be performed per standard interventional technique)

1. Prior to use of the Celero breast biopsy device, inspect the protective packaging and device to verify that neither has been damaged during shipment. If it appears that the packaging has been compromised, do not use the device.
2. Carefully remove the Celero breast biopsy device from its protective packaging.
3. Remove the protective sheath from the Celero breast biopsy device.
4. If it is desired to use the introducer system, remove it from its packaging and remove its protective sheath.
5. Depress and release the black latch on the Celero breast biopsy device to release the actuating lever.
6. Completely depress the lever once to retract the vacuum source.
7. Completely depress the lever a second time to retract the cutting cannula.
8. Completely depress the lever a third time to retract the inner cannula. The lever will automatically latch. Do not proceed with the procedure if the lever does not automatically latch at this point, as the device may not function properly.
9. Check the indicator window. If the window is green, the device is ready for biopsy. If the window is not green at this point, do not proceed with the procedure, as the device may not function properly.

Performing a Biopsy (all steps to be performed per standard interventional technique)

1. Identify targeted area.
2. Prepare the site as required.
3. Insert the tip of the needle and advance to the desired location.
 - a. If using the optional introducer system, insert the introducer with stylet and advance to desired location.
 - b. Remove the stylet, leaving the introducer in place. The stylet may be detached from introducer by depressing the snap mechanism on the introducer hub.
 - c. Advance the Celero breast biopsy device needle to the targeted area by sliding it into the introducer hub and snap into place.

NOTE: If desired, the inner cannula may be fired prior to insertion of the Celero breast biopsy device into the patient.

4. While maintaining the position and the needle orientation of the Celero breast biopsy device, depress the green button to fire the inner cannula (25mm stroke).
5. Depress the blue button to fire the cutting cannula (22mm aperture) and obtain a tissue core. Vacuum will be generated automatically during this step.
6. Remove the Celero breast biopsy device leaving the introducer in place (if applicable). The introducer may be detached from the Celero breast biopsy device by depressing the snap mechanism on the introducer hub.
7. Depress the black latch to release the lever.
8. Completely depress the lever once to retract the vacuum source.
9. Completely depress the lever a second time to retract the cutting cannula and expose the tissue core.
10. Retrieve the tissue core. The lever may be either fully depressed or fully released when retrieving tissue cores.

NOTE: Caution when retrieving tissue cores from the aperture of the Celero breast biopsy device. Do not place fingers near or within the aperture at any time during a biopsy procedure, as the cutting cannula may be armed and could potentially fire, resulting in user injury.

11. To obtain additional cores, completely depress the lever a third time to retract the inner cannula and repeat steps 1-8 of this section.
12. If finished obtaining cores, remove the Celero breast biopsy device and introducer (if applicable) and properly dispose of them.
13. If a biopsy site marker will be placed, the introducer may be left in place.
14. If desired, deploy biopsy site marker according to manufacturer's instructions for use.

Warnings and Precautions

- The Celero breast biopsy device is not recommended for use with MRI.
- The Celero breast biopsy procedure should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.
- The Celero breast biopsy device should be used only by physicians trained in percutaneous biopsy procedures.
- **Rx ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

- Sound professional judgment should be used when the Celero breast biopsy device is used on patients with breast implants
- Avoid operator or instrument contact with the needle portion of the Celero.
- The Celero breast biopsy device should not be tested by “dry” firing the cutting cannula. Damage to the device may occur, dulling the cutting edge.
- In the event of changes in the performance of the device, please dispose of the device and resume the procedure with a new device.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by Hologic, Inc. may not be compatible with the Celero breast biopsy device. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- Do not resterilize and/or reuse the Celero breast biopsy device. Resterilization and/or reuse may compromise the integrity of the instrument. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned and sterilized devices.
- Potential complications are associated with any percutaneous removal/biopsy technique for tissue collection. Potential complications are limited to the region surrounding the biopsy site and include hematoma, hemorrhage, infection, pain and tissue adherence to the biopsy probe while removing it from the breast.

How Supplied

The Celero breast biopsy device is sterilized by gamma radiation and supplied for single patient use. Discard into an appropriate container after use.

As Identified on Labels:

QTY

Number of Devices Enclosed.

YYYY-MM-DD

Expiration date is represented by the following:

YYYY represents the year

MM represents the month

DD represents the day

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.
250 Campus Drive
Marlborough, MA 01752 USA
Phone: 877-371-4372
BreastHealth.Support@hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

EC **REP**

European Representative
Hologic BVBA
Da Vincilaan 5
1930 Zaventem
Belgium
Tel: +32 2 711 46 80

This page is intentionally left blank

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