# SecurMark<sup>®</sup> for celero<sup>®</sup>

# SecurMark<sup>®</sup> ultrasound

Biopsy Site Marker

Titanium Biopsy Site Marker

Stainless Steel Biopsy Site Marker

Instructions For Use

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Mode d'emploi

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**HOLOGIC<sup>®</sup>**

# SecurMark® for Celero® and SecurMark®

## Ultrasound Biopsy Site Marker for use with Ultrasound-guided Breast Biopsy Device

### Instructions for Use (IFU)

Please read all information carefully. Failure to properly follow the instructions may lead to unanticipated surgical consequences.

Important: This package insert is designed to provide Instructions for Use for the SecurMark® for Celero® and SecurMark® Ultrasound biopsy site marker. It is not a reference to surgical techniques.

Upon completion of the breast biopsy procedure, the user will have the option of using the SecurMark for Celero or the SecurMark Ultrasound biopsy site marker.


---

### Indications

The SecurMark biopsy site marker is indicated for the permanent radiographic marking of sites in soft tissue.

### Contraindications

None known.

*\* Available in multiple shapes*

## Device Description

The SecurMark biopsy site marker is a sterile, single patient use device comprised of a single, permanent biocompatible marker surrounded by a bioabsorbable suture-like material and a deployment device. The permanent marker is manufactured from either Titanium or Stainless Steel. Refer to the product label for the marker material.

The deployment device is a hand-held device that delivers the marker from the distal tip. The deployment device consists of a rigid cannula, handle, rigid push rod, and plunger. The marker is located at the distal end of the deployment device.

The SecurMark Ultrasound and SecurMark for Celero biopsy site markers are designed to be inserted into the breast directly or through the Celero Introducer. Other 13G introducers (or coaxial), designed for use with 14G biopsy needles, may be compatible with the SecurMark biopsy site markers. In order for an introducer to be compatible, it must have a minimum inner diameter of 2.15mm and maximum length of 100.05mm.

The titanium or stainless steel marker is classified as magnetic resonance (MRI) conditional at 3.0 Tesla field strength or less. The marker, when present in a patient undergoing an MRI procedure at 3.0 Tesla or less, will not create an additional hazard or risk with respect to magnetic field-related interactions, movement/dislodgement or heating.

## MRI Artifact Considerations

Artifacts for the SecurMark biopsy site marker have been characterized using a 3.0 Tesla MRI system with spin echo and gradient echo pulse sequences. Based on this information, imaging quality may be slightly compromised if the area of interest is in the exact same area as the SecurMark biopsy site marker.

Artifact size is dependent on the type of pulse sequence used for imaging (larger for gradient echo pulse sequences and smaller for spin echo pulse sequences) and the size of the field of view. Image artifact will be smaller for MRI systems with lower static magnetic field strengths using the same imaging parameters as those operating at higher static magnetic field strengths.

## Device Preparation and Use

1. Prior to use of the SecurMark biopsy site marker, inspect the protective packaging to verify that no damage has occurred during shipment. If it appears the packaging has been compromised, do not use the device(s).
2. Remove the devices from their protective packaging using standard interventional technique. Inspect the devices to verify that no damage has occurred during shipment. If it appears a device has been damaged, do not use the device(s).
3. When using a non-Hologic introducer: During the biopsy procedure set-up, insert the SecurMark deployment device through the introducer (or co-axial) to verify compatibility.

4. Perform the biopsy procedure according to the biopsy device's Instructions for Use.

5. Remove the biopsy device from the access site. If using an introducer, hold the introducer steady, release the latching component (if applicable) and pull back on the biopsy device to remove it from the breast.

6. Insert the SecurMark deployment device. If using an introducer, hold the introducer steady, place the distal end of the deployment device into the proximal end of the introducer and carefully advance the deployment device until it engages the introducer hub.

7. Deploy the biopsy site marker by advancing the deployment plunger with your thumb all the way forward until it latches to the hub. Make sure to hold the introducer and/or deployment device steady during deployment of the marker.

\*NOTE: Insert the tip of the deployment device to the distal area of interest and draw the hub of the device back while advancing the deployment plunger to deploy the marker.

8. Verify the deployment and proper position of the marker prior to removal of the device with the appropriate imaging modality.

9. Slowly remove the deployment device, or introducer and deployment device, from the breast as a unit and properly dispose.

## Warnings and Precautions

- A small percentage of patients experience an allergic reaction to stainless steel due to the presence of nickel. It is not recommended to use the stainless steel SecurMark biopsy site marker in patients with known metal allergies.
- The SecurMark deployment devices are not recommended for use within the bore of an MRI magnet.
- The SecurMark biopsy site markers are not recommended for use in patients with breast implants.
- The biopsy site marking procedure should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications, and hazards prior to performing any minimally invasive procedure.
- These devices should be used only by physicians trained in open or percutaneous biopsy procedures.
- The SecurMark biopsy site markers should be deployed into the cavity created during the biopsy procedure. Deployment into tissue outside of the biopsy cavity is not recommended.

- Marker position relative to established landmarks may change under mammography upon subsequent breast compressions.
- The SecurMark biopsy site markers are not intended to be repositioned or recaptured after deployment.
- Users should take care not to unintentionally deploy the marker.
- Excess hematoma within the biopsy cavity and/or introducer can lead to marker adhesion to the deployment device, increasing the risk of marker drag out.
- Care should be taken to avoid damaging the cannula. Avoid operator or instrument contact with the SecurMark biopsy site marker or the distal end of the deployment device.
- The implanted SecurMark biopsy site markers are magnetic resonance imaging (MRI) conditional. The implanted marker presents no additional risk to the patient or operator from magnetic forces, torque, heating, induced voltages, or movement, but it may affect MRI image quality.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by Hologic, Inc. may not be compatible with the SecurMark biopsy site marker device. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- Do not resterilize and/or reuse the SecurMark biopsy site marker device. Resterilization and/or reuse may compromise the integrity of the instrument. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned and sterilized devices.
- **Rx ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- Store the SecurMark biopsy site marker device in a clean and dry area. Avoid storage or handling temperatures above 50°C (122°F).

## How Supplied

The SecurMark biopsy site marker devices are sterilized by gamma radiation and supplied preloaded for single patient use. Discard into an appropriate container after use.

As Identified on Labels:

**QTY**

Number of Devices Enclosed.

YYYY-MM-DD

Expiration date is represented by the following:

YYYY represents the year

MM represents the month

DD represents the day

## For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.  
250 Campus Drive  
Marlborough, MA 01752 USA  
Phone: 877-371-4372  
BreastHealth.Support@hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

**EC** **REP**

European Representative  
Hologic BVBA  
Da Vincilaan 5  
1930 Zaventem  
Belgium  
Tel: +32 2 711 46 80

This page is intentionally left blank

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