celeroMark™

Biopsy Site Identification System
Titanium Biopsy Site Marker

Instructions for Use
Gebruiksaanwijzing
Mode d'emploi
Bedienungsanleitung
Istruzioni per l'uso
Instruções de Utilização
Instrucciones de uso

HOLOGIC®

CeleroMark™ Biopsy Site Marker for Celero® Biopsy Device

Instructions for Use

Please read all information carefully. Failure to properly follow the instructions may lead to unanticipated surgical consequences.

Important: This package insert is designed to provide Instructions for Use for the CeleroMark™ biopsy site marker for the Celero® biopsy device. It is not a reference to surgical techniques.


Upon completion of the breast biopsy procedure, the user will have the option of using the CeleroMark biopsy site marker for the Celero biopsy device.

Indications

The CeleroMark biopsy site marker is indicated to mark an open or percutaneous biopsy site to radiographically mark the location of the biopsy site.

Contraindications

None known.

** Available in multiple shapes*

Device Description

The CeleroMark biopsy site marker is a sterile, single patient use device comprised of a single, biocompatible titanium marker and a deployment device.

The deployment device is a hand-held device that delivers the marker from the distal tip. The deployment device consists of a rigid cannula, handle, rigid push rod, and plunger. The marker is located at the distal end of the deployment device.

Device Preparation and Use

1. Prior to use of the CeleroMark biopsy site marker, inspect the protective packaging and device to verify that neither has been damaged during shipment. If it appears the packaging has been compromised, do not use the device.
2. Remove the deployment device from its protective packaging using standard interventional technique.
3. Remove the biopsy device from the access site. If using the introducer, hold the introducer steady, release the latching component and pull back on the biopsy device to remove it from the breast.
4. Insert the CeleroMark deployment device. If using the introducer, hold the introducer steady, place the distal end of the deployment device into the proximal end of the introducer and carefully advance the deployment device until it snaps to the introducer hub.
5. Deploy the biopsy site marker by advancing the deployment plunger with your thumb all the way forward until it latches to the hub. Make sure to hold the introducer and/or deployment device steady during deployment of the marker.
6. Verify the deployment and proper position of the marker prior to removal of the device with the appropriate imaging modality.
7. Slowly remove the deployment device, or introducer and deployment device, from the breast as a unit and properly dispose.

Warnings and Precautions

- There are possible adverse reactions when an object is implanted in the body. It is the responsibility of the physician to evaluate any risk or benefit prior to the use of this device.
- Potential complications of marker clip placement consist of pain, seroma formation, inflammation, bruising, hematoma, hemorrhage, infections, hypersensitivity or allergic reaction, soft tissue damage, misdiagnosis (due to marker clip migration), perforation or scar tissue.
- The CeleroMark deployment device is not recommended for use within the bore of an MRI magnet.

- The CeleroMark biopsy site marker is not recommended for use in patients with breast implants.
- The biopsy site marking procedure should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications, and hazards prior to performing any minimally invasive procedure.
- These devices should be used only by physicians trained in open or percutaneous biopsy procedures.
- **ONLY** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- The CeleroMark biopsy site marker should be deployed into the cavity created during the biopsy procedure. Deployment into tissue outside of the biopsy cavity is not recommended.
- Marker position relative to established landmarks may change under mammography upon subsequent breast compressions.
- The CeleroMark biopsy site marker is not intended to be repositioned or recaptured after deployment.
- Users should take care not to unintentionally deploy the marker.
- Excess hematoma within the biopsy cavity and/or introducer can lead to marker adhesion to the deployment device, increasing the risk of marker drag out.
- Care should be taken to avoid damaging the cannula. Avoid operator or instrument contact with the CeleroMark biopsy site marker or the distal end of the deployment device.
- The implanted CeleroMark biopsy site marker is magnetic resonance imaging (MRI) conditional. The implanted marker presents no additional risk to the patient or operator from magnetic forces, torque, heating, induced voltages, or movement, but it may affect MRI image quality.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by Hologic, Inc. may not be compatible with the CeleroMark biopsy site marker device. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.

- Do not resterilize and/or reuse the CeleroMark biopsy site marker device. Resterilization and/or reuse may compromise the integrity of the instrument. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned and sterilized devices.

How Supplied

The CeleroMark biopsy site marker device is sterilized by gamma radiation and supplied preloaded for single patient use. Discard into an appropriate container after use.

As Identified on Labels:

QTY Number of Devices Enclosed.

YYYY-MM-DD Expiration date is represented by the following:

YYYY represents the year

MM represents the month

DD represents the day

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.
250 Campus Drive
Marlborough, MA 01752 USA
Phone: 877-371-4372
BreastHealth.Support@hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

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