



User Companion

Selenia[®]

Digital Mammography System

Selenia[®]

Digital Mammography System

Selenia User Companion

Part Number MAN-02254

Revision 001

March 2011

Technical Support:

USA: +1 877 371 4372

Europe: +32 2 711 4690

Asia: +852 37487700

All Other: +1 781 999 7750

HOLOGIC[®]

Corporate Headquarters

35 Crosby Drive,
Bedford, MA 01730-1401 USA
Tel: +1 781 999 7300
Sales: +1 781 999 7453
Fax: +1 781 280 0668
www.hologic.com

Europe (EU Representative)

Hologic NV
Leuvensesteenweg 250A
1800 Vilvoorde, Belgium
Tel: +32 2 711 4680
Fax: +32 2 725 2087



Refer to the corporate website for more facilities worldwide.

© Copyright Hologic 2011. All rights reserved. Printed in USA. This manual was originally written in English.

Hologic and the Hologic Logo are trademarks or registered trademarks of Hologic, Inc. Other trademarks registered or used by Hologic and its divisions and subsidiaries in the United States and other countries include: Affirm, Dimensions, DSM, FAST Paddle, HTC, MIMS plus, M-IV, MultiCare, SecurView, Selenia, Smart Paddle, StereoLoc, and TechMate. Solaris, Sun, Sun Blade, Sun Ultra, and Ultra are trademarks or registered trademarks of Sun Microsystems, Inc. in the United States and other countries. Microsoft and Windows are trademarks or registered trademarks of Microsoft Corporation in the United States and other countries. Any other product and company names mentioned herein are the trademarks or registered trademarks of their respective owners.

Table of Contents

List of Figures	ix
List of Tables	xiii
Chapter 1 - How to Start Up and Shut Down the System	1
1.0 Quick Reminder for System Startup in “Sleep” Mode	1
2.0 Quick Reminder for System Startup when Completely Off	1
3.0 How Should I Turn Off the System at the End of the Day?	1
4.0 How Do I Turn the System Completely Off?	1
5.0 How Do I Start Up the TechMate™ if It Is Turned Off?	1
Chapter 2—Acquisition Workstation Q and A	3
1.0 How Can I Switch Users?	3
2.0 What Are the AEC Modes I Have Available for Use?	3
3.0 What Exactly is TEC (Tissue Exposure Control)?	3
4.0 What Must I Know to Acquire Images on Patients with Implants?	4
4.1 Method One—Imaging Patients with Implants Away From the Nipple	4
4.2 Method Two—Imaging Patients with Implants Near the Nipple	5
4.3 Method Three—Imaging Patients with Dense Tissue Next To the Implants	6
4.4 Method Four—Imaging Implant Displaced Views	6
4.5 How Should I Perform Implant Procedures?	7
4.6 When I Print Implant Images, What Adjustments Might I Need to Make?	8
5.0 Hints for a Localization Procedure	8
6.0 What is the Recommended Exposure Technique for a Specimen?	9
7.0 How Can I See Which Patients Are on the CD in the Drive?	9
8.0 What Must I Do when I Get a Paddle Message?	9
9.0 How Can I Send Images to an Output after the Procedure Is Completed?	10
10.0 How Do I Edit the Patient Information on Selenia, and PACS?	10
10.1 Correct the Information at the Acquisition Workstation	11
10.2 Resend the Corrected Images to the Outputs	11
11.0 What if I Reject an Image, but Then Change My Mind?	12
11.1 Scenario #1	12
11.2 Scenario #2	13
11.3 Scenario #3	13
12.0 What are the Quality Control Tests We Have to Perform?	14
13.0 Quality Control Hints—Acquisition Workstation	14
13.1 What Is the Best Way to Organize the Weekly Quality Control?	14
13.2 How Can I Be Sure My Weekly Quality Control “Patient” Is Not Deleted?	14
13.3 Can I Run the Repeat/Reject Analysis from an Office Computer?	15
13.4 Can I Save Repeat/Reject Analysis as a Shortcut on a Computer’s Desktop	15
13.5 How Should I Empty the Reject Bin, and How Often Should I Do It?	16
13.6 Hints for the Flat Field Calibration Test	17
13.7 Hints for the Artifact Evaluation Test	17
13.8 Hints for the Phantom View and SNR/CNR Tests.	17

Chapter 3—The User Interface **19**

- 1.0 Details about Login, Logout, and Shutdown 19
 - 1.1 Startup Wait Time 19
 - 1.2 System Messages 19
 - 1.3 How to Log Out or Shut Down 20
 - 1.4 Using the Emergency Off Switches 21
- 2.0 The Select a Patient Screen 23
- 3.0 Patient Demographic Database Organization 25
 - 3.1 Available Databases 25
 - 3.2 Available Patient Lists 25
 - 3.3 The Acquisition Workstation Local Database 25
 - 3.4 Acquisition Workstation Local Database and Modality Worklist Provider Interactions 26
- 4.0 Working with the Databases 27
 - 4.1 Using the Bar Code Scanner 27
 - 4.2 Querying the Modality Worklist Provider 27
 - 4.3 Searching Local Exams—Acquisition Workstation Local Database 30
 - 4.4 Searching Patient Data on a System Without a Worklist 32
- 5.0 Choose a Procedure 33
- 6.0 Patient View Screen 34
 - 6.1 Patient View Screen Options 35
 - 6.2 Selecting Exposure Modes 36
 - 6.3 TEC Exposure Mode 36
 - 6.4 AEC Exposure Mode 38
 - 6.5 Select the Output 41
 - 6.6 Implant Present 41
- 7.0 Acquiring Images 41
- 8.0 Working with Images 41
 - 8.1 Preview Image 42
 - 8.2 Comments 43
 - 8.3 Implant Present 43
 - 8.4 Zoom/Pan 44
 - 8.5 ROI 44
 - 8.6 Measurement 45
 - 8.7 Magnification 45
 - 8.8 Window/Level 45
 - 8.9 Cropping 45
 - 8.10 Accept or Reject Images 46
 - 8.11 Accept a Rejected Image 47
- 9.0 Working with Images After They Are Accepted 49
 - 9.1 Send the Images to the Output Devices 49
 - 9.2 Patient View Screen with Thumbnails 49
 - 9.3 Exposure Taken Using an Incorrect Icon 50
 - 9.4 Close a Procedure 55
- 10.0 Additional Actions on the Patient View Screen 56
 - 10.1 Add a View to the Procedure 56
 - 10.2 Remove Unexposed Views Button 57
 - 10.3 Edit Patients and Procedures 57
 - 10.4 Non-Imaging Mode Button 60
 - 10.5 Simulate Capture Button 61
 - 10.6 Review Button 61

11.0 Clinical Procedures	62
11.1 Manual Mode	62
11.2 TEC Mode	63
11.3 AEC Mode	64
12.0 Performing a Screening Procedure	65
12.1 Select Patient, Procedure, Output, and Exposure Techniques	65
12.2 Set Up the Gantry	65
12.3 Position the Patient	66
12.4 Acquire the Image and Complete the Procedure	66
13.0 Working with Specimens	67
14.0 Performing a Localization Procedure	68
14.1 How to Use the Standard Crosshair Device	68
14.2 How to Use the Magnification Crosshair Device	69
15.0 Accepting and Storing the Image	69
16.0 Optional Workflow Operations	69
17.0 Main Menu Functions	70
17.1 Main Menu with No Selected Patient	70
17.2 File Menu>Exit	70
18.0 Edit Menu	71
18.1 Standard Setup Screens	71
18.2 Set Technique Defaults Screen	71
18.3 View Order Editor	72
18.4 User Setup—Add or Edit a User	72
18.5 Set Up Outputs	74
18.6 Data Storage (Using Outputs)	79
19.0 Admin Menu	82
19.1 Image/Spool Management	82
19.2 Protect Patients	86
19.3 Manage Queues	87
19.4 Eject	90
19.5 Import	90
19.6 Retrieve Priors	92
19.7 MPPS Status	94
19.8 Test Patterns	95
19.9 Available Disk Space	95
20.0 Info Menu	95
Chapter 4—System Hardware Details	97
1.0 System Power Controls	97
1.1 Acquisition Workstation Circuit Breaker	97
1.2 Gantry Circuit Breaker	97
1.3 Uninterruptible Power Supplies (UPS)	98
2.0 Controls and Displays	99
2.1 Compression Device Controls and Displays	99
2.2 Dual Function Footswitches	99
2.3 C-Arm Controls	100
3.0 Operator Radiation Shield	104
4.0 Background for the VOS Switch on the Mobile Selenia	104

Chapter 5—The Manager Functions	105
1.0 User Identification	105
1.1 User Assignments	105
2.0 Main Menu Functions	106
3.0 Edit Menu—Mgr Level	108
3.1 View Order Editor	108
4.0 Admin Menu Functions—“mgr” User Level	109
4.1 Spool Management	110
4.2 Delete Patients	114
4.3 Reclaiming	116
4.4 Reject Management	117
4.5 Log Backup	122
4.6 Paddle Configuration	123
5.0 Info Main Menu Function	124
6.0 ALARM Function	124
7.0 Repeat/Reject Reports Options Available from the Browser	124
7.1 Administrative Reporting	124
7.2 Accessing the Report from Any Computer on the Network	125
7.3 Change Password	125
7.4 Reject Analysis Report	126
Chapter 6—Troubleshooting	127
1.0 Alarm Types and System Messages	127
1.1 Category One Alarm Functions	128
1.2 Category Two Alarm Functions	130
1.3 Category Four Alarms	132
2.0 Tables of Alarm Messages, System Messages, and User Actions	133
2.1 Category One Alarms	133
2.2 Category Three Alarms	135
2.3 Category Two Alarms	136
2.4 System Message Information	143
3.0 Software Troubleshooting	144
4.0 Worklist Problems	146
4.1 Unknown RIS Code	146
4.2 Conflict Between Data in the Modality Worklist and Local Database	146
4.3 Resolving Patient Data Conflicts	147
Index	149

List of Figures

Figure 2-1: AEC Imaging Used for Breasts with the Implant Away From the Nipple	4
Figure 2-2: Implants Near the Nipple Require Manual Exposure Techniques	5
Figure 2-3: Image Implant Displaced Views with a Manual AEC Sensor Position	6
Figure 2-4: Compression Release Set at Manual.....	8
Figure 2-5: The Review Button	12
Figure 2-6: Rejected Images.....	12
Figure 2-7: The Info Screen.....	17
Figure 3-1: Sign Out Button	20
Figure 3-2: Exit from Acquisition Station Dialog Box	20
Figure 3-3: Power Off the Computer Dialog Box	20
Figure 3-4: Confirm Exit from Program	21
Figure 3-5: Panel Power Dialog Box	22
Figure 3-6: Waiting for Link.....	22
Figure 3-7: Initial Select a Patient Screen with No Patient Information	23
Figure 3-8: Select a Patient Screen Options	24
Figure 3-9: Fields Used to Filter a Query to the Modality Worklist.....	27
Figure 3-10: Results of the Patient Query.....	28
Figure 3-11: List of Patients From the Today's Worklist Query	28
Figure 3-12: Selected Patients Appear in the Local Database, Local Exams Screen	30
Figure 3-13: Selecting an Existing Patient from a Search List	30
Figure 3-14: Typing a Patient Name in the Search Field.....	31
Figure 3-15: Select a Patient Screen	31
Figure 3-16: List is Too Long.....	32
Figure 3-17: No Patients Matching Query Message	32
Figure 3-18: Search for a Patient in a System Without a Worklist	33
Figure 3-19: Choose a Procedure.....	33
Figure 3-20: Patient View Screen Panes.....	34
Figure 3-21: Patient Views Options	35
Figure 3-22: Exposure Modes	36
Figure 3-23: Accepting TEC Techniques	37
Figure 3-24: TEC with Techniques Accepted	37
Figure 3-25: Override Mode	37
Figure 3-26: AEC Modes—Choices That Must be Made by the Technologist	38
Figure 3-27: AEC Exposure Adjustment Setting on Screen.....	39
Figure 3-28: No Outputs Warning	41
Figure 3-29: Image Preview Screen.....	42
Figure 3-30: The ROI Statistics Information on the Preview Screen	44
Figure 3-31: Changing the Window/Level Settings.....	45
Figure 3-32: The Accept Button	46
Figure 3-33: Review Dialog Box Showing Two Rejected Images.....	47
Figure 3-34: Accept Rejected Image	48
Figure 3-35: Send Accepted Image to Outputs.....	48
Figure 3-36: Set Rejected Button.....	49
Figure 3-37: Send Accepted Image to Outputs.....	49
Figure 3-38: Patient View Showing Acquired Image	50
Figure 3-39: Patient View Screen's Close Procedure button	55
Figure 3-40: A Finished Standard Screening Procedure.....	55
Figure 3-41: Patient Views Screen with Icon(s) for Scheduled Views.....	56

Figure 3-42: Results of Adding Another Procedure	58
Figure 3-43: Selecting Two Procedures in the Choose a Procedure Screen	58
Figure 3-44: Non-Imaging Mode Button.....	60
Figure 3-45: Non-Imaging Screen	60
Figure 3-46: Non-Imaging Dialog Box Warning	61
Figure 3-47: Rotate the Crosshair Device Out of the X-ray Field.....	68
Figure 3-48: Edit Menu	71
Figure 3-49: View Order Editor	72
Figure 3-50: Edit Users Dialog Box	72
Figure 3-51: Add a New User	73
Figure 3-52: Edit User Profile	73
Figure 3-53: Edit Outputs Dialog Box.....	75
Figure 3-54: New Output Name Dialog Box	75
Figure 3-55: New Name in Edit Outputs Dialog Box	75
Figure 3-56: Adding Devices.....	76
Figure 3-57: Configuring the Output Device(s).....	76
Figure 3-58: Choose Output.....	78
Figure 3-59: No Output Selected	78
Figure 3-60: Invalid Media Alert	80
Figure 3-61: CD-R Alert	81
Figure 3-62: Admin Menu.....	82
Figure 3-63: Image Management Resend Options	83
Figure 3-64: Resend Confirmation.....	84
Figure 3-65: Repreview Screen	84
Figure 3-66: Protecting Patients.....	86
Figure 3-67: Protect Confirmation	86
Figure 3-68: Manage Queues Suspend Job Spooler button	88
Figure 3-69: Delete Selected Jobs in a Queue	89
Figure 3-70: The Import Dialog Box.....	90
Figure 3-71: Import Progress Dialog Box.....	91
Figure 3-72: Retrieve Priors Screen	92
Figure 3-73: Results of the Retrieve Priors Query.....	93
Figure 3-74: Image Import Started	94
Figure 3-75: Query Retrieve Import Successful Dialog Box	94
Figure 3-76: Search too Broad Dialog Box	94
Figure 3-77: Query Retrieve Errors dialog Box.....	94
Figure 3-78: Test Patterns Dialog Box	95
Figure 3-79: Image Queued Dialog Box	95
Figure 4-1: Acquisition Workstation Circuit Breaker.....	97
Figure 4-2: Gantry Circuit Breaker	97
Figure 4-3: Compression Device	99
Figure 4-4: The Dual Function Footswitch	99
Figure 4-5: AWS Radiation Shield	104
Figure 4-6: Value Console Radiation Shield	104
Figure 5-1: Edit Menu	108
Figure 5-2: View Order Editor	108
Figure 5-3: Admin Menu for the Mgr Logon	109
Figure 5-4: Choose a New Patient Dialog Box	110
Figure 5-5: Spool Management Window—One Image Selected	110
Figure 5-6: Repreview Window	113
Figure 5-7: Spool Management-Delete Image Function.....	114

Figure 5-8: Delete Patient Window.....	115
Figure 5-9: Patient Delete Error Dialog Box	115
Figure 5-10: Delete Selected Patients Dialog Box	115
Figure 5-11: Hard Drive Reclaimed	116
Figure 5-12: Hard Drive Full Alarm	116
Figure 5-13: Reject Management Window—(Resend is Shown).....	117
Figure 5-14: Repreview Tab on Reject Management Screen	119
Figure 5-15: Reject Spool Full Alarm	120
Figure 5-16: Reject Management—Resend a Rejected Image.....	121
Figure 5-17: Log Backup—First Dialog Box	122
Figure 5-18: Log Backup—Second Dialog Box	122
Figure 5-19: Log Backup—Third Dialog Box	122
Figure 5-20: Log Backup—Fourth Dialog Box.....	122
Figure 5-21: Paddle Configuration Dialog Box.....	123
Figure 5-22: Info Window	124
Figure 5-23: Log on to the Browser	125
Figure 5-24: Administrative Reporting Options	125
Figure 5-25: Change Password	126
Figure 5-26: Criteria for Reject Analysis.....	126
Figure 6-1: Location of Alarms by Type	127
Figure 6-2: Category One Alarm Dialog Box	128
Figure 6-3: AEC Error Message—Below Minimum	129
Figure 6-4: AEC Error Message—Above Maximum	129
Figure 6-5: Alarm Icon.....	130
Figure 6-6: Alarm Dialog Box	131
Figure 6-7: System Messages	132
Figure 6-8: Acquisition Workstation Circuit Breaker	144
Figure 6-9: Unknown RIS Code	146
Figure 6-10: Conflict Between the Worklist and Local Database.....	146
Figure 6-11: Sample Dialog Box: Add a New Procedure.....	147
Figure 6-12: Worklist Patient Conflict Dialog Box.....	147

List of Tables

Table 2-1: AEC Modes.....	3
Table 2-2: Exposure Techniques for Implants Near the Nipple.....	5
Table 2-3: Recommended Exposure Techniques for an Implant Near the Nipple.....	7
Table 3-1: AEC Mode—Available Selections	38
Table 3-2: Incorrect View Selection Changes.....	53
Table 3-3: Main Menu Functions with No Procedure Chosen.....	70
Table 3-4: Setup Choices.....	71
Table 3-5: Set Technique Default Choices	71
Table 5-1: Main Menu Functions by User (Patient/Procedure Not Selected).....	106
Table 5-2: Main Menu Functions by User (Patient/Procedure Selected)	107
Table 5-3: The Edit Menu for mgr Logon	108
Table 5-4: The Admin Menu for Mgr Logon.....	109
Table 6-1: Alarm Messages	133
Table 6-2: Image Processing Error Messages	135
Table 6-3: Category Two Alarms	136
Table 6-4: System Messages	143
Table 6-5: Software Troubleshooting	144

Chapter 1 - How to Start Up and Shut Down the System

1.0 Quick Reminder for System Startup in "Sleep" Mode

1. Move the trackball to wake the system up.
2. In the Login dialog box enter `reboot` for the user name and password. Wait for the system to start up. The Logon dialog box appears.
3. Enter your Operating System user name *tech* or *mgr* (case sensitive).
4. Click the **OK** button or press the Enter key.
5. Enter the Operating System password (case sensitive) in the next dialog box.
6. Click the **OK** button or press the Enter key. Click the **Dismiss** button only if the system was in "Sleep" mode.
7. In the Login to System dialog box, click the down arrow and select your name in the list.
8. Click the Password field and enter your password (case sensitive).
9. Click the **OK** button or press the Enter key.

2.0 Quick Reminder for System Startup when Completely Off

1. Press the **On** button to start the system. The green LED Power On indicator illuminates.
2. Wait for the computer to start up. Do not click the Dismiss button in the Launch dialog box. Go to Step 3 above.

3.0 How Should I Turn Off the System at the End of the Day?

1. Click **File** (upper left of main screen), then select **Exit**.
2. Click the Log out of the computer option, and then click the Yes button.
3. Click the Yes button in the second "Exit from acq. station" dialog box. If a message appears that spool jobs are waiting, wait until the jobs are finished. The Logout Confirmation dialog box appears. It warns your current session will not be saved. Do not be concerned about this message.
4. Click the OK button. The OS Login dialog box appears. The unit enters a "Sleep" mode.

4.0 How Do I Turn the System Completely Off?

1. Click **File>Exit**.
2. Click the Power off the computer option.
3. Click Yes. If a message appears that spool jobs are waiting, wait until they are finished.
4. Click the OK button at the Logout Confirmation dialog box, "Are you sure you would like to power off the computer?" The green LED power indicator turns off.

5.0 How Do I Start Up the TechMate™ if It Is Turned Off?

1. To start up the TechMate, press the power button on the right side of the front cover.
2. Enter `scr` in the User name field, then click the Return button.
3. A "Technologist Workstation is loading" message appears.
4. The User login screen appears after the software loading is complete.

User Companion

Chapter 1 - How to Start Up and Shut Down the System
How Do I Start Up the TechMate™ if It Is Turned Off?

Selenia®
Digital Mammography System

Chapter 2—Acquisition Workstation Q and A

1.0 How Can I Switch Users?

Click the **Sign Out** button. The Login to System dialog box appears.

The new user then selects their name from the Identify Yourself drop down list and enters their password. They are signed in.

2.0 What Are the AEC Modes I Have Available for Use?

There are three AEC modes: Auto-Time, Auto-kV, and Auto-Filter. These modes allow various levels of user control of the Automatic Exposure Control techniques:

Table 2-1: AEC Modes

AEC Mode	Filter Parameter	kV Parameter	mAs Parameter
Auto-Time	User selected	User selected	Determined from a scout exposure
Auto-kV	Filter is not selectable. Moly tube uses Mo filter. Tungsten tube uses Rh filter.	Determined from compression thickness and look-up table	Determined from a scout exposure
Auto-Filter	Determined from compression thickness and table	Determined from compression thickness and look-up table	Determined from a scout exposure

3.0 What Exactly is TEC (Tissue Exposure Control)?

Tissue exposure control (TEC) exposure mode is designed to take advantage of the imaging characteristics of the Selenia® FFDM system, which has significantly higher detective quantum efficiency and wider exposure latitude compared to a conventional screen-film mammography system.

When the TEC exposure mode is enabled, the Selenia system uses two parameters to compute the optimal x-ray exposure technique: (1) the compressed breast thickness, and (2) the estimated breast tissue density. The compressed breast thickness is automatically read by the Selenia compression paddle system. The mammography technologist estimates the breast tissue density, and sets the breast density control to one of the following three choices: FATTY, NORMAL, and DENSE. The computed x-ray exposure techniques (kV, filter, and mAs) are displayed on the Selenia acquisition display. Once the mammography technologist accepts the displayed exposure techniques, the Selenia system will change from STANDBY to READY, and an image can be acquired instantly.

TEC mode is an advanced MANUAL exposure mode. The mammography technologist always has the option to override the computed exposure techniques, simply by manually changing one of the exposure parameters (kV, mAs, or filter). After the computed exposure techniques are accepted or changed by the mammography technologist, the Selenia system no longer automatically re-calculates the exposure parameters. The TEC mode can be enabled again by changing the breast tissue density control, or by releasing and re-applying breast compression.

4.0 What Must I Know to Acquire Images on Patients with Implants?

The wide dynamic range of the digital detector on Selenia systems allows more flexibility in selecting techniques for imaging implants. However, image processing does not always produce optimal results. When all factors work together, digital images of implants contain more diagnostic information and are far superior to screen-film images.

Image processing on the system uses special algorithms to detect an implant on a digital mammogram. A key element contributing to the success of this algorithm is the use of correct exposure techniques. Although the digital detector may forgive mistakes in imaging techniques, the right exposure range is still necessary to acquire optimum images. The following sections provide guidance and recommendations to optimize the task of imaging patients with implants.

Compressing breast tissue that contains an implant is challenging to the technologist and the patient. Compression is important and adequate compression must be used when imaging patients with implants. Large variations of compressed thickness between the implant and the rest of the breast tissue may prevent image processing from correctly identifying the implant region.

**Note...**

Always manually position the AEC Sensor when using AEC to image patients with implants. (If you use Auto AEC position, the system will place an AEC sensor under the densest portion of the breast; in a case with implants this will always be under the implant.)

4.1 Method One—Imaging Patients with Implants Away From the Nipple

The following figure demonstrates a breast with ample tissue in front of the implant. The Selenia AEC system may be used in this case to select the exposure techniques more accurately because there is enough breast tissue for AEC imaging. You can use the Auto-Filter mode and move the AEC sensor manually to the middle of the breast tissue area.




Figure 2-1: AEC Imaging Used for Breasts with the Implant Away From the Nipple

4.2 Method Two—Imaging Patients with Implants Near the Nipple

The next figure shows minimal breast tissue around the implant. When there is a small amount of breast tissue in front of the implant, the Selenia AEC can not be used accurately because there is not enough breast tissue to use AEC. Part of the implant can interfere with the AEC sensor. In this case use the Manual exposure techniques instead of AEC.




Figure 2-2: Implants Near the Nipple Require Manual Exposure Techniques

The two tables shown below provide the exposure techniques which apply to patients with implants that fall under this category. Do not use the exposure techniques in these tables to acquire images of patients without implants or to acquire implant displaced views. These exposure techniques were created only for imaging implants with minimal tissue around the implant.

Table 2-2: Exposure Techniques for Implants Near the Nipple

<i>Selenia with Tungsten X-ray tube</i>			<i>Selenia with Molybdenum X-ray tube</i>		
Thickness (cm)	kV	mAs	Thickness (cm)	kV	mAs
Less than 4 cm	28	100	Less than 4 cm	28	60
4 cm to < 6 cm	28	120	4 cm to < 6 cm	28	100
6 cm to < 8 cm	28	150	6 cm to < 8 cm	28	160
8 cm to 10 cm	28	150	8 cm to 10 cm	28	200
More than 10 cm	28	180	More than 10 cm	28	240

4.3 Method Three—Imaging Patients with Dense Tissue Next To the Implants

When a very dense area is next to the implant, the image processing can detect this dense area as an extension of the implant.

When a dense area next to the implant appears unexposed and part of the implant, the following actions can correct the condition and improve the image quality of the dense area:

- When you acquire the image of the implant with METHOD 1, switch to the Auto-Time AEC mode and increase the kV exposure technique by 2 kV. Adjust the manual position of the AEC sensor, if necessary.
- When you acquire the image of the implant with METHOD 2, keep the same mAs and increase the kV exposure technique by 2 kV.

4.4 Method Four—Imaging Implant Displaced Views

The figure below shows the implant completely or partly displaced from the imaging area. The Implant Displaced views are equivalent to normal views even though a small part of the implant remains visible in the image. Use the Selenia AEC (for example, Auto-Filter) to acquire these views. You should use a manual AEC sensor position to acquire implant displaced views to avoid interference from part of the implant still being present.




Figure 2-3: Image Implant Displaced Views with a Manual AEC Sensor Position

4.5 How Should I Perform Implant Procedures?



Note...

Never use AUTO AEC sensor mode when imaging implants. AUTO AEC mode puts an AEC sensor under the densest part of the breast. In a case with implants, this part is always under the implant. When imaging patients with implants with AEC, always use a manual position for the AEC Sensor location.

1. Select the **Screening Mammogram w/Implants** procedure from the drop-down list. Check the **Implant Present** box above the views. Leave this box checked for all views.
2. Position and compress the breast correctly for an augmented breast.
3. Acquire the full implant view:
 - a. *Method One*—If the implant is away from the nipple, use Auto-Filter and manually move the AEC Sensor to the middle of the breast tissue.
 - b. *Method Two*—If the implant is near the nipple, use a Manual exposure technique for the full implant view. Refer to the charts below for recommended exposure techniques to acquire an image of an implant near the nipple. The **Exposure Index** will be very low, but you do not need to be concerned. The exposure index is measured under the implant.

Table 2-3: Recommended Exposure Techniques for an Implant Near the Nipple

Selenia with Tungsten X-ray Tube			Selenia with Molybdenum X-ray Tube		
Thickness (cm)	kV	mAs	Thickness (cm)	kV	mAs
Less than 4 cm	28	100	Less than 4 cm	28	60
4 cm to < 6 cm	28	120	4 cm to < 6 cm	28	100
6 cm to < 8 cm	28	150	6 cm to < 8 cm	28	160
8 cm to 10 cm	28	150	8 cm to 10 cm	28	200
More than 10 cm	28	180	More than 10 cm	28	240

- c. *Method Three*—If the full implant preview includes a dense area next to the implant, re-acquire the image:
 - If you used Method One, switch to Auto-Time and increase kV by 2 kV. Adjust the manual sensor if necessary.
 - If you used Method Two, keep the same mAs and increase the kV by 2 kV.
4. Perform implant displaced views with AEC.
 - a. If the implant is in the field of view, place the photocell beyond the implant to avoid photo timing through the implant.
 - b. If the implant can be almost completely displaced from the imaging area, use Auto-Filter.
 - c. If a small part of the implant appears in the imaging area, use a manual sensor position.



Note...

*If you forget to check the **Implant Present** box on the Patient View screen, you may check the box on the Preview screen then click the **Reprocess Image** button at the bottom of the Preview screen.*

4.6 When I Print Implant Images, What Adjustments Might I Need to Make?

Printed images of the full implants can appear flat or dark if sent to printer as a normal mammogram. The following adjustments may be required.

1. After the case is finished, select **Admin>Image/Spool Management**.
2. Find the case. Click all of the full implant views.
3. Click the **Resend selected image(s)** button.
4. Select your printer as the Output.
5. You can adjust the **Contrast/Density** for your printer. Change the contrast or density with the radio buttons. This change is only temporary for the images you are printing now.
6. Click the **Resend** button.
7. A message appears which indicates how many images were resent.
8. Click the **OK** button.
9. Click the **Close** button when the printing is completed.

5.0 Hints for a Localization Procedure

The Localization procedure has additional steps not needed in other clinical procedures. Refer to Chapter 3, Section 14.0, page 68 for the procedure.

- If a paddle other than the Localization paddle will be used, you must turn off the Automatic Compression Release. Refer to Chapter 3, Section 18.1, page 71.




Figure 2-4: Compression Release Set at Manual

- After closing the Standard Setup window, confirm that Manual Release has been activated by checking the settings in the lower left hand part of screen. See Figure 2-4.
- Upon completion of the Procedure, change the Compression Release back to automatic if this is the default. Confirm Automatic Release is activated by checking the settings.
- Remove the face shield.
- To add views, click the **Add View** button before you begin the exam. (Be aware, any additional view is automatically highlighted, so make sure you scroll back to the beginning of your procedure and highlight the first view).
- Take the scout image using Auto-Filter.
- Do not use Auto AEC position. Use manual position 2 or 3.
- If doing a magnification needle localization attach the magnification stand and magnification crosshairs.
- Remember to move the crosshair device or dial the crosshairs out of the field before you acquire the scout image.

6.0 What is the Recommended Exposure Technique for a Specimen?

Refer to the user manuals for the procedure.

- Position the specimen as close to the chest wall edge as possible.
- Always use a spot paddle.
- Use the magnification stand if desired.
- Take an exposure using the **Manual technique, 25 kVp, and 26 mAs.**

7.0 How Can I See Which Patients Are on the CD in the Drive?

Select **Admin>Import**, and the patient name is listed in the dialog box.



Note...

The CD-RW is very sensitive. Do Not press the Eject button while the system reads from or writes to the disk.

*Check Jobs in Queue to confirm no jobs are waiting to be sent before selecting **Eject** from the Admin menu.*



Note...

Do not use stick-on labels on a CD.

8.0 What Must I Do when I Get a Paddle Message?



This message appears if you have an incomplete paddle shift.



This message appears if your Smart Paddle™ is in the wrong position for the selected View.

The system will not acquire an image if the paddle is not installed in one of the three positions. The paddle position must match the selected View. Move the paddle to the correct position and the message clears.

9.0 How Can I Send Images to an Output after the Procedure Is Completed?

Refer to the *Instructions for Use* for directions and Chapter 3, Section 9.1, page 49 for additional details.

1. Click **Admin>Image (or Spool) Mgmt** and select the patient. All the thumbnails in that procedure will be displayed.
 - Click an individual thumbnail to see the available information regarding this particular view (displayed in the upper right portion of the window).
1. Click the **Resend** tab and either click the **Resend All Images** button, or select images and click the **Resend Selected Images** button.
2. Choose the Output from the drop down list.
3. Click the **Resend** button.



10.0 How Do I Edit the Patient Information on Selenia, and PACS?

- The best place to change the patient information is on the HIS/RIS before any information is sent to the Selenia. If the images are acquired and sent to the PACS, change the information at the PACS. Change Patient Information at the Selenia application as a last resort.



Note...

Some RIS systems eliminate the ability to change these errors at the Selenia.

- If the error involved exposing two patients under the same name, it may be necessary to read the first patient at the diagnostic review workstation, then delete the patient study before sending the second set of edited images, because of the information in the DICOM header.
- Because you can cause a mismatch of information when the Selenia images are sent to PACS, work with your HIS/RIS administrator before you make any changes to the patient information at the Selenia.
- You can only change the demographics for a procedure on the Acquisition Workstation. Demographics cannot be changed once the procedure is reclaimed.
- Be sure the incorrect images are still on the Acquisition Workstation before you delete them from the PACS.
- If incorrect images are already accepted and sent to Outputs, prior to sending corrected images you must:

At the SV: Sign on to the SVDX as **delete** and delete that patient study.

At the PACS: Contact your PACS Administrator and alert them to delete the incorrect study or advise another corrective measure.

10.1 Correct the Information at the Acquisition Workstation

These are some additional hints to change the information on the Patient View screen.

1. In the Patient View screen under Patient (right side of screen) click the **Edit** button. See figure at right.
2. If the Accession Number is also incorrect, Click the **Edit** button under Procedure (right side middle of figure) and make the correction in that dialog box.
3. Close the procedure.



10.2 Resend the Corrected Images to the Outputs

1. Click **Admin>Image (or Spool) Mgmt.**
2. Display the patient and click the **Repreview** tab. See figure at right.
3. Select and save each of the original images. (Follow the procedure in the *Instructions for Use.*) There should now be eight thumbnails (assuming there were four in the original study).
4. Click all of the last four thumbnails (these are the corrected views) and resend them.

After you confirm that the new images have been sent, you may (at *mgr* level) go back into Spool Mgmt and delete the four images labeled incorrectly.



11.0 What if I Reject an Image, but Then Change My Mind?

11.1 Scenario #1

- The first image is taken and rejected.
 - The second image is taken and rejected (worse than the first).
 - You'd like to un-reject one of the images.
1. Click the **Review** button (see Figure 2-5). A dialog box opens in the upper left of screen with tabs for all the views taken.



Figure 2-5: The Review Button

2. Click the Tab of interest [ex: RCC]. There are all the rejected images for that view.




Figure 2-6: Rejected Images

3. Click the image that you want to Accept.
 4. When the preview image appears, click the **Set Accepted** button.
 5. To send to the selected Output, see instructions in Scenario #3.
 6. Click the **OK** button to close the Review dialog box.
- Only one Accepted image can be housed under one view icon at any one time.

11.2 Scenario #2

- The first image is taken and rejected.
 - The second image is taken and Accepted (sent to outputs when accepted).
 - You'd like to un-reject the first image and add it to the study.
1. Click the **Review** button (right side of screen, underneath Jobs in Queue, Figure 2-5, page 12).
 2. A dialog box opens in the upper left of screen with tabs for all the views taken.
 3. Click the Tab of interest [ex: RCC]. There are 2 images, 1 Rejected, 1 Accepted.
 4. Click the **Accepted** image.
 5. Wait for the preview image to appear, click the **Set Rejected** button.
 6. When the Image Rejection Information dialog box opens, check the appropriate reason.
 7. Click the **Reject** button.
 8. Now there are 2 rejected images. Click the first Rejected image.
 9. When the preview image appears, click the **Set Accepted** button.
 10. To send to the selected Output, see instructions in Scenario #3.
 11. Click the **OK** button to close the Review window.

11.3 Scenario #3



Note...

Any image that was rejected prior to the Close Procedure are not sent to CAD. These images may be sent to CAD using the Scenario #3 instructions.

To use this method you must be signed in at the *mgr* level.

1. Click **Admin**.
2. Click **Reject Management**.
3. Click **Choose a New Patient**.
4. Select **Patient Name** tab, in the Patient Name field type in the patient's name. Click **Search**.
5. Select the patient.
6. Click the rejected images you want to send to an output.
7. Select **Resend**.
8. Click either **Resend All Images** for patient or **Resend Selected Images**.
9. Choose the Output from the drop down box.
10. Click **Resend**.
11. Click **Close**.

12.0 What are the Quality Control Tests We Have to Perform?

Refer to the current Selenia Quality Control Manual for the Quality Control Tests you must perform.

These tests, from the ACR Mammography Quality Control Manual, must still be performed:

- Visual Checklist
- Viewboxes and Viewing Conditions
- Compression Force

13.0 Quality Control Hints—Acquisition Workstation

13.1 What Is the Best Way to Organize the Weekly Quality Control?

Every month you should create a new “patient” for Weekly Quality Control.

1. Under Patient, click the **New** button.
2. For the patient info:
 - For the last name enter: WEEKLY Quality Control
 - For the first name enter: CURRENT MONTH (Jan, Feb, etc....)
 - For the ID, DOB, and Accession Number: Enter the current date (numbers only).



Note...

If you have more than one Acquisition Workstation, you should add the room number to the patient ID field.

3. At the Procedure Description, choose **Tech Quality Control procedure** from the drop down menu.
4. Confirm **Weekly Quality Control** is the default.
5. Click the **Accept** button.

13.2 How Can I Be Sure My Weekly Quality Control “Patient” Is Not Deleted?

Make sure you Protect your Weekly Quality Control patient.

1. Select **Admin>Protect Patient**; search and select the current Weekly Quality Control patient.
2. Click the **Protect** button, and click the **Yes** button in the dialog box.
3. Click the **Close** button to close the Protect Patient dialog box.

13.3 Can I Run the Repeat/Reject Analysis from an Office Computer?

Setting up an office computer for Repeat/Reject Analysis

1. Get the IP Address of the Selenia system from Hologic™ personnel or your IT Group.
2. At the Office Computer, log on to the Internet Explorer or your internet provider.
3. Enter the web address:
Example: `http://121.7.32.4.10:8080`
`http:// [the IP address]:8080`
4. Click the **Return** button.
5. The login screen displays.
6. Login with user name: `admi n` (case sensitive, lower case).
7. Click the **Logon to the System** option. Click **Repeat/Reject**.
8. Select the To/From dates and select the technologist. Select the Type (Report).
9. Click **Submit**. The report is compiled.
10. Once the report finishes, click the **Print** button to print the report using the browser's print button.
11. Exit from the Browser when you are done to maintain security.

13.4 Can I Save Repeat/Reject Analysis as a Shortcut on a Computer's Desktop

1. Right click the login page, then, click **Favorites**.
2. The computer prompts you for a link name.
3. Enter **Selenia Repeat-Reject** (and the Room # if there are multiple rooms).
4. Once the favorite is saved, click **My Favorites**.
5. Highlight the Selenia Repeat/Reject web link you just saved.
6. Right click to display the dialog box.
7. Click **Send To** and another dialog box displays.
8. Click **Desktop (create a shortcut)**.
9. A shortcut to the Selenia Repeat/Reject Analysis Report is now on the computer's Desktop.

13.5 How Should I Empty the Reject Bin, and How Often Should I Do It?

Running the Reject report monthly gives the site better control, but it must be performed at least quarterly.

1. You must be signed on as mgr to perform this task.
2. Click **Admin>Reject Management**.
3. Click **Rejected Image Report**.
4. Click **Search**. If message displays, “There are too many results. Please refine your search,” you must narrow the search:
 - a. Enter the Operator’s (Technologist’s) Name. This is their user name for the system, which may be their initials.
 - b. Click the **Search** button.
5. Click the first item in your Search Results.
6. Click the **Delete** Tab.
7. Click the **Delete All for Patient** button.
8. When Confirm Operation box appears, click the **Yes** button.
9. Click the **Rejected Image Report** button.
10. Click the **Search** button again.
11. Click the first item in your Search Results.
12. Follow step 7 to step 11 until you have deleted all rejected images.
13. If you entered the operator’s name, you must enter the next operator’s name and delete all of their images. Follow these steps until you have deleted all rejected images for all technologists.
14. Click the **Close** button (bottom right) to close the Choose an Operator dialog box.
15. Click the **Close** button (bottom right) to close the Reject Management dialog box.

13.6 Hints for the Flat Field Calibration Test



Note...

Clean the detector and acrylic block before Calibration.



Figure 2-7: The Info Screen

1. Follow Instructions printed on the screen.
2. To check the Detector Temperature Click **Info>About the Acquisition Station**.
3. On the left side of the Info screen, midway down, look for Temperature. This is the current detector temperature.
4. If the current detector temperature is not within normal temperature range, between 20 and 37 degrees Celsius, delay the calibration test until it is within normal temperature.



Note...

If you do this test in the middle of the day, wait at least 30 minutes after the last patient to start the calibration procedure.

5. Click **OK** to close the screen.

13.7 Hints for the Artifact Evaluation Test

- Select the weekly QC patient and for Outputs select QC Artifact Evaluation.
- Put a CD in the CDRW drive if you store your QC on a CD.
- To use the Zoom/Pan, either hold down the right button and move the trackball/mouse or hold down the left button when the cursor is in the thumbnail and move the trackball/mouse.

13.8 Hints for the Phantom View and SNR/CNR Tests.

- For Outputs, select QC Phantom.
- If there is an SNR button, these are calculated for you automatically.
- If there is no SNR button, use the logs and worksheets in the QC manual to record your values and calculate your results, then record them on the QC forms from the Quality Control manual.

User Companion

Chapter 2—Acquisition Workstation Q and A
Quality Control Hints—Acquisition Workstation

Chapter 3—The User Interface

This chapter contains additional information about the User Interface available to all users. The first part is organized in the order you might use the system. The next part has suggestions and hints for clinical procedures. The last part has details about using the menu functions and additional information about functions available on the screens.

1.0 Details about Login, Logout, and Shutdown

When the system is started, the computer performs internal startup tests. If the system detects a fault condition during the internal checks, a message displays and startup is suspended until the problem is remedied. When the startup tests are complete, the Sun Operating System login dialog box is displayed.

Users must log in to the Sun Operating System and to the Selenia Application software separately.

1.1 Startup Wait Time

After login into the Operating system, the application starts and a **Launch** dialog box alerts you that the digital detector is warming up and to wait before starting any clinical procedures. The alert does not restrict use of the system before the warm-up time expires.

The temperature of the Selenia Detector is critical for image quality. To stabilize its temperature, the detector should be on for at least 1 hour (warm up) before patient exposure.

- If the system has been off longer than 30 minutes, wait until the message no longer displays before acquiring images for clinical purposes.
- If the system has not been off longer than 30 minutes, click **Dismiss**.

1.2 System Messages

The System Message area displays system status and communication between the detector and the Acquisition Workstation (AWS). The displays in the message area are very important when restarting the system after using one of the Emergency Stop Switches (discussed later in this chapter.)



Note...

*When the application starts, the system evaluates the last calibration date. If the date exceeds the calibration period, an error message appears and use of the unit is restricted until it is calibrated. No patient procedures can be performed. Refer to the **QC Manual** for Calibration procedures.*

1.3 How to Log Out or Shut Down

You can exit the software application to log out, change users, or shut the system down. The dialog box for changing users is available as a button on the main application screen (see Figure 3-1) or from the File Menu (**File > Exit**).



Figure 3-1: Sign Out Button

This **File>Exit** dialog box (see Figure 3-2) presents options for shutting down or restarting the system, logging out of the Application, or changing users.

Select the option you want (Figure 3-2) and click **Yes**. Click **No** to close the dialog box and remain signed on.

The options are discussed in the sections that follow.

After making your selection, follow the on-screen prompts.




Figure 3-2: Exit from Acquisition Station Dialog Box

1.3.1 Power off the Computer

This option shuts down the system.




Figure 3-3: Power Off the Computer Dialog Box

An additional dialog box asks if you are sure you want to power off the computer. Click **Yes** or **No**. To restart the shutdown system, press the Acquisition Workstation On button.



Caution:

Do Not use a Circuit Breaker or Emergency Off switch to routinely turn off the Selenia. Do Not use the task launcher or bottom Solaris™ bar, except in emergencies.



Caution:

Never switch off the Acquisition Workstation Circuit Breaker except in emergency. This turns off the Uninterruptible Power Supply (UPS) and risks data loss.

1.3.2 Restart the Computer

This option reboots the system. After the reboot you must log on to the OS and Acquisition Workstation software. Use this selection for troubleshooting purposes.

1.3.3 Log Out of the Computer

This option exits the Acquisition Workstation application and the OS.



Figure 3-4: Confirm Exit from Program

A dialog box questions if you want to continue to log out of the Acquisition Workstation. After you log out, you can shut down the computer, or a new user can log in at a different user level. The next user must log on to the OS first, then the Acquisition Workstation software.

Use this at night to allow auto shutdown/startup. In the morning, log in by typing `reboot` for the user name and password. When the system restarts, log in normally.



Note...

Reboot the Acquisition Workstation's computer daily to resync system communication.



Note...

If a message appears that spool jobs are waiting, choose if you want to go ahead and log off or wait until the jobs are finished. (The spool picks up again later.) Hologic strongly suggests allowing the jobs to complete before shutting down.

1.3.4 Sign In as a Different User

This option signs out the current technologist without changing the OS logon. A new technologist signs in with the same OS access level. This works the same way as clicking the Sign Out button.

1.4 Using the Emergency Off Switches

There are three Emergency Off switches on the Selenia, one on either side of the Gantry and one on the Acquisition Workstation. Do not use the Emergency Off switches for routinely turning off Selenia. They only affect the Gantry and should be used only in an emergency.

1.4.1 Use the Emergency Off Switch

1. Press any Emergency Off switch. The power to the Gantry is immediately removed.
2. If the emergency occurred during a Patient Exam and the patient is compressed, use the handwheels to manually release the patient.

1.4.2 Restart the System After Activating the Switch

1. Turn the Emergency Off switch one-quarter turn clockwise.

- When the Power Panel dialog box displays, “**Communication with the generator is lost,**” examine the Emergency Off Switches and the circuit breaker.



Caution:

Do Not attempt to close this dialog box using the ‘x’ in the upper corner.

- Click **Cancel** if the Gantry is already on.
- If the Gantry is off, and the Emergency Off Switches and the circuit breaker are set correctly, click **OK** to restart the Gantry.
- Click **OK** again when the “**Link established with generator**” dialog box appears.




Figure 3-5: Panel Power Dialog Box

If the dialog box reappears:

- Observe the Gantry status before you respond to the Panel Power dialog box.
- If the Gantry displays are *On*, click **Cancel**.
- If they are *Off*, click **OK**.
- After you click **OK**, a new dialog box, “**Waiting for link to be established with generator,**” is displayed.



Caution:

Do Not attempt to close this dialog box using the ‘x’ in the upper corner.




Figure 3-6: Waiting for Link

- Once a link is established, the message changes to “**Link established with generator**” and an OK button is enabled.
- Select **OK** as soon as this is enabled.



Note...

It may be necessary to click Reset Generator after clicking OK if the Generator Problem message isn't cleared.

2.0 The Select a Patient Screen

This screen appears when the application first loads.



Figure 3-7: Initial Select a Patient Screen with No Patient Information

This chart diagrams the choices available on the Select a Patient screen.




Figure 3-8: Select a Patient Screen Options

The screen enables you to select a patient from existing databases or create a new patient. This section covers additional details about working with the system databases.

3.0 Patient Demographic Database Organization

3.1 Available Databases

The Selenia Acquisition Workstation has a local database on the workstation's computer that stores the patients and image data. The Selenia system also can be configured to have a network connection to a Modality Worklist Provider. This provider can be queried for scheduled patients and their procedures. This process creates a Modality Worklist. Selecting patients from the worklist adds them to the local database.

The Selenia can be configured to communicate with other systems over a network or as an independent unit. When connected to a network, Selenia can query a Modality Worklist Provider, such as a RIS (Radiology Information System) or PACS, for a list of scheduled patients, and the Local Exams (the local database stored on the Acquisition Workstation's computer).



Note...

The Modality Worklist is a list of scheduled exams, typically maintained by a RIS (Radiology Information System) or PACS. For the Modality Worklist to be operational, this feature must be installed and configured for your specific Modality Worklist Provider.

When Selenia is not connected to a Modality Worklist Provider, the local database can be queried for local exams. If your system is not connected to any PACS or RIS system, there are not any Worklist or Local Exams tabs.

3.2 Available Patient Lists

Selenia can display the following types of patient lists:

1. A list of patients resulting from an "On Demand" query to a Modality Worklist Provider (RIS or PACS).
2. A list of patients in the local database. This list includes:
 - Local Exams.
 - Patients entered locally using the "New" Patient button. These patients include completed and uncompleted exams.
 - Uncompleted patient exams that were originally part of a query to a Modality Worklist Provider.

3.3 The Acquisition Workstation Local Database

The local database provides an automatic means to identify patients and temporarily store information about acquired studies so that the output from the Acquisition Workstation is as accurate and complete as possible.

Patient identification involves three functions of the Acquisition Workstation:

- Acquiring Patient demographic information
- Storing Patient demographic information in a local database
- Identifying the Patient and associating the correct demographic information with that Patient upon examination

3.4 Acquisition Workstation Local Database and Modality Worklist Provider Interactions

When the system is configured with a Modality Worklist Provider (MWP), which is typically part of a RIS or a PACS, Patient Demographic information may be obtained from it. There are three basic interactions between the Acquisition Workstation and the Modality Worklist (MWL SCP) functionality:

- **Local Worklist:** The purpose of the Local Worklist is to cache the worklist items obtained from the MWL by the query function. It represents a snapshot of the MWL SCP's worklist obtained by the MWL query. The Local Worklist is distinct from the Local Database, and may contain information that is not included in the Acquisition Workstation information for a given Patient.
- **Local Database:** The Acquisition Workstation maintains its local database with the sole purpose of completely and correctly identifying acquired images. The Local Database is updated from a worklist item when a worklist record is selected.
- **Synchronization:** The Acquisition Workstation synchronizes the Local Worklist and the Local Database when you select a worklist item(s). At Synchronization time, if the information from the selected worklist is found to be inconsistent with the local database, you are given the option of resolving the inconsistencies via a dialog box showing both Local and Worklist data. Inconsistencies are not resolved automatically.

For help with Worklist Problems, refer to Chapter 6, Section 4.0, page 146.

4.0 Working with the Databases

4.1 Using the Bar Code Scanner

The Acquisition Workstation supports a single bar code scanner that identifies a Patient already in the Worklist or the local database. Scanning a barcode automatically generates a query to match the Patient ID or Accession Number (whichever is configured). It queries whichever database is selected before the barcode is scanned.

4.2 Querying the Modality Worklist Provider




Figure 3-9: Fields Used to Filter a Query to the Modality Worklist

To query for a Patient from the Modality Worklist Provider:

1. Click the **Worklist** tab.
2. To input the query information:
 - a. Bar Code Scanner—scan the patient ID barcode. The patient's scheduled procedure displays and the patient is added to the local database.
 - b. Keyboard—one or more fields may be used to query the Modality Worklist Provider. Enter the information and click **Query** or press **Enter**.

A sample query for a patient scheduled for today's date (enter into the appropriate field):

- The patient's name:
- Patient ID:
- Accession Number:
- Requested Procedure ID:
- From the Date drop-down menu select **today**.






Figure 3-10: Results of the Patient Query

3. The search results are displayed on the screen.
4. If you changed the date from the 'today' default, you can click **Today's Worklist** tab to display today's appointments.
5. The Patients who are scheduled for today and meet the criteria you specified are displayed on the screen.





Figure 3-11: List of Patients From the Today's Worklist Query

6. Select the patient to open their procedure.
7. If your search for a single match with a patient ID or Accession Number is successful, the patient's Procedure opens automatically.

4.2.1 Additional Query Information and Hints

- Typing partial information for the patient name assumes a wildcard character at the end of the name.



Note...

You can use the asterisk () character as a wildcard character to broaden your search. For example, entering R*ph would display Patients named Randolph and Rudolph. If you are using the wildcard character, make sure you include at least one other character (that is, don't enter only an asterisk) to avoid listing all Patients or Procedures, which might take several minutes to complete.*

- One or multiple fields may be used as search criteria.
- You can enter a partial name and press Enter on the keyboard to initiate the search (no need to click the Query button).
- Limit the query to a date range using the today drop-down field and selecting a date range.
- Most Worklist Providers require exact matching on Patient ID and Accession Number fields and are case sensitive.
- Names are stored in the database as LAST^FIRST^MIDDLE. You must include the caret (^) character between the name parts (for example, DOE^J).

4.2.2 Using the Patient Query feature

The resulting list displays Patients that met the criteria. The number of results is in the Query pane and the data appear in the **Patient List**. To scroll through a list longer than one screen, use the arrows on the right of the screen. Depending on the criteria you selected, patients that are not scheduled for today may be in this list.

4.2.3 Querying Modality Worklist Connections

In sites with a Modality Worklist Provider, select the **Worklist** tab (if it is not already highlighted) to query the MWP.

4.2.4 Using the Today's Worklist Feature

If you did not use the 'today' date query, you can list the appointments for today by clicking on the **Today's Worklist** tab. (If you used the 'today' date query, the results are the same.)

Click the Patient name in the list to add their data to the local database (Local Exams) and open their procedure.

4.3 Searching Local Exams—Acquisition Workstation Local Database




Figure 3-12: Selected Patients Appear in the Local Database, Local Exams Screen

Patients and related information are kept temporarily, storage space availability permitting, in the local database to save having to re-enter information for a recently X-rayed Patient.

4.3.1 How Patients are Added to the Local Database

The Acquisition Workstation Local Database acquires a Patient/exam entry by:

- **Manual Creation:** Creating a new Patient/exam at the Acquisition Workstation.
- **Modality Worklist:** The selection of a worklist item from a list of items obtained by querying a DICOM Modality Worklist Service Class Provider (MWL SCP).

4.3.2 Search the Local Exams




Figure 3-13: Selecting an Existing Patient from a Search List

- To search for a Patient in the Local database:
 - Click the **Patient Name** or **ID** tab. Enter the Name or ID in the Search box and either click **Search** or press **Enter**.
 - Enter part of a name or ID number and click **Search**. If you are not typing the first part of the ID, you must put an asterisk in the beginning of the name or ID to stand for the part that you did not enter.




Figure 3-14: Typing a Patient Name in the Search Field

- The system searches its local database for the Patient or Procedure ID and displays the list of patients matching the search criteria (see Figure 3-15).




Figure 3-15: Select a Patient Screen

- Select the Patient in the list:
 - Their procedure opens in the Patient Views screen.
 - If the Patient has more than one Procedure scheduled, the Choose a Procedure Screen displays with a list of their procedures.

4.3.3 If the List Is Too Long



Figure 3-16: List is Too Long

If your results give a very long list, you may get an error message. Include more of the name or ID number and try searching again. If the list is long, you can search for an exact name or Patient ID.

1. Select the appropriate tab.
2. Enter a name, or partial name or partial ID in the field to the left of the Search button.
3. Click Search. The resulting list contains only those Patients that meet the requested criteria.
4. Select the Patient you want from the shorter list.

4.3.4 Guidelines for Entering Names in a Search

Note the following guidelines when entering names:

- Names are stored in the database as LAST^FIRST^MIDDLE. You must include the caret (^) character between the name parts (for example, DOE^J*).
- Because names are stored in the format LAST^FIRST^MIDDLE, to search for the first name Jane you must use an asterisk (for example, *Jane).
- Search entries are not case sensitive; you can enter upper or lower case letters.

4.3.5 Patient Not Found

If there is no match, the No Patients Matching Query dialog box appears.




Figure 3-17: No Patients Matching Query Message

4.4 Searching Patient Data on a System Without a Worklist

The Acquisition Workstation Local database contains the demographic information for un-reclaimed Patients examined on this system, and selected patients from the Worklist Patient Query. The Worklist/Local Exams tabs do not exist if your system is not configured for a Modality Worklist Provider.

1. Select the **Patient Name** or **Patient ID** tab, enter the information to search for, and click **Search**. The search results are listed. The arrows on the left side of the screen assist in scrolling through the list.




Figure 3-18: Search for a Patient in a System Without a Worklist

2. Select the Patient from the list. The patient's procedure displays in the Patient View screen.

5.0 Choose a Procedure



Note...

A Procedure ID tracks a Procedure (series of images). This ID relates the images to a particular Procedure and Patient.

This option only appears when the patient has more than one scheduled procedure.




Figure 3-19: Choose a Procedure

Legend for Figure 3-19

1. List of Procedures for a Patient with more than one scheduled Procedure.
2. **Select** Button to open the procedures you have selected.
3. **Close** Button to close the Patient without opening any Procedure.
4. **New** Button to schedule a new Procedure.

6.0 Patient View Screen




Figure 3-20: Patient View Screen Panes

Legend for Figure 3-20 and section where discussed

1. Patient View Icons Pane, Implant Present: Section 6.6, page 41
2. Patient Pane, Edit button: Section 10.3.1, page 57
3. Procedure Pane
 - New button: Section 10.3.2, page 57
 - Edit button: Section 10.3.3, page 58
 - Add button: Section 10.3.3, page 58
4. Outputs Pane: Section 6.5, page 41
5. Special Mode Buttons
 - Non-Imaging Mode button: Section 10.4, page 60
 - Simulate Capture button: Section 10.5, page 61
 - Review buttons: Section 10.6, page 61
6. Exposure Information and Techniques Pane: Section 6.2, page 36
7. Generator Pane
 - Generator Status: Section 7.0, page 41
 - Generator Error Message: Chapter 6, Section 1.1.3, page 130
 - Reset Generator button: Chapter 6, Section 1.1.4, page 130
8. System Messages Pane: Chapter 6, Section 2.0, page 133

6.1 Patient View Screen Options

This chart diagrams the choices available on the Patient View screen.




Figure 3-21: Patient Views Options

The next step in the normal workflow is to select an exposure mode for the procedure.

6.2 Selecting Exposure Modes

Selenia functions in these Exposure Modes:




Figure 3-22: Exposure Modes

To choose the mode, click the up/down arrows in the Exposure Mode pane (the lower far-left box). The exposure defaults for that mode appear to the right in the pane. For more information about setting defaults, refer to Section 18.2, page 71.

Use the trackball or the function keys to change Exposure techniques.



Note... Changing the mode resets all Exposure Techniques to the defaults for that mode.

6.3 TEC Exposure Mode

1. The Gantry reads the compressed breast thickness and updates the recommended X-ray exposure parameters for kV, mAs, and filter using the standard reference table. It continues to update and recalculate until you override or accept the techniques.
2. Verify the selections are appropriate and:
 - Either click **Accept**, below the breast thickness (or press **Enter**).
 - Or change the techniques you do not agree with (Override mode).




Figure 3-23: Accepting TEC Techniques

If you agree and click **Accept** in the Breast Density box, the system is in the Accept Mode.

- It stops recalculating the Technique settings.
- It displays “Ready” in the Generator Status message on the screen.
- It is ready to acquire an Image.




Figure 3-24: TEC with Techniques Accepted

If you change one or more of the parameters (kV, mAs, or filter), the system enters Override Mode:

- It stops recalculating the Technique settings.
- It displays the Override message on the screen.
- It is ready to acquire an Image.




Figure 3-25: Override Mode

The Mammography unit remains in standby until one of the two modes appears. When the Override message displays, the system is ready to acquire an image

6.4 AEC Exposure Mode

AEC is calibrated to give a set dose for a standard breast (50/50, 4.5 cm thick). During AEC, the kVp is calculated by the set AEC table as a function of compression thickness. The pre-exposure (50 milliseconds) is used to calculate the mAs. During the pre-exposure, the resulting pixel value is evaluated and the appropriate mAs that produces a targeted pixel value is determined.




Figure 3-26: AEC Modes—Choices That Must be Made by the Technologist

Table 3-1: AEC Mode—Available Selections

AEC Mode	Techniques
Auto-Time	Filter = user selected kV = user selected mAs = calculated from Scout Pulse
Auto-kV	Filter = not selectable (Mo if Mo tube, Rh if W tube) kV = calculated from compression thickness and lookup table mAs = calculated from Scout Pulse
Auto-Filter	Filter = calculated from thickness and lookup table kV = calculated from compression thickness and lookup table mAs = calculated from Scout Pulse

6.4.1 Auto-Time Mode

In Auto-Time mode, select all exposure factors except mAs, and use Auto AEC or select the AEC Sensor Position. Auto-Time terminates the exposure after the AEC circuitry senses a pre-determined amount of radiation. A final mAs value displays after the exposure.

6.4.2 Auto-kV Mode

In Auto-kV mode, the system automatically selects the exposure kV and mAs. Once you select the Focal Spot and the AEC Sensor position or Auto AEC, the system is ready for the exposure. A Final kV, Post-mAs value, and Filter selection displays after the exposure. (The Filter is not selectable; the filter defaults to Mo if the system is configured with a Molybdenum x-ray tube and Rh if configured with the Tungsten x-ray tube.)

6.4.3 Auto-Filter Mode

In Auto-Filter mode, the system automatically selects all the Technique factors. Once you select the AEC Sensor position or Auto AEC, this is a fully automatic operation. The system automatically sets the kV using the lookup table based on the compressed breast thickness and changes the filter type when indicated by the table. A Final kV, Post-mAs value, and Filter selection displays after the exposure.

6.4.4 AEC Exposure Adjustment Settings

If you increase or decrease the targeted pixel value by using the Exposure Adjustment Settings, each step changes the target pixel value by approximately 15% which corresponds to a change in dose of the same magnitude.

The control range is from +4 to -3 to increase or decrease the exposure. Each step corresponds to a 15% change in exposure from the previous step. The calculated mAs increases or decreases with this change and remains until you change it or reboot the system.



Warning:

Increasing the AEC exposure adjustment setting will increase patient dose to higher levels. Decreasing the AEC exposure may result in increased image noise or degraded image quality.



Warning:

Increased Exposure adjustment leads to higher dose to patient. Keep patient dose as low as practical to obtain good image quality.




Figure 3-27: AEC Exposure Adjustment Setting on Screen

6.4.5 AEC Sensor Position Indicator Settings

You can select one of the seven positions to define the area of interest for AEC calculations. You can also select Auto AEC where a larger area is evaluated for the densest areas of the breast, and those are considered for the AEC calculation.

**Note...**

Always be aware of the position of the AEC Sensor when making AEC exposures.

Set the AEC Sensor Position after compressing the breast. There are three ways to change the AEC Position selection:

- The buttons on the compression device.
- The AEC Pos key on the keyboard.
- The AEC position arrows on the Acquisition Workstation display.

The displays on the side of the Compression Device and the lower right corner of the screen indicate the position setting.

The AEC Sensor Position has eight available positions.

- Positions 1 through 7 are manual selections, from the chest wall edge (position 1) to the nipple edge (position 7).
- Position 8 (Auto) automatically positions the sensor after sampling, and adjusts the sensor position for the two 1 x 1 cm regions.

The AEC Sensor uses areas of the detector rather than physical radiation sensors used in screen-film systems. There are two types of AEC Sensors. The Fixed AEC Sensor is used when a selection from 1 to 7, corresponding to the markings on the screening paddles, is selected. Manually selected Sensor positions do not change from your selection unless the collimator size decreases. This happens when:

- A smaller compression paddle changes the collimator.
- Shifting a paddle will change it to Auto AEC Position.
- The collimator override is used to select a smaller size.

Then, if necessary, the AEC sensor position moves to just within the field. (This change does not occur when the collimator size increases.)

The Auto AEC Sensor covers a larger area. A detector region, with a size defined by paddle and shift position, is used. This region is divided into approximately 1 cm² areas. Auto AEC positions one sensor beneath the densest part, and the second beneath the next densest area and averages them.

When viewing an AEC Preview, you can display the AEC Sensor area(s) without regard for which method, Manual or Auto Sensor, was used for the exposure. Refer to Figure 3-27, page 39.

6.5 Select the Output

To change an Output, select a new output from the drop-down list. When you change the Outputs option, a dialog box appears, informing you of the change. Click **OK**. If you select an output with no destination devices, the dialog box asks if you really do not want any outputs. If this was a mistake, reselect the output. Outputs can be a group of outputs, for example, a diagnostic review workstation and CAD, a combination of printers, and an archive system. For more details about outputs, refer to Chapter 2, "What Must I Know to Acquire Images on Patients with Implants?," page 4.




Figure 3-28: No Outputs Warning



Note...

Only accepted Images are sent to the outputs. Rejected Images are sent to a storage area.

6.6 Implant Present

If the patient has an implant, select the Implant Present checkbox above the views. All the images are marked as having an implant and Implant Processing is applied to all views until the box is unchecked or you close the procedure.



Note...

Select the **Implant Present** box for all Implant views and Implant Displaced (ID) views.



Note...

For full Implant views, use either AEC (Auto Filter) or Manual Exposure Modes. For Implant Displaced views, use AEC and manual placement of the AEC sensor. Refer to Chapter 2, "What Must I Know to Acquire Images on Patients with Implants?," page 4.

7.0 Acquiring Images

The system is ready to acquire Images when the Generator Status indicates "Ready" and the X through the X-ray symbol in the Generator Status box is not displayed. The Generator Status Message changes during the X-ray: Boost, Rotor, X-ray On (the high-pitched tone starts), X-ray Off.



Caution:

Do not move the C-arm while the system retrieves the image. Do not change the paddle until after you accept or reject the image.

8.0 Working with Images

You use the Preview Image screen to:

- Evaluate the Exposure Index and the Preview Image
- Add Comments
- Edit the View

- Make enhancements with the Preview tools
- Accept and reject images

8.1 Preview Image




Figure 3-29: Image Preview Screen

The Acquisition Workstation proceeds into the Preview mode as the image is processed. The Preview Image is a non-diagnostic version of the image with processing applied for confirming anatomy and positioning. The Image pane displays how the image looks at the output device, scaled down 4:1. The Patient Name, ID, and View appear on the top bar of the Preview.

You may apply Window/Level adjustments or use other Preview tools located on the right side of the screen. Any changes to the Window/Level will not be saved to Outputs.

8.1.1 Progressive Preview

Generally, the first preview image is available in a short time for a quick evaluation of positioning. This first preview image updates with an improved preview image, which provides better image quality to evaluate motion or image artifacts.

8.1.2 Film Label Area

The small rectangle in the upper left corner is the film label area and is blank in the Preview. The label box is shown to ensure that there is no anatomy underneath the box, since that is where the label appears if the image is printed.

8.1.3 Digital Markers

A Digital Marker next to the label appears automatically. The technologist's initials are displayed and the chosen View selects the appropriate marker. You can click and drag the marker or change it using the Markers buttons on the Tool Bar.

8.1.4 Hanging Options

Images are presented on the screen as chest wall on the left or right (Dorsal or Ventral.) The graphics in this manual are all in **chest wall on the right** (or Dorsal) mode.

These options can affect how images display on your diagnostic review workstation and should be confirmed with the Vendor.

8.2 Comments

The Comments button allows you to enter or modify Image Comments during the Procedure while the Preview screen is open. The exact comment text is copied into the DICOM header of the Acquired Image. The diagnostic review workstation displays any comments. The comments can be printed on film output from the diagnostic review workstation.



Note...

*Clicking **Clear Comments** removes all of the image's comments.*

8.3 Implant Present

A checkbox was available in the Patient View screen to allow you to mark when an implant is present. The Preview screen displays an Implant Present checkbox which you can select or unselect to change the information.

8.3.1 How to Correct Implant Information Before You Accept the Image

If you realize the implant information on an acquired image is incorrect before you Accept the image:

1. Click the **Implant present** checkbox on the Preview Screen.
2. Click the **Reprocess View** button on the bottom of the preview screen. The image is reprocessed and you can then click **Accept**.
- 3.

8.3.2 How to Correct Implant Information After You Accept the Image

If you realize the implant information is incorrect and you have already accepted the image:

1. Repreview the image by clicking the thumbnail.
2. Select or unselect **Implant Present** on the Preview screen.
3. Click the **Reprocess View** button.
4. Click **Save**. A Send Image to Output dialog box asks you to confirm that you want this image sent to the selected outputs.
 - If you click Yes, images are sent to the selected output(s) group.
 - If you click No, the image can be resent using Image/Spool Management.

5. The incorrect image was sent to the chosen outputs. Corrections must be made and then you must use Image or Spool Management and resend the corrected image. Refer to Section 19.1.2, page 84.
 - If you printed the image, collect and discard the incorrect image.
 - If the image was sent to a diagnostic review workstation, delete the incorrectly processed image.
 - If the image was sent to the PACS, contact the PACS administrator and have them delete the incorrectly processed study.

8.4 Zoom/Pan

Two methods of Zoom/Pan are available: Quick Zoom/Pan and Full Zoom Pan.

8.4.1 Quick Zoom/Pan

When Quick Zoom/Pan is active, a magnification of the Preview Image displays using the entire Preview area. The magnification factor depends on the paddle. The 24 x 30 cm paddle has less magnification than the others. The placement defaults on the centerline of the chest wall. Activation and deactivation of this feature does not affect any crosshairs or line measurements.

8.4.2 Full Zoom/Pan

Full Zoom/Pan is available after the final image is prepared. When Full Zoom/Pan is active, a 70-micron magnification of the Preview Image is shown using the entire Preview area. Each pixel of the image display represents a detector pixel. Activation and deactivation of this feature does not affect any crosshairs or line measurements.

This thumbnail view, seen at the upper right part of the screen, indicates where the magnified view is located in relation to the entire image. You can drag this box to get to where you want to be in the image.

8.5 ROI

The ROI Statistics dialog box shows the:

- location of the left, top and the width and height of the ROI box in pixel units.
- minimum and maximum pixel value with in the ROI box.
- mean, standard deviation, and signal-to-noise ratio (SNR) measurements for the ROI. SNR does not compensate for the DC offset and should not be used for QC.




Figure 3-30: The ROI Statistics Information on the Preview Screen

8.6 Measurement



This feature measures distances on the image.

1. To clear the line, click the '+' at either end of the line.
 - Click the **Measurement tool** to turn Measurement off or to hide the numbers.
 - Turning the Measure tool off does not clear the caliper lines. They are automatically cleared when the Preview window is closed.
 - Caliper lines are not saved or recalled.
 - You can draw multiple lines.

8.7 Magnification



Magnification of the image is available with a short delay after the image appears. The button labels appear when this option is available. Using the Magnifier, you can view part of the image at 2X the currently displayed resolution.

A full resolution view of a part of the image is available when the Acquisition Workstation receives the entire image. You can also use the magnifier on that full zoom image after locating an area of interest and clicking the magnification button to further magnify that area.

8.8 Window/Level



Note...

This operation has no effect on the final processed Image.




Figure 3-31: Changing the Window/Level Settings

You can edit the **Window/Level** by clicking in the box you want to change and entering a new number. When you click **Set** the image display changes to the new setting. Double click the image to return to the original setting.

8.9 Cropping

The Acquisition Workstation images are acquired with a collimated X-ray source. When an 18 x 24 cm or smaller paddle is installed with automatic collimation, the system crops the image to 18 x 24 cm. The area outside is not stored. You cannot change cropping. When manual collimation is selected, cropping follows the field.

8.10 Accept or Reject Images

8.10.1 Accept Button Details




Figure 3-32: The Accept Button

When you accept the Preview Image:

1. The full resolution image with all associated attributes is marked as accepted in the pool space. Acceptance “freezes” the Acquired Image, and no subsequent action changes any attribute of an Acquired Image. It can be rejected later.
2. The system transmits the image to the selected destination output devices.

8.10.2 The Image Repetition Information Dialog Box

When you have already accepted an image for a view and then acquire and accept a second image for the view, the Image Repetition Information dialog box appears. Except for the choices Body Habitus and Wire Localization, the information collected here goes in the Repeat/Reject report count for the logged in technologist.

Select the reason for the repeat and close the dialog box.

8.10.3 Reject Button Details

To reject an image:

1. Click **Reject**. An **Image Rejection Information** dialog box appears, listing possible reasons for rejection and providing a text field for typing a reason.
2. Select a reason for rejection.



Note...

This information is used in the Reject Analysis (required for Quality Control.)

3. Click **Reject** to close the dialog box, or click **Return to Preview** if you decide not to reject the image. The system moves the image to the Reject Bin.
4. After you reject an image, the Preview closes and you can:
 - Re-acquire the rejected View.
 - Select a different View.
 - Close the Procedure.

An asterisk is displayed next to the view label of any views with a rejected image.



8.11 Accept a Rejected Image

You can reverse the Rejection process at any time. If you reverse the rejection process before the procedure is closed, the change is reflected in the Reject Analysis Report. If you reverse the rejection process after the procedure is closed, the change is not reflected in the Reject Analysis Report.

If you decide that a Rejected Image is better than the following image:

1. Reject the second image following the steps in Section 8.10, page 46.
2. In the Patient View screen, click **Review** (below the Outputs box).





Figure 3-33: Review Dialog Box Showing Two Rejected Images.

3. In the **Review** dialog box, select the **View** tab for the View. (for example, R MLO).
4. Click the **Rejected Image** thumbnail you now want to accept. A Preview image appears.
5. Click **Comments** if you wish to add or delete comments.
6. Click **OK** to close the **Comments** dialog box.
7. Click **Set Accepted**.
8. Click **OK** to close the **Review** dialog box. The Accepted Image appears as a thumbnail on the Patient View screen.






Figure 3-34: Accept Rejected Image

8.11.1 Sending Accepted Rejected Images to Outputs



Figure 3-35: Send Accepted Image to Outputs

When you accept a rejected Image, you must confirm that you want this Image sent to the selected outputs.

- If you click **Yes**, Images are sent to the selected output(s) group.
- If you click **No**, the image is not sent. It can be resent using Image/Spool Management at a later time.

8.11.2 Reject an Accepted Image

If you decide to reject the Accepted Image:

1. In the **Patient View** screen, click **Review**.
2. In the dialog box, select the **View** tab for the View you are rejecting.
3. Click the Accepted Image (Figure 3-34). The Preview screen displays.
4. Click **Set Rejected**. The **Image Rejection Information** dialog box appears.
5. Choose a rejection reason.
6. Click **Reject**. The system reprocesses it and labels this image rejected.

- Click **OK**. The displayed thumbnail disappears from the Patient View screen.




Figure 3-36: Set Rejected Button



Note...

After you reject an Image using this method, it is **not** removed from the outputs destination queues, for example, diagnostic review workstation, or archive. It can be deleted from the queues. Refer to Section 19.3, page 87. The Rejected Image is added to the reject analysis bin for tracking purposes.

9.0 Working with Images After They Are Accepted

9.1 Send the Images to the Output Devices

After you accept an image, the Acquisition Workstation sends the image to the destination devices in the selected output. The system makes entries in all queues that are on the output list for the image. If you accept a rejected image, confirm that you want this image sent.



Figure 3-37: Send Accepted Image to Outputs

If you ignore the warning and acquire images with no output selected, use the Resend function to send those images to an output. The system saves all Acquired Images on the Acquisition Workstation hard drive regardless of the output selection. The system reclaims space in the Image pool only after the Image is successfully sent to all of its output destinations. Refer to Section 18.6.2, page 79, for details.

9.2 Patient View Screen with Thumbnails

After the Image is accepted, a thumbnail view of the Image replaces the icon on the Patient View screen, and the system transmits the Image to its outputs. Click the thumbnail to repreview the image.




Figure 3-38: Patient View Showing Acquired Image

9.3 Exposure Taken Using an Incorrect Icon

9.3.1 You Have an Edit View Button

If your Preview screen looks like this in the top right corner and an incorrect view is acquired (e.g., RCC icon highlighted, but patient is positioned for LCC):

- If you *have not* accepted the image:
 - A. Select the **Edit View** button.
 - B. Select the correct view.
 - C. Select **OK**. A corrected thumbnail replaces the incorrect icon. Add another view for the icon that was replaced.
- If you *have* accepted the image:
 - A. Repreview the image.
 - B. Make the view corrections. Refer to Chapter 3, Section 19.1, page 82 for Repreview and Resend directions.



9.3.2 You Do Not Have an Edit View Button


If your Preview screen looks like this in the top right corner and an incorrect view is acquired (e.g., RCC icon highlighted, but patient is positioned for LCC) there are two scenarios to consider:



Scenario One

The mistake is discovered *before* the image is accepted:

1. With the Preview screen displayed, click the **Clear All Markers** button (upper right side of screen); the incorrect marker is removed.
2. Click the correct marker tab.
3. Click the correct marker. A new marker appears on the preview image.
4. Use the table on the following pages as a guide to determine if additional changes are required. These changes may include:
 - Drag the new marker to the bottom of the image. (This is done by putting your cursor on the new marker, then holding the left trackball button down while moving the trackball to move the new marker to the desired area.)
 - Change the Image Laterality (middle right side of screen).
 - Change the Patient Orientation Column.
5. Click the **Accept** button, after the changes have been made.



Note

The flipped images do not show on the Acquisition Workstation, but are shown on the SecurView®. They may not show on 3rd party systems.



Note

CCL, CV, RL, RM, TAN, CCID - these views all behave the same as CC. ML, AT, MLOID - these views all behave the same as MLO. LMO, SIO - these views all behave the same as LM.

Scenario Two

The mistake is discovered *after* an image once is accepted:

If the error is not corrected before you Accept the image, the system sends the incorrect image to the chosen outputs.

1. Corrections need to be made at the outputs before you resend the corrected image.
 - Printer—Collect incorrectly printed images and discard.
 - SecurView—Log on, user name-delete, click *Administration*, click the patient to be deleted, click the *Delete* button, then click the *Exit* button to close the window.
 - PACS—Contact your PACS administrator to delete the incorrectly marked study.
2. Corrections must be made in Image (or Spool) Management.
3. Select **Admin>Image (or Spool) Management**.
 - a. Search for the patient.
 - b. Click the thumbnail of the image that needs to be corrected.
4. Repreview the image to make the corrections.
 - a. Click the **Repreview** tab then the **Repreview** button.
 - b. Make the corrections as previously indicated in Scenario One.
 - c. Click the **Save** button. This creates an additional thumbnail showing the corrected image.
5. Click the **Resend** tab.
 - a. De-select the image that you have corrected.
 - b. Click the last thumbnail (corrected image), along with any other images that need to be resent.
 - c. Click the **Resend selected image(s)** button.
 - d. Go to Outputs, and select the correct outputs.
 - e. Click the **Resend** button.
 - f. Click the **OK** button in the Resend Confirmation dialog box.
6. Click the **Close** button (bottom right hand corner).
7. A manager can delete the incorrect image. Refer to Chapter 5, Section 4.1.8, page 114.

**Note**

If there is any doubt of the accuracy of the corrected images, repeat the image.

**Note**

The correction may only be made if the procedure is still on the Acquisition Workstation.

**Note**

Be sure the incorrect images are still on the Acquisition Workstation before you delete them from the PACS.

9.3.3 Changing the View Selection

Table 3-2: Incorrect View Selection Changes

The view you wanted	View icon chosen in error	Change Marker to:	Drag	Image Laterality	Patient Orientation Column
RCC	RMLO	RCC	no	no change	no change
	RLM	RCC	bottom	no change	H
	RFB	RCC	bottom	no change	R
	LCC	RCC	bottom	R	no change
	LMLO	RCC	bottom	R	H
	LLM	RCC	no	R	no change
	LFB	RCC	no	R	L
RMLO	RCC	RMLO	no	no change	no change
	RLM	RMLO	bottom	no change	H
	RFB	RMLO	bottom	no change	R
	LCC	RMLO	bottom	R	no change
	LMLO	RMLO	bottom	R	H
	LLM	RMLO	no	R	no change
	LFB	RMLO	no	R	L
RML	RLM	RML	bottom	no change	H
RLM	RCC	RLM	bottom	no change	R
	RMLO	RLM	bottom	no change	H
	RFB	RLM	no	no change	no change
	LCC	RLM	no	R	L
	LMLO	RLM	no	R	no change
	LLM	RLM	bottom	R	H
	LFB	RLM	bottom	R	R
	RML	RLM	bottom	no change	H
RFB	RCC	RFB	bottom	no change	R
	RMLO	RFB	bottom	no change	H
	RLM	RFB	no	no change	no change
	LCC	RFB	no	R	L
	LMLO	RFB	no	R	no change
	LLM	RFB	bottom	R	H
	LFB	RFB	bottom	R	no change
LCC	LMLO	LCC	no	no change	no change
	LLM	LCC	bottom	no change	H
	LFB	LCC	bottom	no change	L
	RCC	LCC	bottom	L	no change
	RMLO	LCC	bottom	L	H
	RLM	LCC	no	L	no change
	RFB	LCC	no	L	R

Table 3-2: Incorrect View Selection Changes

The view you wanted	View icon chosen in error	Change Marker to:	Drag	Image Laterality	Patient Orientation Column
LMLO	LCC	LMLO	no	no change	no change
	LLM	LMLO	bottom	no change	H
	LFB	LMLO	bottom	no change	L
	RCC	LMLO	bottom	L	no change
	RMLO	LMLO	bottom	L	H
	RLM	LMLO	no	L	no change
	RFB	LMLO	no	L	R
LML	LLM	LML	bottom	no change	H
LLM	LCC	LLM	bottom	no change	L
	LMLO	LLM	bottom	no change	H
	LFB	LLM	no	no change	no change
	RCC	LLM	no	L	R
	RMLO	LLM	no	L	no change
	RLM	LLM	bottom	L	H
	RFB	LLM	bottom	L	L
LML	LLM	LLM	bottom	no	H
LFB	LCC	LFB	bottom	no change	L
	LMLO	LFB	bottom	no change	H
	LLM	LFB	no	no change	no change
	RCC	LFB	no	L	R
	RMLO	LFB	no	L	no change
	RLM	LFB	bottom	L	H
	RFB	LFB	bottom	L	no change

9.4 Close a Procedure

9.4.1 The Close Procedure Button



Figure 3-39: Patient View Screen's Close Procedure button

1. Click the **Close Procedure** button or press the **Close Exam** key on the keyboard when you have finished. The Procedure remains open until you specifically close it. The **Patient View** screen's **Close Procedure** button allows you to end the Procedure any time.




Figure 3-40: A Finished Standard Screening Procedure

2. If the Procedure is not finished, a dialog box appears. Select **Yes** to close. The Select a Patient screen displays. Choose the next Procedure.

9.4.2 Close a Procedure with MPPS Enabled

An optional feature is available for sites with an MPPS Service Class Provider. When this feature is installed, you have three buttons to use to close a procedure:

- Complete button: closes the active procedure and its status is labeled completed.
- Pend button: closes the active procedure and its status is labeled in progress.
- Discontinue button: closes the active procedure and its status is labeled discontinued. Discontinuing a procedure requires a reason.

10.0 Additional Actions on the Patient View Screen




Figure 3-41: Patient Views Screen with Icon(s) for Scheduled Views

Legend for Figure 3-41

1. Add View
2. Add or Edit a Procedure
3. Output selection drop-down list

10.1 Add a View to the Procedure

1. Click **Add View**, (Figure 3-41, #1).
2. In the dialog box, select an additional View(s) from the list of Standard Mammography Views. You can add multiple views by holding the Control key down while clicking.
3. Click **OK** to continue. The new View Icon(s) appears with the other views.
4. Select the new view before acquiring the image.

10.2 Remove Unexposed Views Button

Unused views or procedures are automatically reclaimed, or can be deleted manually with the Remove Unexposed Views button. When you click the button you need to select **All** or **Selected** to indicate which views are to be deleted.

10.3 Edit Patients and Procedures

10.3.1 Edit Patient Data

If information about the Patient has changed and needs to be edited:



Note...

*This procedure only changes the patient information of images in the study not yet taken. Acquired images will have to be resent using **Image Management** and **Resend** after the information is edited. Refer to Section 19.1, page 82.*

Do not edit patient information if using a Modality Worklist. Work with the HIS/RIS administrator to resolve incorrect patient information. Refer to Chapter 3, Section 19.1, page 82.

1. Click **Edit** in the **Patient Box**.
2. Select the field(s) with incorrect or missing information and make changes to the field(s). Use the Tab key or the trackball to move through the fields.
3. Click **Accept**. The dialog box closes and the main screen updates with the new information.



Note...

The Clear button clears last name, first name, DOB, and age.

4. Verify that the correct changes appear in the screen before acquiring a new image.

10.3.2 Add a New Procedure to an Existing Patient

1. Recall the existing Patient.
2. Click **New** in the *Procedure* box and click **Yes to: Creating a New Procedure will close the current Procedure**. The **New Procedure** dialog box appears.
3. Click the **Accession Number** field. Enter a new **Accession** Number.



Note...

The Clear button clears the Accession number so you can re-enter it.

4. Select the **Procedure Description** from the drop down list(s).
5. Click **Accept** and then click **Close** in the Patient Pane to close this screen.



Note...

If the Patient already has more than one Procedure, you can add a new one directly from the Select a Procedure search screen. Otherwise, when you select the Patient, the scheduled Patient Procedure opens and you must close the Procedure before you can add the new Procedure.

10.3.3 Add a New Procedure to an Open Patient

You can have multiple Procedures open at the same time. This allows you to compare other images to a new image. To add a new procedure when you want the current procedure to remain open:

1. Click **Add** in the **Procedure Box**. The **New Procedure** dialog box appears.
2. Select the **Procedure Description** from the drop down list.
3. Click **Accept**. A new **Procedure** tab appears on the **Patient View** screen above the View icons and you can work with either Procedure (Figure 3-43).




Figure 3-42: Results of Adding Another Procedure



Note...

This new procedure must use the same Accession Number as the currently open procedure. Use the **New** button if this is not your intention. If your facility requires different accession numbers for different procedures you need to use the **New** button. Refer to Section 10.3.2, page 57. It is not possible to change the Accession Number if you use the **Add** button.




Figure 3-43: Selecting Two Procedures in the Choose a Procedure Screen

The images can be taken now, or later. If taken later, the two Procedures are treated as separate Procedures in the Choose a Procedure screen. They can both be opened as one by pressing and holding the Ctrl key and clicking each one. Then release the Ctrl key and click **Select**.

10.3.4 Changing Patient Information After Accepting an Image

Once an image has been accepted, it is immediately sent to the selected outputs. If any change to the accompanying information is needed, it must be reviewed, saved, and resent as an additional record and the previous record deleted.

Delete the incorrect image from the following outputs:

- Printer: Collect the incorrect printed images and discard them.
- DRW: Contact the Diagnostic Review Workstation manager to delete the incorrect image(s).
- PACS: Contact the PACS administrator and alert them to delete the incorrect study or advise another corrective measure.

To correct the information at the AWS:

1. At the AWS, click the **Local Exams** tab, then search using the first letters of the last name.
2. When the patient procedure displays on the screen, click **Edit** in the **Patient Pane**.
3. Change the incorrect patient information.
4. Click the **Accept** button.
5. If the **Accession Number** is also incorrect, click **Edit** in the **Procedure Pane** and make the correction.
6. Click the **Close Procedure** button.
7. Select **Admin>Image Mgt**.
8. Search for the patient using the first letters of the last name.
9. Click the **Repreview** tab.
10. Click the first thumbnail.
11. Click the **Repreview** button. The preview image appears.
12. Click the **Save** button.
13. Repeat steps 10 to 12 for the rest of the images in the procedure. There should now be eight thumbnails (assuming there were four in the original procedure.)
14. Click the last four thumbnails (the corrected views).
15. Click **Resend** to all required outputs.
16. After you confirm that the four images have been resent, you may want to go back into Spool Management and Delete the four images that are labeled incorrectly. (You have to be signed in at the manager level to do this.)

10.4 Non-Imaging Mode Button



Figure 3-44: Non-Imaging Mode Button

Non-Imaging mode allows you to make X-ray exposures without acquiring images. The medical physicist may use it during QC tests to measure the X-ray tube performance without having to wait for Images.

This mode:

- Is not available when you have selected a Patient.
- Deactivates the Image Receptor.
- Allows an X-ray prep/expose cycle to occur while overriding the Patient selection and Image Acquisition sequences.
- Is available to all users




Figure 3-45: Non-Imaging Screen

When you click **Non-Imaging Mode**, a clear message appears stating that you are entering this mode.



Caution:

Be sure to protect the Image Receptor from excessive radiation by covering it with lead while using this mode.



Figure 3-46: Non-Imaging Dialog Box Warning

To use this button:

1. Click **Non-Imaging Mode**. A Non-Image dialog box appears.
2. Click **OK**. The dialog box closes and the Non-Imaging window appears.
3. Follow the directions in the QC Manual or Service Manual to use the window.
4. Click **End Non-Imaging Mode** when finished.

10.5 Simulate Capture Button

This function is for use by service personnel only.

10.6 Review Button

This button, see Figure 3-38, Item 1, is used for Reject /Unreject. Select the thumbnail to review an image.


11.0 Clinical Procedures

This section provides a suggested Sequence of Operation to acquire Patient Images, but it is not intended as a substitute for conventional clinical training and methodology.

The following three sections provide flowcharts for typical clinical procedures.


11.1 Manual Mode

In the Manual mode, you set all the Techniques. The defaults appear when you choose the View, and then you make any necessary adjustment for the individual exposure.




11.2 TEC Mode

A second mode available to you is the Tissue Exposure Control (TEC) Mode. It is an enhanced Manual Exposure Control mode.



11.3 AEC Mode

The third mode available to the technologist is the Automatic Exposure Mode (AEC).



12.0 Performing a Screening Procedure

The Selenia is capable of performing routine screening, diagnostic studies, magnification studies, and localization procedures. Each of these procedures requires some different accessories, but all are performed in the same general sequence of steps. The procedure for selecting the patient and view for a Magnification or Spot procedure is the same as for a Standard Screening. A Localization procedure has additional steps, refer to Section 14.0, page 68.

12.1 Select Patient, Procedure, Output, and Exposure Techniques

At the Acquisition Workstation:

1. Select a Patient from the Worklist or local database, or add a new one.
 - If your installation uses Patient or Procedure ID bar codes, use the bar code scanner to scan patient or procedure (accession) records in both the modality worklist and the local database.
 - If your installation does not use bar codes, search the modality worklist.
 - If your installation does not use a network, search the local database records.
 - If the patient is not in the local database, enter patient data using the keyboard.
2. Identify the required Procedure(s).
 - When working with implants, the beginning procedures are the same as a screening except select **Screening Mammogram w/Implants** procedure in **New Procedure** dialog box.
3. Select the Output Destination(s).



Note... *If using the CR-RW drive as an output, refer to Section 18.6.3, page 80.*

4. Select the first View.
5. Set the Exposure Mode and Techniques.

12.2 Set Up the Gantry



Warning: ***All C-arm movement is motorized.***

1. Install the compression paddle for the procedure and any other necessary accessories.
2. Set the C-arm height and rotation angle.



Note... *To help with patient positioning, the system memorizes oblique C-arm angles. After making an oblique exposure, the C-arm angle stores in memory. When you rotate the C-arm for the next exposure, the system automatically stops the C-arm at the opposite oblique angle. This is true for the C-arm Rotation switches, and the Gantry Rotation switches if equipped. You may then rotate the C-arm to an angle other than the memorized oblique angle to best accomplish breast positioning.*

3. Ensure the light field illuminates the correct area for the installed accessories.


12.3 Position the Patient

1. Use the C-arm height controls to raise or lower the breast platform to accommodate each patient. Do not use any heating source, for example, heating pad, on the image receptor.
2. Use the C-arm Rotation-Release controls to move the tubehead to the required clinical projection angle. Displays on the Gantry show the current angle of the C-arm.
3. Move the patient up to the C-arm.
4. Position the patient as required.
5. Put the arm or hand of the patient on the Patient Handle or at the side of the patient. Inform the patient of the locations of all push button controls and give instructions to avoid contact with any controls.



Warning:

Be aware of the position of the patient's hands and arms relative to any of the C-arm controls.



6. Compress the breast.
 - Use the footswitch controls to provide hands-free compression control and C-arm height adjustment.
 - Use the light field lamp as necessary to visualize the X-ray field. Position the patient's breast within the illuminated field to ensure total coverage.
 - Apply compression slowly. Use the handwheels to make final compression adjustments.

12.4 Acquire the Image and Complete the Procedure

1. Confirm X-ray equipment parameters and Image Acquisition settings.
 - For full implant views, either AEC (Auto Filter) or Manual Exposure Modes may be used. Refer to Chapter 2, Section 4.1, page 4 and Section 4.2, page 5.
 - For implant displaced views use AEC and manual placement of the AEC sensor. Refer to Chapter 2, Section 4.4, page 6.
2. The system is Ready for exposure when:
 - The screen has no "X" across the radiation symbol.
 - The generator status says **Ready**.

If the system does not enter a **Ready** status in 60 seconds, verify that the accessories are correctly installed and the paddle is firmly locked into place.



Warning:

Observe all safety precautions when making X-ray exposures.

3. Press and hold the X-ray buttons for the entire exposure.
The Radiation icon on the screen changes from gray to white during this time and the X-ray Indicators by the buttons on the keypads glow yellow.

4. Release the patient.
 - If the compression release is set to manual or you are using a localization paddle, raise the compression device.
 - If the automatic release feature is set, the compression device automatically rises after the exposure.



Warning:

Do not move the C-arm while the system is retrieving the image.

5. Preview the image. Observe the Exposure Index to ensure the exposure is correct.
6. If required, use the Window/Level tool or other Preview options including Image Comments—all available at any time while the Preview screen is active.
7. Accept or Reject the image.
8. Perform the Acquisition cycle as many times as required for the requested Procedure(s).
9. If necessary, add an additional view or procedure.
10. Ensure the patient has safely moved away from the unit when you have taken all the Views.
11. Close the Procedure.

13.0 Working with Specimens

Follow standard procedure with the following considerations:

- Use a Magnification Stand if desired.
- Use a Spot Paddle.
- Bring the specimen as close as possible to the chest wall edge of the detector or stand.
- Use Manual Exposure Mode.

Note that you can use **Admin>Image Management>Resend** to resend to a printer if you need a second image for pathology.



Note...

You can use the manual collimation button on the tubehead to expose the specimen only. This does not reduce the size of the image, but does mask out all unnecessary space. The minimum cropped size of a Selenia image is 18 x 24 cm. You cannot further reduce the image size for specimens.

14.0 Performing a Localization Procedure



Warning:

The system is normally configured to disable Automatic Compression Release when a Localization Paddle is installed. If you use a paddle other than a Localization Paddle during this procedure, you must disable Automatic Compression Release. On the Main Menu, click Edit>Standard Setup and change the Compression Release option from Automatic to Manual. Return this option to Automatic when the procedure is complete.

14.1 How to Use the Standard Crosshair Device

The Localization procedure has additional steps not needed in other clinical procedures.

1. Follow the standard screening procedures to step 3, "Install the compression paddle for the procedure and any other necessary accessories," page 65.
2. Install the Crosshair Device.
3. Change the Compression Release setting to **Manual** if set on Automatic.
 - a. Select **Edit>Standard Setup**.
 - b. Select **Manual** in the Compression Release field. Wait a few seconds before closing the window for the change to take effect.
 - c. Click **Close**.




Figure 3-47: Rotate the Crosshair Device Out of the X-ray Field

1. Rotate the Crosshair Device out of the X-ray field before taking the scout image.



Note...

Auto AEC is not available with localization paddles. Choose AEC position 2 or 3.

2. Take the scout image.
3. Click **Zoom/Pan** in the Preview window to see a magnification using the entire screen.
4. **Pan** the image by dragging with the right trackball button.
5. Click **Crosshair** and then click on the area of interest to display the crosshair lines.
6. Click **Zoom/Pan** button off to see where the crosshairs intersect with the alpha-numerics on the screen. This determines the horizontal and vertical position of the lesion.
7. Move the crosshair device into the field.
8. Turn on the light field lamp.
9. Adjust the crosshair knobs to place the crosshair shadows at the appropriate reference labels on the localization paddle. The intersection of the crosshair shadows on the skin indicates the position of the lesion on the scout image.
10. Note the position on the breast.
11. Move the crosshairs out of the X-ray field.
12. Continue with the Localization Procedure, changing paddles as required.

13. Upon completion of the procedure, change the Compression Release back to automatic if this is the site default.
 - a. Select **Edit>Standard Setup**.
 - b. Select **Automatic** in the **Compression Release** field. Wait a few seconds before closing the window for the change to take effect.
 - c. Click **Close**.
14. Use Standard Screening steps for the Procedure. Refer to Section 12.0, page 65.

14.2 How to Use the Magnification Crosshair Device

Follow the instructions in Section 14.1, page 68 except:

1. Attach the Magnification stand.
2. Attach the Magnification Crosshair Device to the tubehead.
3. Check to make sure the Grid Out and Small Focal Spot are selected.
4. Dial the Crosshairs out of the way before you take the scout shot.
5. Remove the crosshair device after the needle placement when you are ready to continue with the procedure.

15.0 Accepting and Storing the Image

Use the Preview screen to:

1. Inspect for good positioning and proper exposure level.
2. **Accept** or **Reject** the image.
 - If it is acceptable, select **Accept** and wait for the image to appear as a thumbnail on the Patient Procedure screen.
 - If the image has problems (for example, under-exposure, movement) enter **Reject** and the reason it is unacceptable. Repeat the exposure.
 - Note the displayed actual exposure value. If a minimum exposure (250) is not achieved, the image should be rejected.

Accepting an image transfers it to the selected Output device(s).

Rejecting an image stores the image in the **Reject Bin**, and does not send it to the Output destination. All images, Rejected as well as Accepted, appear on the Procedure's **Review** of images. All **Rejected Images** are stored in a log with their reasons for rejection.

For additional information on accepting and storing images, refer to "Accept Button Details," page 46.

16.0 Optional Workflow Operations

In addition to the workflows described above, optional operations are available to:

- Open multiple Procedures within a Patient. Refer to Section 10.3.2, page 57:
- Search by either Patient Name or ID to Open a new Procedure for an existing patient. Refer to Section 10.3.3, page 58.
- Search by either Patient Name or ID to Recall an open Procedure to continue that Procedure. Refer to Section 5.0, page 33.
- Resend previously acquired Images. Refer to Section 19.1, page 82.

17.0 Main Menu Functions

Main Menu functions provide the technologist with easy access to non-routine Acquisition Workstation functions. They are always available while the Acquisition Workstation application is running unless a dialog box or Procedure is open.

This section describes additional details about the menus in the order they appear on the screen. If a menu item is not listed, there are no additional hints or descriptions about the item available beyond the information in *Instructions for Use*.

17.1 Main Menu with No Selected Patient

The table below lists the Main Menu functions for the technologist when a Patient Procedure is *not* selected. Once a Procedure is open, the Main Menu functions are not available.

Table 3-3: Main Menu Functions with No Procedure Chosen

Menu Function	Tech
File	Exit
Edit	Standard Setup
	View Order Editor
	User Setup
	Outputs
Admin	Image Management
	Protect Patients
	Manage Queues
	Eject
	Import
	Retrieve Priors
	PPS Status*
	Calibrate
	Test Patterns
	DR Device Control
Available Disk Space	
Info	About the Acquisition Workstation
*This only appears when the MPPS Service is installed.	

17.2 File Menu>Exit

The **Exit** function provides options for exiting the system if a Patient/Procedure is not displayed. If a Patient/Procedure is open, a dialog box informs you that you cannot Exit while a Patient/Procedure is open. Refer to Section 1.3, page 20.

18.0 Edit Menu



Figure 3-48: Edit Menu

18.1 Standard Setup Screens

On the Main Menu, click **Edit>Standard Setup** to display a list of startup defaults. Table 3-4 lists available options. Hints appear in the lower pane as the cursor moves over the different choices.

Table 3-4: Setup Choices

Choice	Range of choices
Pre-compression Force	10 - 100%
Full Compression Force	10 - 100%
Compression Force Units	Pounds (lb) or Newtons (N)
Compression Release Mode	Automatic or Manual
Compression Force Mode	Pre, Full, or Dual
Collimator	Automatic or Manual
Default Exposure Mode	Auto-Time, Auto-kV, Auto-Filter, TEC, or Manual
Default Magnification Mode	Auto-Time, Auto-kV, TEC, or Manual
Set Technique Defaults	Refer to "Set Technique Defaults Screen"

18.2 Set Technique Defaults Screen

The last button on the Standard Setup screen, **Set Technique Defaults**, displays the Technique Default Setup screen. Select *Manual* Exposure Mode to start. Remember to set the Magnification modes. Each possible default must be set up.

Table 3-5: Set Technique Default Choices

Mode	Choices
<i>Manual</i> (and Mag)	Filter, kVp, mAs, mAs Table, Focal Spot (Large or Small)
<i>Auto-Time</i> (and Mag)	Filter, kVp, mAs Table, Focal Spot, AEC Region position, and Density
<i>Auto-kV</i> (and Mag)	Filter not selectable (default is Mo for Moly Tube, Rh for Tungsten), mAs Table, Focal Spot, AEC Region position, and Density setting
<i>Auto-Filter</i> (and Mag)	mAs Table, AEC Region position and Density setting
<i>TEC</i> (and Mag)	Filter, mAs Table, Focal Spot, and Breast Density (Fatty, Normal or Dense)

18.3 View Order Editor




Figure 3-49: View Order Editor

To set up your individual view order defaults, drag the **View icons** to change their order. Select the procedure from the drop down list. You may set view orders for each Procedure. Click **Close**. Your **View Order** preferences are saved and load each time you log in.

18.4 User Setup—Add or Edit a User

Click **Edit>User Setup**. The **Edit Users** dialog box appears.

A user can only be added or deleted at the same or lower level as the user signed in. Example:

- *Mgr* can add or delete a mgr or tech user.
- *Tech* can add or delete only a tech, not a mgr.



Figure 3-50: Edit Users Dialog Box

18.4.1 Add a New User

To add a new user:

1. Click **New**.
 2. Enter the requested information in the fields.
 3. Assign a password. The the password must contain at least 6 characters.
 4. Click **Accept** to add the user, or **Cancel** to leave without adding the user.
- The user name appears on the printed images.
 - The initials appear on the electronic image.



Figure 3-51: Add a New User

18.4.2 Edit Your Profile

To edit your initials or password:

1. Select your name from the list.
2. Click **Edit**. The Edit Profiles for <name> appears.
3. Change the Initials if necessary.
4. Enter your old password.
5. Enter a new password. The the password must contain at least 6 characters.
6. Reenter the new password to confirm it.
7. Click **Accept** to make the change.



Figure 3-52: Edit User Profile



Note...

You have to enter the password to change the initials.

18.4.3 Delete a User

To delete a user:

1. Select a User in the **Edit Users** dialog box to delete.
2. Click **Delete**.
3. Confirm the deletion in the **Delete?** dialog box by clicking **Yes**, or cancel the deletion by clicking **No**.

18.5 Set Up Outputs

An Output is a pre-defined list of one or more physical output devices, which determines where to send images from the Acquisition Workstation. Outputs are initially created during the installation process. The Output name is used to identify which devices are included in the Output group. Physical devices (for example, Printer archives, diagnostic review workstations, etc.) are supported via Named Virtual Devices. These become the output destinations of the system and can be grouped together or used individually.

The Edit Outputs function adds or removes Outputs. A list of currently available choices appears in the dialog box and you can group destination devices into an output group, add new groups, delete, copy, or edit current ones. The external CD-RW drive and any network connections are setup at installation. They can be combined into additional group outputs.

Outputs are defined using the Edit Outputs Menu. This option enables you to:

- Select output destination(s) and add/subtract them from the output's device list.
- Select devices to be associated with each output destination(s).
- Save the output under a recognizable name.
- Create a new output list using an existing output list.
- Delete an output list.

In the Edit Outputs function, the lists of configured possible output destinations are sorted by device name:

- Archive = PACS and local storage devices
- CDRW = CD-RW
- Physician Display = diagnostic review workstation, technician review workstation, CAD
- Laser Camera = Printer

Other characteristics of output lists are:

- A single named virtual device that appears in multiple output lists has the same configuration in all those lists.
- At least one output list is named "No Outputs" and contains no devices.
- Named virtual devices appear in the Edit Outputs tool's destination list when the service engineer configures the device. An output for the virtual device appears in the output list when the device is added using the Edit Outputs function.



Note...

Only Accepted Images are sent to the selected list of outputs. Rejected Images are sent to their own list of output devices and must be removed from the system after being reviewed.

18.5.1 Add a New Output Group




Figure 3-53: Edit Outputs Dialog Box

1. Select **Edit > Outputs**, and click **New** to add a new Output (Figure 3-53).




Figure 3-54: New Output Name Dialog Box

2. Enter the name you want to use for this Output (for example, LaserPrinter2) in the *Output Name* dialog box (Figure 3-54).
3. Click **OK**. The new name shows up in the *Edit Outputs* dialog box.




Figure 3-55: New Name in Edit Outputs Dialog Box

4. Select the new name (for example, LaserPrinter2) in the *Edit Outputs* dialog box, and click **Edit** (Figure 3-55).

The *Edit Output: <device name>* dialog box appears, with tabs for the four types of outputs. The Device names currently selected are listed in the right column.




Figure 3-56: Adding Devices

5. Click the tab of the device type you need to add. The driver(s) for that type of output are listed.
6. Select the device you want to use (Figure 3-56).
7. Click **Add** in the middle of the box. The device name appears in the right column, *Devices currently selected*.
8. If you want to add an additional output, click the device tab for it.
9. Click the additional device you want to add.
10. Click **Add** to move it to the right column. You can add many devices to a group.




Figure 3-57: Configuring the Output Device(s)

11. Click **Accept** to close the *Edit Output: <device name>* dialog box. The new group name shows up on the *Edit Outputs* dialog box.
12. Click **OK** to close the *Edit Outputs* dialog box.

18.5.2 Edit Outputs and Make New Groups

To recombine devices or delete a device from a group:

1. Select **Edit > Outputs**.
2. Click an Output name(s) that you want to change or recombine.
3. Click **Copy**. A dialog box appears (Fig 2-9).
4. Enter the new Output name you want to use for this Output (for example, LaserPrinter3).
5. Click **OK**. The new name shows up in the *Edit Outputs* dialog box.
6. Select the new name, in this case LaserPrinter3, in the *Edit Outputs* dialog box.
7. Click **Edit** (Fig 2-10). The *Edit Output: <device name>* dialog box appears, with tabs for the four types of outputs. Devices currently part of the group are listed in the right column.
8. Click the tab of the device type to add, for example, laser_cameras. A list of the driver(s) for that type of device appears.
9. Select the device you want to use (Figure 3-56).
10. Click the **Add** button in the middle of the box. The device name appears in the right column, *Devices Currently Selected*.
11. If you want to add an additional device, click the tab for it (for example, **cdrw** for a CD-RW drive).
12. Click the additional device you want to add.
13. Click **Add** to move it to the right column.
14. Click the name of any output in the *Devices Currently selected* column that you do not want in the group and click **Remove Item**.
15. Click **Accept** to close the *Edit Output: <device name>* dialog box.
16. Click **OK** to close the *Edit Outputs* dialog box.

18.5.3 Delete Outputs

To delete an output group name that is no longer needed:

1. Select **Edit > Outputs**.
2. Click the output name(s) that you want to delete.
3. Click **Delete**. The *Delete output?* dialog box appears.
4. Click **Yes** in the *Delete Output?* dialog box. The dialog box closes and the output name is removed from the *Currently available outputs* list.

18.5.4 Output Queues

Actual transmission to a given output device is handled through output queues associated with the device. The following conditions apply to output queues:

- Output queue states are held through shut down cycles.
- Output queues are automatically restarted in the same state as at shutdown.

An Admin Menu function (Section 19.3, page 87) allows you to:

- Cancel jobs from queues (not including active jobs).
- View job status within queues.
- Resend “stopped” jobs. Note that images for problem jobs are not reclaimed.

18.5.5 Using the New Output in the Patient View Window

1. Click the **Outputs** down arrow.
2. Select the new name (for example, LaserPrinter2). A *Setup output devices* dialog box appears informing you that you are changing to the output <name> and listing the virtual device name(s) in the output.
3. Click **OK** to close the *Setup output devices* dialog box.




Figure 3-58: Choose Output

4. Click **OK** to close the *Setup Warning* dialog box.




Figure 3-59: No Output Selected

The “no outputs” option is used by service when they do not need an output. If you choose this option, none of the images you acquire are sent to any other devices.

18.6 Data Storage (Using Outputs)

The Selenia data storage system is comprised of the following subsystem components:

- Acquisition Workstation hard drive (automatic save)—for acquisition, temporary storage, and QA display of the breast images
- CD-RW—for temporary storage of images on removable storage media

The Selenia System can interface with these components:

- Diagnostic Review Workstation—for display of diagnostic quality images
- Film Printer—for printing on high-resolution films
- Network Storage—for storing data, archiving images

The Selenia System, when connected to a network, can store data on many systems and can output to many others. These integrated systems are user selectable and location specific.

If you do not have a connection to a PACS system, you must send the patient information to an approved film printer.

18.6.1 Hard Drive

The Acquisition Workstation computer's hard drive automatically stores the Acquired Images. Select the Main Menu item **Admin>Available Disk Space** to display the available storage space remaining, in units of images, on the system. The capacity of the hard drive, about 50 Images per Gb, depends on the size of the images. A 24 x 29 cm image uses more storage space than an 18 x 24 cm image, so the number of available images is an estimate.

18.6.2 Automatic Reclamation of Hard Drive Space

The system automatically reclaims Hard disk storage space for new Procedures. Reclaiming is a process of deleting existing patient files from the Acquisition Workstation hard drive in order to make room for new patient information. The image reclaimer deletes images from the system when a number of images have accumulated on the system. Once reclamation of images begins, images continue to be reclaimed until the number of images on the system reaches a set number.

Some images are not deleted:

- No Patients are deleted if they have any Images queued to an output device.
- No Images are deleted for Protected Patients.
- No Images queued up on an output device are deleted.
- No Images that are not successfully committed to an archive device are deleted.
- Rejected Images are not automatically deleted.



Caution:

The system does not acquire a new image if the disk does not have the necessary space for a new image. The space must be reclaimed by deleting patients whose images are no longer needed on the hard drive.

This message should appear only after messages alerting you that there are too many protected patients or too many images backed up on an output device, and so may be due to a configuration issue that can be corrected by service.

18.6.3 CD-RW Drive

A Compact Disc ReWritable (CD-RW) drive, mounted on the underside of the Acquisition Workstation, is available for temporary storage of images. If you do not have a PACS network connection, you must select an approved film printer as one of the outputs for patient images.

1. To use the CD-RW drive as an output, in the **Output** box of the **Patient View** window, select an output that includes the CD-RW. A **Setup Output Devices** dialog box appears, informing you that you have selected a new setup <name> and giving you the virtual device name for the CD-RW.
2. Open the drive drawer by pressing the eject button on the drive.
3. Write the patient name, ID # and date on the envelope and/or CD. Insert the approved CD-R disk in the drive prior to starting a procedure.




Figure 3-60: Invalid Media Alert

4. Select the CD-RW as an output, or part of an output group. Acquired images are now sent to the CD as well as to any other output in the group (if any).
5. After the last view is accepted, it takes a few minutes to finish writing the images to the CD (longer for 24 x 30 cm images).
 - Observe the **Jobs in Queue** message in the Outputs pane to tell when all the images are done.
 - You can use **Admin>Manage Queues** to see what is still in the Queues.
6. After the last view is sent, use **Admin>Eject** to eject the CD.

Writing images to a CD takes longer than saving or printing them. You can use the **Image Management>Resend** option and write them after the procedure is finished. This has the advantage of not slowing down the rest of the procedure while you are waiting for the program to write to the CD during Image Acquisition.

To see which patients/views are on a CD, use **Admin>Import**.

The CD-RW drive is also available for importing images from other Selenias recorded on a CD. Use the **Admin>Import** directions for Importing Images.

18.6.4 Change the Disk in the CD-RW Drive

If the disk is near capacity, the message below appears.



Figure 3-61: CD-R Alert

To change the CD:

1. Click **Eject** on the alert message dialog box and the CD-RW drive drawer opens.
2. Remove the full CD and insert a labeled blank one.
3. Close the drawer. Label the full CD and store it per local requirements.
4. A Continue button appears. Click **Continue** to save the remaining Images.

If you click **Cancel** and do not change the CD, the system generates an **Alarm** and does not archive the current Image on the CD. It is saved on the hard drive and you can **Resend** it from the **Image Management** function.

18.6.5 Film Printer

The system can send the Image to an optional high-resolution DICOM film printer. Details on the function of the printer depend on the printer. The printer configuration (as with any new printer) is setup by an authorized service engineer. You must add it to the Outputs list.

18.6.6 Network Storage

The system supports retrieval of the archived Images and Procedures. A separate diagnostic review workstation can retrieve the Images. Details on the function of these outputs depend on your specific configuration.

19.0 Admin Menu




Figure 3-62: Admin Menu

19.1 Image/Spool Management

This Admin Menu item allows you to:

- Query for Patients whose Images are still within the Image spool.
- Resend selected images in the spool to Output Destination devices. (The system resends Images to the selected output exactly as found in the spool.)
- Repreview an Image; add comments, etc. (which generates a new image file).
- Confirm DICOM storage commitment.

Query details:

- Entering an incomplete string returns all matches to the string.
- Queries are not case sensitive.
- “No entry” or “*” is equivalent to everyone in the Image spool.
- If too many matches are found, the application requires more specific search criteria.
- If no matches are found, the message, “No Patients matching query” is displayed.
- When there are Images for the queried Patient(s) the Acquisition Workstation displays Patient information in a table.

When you have selected a patient to manage, an Image Management screen displays. The upper right pane displays detailed information about the view when you have selected a single view. This information includes the view name, acquisition capture timestamp, technologist, comments, and other data. The lower pane displays thumbnails of all the images of the selected Patient in the order of acquisition.

19.1.1 Using the Resend Options



Figure 3-63: Image Management Resend Options

1. Select a thumbnail image to resend.
2. The output selected when the image was acquired appears as the default. Confirm the Output is correct. If not, select the Output you want from the drop-down list. The Output change you make here does not affect the Output selection on the main screen.
3. If the selected Output has a printer, the following **Print Image Processing Options** appear. Click the option to select or deselect it. See Figure 3-63.
 - Adjust the **Contrast** and **Density** for printing. Change by clicking another number to increase or decrease.
 - **Print True Size** (if available at the printer).
 - **Apply PCE** (the Peripheral Contrast Enhancement Function uses a special algorithm to enhance the appearance of the mammograms).
 - Omit Patient Name.
 - Omit Label options.
4. Click **Resend**. A dialog box states the number of images resent.




Figure 3-64: Resend Confirmation



Reminder...

You cannot change any Patient or image object information from the resend function. This includes image, LUT modifications and DICOM information such as Patient name spelling, ID, comments, etc. The only way to send new information is to make a new file using the Repreview function, save it, and Resend this new file.

19.1.2 Using the Repreview Option




Figure 3-65: Repreview Screen

1. Select one Image thumbnail that you want to review. Use the arrows to scroll when more than four images are available (Figure 3-63).
2. Click the **Repreview** tab, then the **Repreview** button.
3. The Repreview screen appears with the same tools as you had in the Preview screen. If you only want to review the Image, click **Cancel** when you are done.

If you select a FOR PRESENTATION image, you cannot change the markers and ROI is not available.

4. If you want to create a new file to resend, make any required changes to the markers, comments, or use the **Edit View** button and then click **Save**.
5. Make any necessary changes to the output selection. When the Send Image to Output dialog box displays, click **Yes** to send the images to the selected output group.
6. If you are reprinting, make any necessary changes in the **Contrast** or **Density**.
7. Click **Save**. The **Repreview** screen closes and the **Image Management** screen appears again. You have created a new copy of the original file.



Note...

Be careful when resending that you send the newly saved image. The original is still there as a possible selection.

19.1.3 Correcting Incorrect View Selections with an Edit Button

One use for the Repreview function is to correct an incorrect view exposure if the error was not corrected before clicking Accept.

1. In the Repreview screen, select the **Edit View** button.
2. Select the correct view.
3. Select **OK**. The incorrect view is replaced with a corrected one.
4. Click **Save** to save the new version with corrected markers and orientation.
5. Select the Resend tab. Deselect the wrong image and select the last new thumbnail.
6. Select the correct output and select the Resend button. Click OK in the confirmation box and click Close to leave the Image Management function.
7. Delete the incorrect image sent to the chosen outputs.
 - If the image was printed, collect and discard the incorrect image.
 - If the image was sent to a diagnostic review workstation, log on and delete the incorrect image.
 - If the image was sent to the PACS, contact the PACS administrator and have them delete the incorrectly marked study.
8. A manager can delete the incorrect image. Refer to Chapter 5, Section 4.1.8, page 114.

If your Preview Screen does not have an Edit button and you have accepted the image, refer to Section 9.3.2, page 51.

19.1.4 Other Options on the Image Management Screen

- Query for another Patient by clicking **Choose a New Patient**.
- Close the **Image Management** screen without taking any actions.

19.2 Protect Patients

This menu item protects or unprotects Patient's records on the Acquisition Workstation hard drive. The system cannot automatically reclaim Protected records.

1. Select **Admin>Protect Patients**.
2. Search for the Patient you want to protect and select the Patient from the search results.




Figure 3-66: Protecting Patients

3. Click **Protect** at the bottom of the screen. A dialog box appears prompting you to confirm you want to protect the Patient. Select **Yes**.




Figure 3-67: Protect Confirmation

To Unprotect a Patient:

1. When you no longer need to protect a Patient, select **Admin>Protect Patients**.
2. Select the name and click **Unprotect**.
3. Repeat for each patient who no longer needs to be protected.
4. Click **Close** to close the dialog box.



Note...

Be sure to **Unprotect Patients** when you no longer need them to avoid filling up the Acquisition Workstation hard drive.

19.3 Manage Queues

Click **Admin>Manage Queues** to find problem jobs or locate a particular job in the queues.

Actual transmission to any output is handled through Output Queues. Each connected output device is associated with an Output Queue. If a physical device does not accept a job or if the job does not complete successfully, the job may be retried or held as a Problem Job, depending on the reason for the error.

A service configurable number of retries are performed before an alarm is generated. A service configurable number of retransmission attempts occurs after the alarm is generated. Once the alarm condition clears, the alarm must be erased to be cleared.

Jobs may fail because there is either a problem with a specific characteristic of the job (for example, configured LUT not available on the physical printer) or because there is a general problem with the physical device (for example, device is off-line). In order to keep Problem Jobs from clogging an otherwise operational queue with retries, the output queues have two levels of priority. The levels are managed automatically and there is no direct user interaction.

You can view the states of the output queue:

- Waiting—job has not made it to the top of the queue yet.
- Running—job is active.
- Sleeping–Error—job has an unresolved alarm that has continued for a configurable number of retries (see below).
- Stopped–Error—job has a non-recoverable alarm and is not being retried.
- Output queue states are held through shutdown cycles.
- Output queues are automatically restarted in the same state they were in at shutdown.



Note... *Images for problem jobs are not reclaimed until the job is removed from the queue.*

19.3.1 Job Retry Options

The Acquisition Workstation supports the following job retry actions. The parameters used in the action processing are service configurable, and the factory defaults are listed below.

- Number of failed retries before an alarm is posted (default = 3).
- Initial delay time between retries that occur after the alarm is posted (default = 90 sec.).
- Number of failed retries using the initial delay time (default = 10).
- Final delay time between retries (default = 900 seconds) Note: This feature helps to prevent failed jobs from “clogging” the queue.
- Total number of retries attempted after an alarm has been posted until the job is marked “Stopped,” and retries are discontinued.

19.3.2 Available Manage Queues Functions

You have the following functions (refer to the numbers to the left of Figure 3-68):

1. *A list of available queues.* Different spoolers may have been configured to handle different communication services.




Figure 3-68: Manage Queues Suspend Job Spooler button

2. Two buttons under the list of available queues
 - A *Suspend Job Spooler* button that suspends images from going to the selected queue. This changes the current state from Enabled to Suspended.
 - A *Refresh* button that refreshes the queued jobs list for the selected device.
3. A *scrolling list of queued jobs*. This list identifies all of the jobs in the queue by:
 - Patient Name
 - Patient ID
 - Study
 - Image
 - Job Status: Running, Waiting, Problem, Sleeping
 - Retry Count
 - Device
4. Selecting an inactive job from the list, activates buttons that allow you to:
 - Delete a job from the queue.
 - Resend a job, if the job is in the stopped status.

You cannot delete currently Active jobs from the queue, and the “Delete Job” button is not active if you select such a job. If a job is stopped, the **Resend Stopped Job** button becomes active.

You may delete all jobs the system is not currently sending from the queue regardless of who originally queued it.

19.3.3 Using Manage Queues

Click **Admin>Manage Queues** to check for problem jobs or locate a particular image in the queues.

19.3.4 Determine the Status of a Job

- Click the queue you want to view in the Devices list.
 - The jobs in that queue are listed in the **Queued Jobs** list.
 - Their status is in the **Job Status** column.
 - A status of “sleeping” means a job has been retried several times in a row and is currently suspended for an hour or so before retrying.
- Click **OK** when you are finished.

Manage Queues						
Devices						
Name of Device	Description				Status	
DicomPrintSpooler	single priority, single server, Dicom print job spooler				ENABLED	
DicomStoreSpooler	single priority, single server, Dicom store job spooler				ENABLED	
DicomCommitSpooler	single priority, single server, Dicom commit job spooler				ENABLED	
CdrwSpooler	single priority, single server, CDRW job spooler				ENABLED	
<input type="button" value="Suspend CdrwSpooler"/> <input type="button" value="Refresh"/>						
Queued Jobs						
Patient Name	Patient ID	Study	Image	Job Status	Retry Count	Device
Patient^One	00581381	Screening-Bilateral Mammography	R CC	RUNNING	0	/dev/rdisk/c1t4d0s2
Patient^One	00581381	Screening-Bilateral Mammography	L CC	WAITING	0	/dev/rdisk/c1t4d0s2
Patient^Two	00581382	Screening-Bilateral Mammography	R CC	WAITING	0	/dev/rdisk/c1t4d0s2
<input type="button" value="Delete From Queue"/> <input type="button" value="Resend Stopped Job"/> <input type="button" value="OK"/>						

Figure 3-69: Delete Selected Jobs in a Queue



Warning:

If you click Delete Job, you cannot undo the deletion. If it still needs to be sent, resend the job to a storage device later.

19.3.5 Delete (Canceling) a Job

- Click the queue you want to manage in the Devices list.
- Click the **Suspend <name> Spooler** button.
- Select a job from the **Queued Jobs** list, or Ctrl+Click to select multiple jobs. Verify the job(s) you want to delete are selected correctly. There is no “undo”.
- Click the **Delete From Queue**—below the list.
- Click **OK** when you are finished.

19.3.6 Resend a Stopped Job

- Click the queue you want to manage in the **Devices** list.
- Select a job or jobs labeled ‘Stopped’ from the **Queued Jobs** list.
- Verify the job(s) are selected correctly.
- Click the **Resend Stopped Job** button—below the list.
- Click **OK** when you are finished.

19.4 Eject



Note...

The drive drawer only opens from the menu when there is a disk in the drive.

Click **Admin>Eject** to open the CD-RW drive drawer when there is a disk in the drive. If there is no media in the CD-RW drive, clicking Eject in the menu results in a CD-R alert window informing you there is no media in the drive. Press the button on the CD-RW drive to open the drawer. To close the drive, in either case, press the CD-RW drive button.

19.5 Import



Figure 3-70: The Import Dialog Box



Note...

*Do Not press the **Eject** button on the drive while the system reads from or writes to the disk. It may cause the drive to stall.*

Imported images can be reviewed and sent to outputs including store-committing them.

1. With the CD containing the images to import in the drive, select **Admin>Import** to import a DICOM Image file from a CD. A dialog box appears and the system reads the CD in the drive. The Patient's name displays in the list.
2. In the dialog box, find the file(s) you wish to import.
 - Clicking on a "+" opens a list of folders and/or files below.
 - Clicking on the empty box selects the file or folder.

- Click the file name and a check appears in front of all files you have selected. Alternately, click a folder to see the list of files in the folder. If you click the box in front of the folder, you select all the files in the folder. All checked files are imported.



Hint...

If you have selected a file that you do not want to import, click the check mark again to clear it. Clicking on a folder selects all subfolders below it. Re-clicking on the folder unselects the folder, but does not unselect all the files below it. You have to clear each one.

3. Click **Import**. A dialog box indicates the progress of the import.
4. Click **Cancel** to close the dialog box without importing anything.




Figure 3-71: Import Progress Dialog Box

5. The **Import Progress** dialog box indicates the files imported and whether the import is successful or not.

The imported Patient data is available in the **Select a Patient** screen.

19.6 Retrieve Priors

All users have the ability to view prior Images on the Acquisition Workstation if the Images are in the Local Database or stored on the PACS system. Locate Images in the Local Database by clicking on the Local Exams tab and searching for the Patient.

Images on the PACS require an additional step. The Selenia system queries and retrieves the prior Images for a Patient from the PACS system. The program asks you to choose the source if more than one source has been configured.

The system imports Images into the Acquisition Workstation Local Database. Images already on the Acquisition Workstation generate an alarm if you try to overwrite them with the PACS version.

The rules for handling both FOR Processing and FOR Presentation retrieves of the same image are:

1. If the image being retrieved is a for processing image and is already present in its for processing form in the database, you are alerted to the fact that the image is already in the database, and therefore, the import of that image did not occur.
2. If the image being retrieved is a for processing image and the for presentation version is already in the database, you are alerted to the fact that the image is already in the database, and therefore, the import of that image did not occur.
3. If the image being retrieved is a for presentation image and it already exists as a for presentation image in the database, you are alerted to the fact that the image is already in the database, and therefore, the import of that image did not occur.
4. If the image being retrieved is a for presentation image and the for processing version is already in the database, then the data for the for processing version should be discarded and replaced with the for presentation version. The log file notes this replacement activity.
5. After the import is completed, the Patient can be opened using the Completed Exams tab and located like any Patient captured on the Acquisition Workstation. Any import errors are displayed as Alarms.

To Retrieve Priors from the PACS system:

1. Click **Admin > Retrieve Priors**. The **Retrieve Priors** screen appears
2. Search for the Patient's name in the **Name** field or ID in the **ID** field (wildcards are allowed).



Figure 3-72: Retrieve Priors Screen

3. After searching the PACS system, the program displays a list of all Patients meeting the query criteria (Figure 3-73 #1). The progress bar on the lower left (Figure 3-73 #2) indicates the system is busy searching.
4. Select the Patient name(s) you want to import. Selecting a Patient activates the **Import** button.
 - You can choose several Patients at once.
 - Clicking a selected Patient deselects them.
 - The title bars allow sorting by field.
5. Click **Retrieve and Import** (Figure 3-73 #3). The PACS system locates the Images and imports them into the Acquisition Workstation Local Database. A dialog box appears and the progress bar indicates the status of the import.
6. Click **OK** to close the dialog box.
7. Click **Close** to close the Retrieve Priors screen and return to the standard acquisition screens while Import operates in the background.




Figure 3-73: Results of the Retrieve Priors Query

Legend for Figure 3-73

1. List of Patients
2. Import Progress Bar
3. Retrieve and Import Button



Figure 3-74: Image Import Started



Figure 3-75: Query Retrieve Import Successful Dialog Box

After the import is completed, you can display the Patient Images using the Local Exams tab in the same method as any Patient Image acquired on the Acquisition Workstation. Alarms display any import errors.

If your query request resulted in more than 50 records, this dialog box appears. Refine your search.




Figure 3-76: Search too Broad Dialog Box

If there are errors in the results of the query, this error dialog box appears.



Figure 3-77: Query Retrieve Errors dialog Box

19.7 MPPS Status

An optional feature available for sites with an MPPS Service Class Provider. When this feature is installed, workflow changes with three close procedure button options instead of one. Adding a procedure and starting a new procedure also have slightly different functionalities.

19.8 Test Patterns




Figure 3-78: Test Patterns Dialog Box

Selecting **Admin>Test Patterns** displays a **Test Patterns** dialog box. Use Test Patterns for the Laser Printer Test (make different selections for other uses). Refer to the QC manual for details about the test.




Figure 3-79: Image Queued Dialog Box

When you send the test pattern to the printer you get an Image Queues dialog box. To see where the image is in the queue, refer to Section 19.3, page 87.

The **Clean Up** button is for use if the system is having problems with disk space or reclamation warnings. Click on it to clean up any files from previous “test pattern sends” before sending the current image.

19.9 Available Disk Space

The available storage space remaining is indicated in units of images. This information is displayed on the window by selecting **Admin>Available Disk Space**. The capacity of the hard drive, about 50 images per Gb, is dependent on the size of the images. A 24 x 29 cm image uses more storage space than an 18 x 24 cm image.


20.0 Info Menu

The Info menu identifies the Acquisition Workstation software version. The screen displays current temperature of the detector, and the temperature, date and time of the last calibration. It also provides the maximum temperature the detector has reached and date and time this occurred.


Chapter 4—System Hardware Details

1.0 System Power Controls

1.1 Acquisition Workstation Circuit Breaker

Location of Circuit Breaker	Detailed Description
 <p data-bbox="467 890 964 953"><i>Figure 4-1: Acquisition Workstation Circuit Breaker</i></p>	<p data-bbox="1015 546 1464 730">The Acquisition Workstation Circuit Breaker, Figure 4-1, is found on the back of the system behind the shield. The Value Console Circuit Breaker (not shown) is found on the back of the Console.</p> <p data-bbox="1015 737 1464 921">Do not normally turn the circuit breaker off. Use the circuit breaker only in an emergency, because the UPS power is removed and its battery begins to discharge. This can result in data loss.</p> <p data-bbox="1015 932 1464 1117">If you use the circuit breaker, wait until the UPS does not make a beep to turn the circuit breaker back On. The system automatically restarts, but the UPS might not respond automatically and then you must reset it.</p> <p data-bbox="1015 1127 1390 1190">Refer to Section 1.3, page 98 for directions.</p>

1.2 Gantry Circuit Breaker

Location of Circuit Breaker	Detailed Description
 <p data-bbox="522 1623 915 1654"><i>Figure 4-2: Gantry Circuit Breaker</i></p>	<p data-bbox="1015 1346 1464 1568">The Gantry circuit breaker, Figure 4-2, is found at the lower right corner of the back panel. It provides overload protection and removes power from the Gantry only. It should be used only to remove power for servicing or in an emergency.</p> <p data-bbox="1015 1579 1321 1610">Up is On and Down is Off.</p> <p data-bbox="1015 1621 1464 1745">When the circuit breaker is used to turn off the Gantry, the system needs to be restarted. If the Acquisition Workstation is off, restart the system normally.</p> <p data-bbox="1015 1755 1464 1879">If a Panel Power dialog box appears on the Acquisition Workstation, follow the procedure in the <i>Instructions for Use</i> to restart the system.</p>

1.3 Uninterruptible Power Supplies (UPS)

When there is no power to the system, the Value Console or the Acquisition Workstation and TechMate automatically begin to shut down. (For example, if the site has a power failure or if a Main power supply circuit breaker is shut off.) When power is lost, an ordered shutdown begins after approximately two minutes. This shutdown must be allowed to complete before the UPS is turned Off. Then both UPS systems chirp repeatedly and can be shut down until the power is restored. Refer to the *Instructions for Use* manual for UPS directions.

If the Acquisition Workstation or the TechMate does not respond to the Power button, its Uninterruptible Power Supply (UPS) may need to be reset.

For example: If power is removed from the Acquisition Workstation (such as power source disconnected or circuit breaker switched off) and then re-applied, the UPS may not automatically respond.

In rare situations the Acquisition Workstation or Value Console computer can become locked up (unresponsive to user input) and must be turned off and restarted. In some cases the UPS can also require a reset.

2.0 Controls and Displays

2.1 Compression Device Controls and Displays

Legend for Figure 4-3

1. Compression Handwheels.
2. AEC Sensor Position controls.
3. AEC Sensor Position Display.
4. Compression Force Display shows the compression force through the range of 10 lb \pm 1 lb to 67.4 lb (44.5 N to 300 N) in 1 lb (4.4 N) increments. Dashes are shown for force less than 10 lb.
5. Compression Thickness Display shows the thickness between 0 and 15 cm above the image receptor in 0.1 cm increments.




Figure 4-3: Compression Device

2.2 Dual Function Footswitches



Warning:

Always position the footswitches so as to eliminate accidental activation by patient, operators, or wheelchairs.

The system incorporates two footswitches for C-arm and Compression Up/Down movement.

Locate the footswitches for easy access from either side of the C-arm.

Legend for Figure 4-4

1. C-arm Down
2. C-arm Up
3. Compression Down
4. Compression Up















Figure 4-4: The Dual Function Footswitch

2.3 C-Arm Controls

Function	Detailed Description
<p>Counterclockwise C-arm Rotation</p> 	<p>When you press a C-arm Rotation switch, rotation commences immediately and ceases when you release the switch.</p> <p>Press this button to begin motorized counterclockwise C-arm rotation. C-arm movement is disabled when compression force of 58 N (13 lb) or greater is applied.</p> <p>The C-arm Angle Rotation Displays on the Gantry reflect the changes in rotation. Rotating the C-arm counterclockwise changes the display to a more negative reading.</p> <p>After making an oblique exposure, the C-arm angle is stored in memory. When you rotate the C-arm past 0 degrees for the next exposure, the system automatically stops the C-arm at the opposite oblique angle.</p>
<p>Clockwise C-arm Rotation</p> 	<p>When you press a Gantry Rotation switch, rotation commences after a short delay and stops after a pre-programmed amount of travel. Service configures the rotation detent position for the Gantry switches from 10-100 degrees of rotation; the default is 45 degrees.</p> <ul style="list-style-type: none"> • Press and hold a Gantry Rotation switch; audible beeps indicate the C-arm is about to move, then movement starts. Movement stops when the pre-programmed rotation angle is reached • Press the center button and the C-arm moves to the zero degree position.
<p>Compression Down</p> 	<p>Pressing a Compression Down button:</p> <ul style="list-style-type: none"> • Engages the compression brake. • Turns on the light field lamp. • Lowers the compression device. <p><i>Note... Once activated, the compression brake remains engaged until compression release is pressed.</i></p> <p>Compression down movement stops:</p> <ul style="list-style-type: none"> • Upon release of the button. • At the compression down force limit. • At the lower compression travel limit.

Function	Detailed Description
<p>Compression Up</p> 	<p>Pressing a Compression Up button:</p> <ul style="list-style-type: none"> • Raises the compression device. • <i>Does not</i> release or activate the compression brake. <p>Motorized compression up movement automatically stops:</p> <ul style="list-style-type: none"> • Upon release of the button. • Upon reaching the upper compression travel limit.
<p>Compression Release</p> 	<p>Pressing the Compression Release button:</p> <ul style="list-style-type: none"> • Releases the compression motor brake. • Raises the compression device approximately 10 cm. <p><i>Note...</i> The system disables all compression release functions when a localization paddle is installed (including automatic compression release, if enabled.) Patient release is available through the compression up functions.</p>
<p>C-arm Up</p> 	<p>Pressing this button moves the C-Arm up. Motorized C-arm movement automatically stops upon release of the button or upon reaching the upper C-arm travel limit. C-arm movement is disabled when compression force of 58 N (13 lb) or greater is applied.</p>

Function	Detailed Description
<p>C-arm Down</p> 	<p>Pressing this button moves the C-Arm down. Motorized C-arm movement automatically stops upon release of the button or upon reaching the lower C-arm travel limit. C-arm movement is disabled when compression force of 58 N (13 lb) or greater is applied.</p>
<p>Light Field Lamp</p> 	<p>Pressing the light field lamp button turns on the light field lamp for approximately 30 seconds, permitting you to visualize the x-ray field. The light field lamp automatically turns on whenever a Compression Down button is activated.</p>
<p>Collimator Override</p> 	<p>The Collimator Override button allows the collimator to cycle through the available calibrated x-ray fields. Choice of the collimator mode, Automatic or Manual, is made in Edit>Standard Setup. Press the light field lamp button to visualize the x-ray field, and then press the Collimator Override button. The collimator device moves to another field size and the light field size changes.</p>

Function	Detailed Description
<p>Smart Paddle System</p> 	<p>Smart Paddles shift ~2.5 cm into a left, center, or right position. While compression is applied, the paddle cannot be unlocked. The collimator is programmed to follow the position of the Smart Paddles.</p> <ol style="list-style-type: none"> 1. Rotate and hold the knob on the side of the paddle to shift the paddle into one of the detent positions. 2. Move the paddle into the correct position. 3. Release the knob to lock the paddle in place.
<p>AEC Position</p> 	<p>Change the AEC Position by:</p> <ul style="list-style-type: none"> • pressing one of the arrows on the compression device top • using a function key at the Acquisition Workstation or Value Console • selecting an up/down arrow on the display <p>The AEC has an Auto setting. Move the AEC position beyond #7 or #1. The Auto position indicator on the compression device illuminates and the Acquisition Workstation display indicates Auto.</p>

3.0 Operator Radiation Shield

The system's radiation shield is permanently mounted on the Acquisition Workstation or the Value Console.




Figure 4-5: AWS Radiation Shield




Figure 4-6: Value Console Radiation Shield

4.0 Background for the VOS Switch on the Mobile Selenia

Switches located in the Gantry limit the C-arm vertical travel to a pre-defined range. When the C-arm reaches its highest or lowest point of normal travel, the switches engage to stop further motion. If the C-arm should travel past the Limit switches (in the event of switch failure or other electrical malfunction), a second switch engages that removes power from the system. A Field Engineer would normally be contacted to repair and restart the Selenia.

In a mobile environment (during transport), it is possible for the C-arm to 'creep' downward and beyond the two Limit switches. If this occurs, the Selenia would be disabled—the lower limit switch would engage, preventing the system from normal start-up. This is not a failure condition that would require a Service Engineer's attention. The Vertical Position Override Switch is used to allow the operator to start the system and raise the C-arm to within its normal working limits.

Chapter 5—The Manager Functions

1.0 User Identification

The Acquisition Workstation first identifies the user at the login process with these Operating System (OS) and User Accounts:

- *tech* - a single X-ray technologist user account, a group log in for all technologist users.
- *mgr* - a manager user account for the Acquisition Workstation system administrator.
- *phys* - a medical physicist user account with special tools for the physicist.

1.1 User Assignments

When logged into the Acquisition Workstation, each individual user views a consistent set of screens saved by the OS for them. The OS user login determines which access level a user has to the Acquisition Workstation functionality. When the OS user is *tech*, *mgr*, or *phys* an Acquisition Workstation **Login to system** dialog box appears.

The manager configures user login information. The Acquisition Workstation maintains both user identity and passwords in a file. The user logs in by selecting from a list of users. The password is typed in and must contain at least six characters. These designators identify the individual who acquired the images.

A user can be assigned more than one OS access level and thus have their application sign-in appear at more than one level. For example, the manager may configure the following user identifications:

- Chief - both *tech* group and *mgr* user
- King - *mgr* group
- Joe - *tech* group
- Mary - both *tech* group and *mgr* user
- Sally - *tech* group
- John - *phys* group

The technologists identified by “Chief” and “Mary” may log in as either the OS user of *tech* or *mgr*. Their IDs appear in the User Identification list in either case. Depending on which OS login they are under, they may perform either normal tech functions or manager functions. The technologists identified as “Joe” and “Sally” only appear when the OS login is *tech*, and they have only normal tech level access to the Acquisition Workstation. The technologist identified as “King” only appears when the OS login is *mgr*, and “King” only has mgr level access to the Acquisition Workstation. The physicist identified as “John” only appears when the OS login level is *phys*.

2.0 Main Menu Functions

The Main Menu functions are always available except when a dialog box is presented. The table below indicates the state of the Main Menu functions for the *tech* and *mgr* OS users when a Patient/Procedure is **not** selected—*italics* indicates that the button is visible but inactive. Your logon level determines the availability of some Main Menu functions.

Table 5-1: Main Menu Functions by User (Patient/Procedure Not Selected)

Menu Function	Tech	Mgr
File	Exit	Exit
Edit	Standard Setup	Standard Setup
	View Order Editor	View Order Editor
	User Setup	User Setup
	Outputs	Outputs
Admin	Image Management: Resend Repreview	Spool Management: Resend Delete Image Repreview
	Protect Patients	Protect Patients
		Delete Patients
		Reject Management
	Manage Queues	Manage Queues
	Eject	Eject
	Import	Import
	Retrieve Priors	Retrieve Priors
		Log Backup
	Calibrate	Calibrate
	Test Patterns	Test Patterns
	DR Device Control	DR Device Control
	Available Disk Space	Available Disk Space
	Paddle Configuration	
Info	About the Acquisition Workstation	About the Acquisition Workstation
Alarms	Alarms	Alarms
Procedure Buttons	Sign Out	Sign Out
	Search	Search
	Patient (New/Edit/ <i>Close</i>)	Patient (New/Edit/ <i>Close</i>)
	Procedure (New/Edit/Add)	Procedure (New/Edit/Add)
	Outputs	Outputs
	Non-Imaging Mode	Non-Imaging Mode
	Simulate Capture	Simulate Capture
	Review	Review

The ***Close*** button in the Patient box, active only when you have the Choose a Procedure window displayed, enables you to close the Choose a Procedure window without opening any procedures.

These Main Menu functions are available when a procedure is open. **Bold** indicates operational. *Italics* indicates visible but inactive

Table 5-2: Main Menu Functions by User (Patient/Procedure Selected)

Menu Function	Tech	Mgr
File	Exit	Exit
Edit	Standard Setup	Standard Setup
	View Order Editor	View Order Editor
	User Setup	User Setup
	Outputs	Outputs
Admin	Image Management	Spool Management
	Protect Patients	Protect Patients
		Delete Patients
		Reject Management
	Manage Queues	Manage Queues
	Eject	Eject
	Import	Import
	Retrieve Priors	Retrieve Priors
		Log Backup
	Calibrate	Calibrate
	Test Patterns	Test Patterns
	DR Device Control	DR Device Control
	Available Disk Space	Available Disk Space
	Paddle Configuration	
Info	About the Acquisition Workstation	About the Acquisition Workstation
Alarms	Alarms	Alarms
Procedure Buttons	<i>Sign Out</i>	<i>Sign Out</i>
	Close Procedure	Close Procedure
	<i>Patient (New/Edit/Close)</i>	<i>Patient (New/Edit/Close)</i>
	Procedure (New/Edit/Add)	Procedure (New/Edit/Add)
	Outputs	Outputs
	Non-Imaging Mode	Non-Imaging Mode
	Simulate Capture	Simulate Capture
	Review	Review
	Add View	Add View
	Remove Unexposed	Remove Unexposed
Implant Present	Implant Present	

3.0 Edit Menu—Mgr Level




Figure 5-1: Edit Menu

Table 5-3: The Edit Menu for mgr Logon

Menu Item	Where item is discussed
Standard Setup	Chapter 3, Section 18.1, page 71
View Order Editor	Section 3.1, page 108
User Setup	Chapter 3, Section 18.4, page 72
Outputs	Chapter 3, Section 18.5, page 74

3.1 View Order Editor

This Edit Menu item is only available to the Mgr login. It sets the site’s default View Order. This is the order the Views appear if they are not reordered by the technologist.

Each technologist can determine their own View order for performing a procedure. If they have set a View order, their setting determines their View order.




Figure 5-2: View Order Editor

4.0 Admin Menu Functions—“mgr” User Level




Figure 5-3: Admin Menu for the Mgr Logon

The Admin Menu items discussed in this section are those only available to the *mgr* login. The table below indicates where each topic on this menu is discussed.

Table 5-4: The Admin Menu for Mgr Logon

Menu Item	Where discussed
Spool Management	Section 4.1, page 110
Protect Patients	Chapter 3, Section 19.2, page 86
Delete Patients	Section 4.2, page 114
Reject Management	Section 4.4, page 117
Manage Queues	Chapter 3, Section 19.3, page 87
Eject	Chapter 3, Section 19.4, page 90
Import	Chapter 3, Section 19.5, page 90
Retrieve Priors	Chapter 3, Section 19.6, page 92
Log Backup	Section 4.5, page 122
Calibrate	QC Manual
Test Patterns	Chapter 3, Section 19.8, page 95
Dr. Device Control	Service Use Only
Available Disk Space	Instructions for Use
Paddle Configuration	Section 4.6, page 123

4.1 Spool Management

This function changes with the OS logon level. The Spool Management function allows users to manage the image spool space.

- **tech:** resend and repreview images to output devices (labeled “Image Management”).
- **mgr:** repreview, resend captured images to output devices, and delete them.




Figure 5-4: Choose a New Patient Dialog Box

To use Spool Management:

1. Select **Admin>Spool Management**. The **Choose a New Patient** dialog box appears.
2. **Search** for a Patient whose images you want. The query works the same as elsewhere.
3. Select the Patient. The **Spool Management** window appears.



Figure 5-5: Spool Management Window—One Image Selected

4. Select the image(s) you want to work with. The thumbnails are displayed in the order of acquisition. Use the arrows to scroll if there are more than four views.

If you select a Patient with Rejected Images, you have a thumbnail picture of each of the patient's Rejected Images and demographics. When you click on a thumbnail, information about that view is displayed. See Figure 5-5, page 110.



Note...

There is an indication in the upper right pane as to whether or not the archive has committed the View to storage if that function is enabled.

5. Select the tab for the function you want to perform. The details for each function are explained below.
6. Click **Close** to exit Spool Management.

4.1.1 Resend Function

The Resend function allows you to query for Patients whose images are still within the Image Spool and to resend selected images in the spool to Output Destination devices. No Patient or Image object information can be changed from the Resend function. This includes Image or LUT modifications and DICOM information such as Patient Name spelling, ID, Comments, etc. Images are resent to the selected Output exactly as found in the spool.

4.1.2 Printer Contrast and Density

The Print Client in Spool Management applies the Lorad Image Processing algorithm (using the breast density/contrast control parameter) to the Image data as it is sent to the printer.

If a printed film needs to be changed, select Spool Management, and Resend the patient or image. Change the Density and Contrast as necessary for the image. You can resend one or many images, and any changes you make to the Contrast and Density apply to all images you send and only at that time. They are not saved with the images.

4.1.3 What Happens When You Select the Resend Tab

- A list of available outputs becomes active. The default output is the one currently selected in the Patient View Window.
- You can choose to resend all images for the Patient or only selected images.
- Different options are available, depending on the Output(s) selected.
- For example, if you have selected a printer, **Apply PCE** (the Peripheral Contrast Enhancement Function which uses a special algorithm to enhance the appearance of the mammograms) becomes an option.
- A “Resend” button becomes active.

At this point the following functions are also available:

- Query for a New Patient.
- Close Resend without taking any actions.

4.1.4 How to Use the Resend Function

To use the Resend Function:

1. Select a Patient.
 2. Click the **Resend** tab.
 3. Select the image to Resend or select **Resend All**.
 4. Select the **Output**. The Output selected applies only to the Resend function and does not affect Outputs chosen for normal acquisition. Use the **Edit Output** option to add or modify the Outputs.
 5. If the selected Output has a printer, the following **Print Image Processing Options** appear. The printing changes are not saved anywhere. Click the option to select or deselect it.
 - Adjust the **Contrast** and **Density** for printing. Change by clicking another number to increase or decrease.
 - **Print True Size** (if available at the printer).
 - **Apply PCE** (the Peripheral Contrast Enhancement Function uses a special algorithm to enhance the appearance of the mammograms).
 - **Omit Patient Name**.
 - **Omit Label options**.
1. Click **Resend**. A dialog box confirms the resend has been performed and the number of images resent. Click **OK**.
 2. Click **Close** when you are finished with the **Spool Management** function.

Writing images to a CD takes longer than saving or printing them. For this reason it is usually best to use the Resend option and write them after the procedure is finished. This has the advantage of not slowing down the rest of the procedure while you are waiting for the program to write to the CD during Image Acquisition.

4.1.5 Image Review Function (Also Known as Repreview)

The Image Review function allows you to recall a Patient Image for review and modification. All modifications made during the review are saved with the new image. Additionally, the DICOM tag, Image Type attribute (0008,0008) is set to “DERIVED\SECONDARY.” A unique DICOM Image instance UID is created for the image. The image is reclaimed using standard reclaim rules.

You have all the capabilities as in the original Preview (other than reject and accept). The differences from the standard Preview function include:

- The image is not sent to any output devices. If you need to send the image to an output device, use the “resend” function.
- Save and Cancel buttons are available.

Save creates a new image for the Patient. The original image is not modified in any manner. When you click Save, you save only the changes in markers, orientation and comments. Printer changes are not saved.

Cancel closes the Preview window and displays the Spool Management window.




Figure 5-6: Repreview Window

To use the Repreview Function:

1. Select a Patient from the **Choose A New Patient** list.
2. Select the image to Review.
3. Click the **Repreview** tab.
4. Click **Repreview**. The **Repreview** window appears.
5. Perform the actions and then click either **Save** or **Cancel**.
6. Click **Close** when done with the **Spool Management**.

4.1.6 Archiving Patients with Acquired Images

1. Select **Spool Management** from the Admin Menu.
2. Click **Search** and select the Patient name from the list. The **Review Images** window appears.
3. Click the images you want to resend, or select **Resend all the images**.
4. In the **Setup** window, select any other choices. Select the **Output** for archive.
5. Click **Resend**.
6. Close the **Spool Management** dialog box.

4.1.7 Correcting Incorrect View Selections

The Repreview function is used correct an incorrect view exposure if the error was not corrected before clicking Accept. Refer to Chapter 3, Section 9.3, page 50.

4.1.8 Delete Image Function

This function allows you to delete all or selected image(s) of a patient from the hard drive. Use this feature to delete incorrect images that are corrected and resent. When you click **Delete**, a dialog box reminds you that the image will be permanently deleted from the Image pool.

To use the **Delete Image** function:

1. Select the image to delete.
2. Select the **Delete** tab.
3. Click **Delete Selected** to choose the Patient Images, or **Delete All for Patient**.
4. Click **Yes** in the **Confirm Operation** dialog box.
5. Click **Close** when done with **Spool Management**.



Figure 5-7: Spool Management-Delete Image Function

4.2 Delete Patients

This option at the *mgr* level allows a manager to delete Patients from the Local Exams and all of the Patients' Images from the Image pool. Successfully exported images and patients are routinely reclaimed, and do not need to be deleted manually.

You can delete a Patient under the following circumstances:

- They have no rejected images in the reject bin.
- All their associated images are successfully exported or deleted from output queues.

To delete a Patient:

1. Select a Patient name from a list, or use the Patient name query to search for the Patient.
2. The **Delete** button activates when a Patient is selected.




Figure 5-8: Delete Patient Window

3. Click **Delete**.



Figure 5-9: Patient Delete Error Dialog Box

4. A dialog box might display reasons why the deletion cannot be performed, or inform you that the patient is protected. If there are protected patients in your selection, a dialog box tells you this and requests a confirmation to delete a protected patient. If the message indicates, for example, the Patient has jobs in the reject bin and/or the output queue), the jobs must first be deleted or allowed to pass through the queue.




Figure 5-10: Delete Selected Patients Dialog Box

Otherwise, A **Delete Selected Patients and Images** dialog box appears.

5. Click **Yes**. The Patient and all related data are deleted.

4.3 Reclaiming



Figure 5-11: Hard Drive Reclaimed

The Acquisition Workstation computer’s hard drive automatically stores the Acquired Images. Information that is successfully sent to other storage is routinely reclaimed from the hard drive in order to make room for new Patient information.

The system automatically reclaims Hard disk storage space for new Procedures using configurable parameters. The image reclaimer deletes images from the system when a number of images (service configurable) have accumulated on the system. Once reclamation of images begins, images continue to be reclaimed until the number of images on the system reaches a lower service configurable number.

The following Reclamation rules use the factory default parameters:

- If there are more than 2000 images on the system, then the 25 oldest images are reclaimed. No Procedure is left with missing images.
- Patients with no images are reclaimed.
- If there are more than 700 Patients, then the 50 Patients (and their images) with the oldest Procedures are reclaimed.

The following is a list of exceptions to these rules:

- No Patients are removed if there are any images queued to an output device.
- No images are deleted for Protected Patients. (Protect a Patient by using the Main Menu function Admin>Protect Patient. Be sure to unprotect them.)
- No images queued up on an output device are reclaimed.
- No images are deleted that have not been successfully printed on film or committed to an archive device, (Archive devices should be set to send commitments).
- Rejected images are not automatically reclaimed. Managers need to periodically remove all Rejected images using **Admin> Reject Management**.



Caution:

The system does not acquire a new image if the disk does not have the necessary space for a new image. This condition can occur if there are too many protected Patients or images backed up on an output device. Delete Patients with images not needed on the hard disk drive.



Figure 5-12: Hard Drive Full Alarm

4.4 Reject Management

This function allows a **mgr** user to manage the Reject Bin. The System Reject Bin provides the following functions:

- Delete Rejected Images
- Display the Reject Bin percentage full
- Select a Patient with Rejected Images and Review the images
- View a report of Rejected Images

A local procedure should be in place for routine examination and removal of Rejected Images, because they are not automatically reclaimed and can fill up the hard drive.



Figure 5-13: Reject Management Window—(Resend is Shown)

4.4.1 View a Patient's Rejected Images

When you select a Patient with Rejected Images, you have a thumbnail picture of each of the patient's Rejected Images and this information in the upper left pane:

- Patient Name
- If they are protected
- Patient ID
- Patient Birth Date
- Number of images in the Procedure
- Procedure number
- Procedure name

When you click on a thumbnail, information about that view is displayed in the upper right pane:

- Image view
- Date and time of capture
- Tech
- All image and technique information including the paddle used
- If it is Committed
- Comments
- Last six digits of UID

These Tabs allow you to administer the selected patient's images:

- Delete—Delete a selected (or all) image(s) permanently from the Reject Bin.
- Resend—Select an output and re-send a selected image to the output.
- Repreview—View the image

The Rejected Images and the comments are stored in the Reject Bin. DICOM attributes associated with the image and the LUT are held in the Reject Bin in the same manner as Accepted Images.

Proper image evaluation should be performed under the **Repreview** function to ensure effective image evaluation management.

4.4.2 Repreview Rejected Images



Figure 5-14: Repreview Tab on Reject Management Screen

To repreview rejected images:

1. Select **Admin>Reject Management**
2. Click **Choose a New Patient**.
3. Search for the patient in the **Choose a New Patient** dialog box. The patient's rejected images are displayed in the lower left pane.
4. Click the image you want to review. Image information is in the top of the box.
5. Click the *Repreview* tab, the **Repreview** option button is active. If the Repreview button does not display, contact service.
6. Click the **Repreview** button. The Repreview screen appears.
7. Click **Save** or **Cancel** to close the Repreview screen.
8. Click on the **Delete** tab to delete them, or click on the **Resend** tab to resend them to an output and then delete them.
9. Click **Close** to close **Reject Management**.

4.4.3 View Rejected Images Report

There are two different Rejected Images Reports. One on screen where you can see the rejected images, and one in the network browser that presents statistics and can be printed if the computer you are using is connected to a printer. For information on using the printable report, refer to Section 7.0, page 124.

When you choose to view the Rejected Images Report, you are able to filter the report to a single technologist, or view all technologists. The report lists the operator name, date, and time for each Rejected Image. Selecting a Rejected Image in the report shows the same information as explained in the previous section, “View a Patient’s Rejected Images,” page 118.

To review the rejected images report:

1. Select Admin>Reject Management
2. Click **Rejected Image Report**.
3. Search for a technologist in the **Choose an Operator** dialog box or click **Search** to search for all technologists with rejected images.
4. Select a technologist from the list. The technologist’s rejects for the procedure on the date and time listed display on the screen as thumbnails.
5. Click on the **Repreview** Tab to review them. Click on the **Delete** tab to delete them, or click on the **Resend** tab to resend them to an output and then delete them.
6. Click **Close** to close **Reject Management**.

4.4.4 Delete Images from the Reject Bin

Rejected Images should be deleted when they are no longer needed. If this alarm appears on the screen, someone needs to “empty” the Reject Bin.



Figure 5-15: Reject Spool Full Alarm

To delete images from the Reject Bin in Reject Management by patient:

1. Select a **Patient**, (you can review and resend the image before you delete it).
2. Click the **Delete** Tab.
3. Select the thumbnails of the images to delete and click **Delete Selected**, or click **Delete All for Patient** to delete all of that Patient’s images.

To delete them by the technologist who acquired them, use directions in Section 4.4.3, p. 120.

4.4.5 Resend Images from the Reject Bin



Figure 5-16: Reject Management—Resend a Rejected Image

To resend the image before you delete it:

1. Select the **Resend** tab. This window functions the same as those in the regular **Resend Menu**.
2. Click **Resend selected images** or **Resend all images for patient**.
3. Select the **Output**. The Output selected applies only to the Resend function and does not affect Outputs chosen for normal acquisition.
4. If the selected Output has a printer, the **Print Image Processing Options** activate. The printing changes are not saved anywhere.
 - Adjust the **Contrast** and **Density** for printing. Select a different number to make a change.
 - **Print True Size** (if available at the printer).
 - **Apply PCE** (the Peripheral Contrast Enhancement Function uses a special algorithm to enhance the appearance of the mammograms).
 - **Omit Patient Name**.
 - **Omit Label options**.
1. If you are done with the image and are ready to delete it, click **Delete Tab** and select the images to delete if the selection has changed.
2. Click **Delete** or **Delete All for Patient**. This only deletes images in the reject bin.
3. Click **Close** to close **Reject Management**.

4.5 Log Backup

When Service requests, Log Backup saves copies of the system log files to a CD-R.

To back up Log Files:

1. Select **Log Backup** from the **Admin** menu. The Backup dialog box appears.




Figure 5-17: Log Backup—First Dialog Box

2. Select **Backup to CD-R**. Click **Continue**.



Figure 5-18: Log Backup—Second Dialog Box

3. Insert a blank CD-R into the CD-RW device.



Figure 5-19: Log Backup—Third Dialog Box

4. Click the **Test Blank CD-R** button. The message **Blank CD-R found** appears. Click **Continue**.




Figure 5-20: Log Backup—Fourth Dialog Box

5. Click **Done** after the files are backed up. The disk auto-ejects.
6. Label it per your site procedures and store it safely or send to Hologic if requested.

4.6 Paddle Configuration

This menu option sets up the paddles available to your site. If these options are not set up correctly, an “Incorrect Paddle Position” alarm message appears.




Figure 5-21: Paddle Configuration Dialog Box

To check a Paddle Configuration or add/change an available View for a paddle:

1. Click the arrow for the drop-down list in the center.
2. Select the paddle you want to View/Configure. Your selection appears on the **Paddle ID info** line.
3. Select **Left, Center, or Right** paddle position for all shifting paddles.
4. Check (or clear) the **Allowable Views** for that paddle.
5. Repeat steps 1-3 for each paddle you want to configure and for each paddle position.
6. Click **Done** when finished.

5.0 Info Main Menu Function

The **Info** function identifies the Acquisition Workstation software versions, the current temperature of the array, and the date, time, and temperature of the last calibration.



Figure 5-22: Info Window

6.0 ALARM Function

The Alarm function manages alarms that have occurred due to failed jobs to destination devices. For more information on Alarms, refer to Chapter 6, page 127.

7.0 Repeat/Reject Reports Options Available from the Browser

7.1 Administrative Reporting

The Reject Analysis Report can be printed for administrative review. It is available from a web browser, not from the Selenia application. Any network connected PC can access them.

7.2 Accessing the Report from Any Computer on the Network

To access these reports:

1. Run the Browser program on your machine.
2. In the Browser Address Bar, type: `http://<the system name you want to access>` or the `<IP Address>:8080` and press Enter. A Password dialog box appears. (The system's IP address is recorded in an Appendix of the Installation Manual.) For additional details, refer to Chapter 2, Section 13.3, page 15.
3. Exit from the Browser when you are done to maintain security.

7.2.1 Log In Using the Admin Password.

Enter Logon Information For This Site

Username	<input type="text" value="admin"/>
Password	<input type="password" value="*****"/>
<input type="button" value="Logon To The System"/>	

Figure 5-23: Log on to the Browser

Administration Tools LOGRAD *Selenia*

Options

[Change Password](#)
[Reject](#)

Done

Figure 5-24: Administrative Reporting Options

You can choose to change the browser's password or open Repeat/Reject Analysis.

7.3 Change Password

This changes the browser's password for a Selenia system. Do not change this without system wide agreement.




Figure 5-25: Change Password

7.4 Reject Analysis Report

1. Click the Reject option.
2. Choose the Report Parameters to view.
3. Click Submit. The request is sent to the computer and the information is collected and formatted into a report.
4. Set the page to Landscape, and print the report using the browser’s print button.

Administration Tools LORAD Selenia

Reject

FROM: Day: Month: Year:

TO: Day: Month: Year:

Technologist: Type:

Figure 5-26: Criteria for Reject Analysis

To Save the Reject Report as a Shortcut on Your Desktop:

1. With the Reject Report Parameters page showing, click **Add to Favorites** in the browser.
2. In the dialog box that appears, type Selenia Reject and the room number if there are multiple units.
3. Click the appropriate button to save the favorite.
4. Once the favorite is saved, click on Favorites.
5. Highlight the link you just saved.
6. Right-click to display a menu of options.
7. Click Send To.
8. Click Desktop (create a shortcut). The shortcut is now on the desktop.

Chapter 6—Troubleshooting

1.0 Alarm Types and System Messages

The goal of an Alarm is to provide feedback about abnormal events associated with either the Image Acquisition on the Acquisition Workstation or the export of Acquired Images to destination devices.




Figure 6-1: Location of Alarms by Type

Legend Figure 6-1

1. Category Two Alarm Icon
2. Category One Alarm Message
3. Reset Generator button
4. System Information and Messages

There are four basic categories of alarms/messages:

1. A condition in which the Acquisition Workstation cannot continue (Category One).
2. A condition in which the operation of the Acquisition Workstation can continue (Category Two—normally associated with a warning or a failure in an output device).
3. A condition in which the image processing has failed, displayed on the Preview screen (Category Three).
4. A condition in which the system produces messages explaining why the system is not going into Ready Mode (System Messages).

1.1 Category One Alarm Functions

Category One Alarms appear on the screen in the Alarm Message box. A dialog box appears for Alarms requiring your immediate input. You must close the dialog box before you can acquire an exposure.

For *Category One* alarms, a dialog box appears immediately to indicate a problem that you are required to acknowledge. Depending on the specific condition, a dialog box asks whether to continue or to abort the current action. Some image receptor, generator, and internal alarms are Category One.




Figure 6-2: Category One Alarm Dialog Box

1.1.1 Sample Category One Alarm Messages

The Operator Message for each alarm contains the alarm number, the associated message as listed below, and the action required (for example, “Press Reset to continue”). If Reset is required, click Reset Generator.



Note...

Some conditions, such as Invalid Paddle Position, remain active until the condition is corrected. The Reset command is not required; the Acquisition Workstation message clears automatically once you correct the condition.

Sample Alarm Messages

Table 6-1 describes Category One Alarm messages.

- “Calculated exposure time exceeds max. Exposure Aborted. Press Reset to Continue.”
- “Premature release of x-ray switch. Exposure Aborted. Press Reset to Continue.”

1.1.2 AEC Alarm Messages

When the system is unable to use the selected AEC mode to acquire an Image, one of the following messages is displayed (Figure 6-3 and Figure 6-4).

When you have an Error during exposure:

1. Click **OK** in the dialog box.
2. Click **Reset Generator**.

If the system is unable to use Auto-Filter or Auto-kV, try changing to Auto-Time and selecting a higher, or lower kV as appropriate



Note...

The operator should always be aware of the position of the AEC Sensor when making AEC exposures.

Calculated Exposure Time Less Than Min




Figure 6-3: AEC Error Message—Below Minimum

When the AEC Sensor is over an unattenuated or very minimal breast tissue area, the alarm message, “Calculated exposure time less than min,” displays in the Generator Status pane. The AEC Sensor should be used in a position under breast tissue. Alternately, use AutoTime and a lower kV.

Calculated Exposure Time Exceeds Max




Figure 6-4: AEC Error Message—Above Maximum

When very dense breast tissue or any object such as an implant, pacemaker, or other anatomy influences the AEC Sensor, the alarm message “Calculated exposure time exceeds max,” is displayed in the Generator Status pane. The manual AEC sensor should be used in a position where it is not obstructed. Alternately, use AutoTime and a higher kV.

1.1.3 Generator Related Alarms

If the Gantry generator sends a message classified as an alarm, it displays as a Category One alarm. See Figure 6-2.

1.1.4 Digital Image Receptor Related Alarms

There are four basic classes of alarms associated with the Digital Image Receptor Interface. All four are Category One alarms.

<i>Start-Up:</i>	<i>The array is not responding at start-up.</i>
<i>Communication:</i>	<i>Communications were lost or incorrect mid-operation.</i>
<i>Image Receptor Error:</i>	<i>Error reported from Image Receptor</i>
<i>Scan ID Error:</i>	<i>The Scan ID does not match the one returned from the array.</i>

1.2 Category Two Alarm Functions

Category Two alarms are generally output problems that occur due to jobs that fail at the destination devices. When a new Category Two alarm occurs, the Alarm icon indicates that an alarm is present. Viewing the new alarm automatically acknowledges the alarm condition. The Alarm icon highlight remains if alarms are still present. Erase the alarm to clear it.

Output queues handle actual output from the Acquisition Workstation. If a physical device does not accept an exported job, or if it does not complete successfully, it is retried or held as a problem job, depending on the reason for the error (Figure 6-6). The system performs a service configurable number of retries before it generates an alarm.

Jobs can fail because there is either a problem with a specific part of the job (for example, printing requirements that a printer doesn't support) or because there is a general problem with the physical device (for example, the printer is off-line or not turned on). The system manages queues automatically to prevent problem jobs from clogging the queues. You can view the status of Images that have not left the queue using the Admin>Manage Queues menu function.

An alarm condition is either a warning or a notice. When a message such as a low-film condition is received, an alarm of type "warning" is posted to indicate this. However, this alarm does not get automatically removed when the low-film condition no longer exists; instead, another alarm of type "notice" is posted indicating that the printer status is now normal. The "From" column is not the name of the device posting the alarm, but the name of the spooler to which the alarm was sent.

The Category Two Alarm provides the following functions:

1. When you click the Alarm icon, the Alarm dialog box becomes active.
2. The Alarm dialog box allows you to:
 - View alarms
 - Acknowledge alarms
 - Erase alarms




Figure 6-5: Alarm Icon

The Alarm dialog box (Figure 6-6, page 131) contains the following:

1. Information about the alarm condition:
 - From: Name of the device with the alarm
 - Type: Warning, Notice
 - Subject: Additional system information
 - Time Sent: Time the alarm was posted to the spooler
2. An area that displays the message contents.
3. A list that shows all of the current alarms sorted in reverse chronological order.
4. A Forward/Backward button set that allows you to quickly navigate the alarm list.
5. An Erase and Erase All button that allows you to erase the alarm.




Figure 6-6: Alarm Dialog Box

Legend for Figure 6-6

- | | |
|-----------------------------|----------------------------|
| 1. Alarm Information | 4. List Navigation Buttons |
| 2. Message Display Area | 5. Erase and OK Buttons |
| 3. Scrolling List of Alarms | |

Erase an alarm by selecting the alarm and then clicking **Erase**. Clicking **Erase All** returns the Alarm icon to normal. Erased alarms are completely deleted. Therefore, they may recur, in which case they are treated as new alarms.



Note...

If you accidentally click the main screen when an Alarm dialog box appears, re-click the Alarm icon and the dialog box re-appears. If a software problem occurs and a screen stops responding, right-click the top bar of the problem screen and close it.

1.3 Category Four Alarms

System Messages appear in the lower right System Message box. They indicate why the system is not going into the Ready Mode. Table 6-4, page 143 lists required user responses to these messages. The dashboard panel displays text and changes color (to white) when a system problem occurs. Loss of PXCM to Generator link, receptor power loss, and “Paddle not locked”, are examples of messages that appear in the area.

An example of a System Message is:

- “Incorrect Paddle Position.”




Figure 6-7: System Messages

2.0 Tables of Alarm Messages, System Messages, and User Actions

2.1 Category One Alarms

Table 6-1: Alarm Messages

Error #	Cause	Action
2	Calculated exposure time exceeds max	reset / retake in AutoTime mode
3	X-ray switch not released after exposure	reset required
4	Stuck x-ray switch on start-up	reset required
5	Calculated exposure time less than min	reset / retake in AutoTime mode
9	Invalid Cal Data (service call required)	"non-resettable error"
20	Rotor error	reset / re-expose
21	No x-ray switch after start from host	reset / re-expose
22	ARC-kV control error	reset / re-expose
23	TUBE OVERCURRENT	reset / re-expose
24	TUBE OVERVOLTAGE	reset / re-expose
25	INVERTER OVERCURRENT	reset / re-expose
26	X-RAY INTERLOCK FAULT	reset / re-expose
27	FILAMENT OVERCURRENT	reset / re-expose
28	FILAMENT OVERVOLTAGE	reset / re-expose
29	FILAMENT GRID FAULT	reset / re-expose
30	HOST TO GEN. COMM FAULT	reset / re-expose
31	HARDWARE BACKUP TIMER TIMEOUT	reset / re-expose
32	FIRMWARE BACKUP TIMER TIMEOUT	reset / re-expose
33	HV GENERATOR TIMING FAULT	reset / re-expose
34	Spot relay error (unable to switch focal spot)	reset / re-expose
35	Tube Undercurrent Fault	reset / re-expose
51	Bucky error: Bucky not responding	reset / re-expose
52	Brick communication error during exposure (message block 1 request not received)	reset / re-expose
53	No DETECTOR_READY (at start of start exposure, or premature de-activation during x-ray)	reset / re-expose
54	READY signal error (premature de-activation during x-ray)	reset / re-expose
55	Invalid parameter field received in Block 2 or 3 command	reset / re-expose
56	Invalid value field received in Block 2 or 3 command	reset / re-expose
57	Post Exposure Data Error (Selenia did not successfully process post exposure data)	reset / re-expose
58	ERROR signal activation (premature activation during x-ray)	reset / re-expose

Table 6-1: Alarm Messages (Continued)

Error #	Cause	Action
59	No mAs received from Brick (New mAs not received during AEC exposure)	reset / re-expose
60	CCW C-ARM Rotation Switch Stuck (stuck switch detected on start-up)	inspect switch
61	CW C-ARM Rotation Switch Stuck (stuck switch detected on start-up)	inspect switch
62	Lamp Switch Stuck (stuck switch detected on start-up)	inspect switch
63	Comp Up Switch stuck (stuck switch detected on start-up)	inspect switch
64	Comp Dn Switch Stuck (stuck switch detected on start-up)	inspect switch
65	Comp Rel Switch Stuck (stuck switch detected on start-up)	inspect switch
66	C-Arm Up Switch Stuck (stuck switch detected on start-up)	inspect switch
67	C-Arm Dn Switch Stuck (stuck switch detected on start-up)	inspect switch
70	No comm from Tubehead	reset
71	filter/mirror error	reset
72	filter error (AF mode)	reset
75	Invalid Paddle Position (Compression Paddle shifted to an "in-between" position.)	*self-correcting
76	Unrecognized Paddle	*self-correcting
79	Compression Device Interference	*self-correcting

*Self correcting errors: Error remains active until the error condition is corrected, and clears once the error condition is corrected—it does not require a "reset." The alarm message on the screen clears once the error clears.

**Undefined error codes: If any error codes not defined above are reported, the alarm message displays "Generator Error#" where # is the error code. Assume that a reset is required.

2.2 Category Three Alarms

Category Three alarms are image processing alarms. An error code in the Exposure Index box on the Preview Screen informs you of problems that have occurred. These alarms appear on the Preview Screen when image processing fails.

Table 6-2: Image Processing Error Messages

Error #	Message
1000	Image processing is not licensed.
1001	Image is not FOR PROCESSING.
1002	Image orientation flag is missing or not supported.
1003	Paddle ID error.
1004	Internal configuration does not allow for any outputs!
1005	Out of range pixel values.
1006	Calling application did not allocate memory for one or more image outputs.
1007	Invalid data in input image. Data may not be of a mammographic breast image.
1008	Memory allocation error.
1009	Error opening configuration file.
1010	Error opening log data file.
1011	Data bit depth not supported.
1012	Missing required input parameter.
1013	Collimation configuration error.
1014	Laterality error.
2000	Density/contrast input parameters out of range
2001	Pixel manipulation input parameter out of range
2002	Histogram input parameter out of range
2003	Unsharp mask input parameters out of range
2004	PCE input parameters out of range.
2005	P-Value input parameters out of range
2006	Print boundary parameters out of range
2007	Diagnostic density/contrast parameters out of range
2008	Implant parameters out of range
2009	Image DC Offset out of range
2010	MCE parameters out of range
2011	Diagnostic logistic parameters out of range
2012	Screening logistic conversion parameters out of range.

2.3 Category Two Alarms

Only the key part of the message is included in this table.

Table 6-3: Category Two Alarms

Message From	Message	Action Required
Printer or Storage Task	Print request for ID: ..., Name: ... Study: ..., Image: ... to device ... has failed and will be retried. Reason: Failed to communicate with the printer pdev :pdevprint_...; will retry Please contact Hologic service.	Please contact Hologic service.
	Print request for ID: ..., Name: ... Study: ..., Image: ... to device ... has failed after ...attempts, and will not be retried. Reason: Failed to communicate with the printer pdev :pdevprint_...; will retry Please contact Hologic service.	Please contact Hologic service.
	Store request for ID: ..., Name: ... Study: ..., Image: ... to device ... has failed after ...attempts, and will not be retried. Reason: Failed to communicate with the remstore pdev :pdevremstore_...; will retry Please contact Hologic service.	Please contact Hologic service.
	Store request for ID: ..., Name: ... Study: ..., Image: ... to device ... has failed and will be retried. Reason: Failed to communicate with the remstore pdev :pdevremstore_...; will retry Please contact Hologic service.	Please contact Hologic service.
	Error happened when trying to initialize CDRW	Try another CD.
	CD-R space full, and user selected not to replace with a new CD-R	Resend the image later if you want it on a CD.
	CD-R was not recordable, unknown error happened	Try another CD.
Output devices	Messages generated from the outputs or workstation and passed through to the user (for example, printer low on film)	No. See device documentation for action.
Database Peer Manager	Version mismatch detected: Local (xxx) Remote (yyy)	Call Hologic for service.
	Unable to connect to peer - please consult log file	Call Hologic for service.
MWL Updater	Error updating local worklist	Contact Network Administrator
	MWL updater query received too many results. Not updating local worklist.	Refine query.
	MWL periodic updater timed out querying SCP	Refine query.
	Worklist updater encountered error querying SCP	Contact Network Administrator.
	Worklist updater: Unexpected error encountered.	Call Hologic for service.

Table 6-3: Category Two Alarms (Continued)

Message From	Message	Action Required
DICOM Print alarms	Reading or processing data (image, lut, or marker) failed locally. Please contact Hologic service.	Please contact Hologic service.
	A DICOM dimse error was returned in a response <msg from DCF> Please contact Hologic service.	Please contact Hologic service.
	Cannot open socket to specified host/port, a network error has occurred or received abort from remote system. Please check destination host setup. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	The remote device did not accept an association. Please check destination host setup. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	The remote device is currently unable to accept an association. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	Print job was not completed in the allotted time, or time-out waiting for response. Status of print job is unknown. Please contact Hologic service.	Please contact Hologic service.
	Requested DICOM SOP class rejected, or an unspecified error occurred. Please contact Hologic service.	Please contact Hologic service.
	(unknown error) Please contact Hologic service.	Please contact Hologic service.
	<Error or warning message from printer such as a film supply problem, processor problem, or other printer malfunction>	Please contact Hologic service. Also, please consult printer documentation.
	Device address not in form <hostname-or-ipaddr>:<port> Please check destination host setup. Please contact Hologic service.	Please contact Hologic service.
	Unable to ping destination host: <hostname> Please check network connections and destination host. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	An unspecified error occurred. Print failed and will not be retried. Please contact Hologic service.	Please contact Hologic service.
	Unexpected CORBA error. Print failed and will not be retried. Please contact Hologic service.	Please contact Hologic service.
	Unexpected local processing error. Print failed and will not be retried. Please contact Hologic service.	Please contact Hologic service.

Table 6-3: Category Two Alarms (Continued)

Message From	Message	Action Required
DICOM Store Commit	Reading or processing image data failed locally. Please contact Hologic service.	Please contact Hologic service.
	A DICOM dimse error was returned in a response. Please contact Hologic service.	Please contact Hologic service.
	Cannot open socket to specified host/port, a network error has occurred or received abort from remote system. Please check destination host setup. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	The remote device did not accept an association. Please check destination host setup. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	The remote device is currently unable to accept an association. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	Storage Commit job was not acknowledged in the allotted time. Status of commit job is unknown. Please contact Hologic service.	Please contact Hologic service.
	An unspecified error occurred. Please contact Hologic service.	Please contact Hologic service.
	(unknown error) Please contact Hologic service.	Please contact Hologic service.
	Unable to ping destination host: <hostname> Please check network connections and destination host. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	Device address not in form <hostname-or-ipaddr>:<port> Please check destination host setup. Please contact Hologic service.	Please contact Hologic service.
	StorageCommit failed and will not be retried. <error message from database or from DCF> Please contact Hologic service.	Please contact Hologic service.
	StorageCommit failed and will not be retried. Unexpected CORBA error. Please contact Hologic service.	Please contact Hologic service.
	StorageCommit failed and will not be retried. Unexpected local processing error. Please contact Hologic service.	Please contact Hologic service.

Table 6-3: Category Two Alarms (Continued)

Message From	Message	Action Required
DICOM Store Alarms	The remote device is currently unable to accept an association. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	The remote device did not accept an association. Please check destination host setup. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	A DICOM dimse error was returned in a response: <msg from DCF> Please contact Hologic service.	Please contact Hologic service.
	Time-out waiting for response. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	Unable to connect to destination host/port. Please check destination host setup and network connections. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	A network error has occurred. Please check network connections. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	Reading or processing image data (image, lut or marker) failed locally. Please also check local device configuration files. Please contact Hologic service.	Please contact Hologic service.
	Error creating a store job. Please check local device configuration files. Please contact Hologic service.	Please contact Hologic service.
	Unable to ping destination host: <hostname> Please check network connections and destination host. Please contact your PACS or network administrator.	Please contact your PACS or network administrator.
	Device address not in form <hostname-or-ipaddr>:<port> Please check destination host setup. Please contact Hologic service.	Please contact Hologic service.
	An unspecified error occurred. Store failed and will not be retried. Please contact Hologic service.	Please contact Hologic service.
	Unexpected CORBA error. Store failed and will not be retried. Please contact Hologic service.	Please contact Hologic service.
	Unexpected local processing error. Store failed and will not be retried. Please contact Hologic service.	Please contact Hologic service.
	<Error message from DCF> Please contact Hologic service.	Please contact Hologic service.

Table 6-3: Category Two Alarms (Continued)

Message From	Message	Action Required
CDRW Alarms	Part10Maker:checkDiskCompliant(). CD-RW Drive was busy, will retry later.	No.
	Part10Maker:checkDiskCompliant(). Exception caught when trying to check file system ISO 9660 Level 1 conformance with <disk_status.mount_point>.	User has to resend the image at a later time if they want it on a CD.
	Part10Maker:checkDiskCompliant(). Exception caught when trying to guess filename from <disk_status.mount_point>.	Resend the image later if you want it on a CD.
	Part10Maker:checkDiskCompliant(). Error happened when trying to initialize CdrwSpaceFullDialogObject	Resend the image later if you want it on a CD.
	Part10Maker:checkDiskCompliant(). CD-R does not conform to ISO 9660 level 1 Standard, and user selected not to replace with a new one.	Try another CD-R. Resend the image later if you want it on a CD.
	Part10Maker:checkDiskCompliant(). CD-R does not conform to ISO 9660 level 1 Standard, unknown error happened in CdrwSpaceFullDialog Object.	Try another CD-R. Resend the image later if you want it on a CD.
	Part10Maker:checkDiskCompliant(). Failed to find any valid media from CD-RW drive. and user selected not to insert with a valid one.	Try another CD-R. Resend the image later if you want it on a CD.
	Part10Maker:checkDiskCompliant(). Failed to find any valid media from CD-RW Drive. unknown error happened in CDrwSpaceFullDialog Object.	Try another CD-R. Resend the image later if you want it on a CD.
	Part10Maker:doit(). Exception caught while trying to reset the temp dir:<TMPDIR>	Contact Hologic Service.
	Part10Maker:doit(). Exception caught while trying to decrypt image file or convert image to part-10 format.	Resend the image later if you want it on a CD.
	Part10Maker:doit(). A matching image record already exists under that series.	No. The image already exists in the current CD-R.
	Part10Maker:doit(). Exception caught while trying to create DICOMDIR file.	Resend the image later if you want it on a CD.
	Part10Maker:doit(). Failed to create temp dicomdir hide file:<dirHideFile>.	Resend the image later if you want it on a CD.
	CdrwJob_j failed to created CdrwOutput Object to handle output task. device resource name=<m_description.resource_name>	Resend the image later if you want it on a CD.
	CD-RW Fail: error detecting burning speed. Check CD-RW device path, and make sure the CD is recordable.	Retry. If condition persists with different media, contact Hologic service.

Table 6-3: Category Two Alarms (Continued)

Message From	Message	Action Required
CDRW Alarms (continued)	CD-RW Fail: error detecting burning speed, CD-RW or CD-R not compatible.	Try another CD-R. Be sure you are using approved CDs.
	CD-RW Fail: error reading CD-RW or CD-R, make sure the CD-R is recordable.	Try another CD-R. Be sure you are using approved CDs.
	CD-RW Fail: error finding session address, make sure the CD-R is recordable.	Try another CD-R. Be sure you are using approved CDs.
	CD-RW Fail: unable to detect the next available recordable address, CD-R probably unrecordable.	Try another CD-R. Be sure you are using approved CDs.
	CD-RW Fail: unable to decide the space needed of the CD-R	Retry. If condition persists with different media, contact Hologic service.
	CD-RW Fail: not enough space on the CD-R. CD-R space full, and user selected not to replace with a new CD-R.	Resend the image later if you want it on a CD.
	CD-RW Fail: not enough space on the CD-R. Error occurred when trying to initialize CdrwSpaceFullDialogObject.	Try another CD-R. Resend the image later if you want it on a CD.
	CD-RW Fail: error reading Cd-RW or CD-R, make sure the CD-R recordable. Error happened when trying to initialize CdrwSpaceFullDialogObject	Try another CD-R. Be sure you are using approved CDs. Resend the image later if you want it on a CD.
	CD-RW Fail: error reading Cd-RW or CD-R, make sure the CD-R recordable. CD-R was not recordable, and user selected not to replace with a recordable CD-R.	Try another CD-R. Be sure you are using approved CDs. Resend the image later if you want it on a CD.
	CD-R Partition being used by other applications (usually means File Manager or terminals are using CD-R). cannot continue, before retry, please exit those applications first.	Exit any other applications and retry.
	CD-RW Drive busy, will retry later.	No
	CD-R Warning: After the record process, CD-R could not be successfully remounted, Please unmount and mount CD-R manually to access the image on CD-R.	Remount CD-R manually.
	CD-RW Fail: error Stopping Volume Management.	Reboot the system. If condition persists, contact Hologic Service.
CD-RW Fail: error Restarting Volume Management. Please start it manually.	Reboot the system. If condition persists, contact Hologic Service.	
Reclamation Alarms	The system is unable to reclaim space for new patient records. Please ensure that there are not an excessive number of protected patients. Call service if problem persists.	Delete unnecessary protected patients.

Table 6-3: Category Two Alarms (Continued)

Message From	Message	Action Required
Repeat Spool Alarms	The repeat bin is full! Please have a manager enter "Repeat Management" through the toolbar ADMIN function to delete repeated images and create more space. Note... The above alarm is for tech logins only.	Delete rejected images in spool.
	The repeat bin is full! Enter "Repeat Management" through the toolbar ADMIN function to delete repeated images and create more space. Note... The above alarm is for all logins but tech logins.	Delete rejected images in spool.
	The repeat bin is nearly full! There is only enough space for <space> more repeated image(s). Please have a manager enter "Repeat Management" through the toolbar ADMIN function to delete repeated images and create more space. Note... The above alarm is for tech logins only	No. Delete rejected images in spool.
	The repeat bin is nearly full! There is only enough space for <space> more repeated image(s). Enter "Repeat Management" through the toolbar ADMIN function to delete repeated images and create more space. Note... The above alarm is for all logins but tech logins.	No. Delete rejected images in spool.
Sheet Alarms	<print job description>	Review device documentation.
Setup Alarms	Error constructing setup:<setup>	Contact Hologic Service.
DICOM Query/Retrieve	There were error(s) when pulling back the image(s).	Contact Network Administrator.
	(number) image(s) moved successfully. There were error(s) when importing the image(s).	If the message in alarm body doesn't help. Contact Hologic Service.
	Unknown error encountered when pulling back the image(s).	Contact Hologic Service.
	(number) image(s) moved and imported successfully.	No.

2.4 System Message Information

Table 6-4: System Messages

Message	Meaning	Action
Incorrect paddle position	This indicates that the paddle position is not valid for the selected patient view	Correct the paddle position. The Selenia system only allows the use of certain paddles with specific views. If this message appears, you may need to change the position of the paddle or use a different paddle.
Paddle not locked	This indicates that the paddle is between positions (left, center, right)	Shift the paddle until it clicks into position.
Remove paddle for calibration	This indicates that a paddle is present and a calibration scan is being requested.	Remove the paddle for calibration.
Detector power	The PXCM is reporting that the power to the detector has been removed	Make sure the gantry has power. Examine the emergency stops switches (three) to be sure the switches are off, and the circuit breaker to be sure it is on (up). Reboot the system. If the condition persists, contact Service.
Detector service mode	This indicates that the detector is in service mode.	Reboot the system. If condition persists, contact Service.
Calibration needed	This indicates that a calibration is needed and could be the result of a calibration file that is beyond the calibration life span	Perform a calibration.
Generator problem	The PXCM has lost its link with the generator	Reboot the system. If condition persists, contact Service.
Detector problem	The PXCM has lost its link with the detector	The imaging subsystem has reported that it lost the link with the detector. Reboot the system. If condition persists, contact Service.
Detector hi temp	The PXCM is reporting that the detector hi temperature has been reached and is in danger of being turned off	The detector reached the high temperature. Shut down the system. Contact Service.
Detector not present	This is a reminder that the acquisition workstation was started without a detector present	Not Applicable
Non Acquisition Mode	This is displayed during review and test pattern generation and indicates that the system cannot be used to take patient exposures	Not Applicable
Detector shutdown in (n minutes)	This indicates that the detector is about to be shutdown in n minutes	Click Cancel Shutdown if you are still working. The detector does not shutdown until the backup shutdown time.

3.0 Software Troubleshooting

Table 6-5: Software Troubleshooting



Symptom	Probable Cause	User Directions
<p>A message appears saying: " The Acquisition Station database is corrupt. The system may have been shutdown incorrectly." It will then give the options of:</p> <ul style="list-style-type: none"> • "Clear database and restart" • "Repair database and restart" • "Exit the application" 	<p>The system detects a database corruption</p>	<p> Caution: <i>If you see this message, you should first attempt to repair the database. If after performing this operation, and rebooting, you still receive this message, contact Service.</i> <i>The "Clear the database and restart" should only be used under direction of a service engineer.</i></p>
<p>White screen with message: program terminated ok</p>	<p>Power problems</p>	<ol style="list-style-type: none"> 1. Turn off the circuit breaker on back of Acquisition Workstation (Figure 6-8). 2. Wait for UPS to stop beeping (a minimum of 10 minutes). 3. Turn the circuit breaker switch back on. 4. Press the AWS On button. The computer restarts. 5. Log on to the OS.  <p><i>Figure 6-8: Acquisition Workstation Circuit Breaker</i></p>
<p>Forgot user password</p>		<p>Delete user account and re-enter.</p>
<p>Can't make a "tech" into a "manager"</p>	<p>Logged in as tech</p>	<p>This must be done by a Manager. Delete user account and re-enter it when logged on as a Manager.</p>
<p>Panel Power dialog box continues to reappear after using the Emergency Off Switch.</p>	<p>System may need to be reset or restarted.</p>	<p>Select the Reset Generator button after selecting OK. If that doesn't work, System may need to be restarted. Shut down the System following the procedure Chapter 1, "How Do I Turn the System Completely Off?," page 1. After waiting a few minutes, press the power button to restart the system.</p>

Table 6-5: Software Troubleshooting (Continued)

Symptom	Probable Cause	User Directions
No screen function, cursor still moves.	Dialog box hidden behind main screen	<p>If you accidentally click the Main screen when an alarm dialog box is displayed, re-click the alarm icon and the dialog box is re-displayed.</p> <ol style="list-style-type: none"> If the main screen appears to be non-functional, click the upper right corner dot of the Main screen bar, and the screen minimizes to an icon that says Lorad. Close any open application dialog boxes (not the little icon boxes on the left). Double Click the icon Application to return it to a full screen. <p>If the Task Launcher screen is accidentally closed, restart the application.</p> <ol style="list-style-type: none"> If a screen does not function—and, there is nothing behind it, after you have minimized the program—click the exit button on the OS task bar at the bottom of the screen. Log out of the session. Click OK on the logout confirmation dialog box. <p>Call field service if these solutions do not work. <i>Do not</i> turn off the Acquisition Workstation circuit breaker with the program running.</p>
Lorad wallpaper appears	Open key pressed.	Press the Open key again or click the application icon on the left and program should reappear.
Incorrect paddle position message doesn't go away	Changed or moved paddle while Preview screen was displayed Paddle not recognized.	<ul style="list-style-type: none"> Restore the previous paddle or position. You may need to also click Reset Generator. Re-seat paddle to confirm good connection. Examine Paddle Configuration table.
Collimation on spot paddles is wrong	Automatic Collimation was turned off in setup	Turn Automatic Collimation back on and leave on.
Error 100	Error transferring preview Image	Restart the system. Be sure to restart the computer every morning.
RPC3-Write error	Faulty generator communications, brick, cabling, etc.	Call Service.
Cannot delete patients	Logged in as tech	Log in as Manager.
Image Washed out		Evaluate exposure bar, should be above 250.
CD-RW will not respond		<ul style="list-style-type: none"> Wait five minutes before trying again. Look at CDRW spooler in Manage Queues to see if Images are still being transferred. Wait until Image transfer has finished before retrying. Try a different media.

4.0 Worklist Problems

4.1 Unknown RIS Code




Figure 6-9: Unknown RIS Code

This dialog box appears if the RIS code returned by the Modality Worklist Provider (Radiology Information System or PACS) does not map to an existing procedure on the Acquisition Workstation. Your options at this point are:

- Click **Cancel** because the previous selection was incorrect.
- Select the correct Procedure from the drop-down list and click **OK**.
- After the Patient is completed, contact your local IS department to ensure the RIS code matches your expectation. Then, if it is still a problem, contact the Selenia help desk for assistance in using the Procedure Editor to add the new RIS code.

4.2 Conflict Between Data in the Modality Worklist and Local Database

When a conflict between the worklist and local database does occur, work with the HIS/RIS administrator to resolve this issue.




Figure 6-10: Conflict Between the Worklist and Local Database

This dialog box indicates there is a conflict between the worklist data and an existing Patient in the local database. Verify that you have selected the proper record and that a mistake hasn't been made during scheduling and identifying the patient and procedure.

If you are to assign a new Accession Number, follow this procedure.

1. Click **OK** in the MWL error data box (Figure 6-10). A New Procedure dialog box appears.

The screenshot shows a dialog box titled "New Procedure" for "Patient, Six". It has two main input fields: "*Accession Number:" followed by a text box, and "*Procedure Description:" followed by a dropdown menu currently showing "Mammography Screening-Bilateral Mammography". At the bottom, there is a note: "*NOTE: Fields with an asterisk '*' are mandatory and require input." and three buttons: "Clear", "Cancel", and "Accept".

Figure 6-11: Sample Dialog Box: Add a New Procedure

2. Type a new **Accession number** and select the procedure.
3. Click **Accept**. A new Procedure is added for the Patient.

Another workaround to consider would be to have a manager delete the old record and then reselect the worklist item.

4.3 Resolving Patient Data Conflicts

Depending on which field doesn't match, a dialog box may be displayed allowing you to enter the Local Database information. These changes only apply to the Local Database on the Acquisition Workstation, and do not correct mistakes on the Modality Worklist Server. Inform your HIS/RIS administrator of this issue and work with them to resolve it so that images do not get lost. When you have finished correcting the discrepancies, click OK. The Patient's screen appears.

The screenshot shows a dialog box titled "Patient Data Conflict Resolution" with the subtitle "Resolve Patient". It is split into two columns: "Local" and "Network".
 Local fields: *Last Name: Patient, *First Name: One, Middle Name: (empty), *Patient ID: 00581381, *DOB: (MMDDYYYY) 01 / 01 / 1950, *Age of Patient: 054Y, *Gender: F.
 Network fields: Patient, One, (empty), 00581381, DIFF 04 / 01 / 1948, DIFF 055Y, F.
 At the bottom, there is a note: "*NOTE: Fields with an asterisk '*' are mandatory and require input." and an "OK" button.

Figure 6-12: Worklist Patient Conflict Dialog Box

Index

A

accept

- image 46
- rejected image 12, 47

Accept button 46

accepted image

- destination devices 74
- patient view screen thumbnail 49
- review dialog box thumbnails 47

access levels 105

Acquisition Workstation

- circuit breaker 97
- local database (Local Exams) 25

add

- additional procedure 33, 58
- additional view to procedure 56
- new procedure
 - to Open Patient 58
 - with new accession number 57
- new user 73
- patient 30

Admin menu

- mgr 109
- tech 82

administrative browser, change password 125

Administrative Reporting 124

AEC

- alarm messages 128
- Auto sensor position 40, 103
- exposure adjustment settings 39
- modes 3, 38–39
- sensor position 99, 103
- sensor position settings 40
- technique defaults 71
- using 3, 38

alarms 127–135

- Category One 128, 133–134
- Category Three 135
- Category Two
 - alarm icon and dialog box 130
 - alarm list 136
 - CDRW alarms 140
 - DICOM alarms 138
 - number of alarms 130
 - Repeat Bin alarms 142
- CD 81
- functions 124
- generator 128

- image receptor 128, 130
- image spools 130
- incorrect paddle position 9, 143
- managing 124
- output queues 130
- paddle not locked 9
- re-display 132

allowable views, paddles 123

AUTO AEC sensor position 40

B

bar code scanner, use 27

bar codes 65

beeping UPS 97

button (software)

- Accept 46, 57
- Add (Procedure) 58
- Add View 8, 56
- Brightness/Contrast (Window/Level) 45
- button locations on Patient View screen 34
- Cancel/Clear 57
- Clean Up 95
- Close (in Patient box) 106
- Close Procedure 34, 55
- Comments 43
- Edit 106
- Edit Patient 34
- Edit Procedure 34
- Eject 9
- Exit 20
- Full Zoom/Pan 44
- Import Priors 93
- Magnification 45
- Measurement 45
- New Procedure 33, 57
- Non-Imaging Mode 60
- Quick Zoom/Pan 44
- Reject 46
- Remove Unexposed Views 57
- Retrieve and Import 93
- Review 61
- ROI (Region of Interest) 44
- Search 31–32
- Set Accepted 47
- Set Rejected 48
- Set Technique Defaults 71
- Simulate Capture 61
- Window/Level (Brightness/Contrast) 45

buttons

- C-arm rotation and travel 100, 102
- compression release 101

C**calibration** 19**captured images, resending** 110**C-arm**

- alarms 134
- controls 66
- movement 66
 - disabled 100
- rotation displays, memory, and switches 100

CD, CD-RW

- alarms 81
 - initialization error 136
 - not recordable 136
- eject disks 90
- import images from 80, 90
- Outputs 79
- patients on a 9
- space full alarm 136
- use 80–81

CD-RW drive 80

- as output 74
- virtual device name 80

Center/Width, set 45**change Patient Information, Accepted Image** 59**change users** 3**Choose a Procedure, Close button** 106**circuit breaker**

- Acquisition Workstation 97
- Gantry 97

Clean Up button 95**clinical procedures** 62

- AEC mode 64
- manual mode 62
- TEC mode 63

Close

- Exam key 55
- Patient button 57
- Procedure button 34, 55

collimator

- override 102
- startup defaults 71

comments

- enter or modify 43
- required 46

Communication with the generator is lost 22**compression**

- controls and displays 99
- release buttons 101
- release modes 71

compression brake 100**compression force** 71

- display and range 99
- startup defaults 71

compression thickness display 99**Compression Up button** 101**contrast for printer, adjust in**

- Image Management 83, 112
- Reject Management 121
- Resend 8, 85

controls 97

- AEC sensor position 99
- automatic compression release 101
- C-arm 99–102
- collimator override 102
- compression 99–102
- compression release buttons 101
- handwheels 99
- light field lamp 102
- location, function and appearance 100

cropping, automatic 45**Crosshairs, use** 68–69**D****daily restart** 21**data storage (Outputs)** 79**database**

- access modality worklist 27
- available 25
- conflict 146
- corrupt 144
- local 26, 30
 - Local Exams tab 30
- organization 25
- synchronization 26

default view order 108**defaults**

- exposure modes 36
- set system 71

delete

- images 114
- images from Reject Bin 118, 120
- job in queue 89
- output 77
- patients 114

Delete tab 114
 in Repeat Management 121
density, adjusting for reprinting in Spool Management, Resend 111
desktop shortcut, Reject Analysis report 15
determine job status 89
dialog box, lost or missing 132
DICOM storage commitment 82
digital markers 43
disk
 eject CD from drive 81, 90
 record images on CD 80
Dorsal/Ventral hanging 43

E

Edit menu 71, 108
Edit Outputs
 new groups 77
 tab 76
Edit Patient button 106
edit patient data 57
 on an accepted image 59
 with MWL and PACS 10
Edit Procedure button 106
 location 34
Edit Users 72
Edit View 50
Eject CD 90
Emergency Off switch 21–22
 restart system after use 21
 use 21
error message, Invalid Paddle 123
error messages see alarms
Exit 70
exposure
 AEC adjustments settings 39
 information pane 34
 status 41
 with incorrect icon and edit view button
 corrected after accepted 50
 corrected before accepted 50
 with incorrect icon and no Edit View button
 correct after accepted 52
 correct before accepted 51
exposure modes
 AEC 38
 change resets exposure techniques 36
 default 71
 diagram 36
 select 36
 TEC 36

exposure techniques 34
 default settings 71
 reset to system defaults 36
 select 36

F

File menu 20, 70
film label area 42
flowchart
 AEC Modes 64
 Manual Mode 62
 TEC mode 63
footswitches 99
Full Compression Force, Startup defaults 71
Full Zoom/Pan 44
function
 C-arm rotation switches 100

G

Gantry
 C-arm rotation displays 100
 circuit breaker 97
 set up for screening 65
Gantry Rotation switch 100
generator
 message pane 34
 problem 22
 related alarms 130
 Reset button 34
 status 41, 66

H

Hanging options 43
Hard Drive Full alarm 116
hard drive, AWS computer 79
hints
 Artifact Evaluation test 17
 Flat Field Calibration test 17
 Phantom View and SNR/CNR tests 17

I

icons, system ready for exposure 66
Identify Yourself 3
Image Management
 Repreview 84
 Resend 83
 printing options 83
 use 82–85
Image Pane, film label area 42
Image Preview screen 42

Image Repetition Information dialog box 46

Image Review Function 112

image specimen 9

image spool 120

Image Status, alarms in queues 130

images

accept or reject 46, 69

acquire 41, 65–66

cropping, automatic 45

delete 114

from Reject Bin 120

rejected 118

DICOM attributes 118

import 90

files from CD 80, 90

manage 82

outputs see Outputs

preview 42, 69

printing 111

problem jobs 77, 87

query spool 82, 111

recall 112

reclaimed 79

rejected 46, 69, 74, 117

log 69

outputs 41

repreview 84

reprint 82

resend 82, 110–111, 113

post procedure 10

rejected 118

to new output 49

retrieve on review workstations 81

rules for automatic reclamation 79

send accepted to outputs 74

store 69

technologist initials 105

temporary storage 79

thumbnails 49–50

view full resolution 45

view rejected 118

as report 117

in Review screen 69

work with 41–55

Implant exposure techniques 7

implant present

button 41

correct information

after accepting 43

before accepting 43

implants

displaced views 7

general information 4

print images 8

procedures 7

import

image files 90

priors in background 93

incorrect information, places to remove 59

incorrect view selections

correcting with an Edit button 85

correcting without an Edit button 51

how to correct 50

Info menu 95, 124

initials, where used 73

invalid paddle, reconfiguring 123

J

job status

determine 89

viewing 77

L

laser printer 81

Launch, see Wait Time

light field lamp 102

activated by compression down button 100

use 66

Local Exams

delete Patients from 114

description 25

names 32

results too long 32

search 30

sources of data 25–26

Local Exams tab 25

localization

hints for procedure 8

paddle disables compression release 101

procedures 68

log backups, mgr function 122

log in and out of application 19–20

Log Out 21

Log out of the computer 1

Logout Confirmation dialog box 1

M

Magnification 45

Magnification modes, defaults 71

Main menu functions 70

by user 106–107

maintenance, system 105

Manage Queues 87–88

delete (cancel) job status 89

determine job status 89

resend stopped job status 89

Manual exposure mode, technique default options 71

Markers, digital 43

Measurement tool 45

memory, rotation switches 65

menus

Admin 82

Edit 71

File 70

Info 95

messages see alarms

Modality Worklist 25, 27, 29

MPPS, close a procedure 56

N

network storage 81

New

Procedure 57–58

User 73

No Outputs warning 41

Non-Imaging Mode

button 61

use 60–61

O

OS logon level

menu options 110

output images to CD 80

outputs 74–81

accepted images 46

add new 75

alarms for queues 130

change Resend printer options 83, 112, 121

change/select destinations button 56

delete 74, 77

described 74

edit 77

none 41, 49

on network 79

pane location 34

printer 81

queues 77, 87

rejected images 41

remove incorrect information from 59

retry job options 87

select 41

send

accepted rejected images to 48

images to destination devices 49

spools 49

use after procedure closed 10

use new 78

Override mode (TEC) 37

P

paddle alarms 9

paddle configuration, add or change views 123

Panel Power dialog box 22

password

assign 73

change 73

for Administrative reports 125

patient

archiving many 113

correcting wrong information 11

data conflicts, resolve 147

data, how to enter 30

deleting 114

demographics 25

edit data 57

with MWL and PACS 10

find on a CD 9

lists 25

Name tab for local database 31

pane on Patient View screen 34

protect from reclamation 86, 116

query function 29

search local database 32

with rejected Images 117

Patient Close button 106

Patient Query 29

how to enter patient names 32

Patient View screen 55

options chart 35

thumbnails 49

position AEC sensor control 103

position patient 66

post start-up tests 100

Power off the computer 20

PPS status 94

menu item 70

Pre-compression force, startup defaults 71

Preview Image 42, 69

Print Image Processing Options 83

printer 85

- adjust contrast and brightness 83, 112, 121
- alarms 130
- change resend options 83, 112, 121
- configuration 81
- copies of images 83
- test patterns 95

problem jobs

- in output queues 87
- not reclaimed 77

procedure

- add another
 - to open patient 58
 - with new accession number 57
 - with same accession number 58
- Add or Edit buttons 56
- choose 33
- saving incomplete 55

program terminated 144**Progressive Preview** 42**Protect Patients** 79, 86**protect Weekly QC** 14**Q****QC**

- hints 17
- protect Weekly QC patient 14

query 29–33

- hints for use 29
- images in spools 82
- Modality Worklist Provider 27, 30
- patient tab 29
- retrieve priors results 93

queues 77, 87

- manage 87
- number of Jobs in Queue 34
- output 77

Quick Zoom/Pan 44**R****radiation icon**

- changes color during x-ray exposure 66
- system ready indicator 66

Ready 41

- status 66

rebooting the computer 21**reclaimed, problem jobs not** 77**reclamation**

- Image spool 49
- rules for automatic 79, 116

Region of Interest (ROI), button 44**reject accepted image** 48**Reject Analysis**

- Office Computer 15
- Report 126
- report criteria 126
- Shortcut on the Desktop 15

Reject Bin 114–115, 117, 120

- empty 16
- managing 117
- see Reject Management

reject image 46**Reject Management** 117

- deleting images 120
- reject bin percentage full 117
- resending images 121
- view patient's rejected images 118
- view rejected images report 120

Reject Spool Full Alarm 120**Rejected Images** 117**removable storage** 79**Remove Unexposed Views button** 57**Repreview**

- images 82
- tab 112
- use 84
- window 113

Reprint image 111**Reprocess View** 43**Resend** 82–84, 111–112

- images 8, 10, 83, 110
- Images from Reject Bin 121
- mgr functions 111
- tech functions 82

resend stopped job 89**reset generator** 22, 128**reset UPS** 98**resolve patient data conflicts** 147**restart computer** 21**restart system after Emergency** 21**Retrieve Priors** 92**Review** 61

- use to reject Accepted image 48

Review dialog box

- use 47–48
- view tab 47

RIS code, unknown 146**rotation switch memory** 65, 100**Ruler tool** see **Measurement tool****Rules for automatically deleting images** 116

S

screening procedure 65

Mammogram w/Implants 7

screens

locked up or missing 132

non-functional, but cursor moves 145

Patient View 49, 55

Preview Image 42

Review dialog box 47

Select a Patient

functions available on 23

Standard Setup 71

search 27–33

for Patients 30

local database 30

Local Exams 30

Modality Worklist 27, 29–30

retrieve priors 92

select

AEC sensor position 40

exposure modes 36

exposure techniques 36

more than one procedure 59

patients in Today's Worklist 29

screening options 65

Select a Patient screen 23

chart of options 24

send

accepted rejected images to outputs 48

images to outputs post procedure 10, 83

sensor position, AEC 40

sequence for screening 65

set Center/Width 45

set exposure techniques 36

shutdown, how to 20

sign in as a different user 21

Sign Out 3

Simulate Capture, button 61

Sleep mode 1

SMPTE pattern 95

software

troubleshoot 144

version 95

specimens 67

exposure technique 9

spool

alarms for image status 130

deleting images from 114

jobs waiting when you exit 1, 21

Spool Management 110

Delete tab 114

manage space in image spool 110

menu item 110

mgr 110

Resend 111

density, adjusting for reprinting 111

using 110–111

Standard Setup screens 71

Start Up and Shut Down reference page 1

start up, TechMate 1

Start-up Wait Time 19

store image on removable media 79

switch users 3, 20

switches

circuit breakers 97

Emergency Off 21

Synchronization 26

system

controls and functional tests 100

default setup choices 71

message area 34

messages 19, 133, 143

paddle messages 9

power controls 97

ready indicator 66

start up and turn off 1

system log files, backing up to CD-R 122

T

tab, Delete 114

TEC

(Tissue Exposure Control), described 3

Accept mode 37

Exposure mode 63

flowchart 63

technique default options 71

using 37

tech, access level 105

TechMate, start up 1

techniques, exposure 34

temperature (IR) 95

test patterns 95

thumbnails 47, 49

Today's Worklist tab 29

troubleshooting 144

U

Unexposed Views, deleting 57

unreject an image 12, 47

UPS 98

beeping 97

power 97

reset 98

user interface 19, 70

User Setup function 72

users

access levels 105

add, edit and delete 72–73

identification, profiles, rights 105

log in different user 3

sign in different 21

V

View

add 56

icon 49

select 49

view

image status in output queues 87, 130

paddle configuration 123

prior images on AWS 92

rejected images 46, 118

Rejected Images Report 120

View Order Editor 72, 108

View Selection, change 53

W

Wait Time 19

Weekly Quality Control 14

protect "patient" 14

Window/Level 45

workflow

clinical procedures 65

optional operations 69

worklist 29–147

local 26

Modality 30

problems 146

resolve patient data conflicts 147

tab 27

write images to CD 80

X

X-ray

acquire image 66

indicators 66

change when in use 41

x-ray fields, collimated 102

Z

Zoom/Pan 44