



Instructions For Use

Selenia[®]

Digital Mammography System

Selenia[®]

Digital Mammography System

Instructions for Use

Part Number MAN-01384

Revision 002

March 2011

Technical Support:

USA: +1 877 371 4372

Europe: +32 2 711 4690

Asia: +852 37487700

All Other: +1 781 999 7750

HOLOGIC[®]

Corporate Headquarters

35 Crosby Drive,
Bedford, MA 01730-1401 USA

Tel: +1 781 999 7300

Sales: +1 781 999 7453

Fax: +1 781 280 0668

www.hologic.com

Europe

(EU Representative)

Hologic NV
Leuvensesteenweg 250A

1800 Vilvoorde, Belgium

Tel: +32 2 711 4680

Fax: +32 2 725 2087



Refer to the corporate website for more facilities worldwide.

© Copyright Hologic 2011. All rights reserved. Printed in USA. This manual was originally written in English.

Hologic and the Hologic Logo are trademarks or registered trademarks of Hologic, Inc. Other trademarks registered or used by Hologic and its divisions and subsidiaries in the United States and other countries include: Affirm, Dimensions, DSM, FAST Paddle, HTC, MIMS plus, M-IV, MultiCare, SecurView, Selenia, Smart Paddle, StereoLoc, and TechMate. Solaris, Sun, Sun Blade, Sun Ultra, and Ultra are trademarks or registered trademarks of Sun Microsystems, Inc. in the United States and other countries. Microsoft and Windows are trademarks or registered trademarks of Microsoft Corporation in the United States and other countries. Any other product and company names mentioned herein are the trademarks or registered trademarks of their respective owners.

Table of Contents

List of Figures	xi
List of Tables	xiii
Preface	xv
1.0 Intended Uses	xv
1.1 Indications For Use	xvi
1.2 Contraindications	xvi
1.3 Potential Adverse Effects	xvi
1.4 Product Complaints	xvi
1.5 Summary of Clinical Studies	xvi
2.0 Quality Control	xvi
3.0 User Training	xvi
4.0 Skills Needed for System Use	xvi
5.0 Terms and Definitions	xvii
6.0 Warnings, Cautions, and Notes	xix
7.0 International Symbols	xix
8.0 How to Use this Document	xx
Chapter 1 - General Information	1
1.0 Warnings and Precautions	1
2.0 Radiation Safety	5
2.1 Exposure Duration	5
2.2 Operator Radiation Shield	5
2.3 Patient Face Shield	5
2.4 Exam Room Door Indicators	5
3.0 Mechanical Safety	5
4.0 Emergency Off Switches	5
5.0 Interlocks	6
6.0 Compliance	6
6.1 Compliance Requirements	6
6.2 Compliance Statements	7
7.0 Label Locations	8
Chapter 2 - System Components and Controls	11
1.0 System Description	11
2.0 Film Printer	11
3.0 System Power Controls	12
3.1 Input Power Circuit Breakers	12
3.2 Emergency Off Switches	13
4.0 Acquisition Workstation Controls and Display	14
4.1 Acquisition Workstation	14
4.2 Keypads	15
4.3 Keyboard and Trackball or Mouse	15
4.4 Keyboard and Screen Controls	16
4.5 Bar Code Scanner	16
4.6 The Selenia Display	17

5.0 TechMate Controls and Displays 18

6.0 Tubestand Controls and Displays 19

 6.1 Gantry Controls and Displays 19

 6.2 C-arm Controls and Displays 20

 6.3 Dual Function Footswitches 21

Chapter 3 - Startup, Functional Tests, and Shutdown **23**

1.0 Procedure for Startup 23

2.0 Controls and Functional Tests 25

3.0 Monthly System Tests 29

4.0 How to Restart the System 29

 4.1 Restart From Sleep Mode 29

 4.2 Restart After the Emergency Off Switch was Activated 29

 4.3 How to Turn On or Reset the Uninterruptible Power Supply (UPS) 30

 4.4 How to Turn Off the UPS 30

 4.5 UPS Operation for TechMate 31

5.0 How to Log Out 32

6.0 How to Shut Down the System 32

Chapter 4 - The Selenia Accessories **33**

1.0 Introduction 33

2.0 Patient Face Shield 34

 2.1 How to Install the Face Shield on the Tubehead 34

 2.2 How to Remove the Face Shield 34

3.0 Compression Paddles 35

 3.1 How to Install the Compression Paddles 35

 3.2 How to Remove a Compression Paddle 35

 3.3 Maintenance and Cleaning 35

 3.4 Smart Paddle System 36

 3.5 FAST Paddle Use 36

 3.6 How to Realign the Paddle Front Edge 37

 3.7 Identification of Compression Paddles 38

4.0 Magnification Platform 40

 4.1 How to Install the Magnification Platform 40

 4.2 How to Remove the Magnification Platform 40

5.0 The Localization Crosshair Device 41

 5.1 How to Install the Crosshair Device 41

 5.2 How to Remove the Crosshair Device 41

 5.3 Localization Crosshair Device to the Localization Paddle Alignment 42

6.0 The Magnification Crosshair Device 42

Chapter 5 - The User Interface **43**

1.0 The Three Main User Interface Screens 43

2.0 A Screening Exam Sequence Suggestion 44

Chapter 6 - The Select a Patient Screen **47**

1.0 Introduction 47

2.0 Select a Patient 48

 2.1 How to Query the Modality Worklist Provider 48

 2.2 How to Search the Acquisition Workstation Database (Local Exams) 48

 2.3 More Information about Searches 48

3.0 How to Add a New Patient	49
4.0 How to Edit the Existing Patient Data	50
4.1 About Patient Information Edits with a Modality Worklist and PACS	50
4.2 How to Change the Patient Information After an Image Is Accepted	50
5.0 New, Edit, Add a Procedure	51
5.1 How to Add a New Procedure	51
5.2 How to Edit a Procedure	52
5.3 How to Add a Procedure	52
6.0 Image Output	53
7.0 Exposure Information	54
7.1 The Exposure Pane	54
7.2 About Exposure Modes	55
8.0 Menu Bar	59
8.1 Edit Menu Function	60
8.2 Admin Menu	65
8.3 Info Menu	72
9.0 Special Modes Buttons	72
9.1 Non-imaging Mode Button	72
9.2 The Review Button	72
9.3 The Simulate Capture Button	73
10.0 Generator Status	73
11.0 System Messages	73
12.0 Alarms	73
Chapter 7 - The Patient Views Screen	75
1.0 Work with Views	76
2.0 Acquire Images	77
Chapter 8 - The Preview Image Screen	79
1.0 Introduction	79
2.0 The Preview Image	79
2.1 Progressive Preview	79
2.2 Film Label Area	80
2.3 Hanging Options	80
3.0 Tools for Image Enhancement and Annotation	81
3.1 Comments.....	82
3.2 The Edit View Button and Markers	83
3.3 Implant Present	84
3.4 Quick Zoom/Pan	84
3.5 Full Zoom/Pan	84
3.6 ROI	85
3.7 Crosshairs	86
3.8 Measurement	87
3.9 Magnification	88
3.10 Window/Level	89
3.11 Display AEC Regions	90
3.12 Exposure Index	91
3.13 Exposure and Dose/Entrance Surface Exposure Information	91
4.0 Accept or Reject	91
4.1 Accept	91
4.2 Reject	92

5.0 Additional information about Accepted/Rejected Images	93
5.1 The Image Repetition Information Dialog Box	93
5.2 How to Accept a Rejected Image	93
5.3 How to Reject an Accepted Image	94
Chapter 9 - Maintenance and Cleaning	95
1.0 General Information	95
1.1 For General Cleaning	95
1.2 To Prevent Possible Injury or Equipment Damage	95
2.0 Care and Cleaning—Acquisition Workstation and Value Console	96
2.1 Display	96
2.2 Keyboard and Trackball or Mouse	96
3.0 Preventive Maintenance	97
Appendix A—System Specifications	99
1.0 Dimensional Information	99
1.1 Tubestand (Gantry with C-arm)	99
1.2 Acquisition Workstation	99
1.3 Acquisition Workstation with Dual Swivel Arms	100
1.4 Selenia Value Console	100
2.0 Operating Environment	101
2.1 General Operating Conditions	101
3.0 Storage Environment	101
3.1 Tubestand	101
3.2 Image Receptor	101
4.0 Electrical Input	101
4.1 Tubestand	101
4.2 Acquisition Workstation	101
5.0 Acquisition Workstation Technical Information	102
6.0 Tubestand Technical Information	103
6.1 C-arm	103
6.2 Compression	103
6.3 X-ray Tube: Molybdenum	104
6.4 X-ray Tube: Tungsten	106
6.5 X-ray Collimation	107
6.6 Light Field Indication	107
6.7 X-ray Generator	107
7.0 Imaging System Technical Information	110
7.1 Image Receptor	110
8.0 Tissue Exposure Control (TEC) Mode—Enhanced Manual Mode	111
9.0 Automatic Exposure Control (AEC)	111
10.0 Hardcopy Film Printing Devices	113
10.1 Image Engine	113
10.2 Interface	113
10.3 Printed Film	113
11.0 TechMate	113

Appendix B—The Mobile Selenia	115
1.0 General Information	115
1.1 The Vertical Position Override Switch (VPOS)	115
2.0 Safety Conditions and Other Precautions	116
3.0 Mobile Specifications	117
3.1 Shock and Vibration Limits	117
3.2 Coach Environment	117
3.3 Electrical Input	117
4.0 How to Prepare the System for Transport	118
5.0 How to Test the System Integrity After Transport	118
Index	119

List of Figures

Figure 1-1: Selenia Label Locations	8
Figure 1-2: Selenia Value Console Label Locations	9
Figure 2-1: Selenia System Description	11
Figure 2-2: Gantry Circuit Breaker	12
Figure 2-3: Acquisition Workstation Circuit Breaker	12
Figure 2-4: Emergency Off Switch	13
Figure 2-5: Acquisition Workstation Emergency Off Switch	13
Figure 2-6: Gantry Emergency Off Switches	13
Figure 2-7: Value Console Emergency Off Switch	13
Figure 2-8: The Acquisition Workstation Controls	14
Figure 2-9: The Selenia Value Console Controls	14
Figure 2-10: Keypad Controls	15
Figure 2-11: Selenia Keyboard	15
Figure 2-12: Selenia Function Keys on the Keyboard	16
Figure 2-13: Totoku Display	17
Figure 2-14: Barco Display	17
Figure 2-15: TechMate Controls and Indicators	18
Figure 2-16: Tubestand	19
Figure 2-17: Gantry Controls	19
Figure 2-18: Compression Device	21
Figure 2-19: The Dual Function Footswitch	21
Figure 3-1: Emergency Off Switch	23
Figure 3-2: Gantry Circuit Breakers	23
Figure 3-3: Acquisition Workstation Circuit Breaker	23
Figure 3-4: Power On Buttons	23
Figure 3-5: Operating System Login Screen	29
Figure 3-6: The Acquisition Workstation UPS Switch Lever	30
Figure 3-7: Value Console UPS Switch	30
Figure 3-8: Reset the TechMate UPS	31
Figure 3-9: The Sign Out Button	32
Figure 3-10: Exit from Acquisition Workstation Dialog Box	32
Figure 4-1: Installation of the Face Shield	34
Figure 4-2: Mount the Compression Paddle	35
Figure 4-3: Shifting Compression Paddle	36
Figure 4-4: Alignment of the Paddle	37
Figure 4-5: Alignment of Adjustment Screws	37
Figure 4-6: Magnification Platform Installation	40
Figure 4-7: Installation of the Crosshair Device	41
Figure 4-8: Crosshair Locking Levers	41
Figure 4-9: Adjustment Lock Screw	42
Figure 4-10: Installation of the Magnification Crosshair Device	42
Figure 6-1: The Select a Patient Screen	47
Figure 6-2: New Patient Entry Form	49
Figure 6-3: The Edit Patient Dialog Box	50
Figure 6-4: New Procedure, Edit a Procedure, and Add a Procedure Buttons	51
Figure 6-5: New Procedure Dialog Box	52
Figure 6-6: Add a New Procedure Dialog Box	52
Figure 6-7: Select the Output	53

Figure 6-8: Exposure Techniques	54
Figure 6-9: Breast Density Options.....	56
Figure 6-10: Accept TEC Exposure Techniques	56
Figure 6-11: The Override Mode.....	57
Figure 6-12: AEC Exposure and Sensor Position Adjustment On-Screen Settings	58
Figure 6-13: The Setup Screen	60
Figure 6-14: View Order Editor	61
Figure 6-15: The Edit Users Dialog Box.....	62
Figure 6-16: Edit Output Device Dialog Box	63
Figure 6-17: Choose a New Patient Dialog Box	65
Figure 6-18: Image Management Resend Options	66
Figure 6-19: Manage Queues.....	68
Figure 6-20: Test Pattern Dialog Box.....	71
Figure 6-21: Non-Imaging Mode Button.....	72
Figure 6-22: The Alarm Icon	73
Figure 7-1: Patient View Screen	75
Figure 7-2: Buttons Used on the View Screen	76
Figure 7-3: Ready for Exposure Indicators	77
Figure 8-1: The Preview Image Screen	79
Figure 8-2: Hanging Options.....	80
Figure 8-3: The Preview Screen Tools	81
Figure 8-4: Image Comments	82
Figure 8-5: The Markers in the Preview Pane	83
Figure 8-6: Full Zoom/Pan Image	84
Figure 8-7: ROI Size Drop-down Menu.....	85
Figure 8-8: Preview with Crosshairs	86
Figure 8-9: A Measurement on the Preview Screen	87
Figure 8-10: Preview with Magnification.....	88
Figure 8-11: The Window/Level Settings	89
Figure 8-12: AEC Regions	90
Figure 8-13: Exposure Information	91
Figure 8-14: Exposure Information	91
Figure 8-15: The Reject Reasons List.....	92
Figure 8-16: A Rejected Image	92
Figure 8-17: Set Accepted Button on Review Screen	93
Figure 8-18: Accept Rejected Image.....	93
Figure 8-19: Set Rejected Button	94
Figure A-1: Tubestand Dimensions	99
Figure A-2: Acquisition Workstation Dimensions	99
Figure A-3: Acquisition Workstation with Dual Swivel Arms	100
Figure A-4: Value Console Dimensions.....	100
Figure B-1: Location of Vertical Position Override Switch (VPOS).....	115

List of Tables

Table P-1: International Symbols	xix
Table 1-1: Factors That Limit Exposure Duration	5
Table 2-1: Selenia Function Keys	16
Table 3-1: System Startup Procedures	23
Table 3-2: Monthly Control Function Tests	25
Table 3-3: The Log Out Methods	32
Table 4-1: Available Accessories	33
Table 5-1: Workflow Suggestions to Select a Patient and Acquire an Image.....	44
Table 6-1: Legend for Figure 6-1	47
Table 6-2: Exposure Modes	55
Table 6-3: AEC Alarm Messages	58
Table 6-4: Menu Bar Options and Functions	59
Table 6-5: User Setup Options.....	62
Table 6-6: How to Use Manage Queues.....	68
Table 7-1: Actions on the Patient Views screen	75
Table 7-2: View Options	76
Table 9-1: User Preventive Maintenance	97
Table A-1: Collimation Settings	107
Table A-2: mA setting as a function of kV	108
Table A-3: mA setting as a function of kV	108
Table A-4: mA Factor as a function of mAs.....	108
Table B-1: Integrity Checklist.....	118

Preface

The Selenia[®], based on the M-IV[™] Mammography X-ray System, uses an Image Acquisition system that includes a digital image receptor. This receptor, which covers an area of 24 cm x 29 cm, is a direct-capture detector using an amorphous selenium photoconductor. At the Acquisition Workstation the user selects X-ray exposure technique factors, adds patient identification data, acquires, processes, and displays the digital images. The images are then either processed for printing and transmitted to peripheral hard copy laser film printer or sent to a diagnostic review workstation or both. Contrast and brightness are set automatically and can be adjusted by the user prior to printing or viewing. Hard-copy printers and a diagnostic review workstation are available as options.

The patient is prepared for the procedure in the same manner as for a film-based system. The Acquisition Workstation provides direct digital acquisition system advantages in system efficiency and productivity including:

- Patient demographics are associated with the electronic image.
- Technologist comments can be associated with the image.
- The technologist is assisted through the steps of image acquisition.
- Possible image destinations are provided via configurable settings.
- The technologist may review past acquisitions (including priors), to see previous captures, and then resend, and/or remark them.

1.0 Intended Uses

RxOnly United States Federal Law restricts this device to use by or on the order of a physician.

The Selenia is intended to produce radiographic images of the breast. Its intended use is for the production, storage, and diagnostic review of digital screening, diagnostic, and needle localization mammography.

The system is to be used in a radiology or clinic exam room environment in a hospital, outpatient clinic, or a breast imaging center. Mammography technologists operate the system for the production and storage of digital mammograms. The system may also be used for quality control purposes and other clinical or research related activities by medical physicists and radiologists certified in accordance with MQSA standards.

- The Selenia Acquisition Workstation display is not approved for final interpretation of examinations. Final interpretations should be done from either films or with the diagnostic review workstation. Images shown on the Acquisition Workstation display are for quality assurance or confirmation purposes only.
- Only images produced by recommended laser printers, or an approved diagnostic review workstation, should be used for final interpretation of examinations. For compatible printers, see the latest product data sheets for the system, which can be obtained from Hologic[®] or your sales representative.

1.1 Indications For Use

The Selenia Full Field Digital Mammography System generates digital mammographic images that can be used for screening and diagnosis of breast cancer. The Selenia Full Field Digital Mammography System is intended for use in the same clinical applications as traditional screen-film mammographic systems. Mammographic images can be interpreted on either hard copy film or the diagnostic review workstation.

1.2 Contraindications

There are no known contraindications.

1.3 Potential Adverse Effects

The following is a list of potential adverse effects that apply to mammography and are also applicable to digital mammography using the Selenia.

- Excessive breast compression
- Excessive X-ray exposure
- Electric shock
- Infection
- Skin irritation, abrasions, or puncture wounds

1.4 Product Complaints

Report any complaint or dissatisfaction in the quality, durability, reliability, safety, effectiveness, and/or performance of this product to Hologic. If the device has caused or contributed to a serious injury of a patient, immediately report the incident to Hologic. (See the Title Page for contact information.)

1.5 Summary of Clinical Studies

This information is available in the document MAN-02233.

2.0 Quality Control

Perform all Quality control tests within the correct time frame.

3.0 User Training

Hologic does not accept the responsibility for injury or damage caused by incorrect system operation. Hologic can arrange for training by a clinical application specialist. Refer to this *Instructions for Use* manual for directions on how to use the Selenia.

4.0 Skills Needed for System Use

You must understand how to use a personal computer system, internal and external storage media (DVDs, CDs, diskettes, etc.), and external computer peripherals (keyboard, trackball, printers).

5.0 Terms and Definitions

Accession Number	A DICOM term that refers to a RIS-created number that uniquely identifies a visit to a site by a patient.
ACR	American College of Radiology
Ag	Silver
AEC	Automatic Exposure Control. A method that limits the amount of radiation a patient receives.
collimator	Device at the X-ray tube that limits the area of the receptor that is exposed.
dialog box	Device at the X-ray tube that limits the area of the receptor that is exposed.
DICOM	Digital Imaging and Communications in Medicine. An industry standard specification for communication between medical imaging equipment.
digital marker	A software mechanism that marks an image to indicate some information (normally orientation).
FAST Paddle™	Fully Automated Self-Adjusting Tilt Paddle system
Grid	Element within the Digital Image Receptor that reduces the scatter radiation during exposure
HIS/RIS	Hospital Information System/Radiology Information System. Generic term for non-PACS systems which track the Patient demographics and ordered radiological studies
HTC™	High Transmission Cellular Grid
kV	Kilovolts: One of the x-ray exposure settings.
login/logout	The process of logging into and out of the Operating System of the Acquisition Workstation
LUT	Look-Up Table (LUT). An image processing function that replaces one image pixel value with a different image pixel value.
Mag	Magnification
mA	Milliamperes. One of the x-ray exposure settings.
mAs	Milliampere-seconds. An electrical term used in x-ray exposure settings
Mo	Molybdenum
Modality Worklist (MWL)	A list of scheduled procedures normally kept by a RIS or PACS.
MPPS	Modality Performed Procedure Step. A DICOM service to allow the RIS (or another device) to know about work performed on the AWS.
MQSA	Mammography Quality Standards Act
Operating System (OS)	The software control system which runs all functions of a computer.
Outputs	A list of devices to which the accepted image is sent. The devices can be a combination of archives, workstations and film printers.

PACS	Picture Archiving and Communications System. A computer and network system for the transfer and archive of digital medical images
PPS Status	The condition of a Performed Procedure Step being “in progress”, “completed”, or “discontinued”.
Procedure	A generic medical protocol which contains a set of images (Views) which are acquired under certain conditions, and are performed together for a one purpose (for example standard screening). There is no Procedure instance UID item because a Procedure is a generic item in DICOM. DICOM supports the identification of requested Procedures.
Rh	Rhodium
ROI	Region of Interest
RIS	Radiology Information System
Series	A set of images acquired by a single tech for a single Patient and Procedure on a particular modality with a fixed body part, laterality and view position. DICOM uniquely identifies the series with a globally unique instance UID. Based on this description, each individually acquired DR image on the Acquisition Workstation is also an individual series.
Sign-in	The process of user identification to the Acquisition Workstation application.
Sign-out	The process in which a user exits the Acquisition Workstation application, but the user does not logout of the OS.
Smart Paddle™	A paddle with a release knob on each side which allows the paddle to move from one side of the image detector to the other. The system identifies a Smart Paddle installed on the compression device.
Solaris™	The company Sun™ Microsystems’ version of the UNIX Operating System.
TEC	Tissue Exposure Control mode, an enhanced Manual Exposure Control (MEC) mode.
TechMate™	The Hologic SecurView _{RT} technologist review workstation accessory, as packaged with the Selenia to reduce space in the examination room.
technique	Combination of X-ray parameters (kV, mA, etc.) for a specified view in a procedure.
UPS	Uninterruptible Power Supply.
View	The combination of a single X-ray image and a specified set of conditions under which the image was acquired. The View is not part of DICOM nomenclature, but in the context of DR, is approximately synonymous with a DICOM image object.
W	Tungsten
Wildcard Character	A keyboard character that represents one or many characters, for example the * or ?.

6.0 Warnings, Cautions, and Notes

Descriptions of Warnings, Cautions, and Notes used in this manual:



WARNING!

Procedures that you must follow accurately to prevent possible dangerous or fatal injury.



Warning:

Procedures that you must follow accurately to prevent injury.



Caution:

Procedures that you must follow accurately to prevent the damage to equipment, loss of data, or damage of files in software applications.














Note...

Notes indicate additional information.

7.0 International Symbols

This section explains the International Symbols used on this system.

Table P-1: International Symbols

	Potential Equalization terminal	Connection for a conductor other than the Protective Earth for a direct connection between two or more pieces of electrical equipment.
	Protective Earth terminal	Used for the connection of the ground of the line cord or ground cable of the equipment and no other purpose.
	Off	Power disconnection from the mains
	On	Power connection to the mains
	Off	Off, only for a part of equipment
	On	On, only for a part of equipment
	WEEE	Symbol indicating separate collection for electrical and electronic equipment
	Dangerous Voltage	Identifies area of potentially lethal voltage
	Manufacturer	
	Date of Manufacture	
	X-ray Radiation	Caution—Radiation

8.0 How to Use this Document

- The Preface describes intended system use and identifies the skills a system user needs.
- Chapter 1 provides general safety information.
- Chapters 2, 3, and 4 describe system hardware, system startup/shutdown and functional tests.
- Chapter 5 steps you through a screening exam and the user interface.
- Chapters 6, 7, and 8 describe the options on the user interface screens while logged on as a “tech”.
- Chapter 9 describes system Maintenance and Cleaning procedures.
- Appendix A states specifications for the system.
- Appendix B provides information about Selenia in a mobile environment.

Chapter 1 - General Information

1.0 Warnings and Precautions

Read and understand this manual before you use the system. Always follow all the instructions in this manual.

This system is for use by qualified Operators only. The Operators must have experience in the principles of mammography before following the instructions contained in this manual.



The Selenia system is classified as CLASS I, TYPE B permanently connected equipment as per IEC 60601-1. There are no special provisions to protect the system from flammable anesthetics or ingress of liquids.



WARNING!

Do not open any of the panels. Lethal voltages are present within the interior of this unit.



WARNING!

Do not use the electrical equipment near flammable anesthetics.



WARNING!

The user or the service personnel must correct problems before the system is used.



WARNING!

The user must arrange for preventive maintenance by an authorized service representative.



WARNING:

This X-ray unit can be dangerous to patient and the user. Always follow the safety precautions for x-ray exposures.



WARNING!

After a power failure, remove the patient from the system before you apply power.



WARNING!



To keep the isolation quality for the system, attach only approved accessories or options to the system. Only the authorized personnel can make changes to the connections.



WARNING!

Keep a 1.5 meter safe distance between the patient and any non-patient devices.

Non-patient system components (like the Workflow Manager, the diagnostic review workstation, or the hardcopy printer) must not be installed in the Patient Area



WARNING!

Per North American electrical safety requirements, grounding reliability can only be achieved when the Acquisition Workstation is connected to a receptacle marked Hospital Grade.



WARNING!

The disk drives installed in this system are a Class I Laser Product. Prevent direct exposure to the beam. Hidden laser radiation exists if the case to a disk drive is open.



Warning:

This unit is for use only by qualified personnel.



Warning:

Control the access to the equipment according to local regulations for radiation protection.



Warning:

You must keep your complete body behind the radiation shield for the time of the exposure for maximum protection from x-ray exposure.



Warning:

This device contains dangerous material. Return to Hologic all material removed from service.



Warning:

Motorized equipment. Use care when adjusting for patient use. Observe equipment and patient at all times during set up. If a chair is necessary use an adjustable chair set above its minimum height.



Warning:

Put the footswitches where a patient user, or wheelchair cannot accidentally activate a switch.



Warning:

Never leave the patient during the procedure if in contact with the mammography unit.



Warning:

Keep the hands of the patient away from all buttons and switches at all times.



Warning:

You increase the patient dose to high levels if you increase the AEC exposure adjustment setting. You increase the image noise or decrease image quality if you decrease the AEC exposure adjustment setting.



Warning:

Increased Exposure adjustment leads to higher dose to the patient. Keep the patient dose as low as practical to get good image quality.



Warning:

If a paddle touches possible infectious materials, call your Infection Control Representative for decontamination instructions.



Warning:

The bar code scanner installed in this system is a Class II Laser Product. Prevent direct exposure to the beam. Hidden laser radiation exists if the cover is opened.



Warning:

The Face Shield does not protect from radiation.



Caution:

Do not use any heat source on the image receptor.

Instructions for Use

Chapter 1 - General Information
Warnings and Precautions



Caution:

Never turn off the Acquisition Workstation Circuit Breaker except in emergency. The circuit breaker can turn off the Uninterruptible Power Supply (UPS) and risk data loss.



Caution:

To minimize possible damage from thermal shock to the Digital Image Receptor, follow the recommended procedure to turn off the equipment.



Caution:

Do not put any magnetic media near or on devices that create any magnetic fields, because stored data can be lost.



Caution:

This system is a medical device and not a normal computer. Do not make changes to the hardware or software that are not authorized. Install this device behind a firewall for network security. The computer virus protection or network security for this medical device is provided (for example, a computer firewall). The network security and anti-virus provisions are the responsibility of the user.



Caution:

Software other than those provided specifically for use with this system must not be loaded onto the system.



Caution:

Only use the approved accessories with this equipment. The failure to follow this caution can cause errors and possible data loss.

2.0 Radiation Safety

2.1 Exposure Duration

Table 1-1: Factors That Limit Exposure Duration

Normal Factors	Abnormal Factors
Manual mAs Timer	<ul style="list-style-type: none"> • Release of the x-ray button before the exposure ends. • Exposure longer than the preset “back-up time” • The independent safety hardware back-up timer. • Detection of a generator fault

2.2 Operator Radiation Shield

The radiation shield meets the requirements of 21CFR 1020 and IEC 60601-2-45.2001. The shield has a rating of a 0.5 mm Pb (lead) equivalent and is permanently installed.



Note... Verify the shield integrity every day before use.

2.3 Patient Face Shield

The Patient Face Shield keeps the face of the patient out of the radiation field.

2.4 Exam Room Door Indicators

- The system provides for an external Power-On indicator.
- The system provides for an external X-ray Exposure In Progress indicator.

3.0 Mechanical Safety

The equipment meets the requirements of IEC 60601-1, UL 60601, and CSA 22.2 No. 601.1 by these safety features:

- The C-arm rotation braking is ensured upon loss of power.
- Compression release is disabled when a localization paddle is installed.
- Motorized compression force has a limit of 200 N (45 lb).

4.0 Emergency Off Switches

The Emergency Off switches remove the power from the Gantry. Do not normally use the Emergency Off switches to turn off the system. Refer to Chapter 2, page 13.

5.0 Interlocks

- Display of 58 Newtons (13 lb) or greater of compression force disables the C-arm vertical drive and rotation.
- Installation of a Localization Paddle disables the automatic compression release functions.
- Release of an X-ray button before the end of the exposure ends the exposure.
- Misalignment of the Light Field Mirror or the Filter prevents an X-ray exposure.
- Installation of an 18 x 24 cm shifting paddle prevents the X-ray exposure if the detected paddle position does not match the View.

6.0 Compliance

This section describes the mammography system compliance requirements and the manufacturer responsibilities.

6.1 Compliance Requirements

The manufacturer is responsible for the effects of safety, reliability, and performance of this equipment with the following provisions:

- The electrical installation of the room complies with the appropriate requirements.
- The equipment is used in accordance with instructions for use.
- Assembly operations, extensions, re-adjustments, modifications, or repairs are performed by authorized persons only.
- The installed network and communications equipment must comply with an IEC Standard, and the complete system (network/communications equipment and Selenia Mammography System) must be installed to comply with IEC 60601-1 and IEC 60601-1-1.



Caution:

Medical Electrical Equipment needs special precautions about EMC and must be installed, put into service and used according to the EMC information provided.



Caution:

Portable and Mobile RF communications can affect Medical Electrical Equipment.



Caution:

The use of unauthorized accessories and cables can result in increased emissions or decreased immunity. To keep the isolation quality for the system, attach only approved accessories or options to the system.



Caution:

The Medical Electrical (ME) Equipment or ME System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the ME Equipment or ME System should be observed to verify normal operation in the configuration in which it is used.



Caution:

This system is intended for use by healthcare professionals only. This system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the equipment or shielding the location.



Caution:

Changes or modifications not expressly approved by Hologic could void your authority to operate the equipment.

6.2 Compliance Statements

The manufacturer states that this device is manufactured/conforms to:

- EN ISO 13485:2007
- FDA, 21 CFR [Parts 820, 900 and 1020]
- IEC 60601-1:1988 +A1+A2:1995 +A13:1996 Medical Electrical Equipment—General requirements for safety
- IEC 60601-1-1:2000-12 Medical Electrical Equipment—Collateral Standard: Safety requirements for medical electrical systems
- IEC 60601-1-2:2001 Medical Electrical Equipment—Collateral Standard: Electromagnetic compatibility for medical electric systems
- IEC 60601-1-3:1994 Medical Electrical Equipment—Collateral Standard: Requirements for radiation protection in diagnostic x-ray equipment
- IEC 60601-1-4:1996 +A1:1999 Medical Electrical Equipment—Collateral Standard: Programmable electrical medical systems
- IEC 60601-1-6 Medical Electrical Equipment—Collateral Standard: Usability
- IEC 60601-2-28:1993-03 Medical Electrical Equipment—Particular requirements for the safety of x-ray source assemblies and x-ray tube assemblies for medical diagnosis
- IEC 60601-2-32:1994 Medical Electrical Equipment—Particular requirements for the safety of associated equipment of x-ray equipment
- IEC 60601-2-45:2001 Medical Electrical Equipment—Particular requirements for the safety of mammographic x-ray equipment and mammographic stereotactic devices
- IEC 62304:2006 Medical Device Software—Software life cycle processes
- UL 60601-1: Medical Electrical Equipment, Part 1—General Requirements for Safety
- CSA: Medical Electrical Equipment Part 1: C22.2 No. 601.1-M90—General Requirements for Safety
- ISO 14971

7.0 Label Locations




Figure 1-1: Selenia Label Locations




Figure 1-2: Selenia Value Console Label Locations

Instructions for Use

Chapter 1 - General Information
Label Locations

Selenia®
Digital Mammography System

Chapter 2 - System Components and Controls

1.0 System Description

The Selenia is available with an Acquisition Workstation or a Value Console. The Acquisition Workstation can also include a technologist review workstation accessory (the Hologic SecurView_{RT} TechMate™) packaged with the Selenia to reduce space requirements.




Figure 2-1: Selenia System Description

Legend for Figure 2-1



- The **Selenia Acquisition Workstation** (1) and the **Value Console** (2) contain the image processing electronics and provide the user interface.
- The **SecurView_{RT} TechMate** (3) provides a Technologist review workstation accessory.
- The **Tubestand** is the **C-arm** (4) and the **Gantry**(5).
The **C-arm** provides a platform for the X-ray tube, compression systems, and the digital image receptor. A pivot mechanism connects it to the Gantry.
The **Gantry** contains the electrical and mechanical subsystems for the Selenia.
- Dual-function **Footswitches** (6) allow hands-free C-arm vertical travel and compression movements.

2.0 Film Printer

Selenia accepts film printers. Printed films display the patient information, exposure techniques, projection, and facility information.

3.0 System Power Controls

3.1 Input Power Circuit Breakers

Location of Circuit Breaker	Description
 <p data-bbox="483 751 875 785"><i>Figure 2-2: Gantry Circuit Breaker</i></p>	<ul style="list-style-type: none"> • At the lower right corner of the back panel • Removes the power for service or in an emergency • Up is turned On. Down is turned Off. • When turned off, a system restart is required. Restart the system following Chapter 3, Section 1.0, page 23.
 <p data-bbox="388 1186 971 1220"><i>Figure 2-3: Acquisition Workstation Circuit Breaker</i></p>	<ul style="list-style-type: none"> • On the back of the Acquisition Workstation behind the shield. The Value Console Circuit Breaker (not shown) is on the back of the Value Console. • Use only in an emergency. • If the circuit breaker was used, wait until the beeps stop, then turn on the circuit breaker. The system automatically restarts. If the UPS does not respond automatically, reset the UPS. See Chapter 3, Section 4.3, page 30.

3.2 Emergency Off Switches

- Press any Emergency Off Switch to remove Gantry power.
- Rotate the switch one-quarter turn to reset.
- Restart the system after any Emergency Off Switch is pressed. Refer to procedure in Chapter 3, Section 4.2, page 29.




Figure 2-4: Emergency Off Switch

Location of Emergency Off Switches




Figure 2-5: Acquisition Workstation Emergency Off Switch




Figure 2-7: Value Console Emergency Off Switch




Figure 2-6: Gantry Emergency Off Switches

4.0 Acquisition Workstation Controls and Display

4.1 Acquisition Workstation




Figure 2-8: The Acquisition Workstation Controls

Legend for Figure 2-8

1. Circuit Breaker
2. Emergency Off Switch
3. Selenia Display
4. Trackball
5. CD-RW Drive
6. Left and Right Keypads
7. Selenia Acquisition Workstation UPS (Uninterruptible Power Supply) Reset Button (on lower right back)
8. Bar Code Scanner
9. Selenia Keyboard




Figure 2-9: The Selenia Value Console Controls

Legend for Figure 2-9

1. Computer
2. Bar Code Scanner
3. Left and Right Keypads
4. Keyboard
5. Display
6. Emergency Off Switch
7. Mouse
8. Circuit Breaker
9. CD-RW Drive
10. UPS

4.2 Keypads

The Keypads provide **Power On**, **X-ray**, and **Compression Release** functions.




Figure 2-10: Keypad Controls

Legend for Figure 2-10

1. Power On button.
2. Power On indicator light.
3. X-ray Indicator lights on both panels.
4. X-ray buttons, one on each side of the keyboard. Press both at the same time to begin an x-ray exposure and hold until the tone stops.
5. Compression Release button.

4.3 Keyboard and Trackball or Mouse




Figure 2-11: Selenia Keyboard

Legend for Figure 2-11

1. Selenia Function keys
2. Trackball (Acquisition Workstation)
3. Future Use
4. Standard keyboard keys
5. Standard arrow keys
6. Standard number keypad
7. Mouse (Value Console)

4.4 Keyboard and Screen Controls

Use the trackball or the function keys to change Exposure techniques.




Figure 2-12: Selenia Function Keys on the Keyboard



Table 2-1: Selenia Function Keys

Key	Name	Function
	Exposure Mode	Changes AEC, TEC, and Manual exposure modes.
	Filter	Changes the Filters. The Filter options depend on the X-ray tube. The Molybdenum system has Mo and Rh. The Tungsten system has Rh and Ag.
	kV Up	Increases kV.
	kV Down	Decreases kV.
	mAs Up	Increases mAs.
	mAs Down	Decreases mAs.
	Grid	Changes the Grid Position (In or Out).
	Focal Spot	Changes the Focal Spot size (Small or Large).
	Breast Density	Changes the TEC Breast Density settings.
	AEC position	Moves the AEC Sensor position.
	New Entry	Opens New Patient Entry dialog box.
	Close Procedure	Closes the Procedure. (Disabled when MPPS is installed).


4.5 Bar Code Scanner

The Bar Code Scanner speeds the entry of data from records with bar codes.

4.6 The Selenia Display

Features	Display
<p>Totoku Display</p> <ul style="list-style-type: none"> • Has a Power button (2)/On-screen Display Adjustment dial(1). • Green indicator lights when Display is turned on. • Has the ability to tilt. 	 <p>Figure 2-13: Totoku Display</p>
<p>Barco Display</p> <ul style="list-style-type: none"> • Meets the DICOM requirements. • Press the recessed dial (on the right side) to display a menu. • Press and hold the recessed dial to turn the Display On/Off. • Green indicator lights when Display is turned on. • Has the ability to tilt. • The Value Console uses only the Barco display. 	 <p>Figure 2-14: Barco Display</p>

5.0 TechMate Controls and Displays



Legend for Figure 2-15

1. TechMate Display
2. TechMate Keyboard
3. TechMate Trackball
4. Power Button for TechMate
5. TechMate UPS Reset

Figure 2-15: TechMate Controls and Indicators

6.0 Tubestand Controls and Displays

Tubestand Components

1. The Gantry
2. The C-arm




Figure 2-16: Tubestand

6.1 Gantry Controls and Displays

Gantry Controls

1. The *C-arm Rotation Angle Displays*, on both sides of the Gantry, show the C-arm angle.
2. The *Emergency Off Switches*, on each side of the Gantry, remove power from the Tubestand.
3. *Gantry Rotation Switches* (if installed), on each side of the Gantry, move the C-arm to a programmed position.
4. The *Input Power Circuit Breaker*, at the rear of the Gantry, provides the overload protection.
5. *Dual-Function Footswitches* activate C-arm and Compression Device movement.





Figure 2-17: Gantry Controls


Instructions for Use

6.2 C-arm Controls and Displays


1. Tubehead Controls
2. Compression Device
3. Rear of the C-arm Side Rails
4. Rear of the Image Receptor




The *Tubehead Controls*, on either side of the X-ray tube side covers, provide the Collimator Override, C-arm Rotation, and Light Field functions.



The *Compression Device* contains the AEC Sensor control and position LEDs, and displays compressed breast thickness and force. *Handwheels* allow manual compression control.



The rear edges of the *C-arm Side Rail* have recessed push-button controls for motorized compression and C-arm movement. An additional C-arm Rotation button is available on the rear of the *Image Receptor*.



6.2.1 Compression Device Controls and Displays

Legend for Figure 2-18

1. Compression Handwheels
2. AEC Sensor Position controls
3. AEC Sensor Position Display
4. Compression Force Display
5. Compression Thickness Display




Figure 2-18: Compression Device

6.3 Dual Function Footswitches

1. Press the footswitch to activate.
2. Release the footswitch to stop motorized movement.

Legend for Figure 2-19

1. C-arm Down
2. C-arm Up
3. Compression Down
4. Compression Up




Figure 2-19: The Dual Function Footswitch



Warning:

Put the footswitches where a patient user, or wheelchair cannot accidentally activate a switch.

Instructions for Use





Chapter 2 - System Components and Controls
Tubestand Controls and Displays

Selenia®
Digital Mammography System

Chapter 3 - Startup, Functional Tests, and Shutdown

1.0 Procedure for Startup

Table 3-1: System Startup Procedures

Step	Description
1. Reset the Emergency Off Switches. <ul style="list-style-type: none"> • Two on the Gantry • One on the Acquisition Workstation (or Value Console) 	 <p>Figure 3-1: Emergency Off Switch</p>
2. Make sure that the circuit breakers are turned On.	 <p>Figure 3-2: Gantry Circuit Breakers</p>  <p>Figure 3-3: Acquisition Workstation Circuit Breaker</p>
3. Perform Pre-Startup Checks.	<ol style="list-style-type: none"> 1. Look for open or loose panels, missing hardware, and indications of damage. 2. Inspect the radiation shield for chips, cracks, breaks, and for tight attachment. 3. Inspect the paddles for small cracks. 4. Remove any obstructions to the Operator view and C-arm movement.
4. Turn On the Selenia. Turn on the TechMate. (If the Acquisition Workstation or TechMate does not start, reset the UPS for that unit. See Section 4.3, page 30.)	 <p>Figure 3-4: Power On Buttons</p> <ol style="list-style-type: none"> 1. Press the Power On button (No. 1). The green LED (No. 2) illuminates. 2. Press the Power On button (No. 3) and the TechMate turns on. 3. Allow the time for initialization and diagnostic tests to complete.

Instructions for Use

Chapter 3 - Startup, Functional Tests, and Shutdown Procedure for Startup

Table 3-1: System Startup Procedures

Step	Description
5. Log on to the Operating System.	When the Logon dialog box appears: <ol style="list-style-type: none">1. Enter your user name (case sensitive) for the Operating System.2. Select OK or press the Enter key.3. Enter your Operating System password. (No characters appear in the field.)4. Select OK or press the Enter key.
6. Log in to the Application.	When the Login to System dialog box appears: <ol style="list-style-type: none">1. Select your ID from the drop-down list.2. Enter your password (case sensitive). (Asterisks appear in the field.)3. Select OK or press the Enter key.



Note...

If the system remains on over the night, reboot the system daily to resync system communications.



Note...

If the system was turned off for 30 minutes or longer, the detector temperature needs time to adjust. Allow a minimum of 1 hour before you acquire images on a patient. The Launch dialog box disappears when the wait time is finished.

If it has been less than 30 minutes, click Dismiss in the Launch dialog box.



Note...

If during the internal checks, the system detects a fault condition, a message appears, and startup is suspended until the problem is remedied.

2.0 Controls and Functional Tests





Perform these functional tests monthly.



Note...

The C-arm movement is disabled when compression force of 58 N (13 lb.) or greater is applied.

Table 3-2: Monthly Control Function Tests

Function	Control	Test
<p>Counterclockwise C-arm Rotation</p> 	<p>C-arm and Gantry Rotation Switches</p> 	<p>Press and hold the C-arm button and confirm the following actions:</p> <ul style="list-style-type: none"> • C-arm rotates • Angle Display changes • Rotation stops when you release the button. <p>Press and hold a Gentry Counterclockwise button and confirm the following actions:</p> <ul style="list-style-type: none"> • Beeps before the C-arm rotates • C-arm counterclockwise movement stops when the pre-programmed rotation angle is reached. • Press and hold a Center button and confirm C-arm moves to zero degree position.
<p>Clockwise C-arm Rotation</p> 	<p>C-arm and Gantry Rotation Switches</p> 	<p>Press the button and confirm the following actions:</p> <ul style="list-style-type: none"> • C-arm rotates • Angle Display changes • Rotation stops when you release the button. <ul style="list-style-type: none"> • Press and hold a Gentry Clockwise button and confirm the following actions: <ol style="list-style-type: none"> a. Beeps before the C-arm moves. b. C-arm clockwise movement stops when the pre-programmed rotation angle is reached. • Press and hold a Center button and confirm C-arm moves to zero degree position.

Instructions for Use

Chapter 3 - Startup, Functional Tests, and Shutdown
Controls and Functional Tests

Table 3-2: Monthly Control Function Tests













Function	Control	Test
<p>Compression Down</p> 		<ol style="list-style-type: none"> 1. Press the button. 2. Confirm the compression brake engages and the light field lamp illuminates. 3. Confirm the Compression Down movement stops: <ul style="list-style-type: none"> • When you release the button. • At the compression down force limit. • At the lower compression travel limit.
<p>Compression Up</p> 		<ol style="list-style-type: none"> 1. Press the button. 2. Confirm the compression brake does not release. 3. Confirm the Compression Up movement stops: <ul style="list-style-type: none"> • When the Compression device reaches the upper compression travel limit. • On release of the button.
<p>Compression Release</p> 		<ol style="list-style-type: none"> 1. Press the button. 2. Confirm the following actions: <ul style="list-style-type: none"> • The compression motor brake releases. • The compression device lifts approximately 10 cm. <p><i>Note... The system disables all compression release functions when a localization paddle is installed.</i></p>







Table 3-2: Monthly Control Function Tests

Function	Control	Test
<p>C-arm Up</p> 		<ol style="list-style-type: none"> 1. Press the button. 2. Confirm the C-arm Up movement stops: <ol style="list-style-type: none"> a. When you release the button. b. When the C-arm reaches the upper travel limit. 3. Apply more than 58 N (13 lb) compression and confirm the C-arm does not move.
<p>C-arm Down</p> 		<ol style="list-style-type: none"> 1. Press the button. 2. Confirm the C-arm Down movement stops: <ul style="list-style-type: none"> • When you release the button. • When the C-arm reaches the lower travel limit.
<p>Light Field Lamp</p> 		<ol style="list-style-type: none"> 1. Press the button. Confirm the light turns on for a short time. 2. Confirm the light field lamp turns on when a Compression Down button is activated.

Instructions for Use

Chapter 3 - Startup, Functional Tests, and Shutdown
Controls and Functional Tests

Table 3-2: Monthly Control Function Tests

Function	Control	Test
<p>Collimator Override</p> 		<ol style="list-style-type: none"> 1. Press the light field button then the collimator button. 2. Confirm the Collimator moves to the next field size. 3. Repeat the steps 1 and 2 until you move through all Collimator field sizes.
<p>Smart Paddle</p> 		<ol style="list-style-type: none"> 1. Install a Smart Paddle. 2. Rotate and hold the knob on the side of the paddle to move the paddle into a detent positions. 3. Release the knob to lock the paddle in position. 4. Turn on the light field lamp. 5. Confirm the collimator position matches the paddle position. 6. Repeat this procedure for the other two paddle positions. 7. Apply compression and confirm the paddle does not unlock.
<p>AEC Position</p> 		<ol style="list-style-type: none"> 1. Select an AEC mode. 2. Change the AEC position with the compression device arrows. 3. Confirm the indicator positions on the sides of the device and Acquisition Workstation display match the current position. 4. Move past position 7 or 1. 5. Confirm the Auto position indicator on the compression device illuminates and the display indicates Auto. 6. Select the AEC arrows on the display and repeat steps 3 to 5. 7. Use the keyboard function key and repeat steps 3 to 5.

3.0 Monthly System Tests

- Test all Emergency Stop Switches. Follow the on-screen reset procedure.
- Change the AEC Sensor Position switch, and make sure the compression device indicator and the Acquisition Workstation display change.
- Test all buttons on all keypads and footswitches. Include the compression release on the Acquisition Workstation.
- During Compression Device tests make sure the thickness and force displays change.
- During C-arm rotation tests, confirm that the LED Angle displays change.

4.0 How to Restart the System

4.1 Restart From Sleep Mode

1. To activate the system, move the trackball or mouse.




Figure 3-5: Operating System Login Screen

2. In the Login dialog box, enter reboot for the user name and password.
3. Wait for the system to restart.
4. Log on to the Operating System, then Log on to the Application software.
5. If the system was off for less than 30 minutes, select **Dismiss** in the **Launch** dialog box.

4.2 Restart After the Emergency Off Switch was Activated

1. Turn the Emergency Off switch by one-quarter turn clockwise to reset the switch.
2. When the Power Panel dialog box displays, "**Communication with the generator is lost,**" examine the Emergency Off Switches and the circuit breaker:



Caution:

Do not try to close this dialog box using the 'x' in the upper corner.

- If the Gantry power is turned on, select the **Cancel** button.
 - If the Gantry power is turned off, and the Emergency Off Switches and the circuit breaker are set correctly, select the **OK** button.
3. When the "**Link established with generator**" dialog box appears, select **OK**.

4.3 How to Turn On or Reset the Uninterruptible Power Supply (UPS)



Figure 3-6: The Acquisition Workstation UPS Switch Lever

1. Find the Acquisition Workstation UPS switch lever on the lower, right of the back panel.
2. Press the switch lever down and hold for two seconds or until the UPS beeps.
3. Start the Acquisition Workstation normally.




Figure 3-7: Value Console UPS Switch

1. Find the Value Console UPS switch, on the front of the UPS.
2. Press and hold the UPS switch for two seconds or until the UPS beeps.
3. Start the Console normally.

4.4 How to Turn Off the UPS

1. Make sure the Acquisition Workstation and TechMate or the Value Console are off. The UPS system begin to chirp repeatedly and can now be shut down.
2. Find the UPS switch. (See Figure 3-6 or Figure 3-7).
3. For the Acquisition Workstation, press and hold the switch lever down for two seconds or until the UPS beeps. For the Value Console, press and hold the UPS switch for two seconds or until the UPS beeps.
4. For the TechMate, perform steps 1 and 2 in Section 4.5, page 31.
5. When the power is restored, press the UPS switch one time to turn on the UPS.

4.5 UPS Operation for TechMate

To reset the TechMate UPS:




Figure 3-8: Reset the TechMate UPS



1. Find the UPS access hole in the middle of the left side of the front cover. This hole aligns with the Power button on the UPS.
2. Put the eraser-end of a pencil into the hole to press the UPS button until the UPS beeps.
3. Press the button again and listen for a second beep.
4. Start the TechMate normally.

To turn off the TechMate UPS, follow steps 1 and 2.

5.0 How to Log Out

There are two methods to log out of the application. This leaves the system on and another user can log in to use it.

Table 3-3: The Log Out Methods

Method 1	Method 2
<ol style="list-style-type: none"> 1. Select the Sign Out button.  <p>Figure 3-9: The Sign Out Button</p> <ol style="list-style-type: none"> 2. Select the Exit button. 3. Select Yes. 	<ol style="list-style-type: none"> 1. Select File>Exit. 2. Select Log out of the computer?.  <p>Figure 3-10: Exit from Acquisition Workstation Dialog Box</p> <ol style="list-style-type: none"> 3. Select Yes.

6.0 How to Shut Down the System

1. Select **File>Exit**.
2. Select **Power off the computer?** from the Exit from acq. station dialog box. See Figure 3-10.
3. Select **Yes**.



Note...

*If a message appears and indicates that there are spool jobs, select **log off** or wait until the jobs are finished. (The spool continues again later.) Hologic recommends that you allow the jobs to complete before shut down.*



Caution:

Do not use a Circuit Breaker or Emergency Off switch as a routine method to turn off the Selenia.



Caution:

Never turn off the Acquisition Workstation or Value Console Circuit Breaker except in emergency. The circuit breaker turns off the Uninterruptible Power Supply (UPS) and risks data loss.

Chapter 4 - The Selenia Accessories

1.0 Introduction

Available accessories depend on your system configuration. Selenia performs screening and diagnostic applications. The Selenia S system performs only screening applications.

Table 4-1: Available Accessories

Accessory		Selenia	Selenia S
Routine Screening Paddles	18 x 24 cm	*	*
	18 x 24 cm FAST	*	*
	24 x 30 cm	*	*
	24 x 30 cm FAST	*	*
	Small Breast	*	*
Contact and Spot Compression Paddles	10 cm Contact	*	
	15 cm Contact	*	
	Ultrasound	*	
	7.5 cm Spot Contact	*	See Note
	Frameless Spot	*	
Magnification Paddles	7.5 cm Spot	*	
	10 cm	*	
	15 cm	*	
	10 cm Mag Localization	*	
	10 cm Mag Perforated Localization	*	
Localization Paddles	10 cm Perforated	*	
	15 cm Perforated	*	
	10 cm Rectangular Open	*	
	15 cm Rectangular Open	*	
Patient Face Shield		*	*
Magnification Platform		*	
Localization Crosshair Device		*	
Magnification Crosshair Device		*	



Note...

On the Selenia S system, only use the 7.5 cm Spot Contact Paddle for compression thickness calibration.

2.0 Patient Face Shield

Remember to examine the Face Shield condition before use each day.

2.1 How to Install the Face Shield on the Tubehead

1. Put the tab ends on the open end of the Face Shield into the slots on the tubehead mount.
2. Slide the Face Shield into the tubehead mount until the Face Shield locks in position.

2.2 How to Remove the Face Shield

1. Pull the rear sides of the shield away from the tubehead.
2. Slide the shield off the mount.




Figure 4-1: Installation of the Face Shield



Warning:

Except for magnification case studies, always use the Face Shield.



Warning:

The Face Shield does not protect from radiation.

3.0 Compression Paddles

There are 19 compression paddles for Screening and Diagnostic Procedures. 12 of the 19 paddles are Smart Paddles and have mechanisms that allow manual lateral movement (both left and right). The x-ray collimation moves with the Shifting Compression Paddle position. The other seven paddles include five magnification paddles and two large, 24 x 30 cm paddles.

3.1 How to Install the Compression Paddles

1. Align the mounting pins on the rear frame of the paddle with the mounting holes in the Compression Device.
2. Push the compression paddle into the Compression Device until the paddle stops.
3. Carefully pull the paddle out approximately 3 mm until the paddle clicks into position.

3.2 How to Remove a Compression Paddle

1. Hold the paddle by the metal sides of the frame.
2. Pull the paddle in a straight line from the mounting holes.




Figure 4-2: Mount the Compression Paddle

3.3 Maintenance and Cleaning

Clean the paddles after each use. Refer to Chapter 9, page 95, for cleaning instructions.

3.4 Smart Paddle System

A release knob on each side of the Smart Paddle allows it to move from one side of the image detector to the other. This enables the Technologist to optimize the breast position and imaging within the active area of the Image Receptor, especially when Lateral Views of a small breast are required.

To shift a Smart paddle:

1. Select a View on the Acquisition Workstation.
2. Rotate the release knob on the side of the paddle from the lock to the shift position.
3. While you hold the knob in the shift position, move the paddle to match the View. The paddle locks into position on a detent.
4. Release the knob to lock it in position.
5. Confirm the paddle is locked into position.
6. Activate the light field lamp and confirm the collimator matches the paddle position.




Figure 4-3: Shifting Compression Paddle

3.5 FAST Paddle Use

The Fully Automatic Self-adjusting Tilt (FAST) Paddle is for use when the composition of the breast tissue does not allow uniform compression across the complete breast with a flat compression paddle. For these patients, not enough compression can cause an image to appear to be out of focus at the anterior region from both involuntary motion and not enough compression.

The FAST paddle use with this type of breast provides these features:

- Reduced motion artifacts, because the compression is more effective.
- The compression is more uniform from the chest wall to the nipple.
- Maximum patient comfort, because over compression at the chest wall is prevented.

The FAST paddle automatically tilts when the compression is applied. The paddle is in a flat position until a compression force of approximately 20 pounds (88 Newtons) is applied. The paddle then tilts until its maximum angle is reached at a force of approximately 30 pounds (132 Newtons).

The FAST paddle does not require excessive compression, but you must use enough compression to prevent the movement of the breast. You should use a consistent amount of compression, especially for related left and right views.

The FAST paddle may not be best for breasts that are equal or symmetrical in thickness from the chest wall to the anterior area of the breast.

3.6 How to Realign the Paddle Front Edge

1. Install the paddle and apply 30 lb (132 N) of compression force.
2. Make sure that the front edge is aligned to the front edge of the Image Receptor. See Figure 4-4.

To align the edge when it is not aligned as shown in Figure 4-4:

1. Release the compression and remove the paddle.
2. Loosen the hardware that holds the paddle to the frame. On a standard screening paddle, or a frameless spot paddle, turn the paddle upside down and loosen the hardware one complete turn. On a FAST paddle, use a hex wrench to loosen the paddle.
3. Turn the paddle over (if necessary) and loosen the two inside screws. See Figure 4-5.
4. Install the paddle and apply 30 lb (132 N) of compression force.
5. Move the paddle to the correct position. (Do not release the compression.)
6. Tighten both inside screws. See Figure 4-5.
7. Release the compression and remove the paddle.
8. Turn the standard screening or the frameless spot paddle upside down and tighten the hardware that fastens the plastic paddle to the frame.
9. Install the paddle and apply 30 lb (132 N) of compression force.
10. Confirm that the front edge of the paddle is aligned to the front edge of the Image Receptor.




Figure 4-4: Alignment of the Paddle






Figure 4-5: Alignment of Adjustment Screws

3.7 Identification of Compression Paddles


3.7.1 Routine Screening Paddles




18 x 24 cm Screening Paddle




18 x 24 cm FAST Paddle



24 x 30 cm Screening Paddle




24 x 30 cm FAST Paddle




Small Breast Paddle


3.7.2 Contact and Spot Compression Paddles




10 cm Contact Paddle




15 cm Contact Paddle



Ultrasound Paddle




7.5 cm Spot Contact Paddle




Frameless Spot Paddle


3.7.3 Magnification Paddles




7.5 cm Spot Magnification Paddle




10 cm Magnification Paddle



15 cm Magnification Paddle




10 cm Magnification Localization Paddle




10 cm Magnification Perforated Localization Paddle

Use the Magnification paddles when the Magnification Platform is installed.


3.7.4 Localization Paddles




10 cm Perforated Localization Paddle



10 cm Rectangular Open Localization Paddle



15 cm Perforated Localization Paddle



15 cm Rectangular Open Localization Paddle

3.7.5 About the Localization Paddles

The Installation method is the same for all Paddles. The system locks out Auto Compression Release when a Localization Paddle is installed. The AEC Sensor may require repositioning.

4.0 Magnification Platform

When the Magnification Platform is installed, the HTC grid automatically retracts and the default values for Magnification x-ray exposure techniques are set.

4.1 How to Install the Magnification Platform

1. Remove the Face Shield (refer to Section 2.0, page 34).
2. Remove the compression paddle.
3. Move the Compression Device above the slots in the C-arm.
4. Hold the Magnification Platform by the support brackets and align the hooks of the bracket with the mounting slots on the C-arm.
5. Put the hooks into the C-arm slots.
6. Push the frame down so that the clips lock the Platform in position.

When the platform is installed, you can use only the Magnification paddles. (See Section 3.7.3, page 39)

4.2 How to Remove the Magnification Platform

1. Remove the Magnification paddle.
2. Press the locking clips on the Magnification Platform down to release the hooks.
3. Lift the Magnification Platform up and out by the support bracket.
4. Reinstall the Face Shield.




Figure 4-6: Magnification Platform Installation

5.0 The Localization Crosshair Device

The Localization Crosshair device, used with a localization paddle enables the location of a specified spot on the breast.

5.1 How to Install the Crosshair Device




Figure 4-7: Installation of the Crosshair Device

1. Remove the Face Shield from the tubehead.
2. Move the Compression Device below the C-arm mounting slots.
3. Hold the Crosshair device by the support brackets and slide the hooks into the C-arm slots.
4. Push the device down to the locked position.
5. Install the Localization paddle on the Compression Device.

5.2 How to Remove the Crosshair Device

1. Rotate the Crosshair device to the left or right.
2. Press the two locking levers inside each mounting arm. See Figure 4-8.
3. Lift the device up and from the C-arm. Be careful. The device is top heavy.
4. Install the Face Shield.




Figure 4-8: Crosshair Locking Levers

5.3 Localization Crosshair Device to the Localization Paddle Alignment

**Note...**

Before you perform the following adjustment, make sure the Localization paddle is aligned to the edge of the image receptor. Refer to "How to Realign the Paddle Front Edge," page 37.

1. Install a rectangular localization paddle.
2. Loosen the adjustment lock screw on the bottom of the Crosshair Device. See Figure 4-9 No. 1.
3. Put a piece of white paper on the breast tray so that you can see the shadows of the crosshairs.
4. Move the Localization paddle approximately 6 cm above the image receptor.
5. Turn on the light field.
6. Move the Crosshair Device until the rectangle of light aligns with the opening in the localization paddle.
7. Tighten the adjustment screw.




Figure 4-9: Adjustment Lock Screw

6.0 The Magnification Crosshair Device




Figure 4-10: Installation of the Magnification Crosshair Device

1. Remove the face shield from the tubehead.
2. Align the Magnification Crosshair Device with the grooves of the face shield rails.
3. Push forward until the device locks into position.
4. To remove the Magnification Crosshair Device, pull and slide the assembly toward you.

Chapter 5 - The User Interface

1.0 The Three Main User Interface Screens

The process of image acquisition, enhancement, and output is done through three screens:


- Select a Patient screen
- Patient Views screen
- Preview Images screen

When you select an option on a screen, other screens may appear for information entry or selection of additional options. This manual covers the screens that appear while you are logged in as a “tech”.


- Select a Patient



- Patient Views



- Preview Image



2.0 A Screening Exam Sequence Suggestion

Table 5-1: Workflow Suggestions to Select a Patient and Acquire an Image



Action	What it looks like												
1. Start up the system and log in to Operating System and Application. Refer to Chapter 3, page 24	 <p>The diagram shows a workstation with a monitor, keyboard, and mouse. Callout 1 points to the monitor area, and callout 2 points to the keyboard area. Below the diagram is a screenshot of the Selenia login screen. The screen has a blue header with the Selenia logo and the text 'Welcome to selenia'. Below this is a text input field labeled 'Please enter your user name'. At the bottom of the screen are four buttons: 'Help', 'Options', 'Start Over', and 'OK'.</p> <p>Technician, Manager, or Physicist</p>												
2. Select the patient. <ul style="list-style-type: none"> • Query the MWL <ol style="list-style-type: none"> a. Select Worklist tab in upper left of screen. b. Use the Bar Code Scanner to input query data, then select Query, or enter the query data then select the Enter key. c. Select the patient from the list. • Search the local database <ol style="list-style-type: none"> a. Select the Local Exams tab in the upper left of the screen. b. Select the Patient Name tab (or ID tab) then enter the data. c. Select the Search button. See Chapter 6, Section 2.3, page 48 for information about wildcard searches. d. Select the patient from the list. • Add a New Patient <ol style="list-style-type: none"> a. Select New from the Patient area of the screen. b. Enter the required information. c. Select the Accept button. 	 <p>The screenshots show the software interface. The top screenshot shows two tabs: 'Worklist' and 'Local Exams'. The middle screenshot shows a search interface with tabs for 'Patient Name', 'Local Exams', and 'Patient, One', and a 'Search' button. Below this is a table of search results:</p> <table border="1"> <thead> <tr> <th>↑</th> <th>Patient, One</th> <th>ID</th> <th>Sex</th> <th>DOB</th> <th>DOB</th> </tr> </thead> <tbody> <tr> <td></td> <td>Patient, One</td> <td>00581581</td> <td>F</td> <td>04/01/1948</td> <td>2008</td> </tr> </tbody> </table> <p>The bottom screenshot shows a 'Patient' form with fields for Name, DOB, ID, and Gender, and buttons for 'New', 'Edit', and 'Close'.</p>	↑	Patient, One	ID	Sex	DOB	DOB		Patient, One	00581581	F	04/01/1948	2008
↑	Patient, One	ID	Sex	DOB	DOB								
	Patient, One	00581581	F	04/01/1948	2008								

Table 5-1: Workflow Suggestions to Select a Patient and Acquire an Image








Action	What it looks like
<p>3. Select the Procedure, when more than one is scheduled.</p> <ul style="list-style-type: none"> • Select from the scheduled procedures list. <ol style="list-style-type: none"> a. Find the procedure in the list. b. Double-click the Select button. • Add a procedure <ol style="list-style-type: none"> a. Select New in the Procedure area. b. Select the procedure from the list. c. Select the Accept button. 	
<p>4. Edit the patient Information, if necessary</p> <ol style="list-style-type: none"> a. Select the patient from the patient list. b. Select Edit from the Patient area. c. Make the changes in the Edit Patient screen. d. Select the Accept button. 	
<p>5. Install the accessories. See Chapter 4, page 33 for installation, removal, and alignment procedures for each device.</p>	<ul style="list-style-type: none"> • Face Shield • Compression Paddle • Magnification Platform • Crosshair
<p>6. Select the view.</p> <ul style="list-style-type: none"> • Use the default view sequence. <ol style="list-style-type: none"> a. Acquire the selected view. b. Acquire an image for each of the remaining views. • Select another view from the displayed views. <ol style="list-style-type: none"> a. Select any view displayed on the Patient View screen. b. Acquire the selected view. c. Repeat for the remaining views. • Add a view. <ol style="list-style-type: none"> a. Click Add View in the Patient View screen. b. In the dialog box, select a view from the list of Standard Mammography Views. To add many views, press and hold the Control key and select the views. c. Select the OK button. 	

Table 5-1: Workflow Suggestions to Select a Patient and Acquire an Image

Action	What it looks like
<p>7. Select the Exposure Techniques.</p> <ul style="list-style-type: none"> • Use default values for selected exposure mode. <ol style="list-style-type: none"> a. Select the Exposure Mode. b. Make sure the default settings are correct for the exam/patient. • Change default values for the selected exposure mode. <ol style="list-style-type: none"> a. Select the Exposure Mode. b. Adjust the default values. 	
<p>8. Position the system and Patient.</p> <ol style="list-style-type: none"> a. Adjust the vertical height of the Breast Platform. b. Set the C-arm angle. c. Put the arm and hand of the patient in a safe location. d. Compress the breast. 	
<p>9. Acquire the image.</p> <ol style="list-style-type: none"> a. Make sure the message area of the screen displays “Ready” (1) and there is no “X” through the x-ray indicator (2). b. Press and hold the x-ray buttons. c. Inspect the image on the Preview Screen then select the Accept or Reject button. 	
<p>10. Close the procedure.</p> <ol style="list-style-type: none"> a. Select the Close Procedure button. b. Select the Yes button. 	

Chapter 6 - The Select a Patient Screen

1.0 Introduction

The Select a Patient screen appears when the system application loads. This screen allows you to locate a patient, and access non-routine Acquisition Workstation options with no patient selected.




Figure 6-1: The Select a Patient Screen

Table 6-1: Legend for Figure 6-1

1. Menu Bar	Section 8.0, page 59.
2. Search for a Patient	Section 2.0, page 48.
3. Patient Query	Section 2.2, page 48.
4. Today's Worklist	Section 2.3, page 48
5. Exposure Techniques	Section 7.0, page 54.
6. Add/Edit the Patient Information	Section 3.0, page 49.
7. Add/Edit a Procedure	Section 5.0, page 51.
8. Outputs	Section 6.0, page 53.
9. Special Modes	Section 9.0, page 72.
10. Generator Status	Section 10.0, page 73.
11. System Message	Section 11.0, page 73.

2.0 Select a Patient

A patient can be added to the system or already exist and be selected from a list. There are two places to search for an existing patient, the worklist which is a list of scheduled patients from a Modality Worklist Provider and the Local Exams which are the patients on the computer in the workstation.


2.1 How to Query the Modality Worklist Provider

1. Select the **Worklist** tab in the upper, left area of the Select a Patient screen.
2. Input the query information with the Bar Code Scanner or the Keyboard.
3. Select the **Query** button or press the **Enter** key.



Note... *This query does not search the local Acquisition Workstation database, Local Exams. An Asterisk is not required for this search.*

4. If you changed the date from the 'today' setting, select the **Today's Worklist** tab.
5. Select the patient from the list to open their procedure.



If your search for one match with a Patient ID or Accession Number has results, the patient Procedure opens automatically. Selecting a patient adds them to the Local Exams database.

2.2 How to Search the Acquisition Workstation Database (Local Exams)

1. Select the **Local Exams** tab in the upper, left area of the Select a Patient screen.
2. Select the **Patient Name** or **ID** tab.
3. Enter the Patient Name or ID in the Search box. Or, enter part of a Patient Name or part of an ID number and use an asterisk for the missing characters.
4. Select the **Search** button or press the **Enter** key.

2.3 More Information about Searches

- When part of the patient name is entered, the system applies a wildcard character at the end of the name.



Note... *You can use the asterisk (*) as a wildcard character to increase your search. For example, when you enter R*ph, Patients with the names Randolph and Rudolph will display. If you use a wildcard character, make sure that you include a minimum of one other character.*

- You can use one or more of the search fields.
- You can limit the query to a date range.
- Most worklist providers require information in the Patient ID and Accession Number fields to match. Data in these fields is case sensitive.
- You must include the caret (^) character between the name parts (for example, DOE^).

3.0 How to Add a New Patient

The Patient area of the screen allows you to Add a New Patient and Edit the existing patient information.

1. Select **New** in the Patient pane to display the New Patient Entry dialog box.




Figure 6-2: New Patient Entry Form

2. Enter the Patient demographic information. Use the Tab key or the trackball to move through the fields. You must enter information in the Fields that have an asterisk.
 - The **Patient Name** fields must contain only letters or numbers.
 - The **Clear** button clears the information you entered.
 - The **Cancel** button closes the form without your changes.
3. Select the **Accept** button to enter the Patient in the Acquisition Workstation Local Database. (This new information is not sent to the PACS system.) The Patient Views screen displays.



Warning:

Verify the Patient Demographic Application information before you acquire an image.

4.0 How to Edit the Existing Patient Data



Figure 6-3: The Edit Patient Dialog Box



Note...

This procedure only changes the patient information for images in the case study not yet taken. The Acquired images must be resent using Image Management and Resend, after the information is edited. Refer to Section 8.2.1, page 65.



Caution:

Do not edit the patient information if you use a Modality Worklist. See Section 4.1, page 50

To edit the Patient information:

1. Select **Edit** in the **Patient Box** to display the **Edit Patient** dialog box.
2. Make the changes to the fields that have wrong or missing information. Use the Tab key or the trackball to move through the fields.



Note...

The **Cancel** button closes the dialog box without any changes. The **Clear** button clears last name, first name, DOB, and age.

3. Select the **Accept** button.
4. Verify that the correct changes appear in the screen before you acquire a new image.

4.1 About Patient Information Edits with a Modality Worklist and PACS

- Make the patient information changes on the HIS/RIS *before* any information is sent to the Selenia.
- Make the patient information changes on the PACS if the images were sent to the Selenia.
- Any patient information changes on the Selenia must be a last option. Some RIS systems do not have the ability to change this information at the Selenia.
- The Modality Worklists transfer the Patient Information from a HIS/RIS system to the Selenia. Any change on the Selenia to this information can cause a mismatch of information when the Selenia images are sent to PACS. Work with your HIS/RIS administrator before you make any changes on the Selenia to the patient information.

4.2 How to Change the Patient Information After an Image Is Accepted

When an image is accepted, the image is immediately sent to the output. If the related information requires any change, then repreview, save, and resend the image as an additional record and delete the previous record.

Remove any wrong images from:

- Printer
- Diagnostic Review Workstation
- PACS

Correct the information at the AWS:

1. Select the **Local Exams** tab and search for the patient.
 - If patient information is wrong, select **Edit** in the Patient Pane.
 - If the accession number is wrong, select **Edit** in the Procedure pane.
2. Change the wrong patient information.
3. Select the **Accept** button.
4. Select the **Close Procedure** button.
5. Select **Admin>Image Mgt.**
6. Search for the patient.
7. Select the **Repreview** tab.
8. Select the first thumbnail image.
9. Select the **Repreview** button.
10. Select the **Save** button.
11. Repeat the steps 10 to 12 for the remaining images in the procedure. Eight thumbnail images display (if there were four in the original procedure.)
12. Select the last four thumbnail images (the corrected views).
13. Select **Resend** to all required outputs.
14. After you confirm that the four images were resent, have a manager delete the four images with the wrong labels.

5.0 New, Edit, Add a Procedure

5.1 How to Add a New Procedure

1. Retrieve the existing Patient.
2. Select **New** in the **Procedure** pane.



Figure 6-4: New Procedure, Edit a Procedure, and Add a Procedure Buttons

3. Select **Yes** in the **Creating a New Procedure will close the current Procedure** dialog box. Any open procedure will close.
4. When the **New Procedure** dialog box appears, enter a new **Accession** Number.

Figure 6-5: New Procedure Dialog Box



Note... The **Cancel** button closes the dialog box without any changes. The **Clear** button clears the Accession number.

5. Select the **Procedure Description** from the drop down list.
6. Select the **Accept** button.
7. Select **Close** in the Patient area.

5.2 How to Edit a Procedure

1. Select the **Edit** button in the Procedure area of the screen.
2. Change the Accession Number. You cannot change the Procedure selection.
3. Select the **Accept** button.

5.3 How to Add a Procedure

The Add a Procedure function allows you to have many Procedures open at the same time.

1. Select the **Add** button in the Procedure pane.

Figure 6-6: Add a New Procedure Dialog Box



Note... The added procedure must use the same Accession Number as the open procedure. If you do not plan to use the same Accession Number, select the **New** button and follow Section 5.1, page 51.

2. From the drop-down list, select the Procedure to add.
3. Select the **Accept** button.

6.0 Image Output




Figure 6-7: Select the Output

Before acquiring an Image, verify that you selected the correct output. All real patient images must be printed or committed to a PACS. Select the Output group to use for this procedure from the drop down list. To edit existing or create new output groups, refer to Section 8.1.5, page 63.

There is also a statement of the current number of Jobs in the Queue.

7.0 Exposure Information

7.1 The Exposure Pane




Figure 6-8: Exposure Techniques






The top, left side displays the status of these options:

- Compression Force Mode
- Compression Release
- Compression Thickness
- Collimation Mode
- Collimation Size
- Image Size
- Paddle ID (type and position)

Select the Exposure Techniques in the bottom area. Use the up/down arrows in the Exposure Mode area (the lower, far-left box) to select the Exposure Mode. The exposure default values for that mode appear to the right.

7.2 About Exposure Modes

Table 6-2: Exposure Modes

Exposure Modes (#1)	Operator Adjustments	Where the function is in the Techniques Pane
Auto-Filter	Select the AEC Exposure (#2); AEC Sensor (#3) (Positions 1-7 or Auto).	
Auto-kV	Select AEC Exposure; AEC Sensor.	
Auto-Time	Select kV; Filter; AEC Exposure; AEC Sensor.	
TEC	Select Breast Density (#4), then select Accept or change the Exposure Techniques, which puts you in Override (Manual) Mode.	
Manual	Manually calculate and select all x-ray techniques and Exposure Factors (kV, Filter, and mAs).	



Note...

When you change the mode, all Exposure Techniques reset to the default values for that mode.

7.2.1 About the Manual Exposure Mode

In the Manual mode the Operator sets all Exposure Techniques. The default values appear when you select the View, then you make any necessary adjustment.

- To change an Exposure technique with the trackball, scroll to the value with the up/down arrows, or select the option box.
- To change an Exposure technique with the keyboard, press the Function Key for the item at the top of the keyboard.

7.2.2 About the TEC Exposure Mode

The Tissue Exposure Control (TEC) Mode is an enhanced Manual Exposure Control mode.

1. Before you position and compress the breast, select the **Breast Density** type on the Acquisition Workstation screen or press the **Breast Density** key.




Figure 6-9: Breast Density Options

2. Compress the breast. The X-ray Exposure Techniques for kV, mAs, and filter update.
3. Select the **Accept** button in the Breast Density area (or press the Enter key). When "Ready" displays in the Generator Status message area, the system is ready to acquire an image.




Figure 6-10: Accept TEC Exposure Techniques

If you change one or more of the parameters (kV, mAs, or filter), the system enters the Override Mode. When the Override message displays, the system is ready to acquire an image.




Figure 6-11: The Override Mode

There are three methods to exit the Override Mode and return to TEC mode:

- Change the Breast Density setting.
- Release the compression then compress the breast.
- Select TEC mode a second time.

7.2.3 AEC Exposure Mode

The Automatic Exposure Mode (AEC) uses the pre-exposure Scout pulse, and for some modes the compression thickness to calculate the exposure.

The AEC options are Auto-Filter, Auto-kV, and Auto-Time. See Section 7.2, page 55 for the exposure technique selections available with each mode.

7.2.4 AEC Exposure Adjustment Settings

You can use the Exposure Adjustment Settings to increase or decrease the target pixel value. Each step changes the target pixel value by approximately 15% and makes a related change in dose. See Figure 6-12, page 58.

Set this control in a range from +4 to -3 to increase or decrease the exposure. This adjustment setting remains until you change the setting or restart the system.



Warning:

You increase the patient dose to high levels if you increase the AEC exposure adjustment setting. You increase the image noise or decrease image quality if you decrease the AEC exposure adjustment setting.



Warning:

Increased Exposure adjustment leads to higher dose to the patient. Keep the patient dose as low as practical to get good image quality.




Figure 6-12: AEC Exposure and Sensor Position Adjustment On-Screen Settings

7.2.5 AEC Sensor Position Indicator Settings

Select one of the seven positions to indicate the area of interest for AEC calculations. Or, select Auto AEC where a larger area is tested.



Note... Always confirm the position of the AEC Sensor before an AEC exposure.

The AEC Sensor Position has eight available positions.

- The Positions 1 through 7 are manually selected, from the chest wall edge (position 1) to the nipple edge (position 7).
- Position 8 (Auto) automatically positions the sensor.

To change the AEC Position selection, use:

- The buttons on the compression device.
- The AEC Pos key on the keyboard.
- The AEC position arrows on the Acquisition Workstation display.

Set the AEC Sensor Position after you compress the breast. The displays on the side of the Compression Device and the lower right corner of the screen indicate the position setting.

7.2.6 AEC Alarm Messages

When the system cannot use the selected AEC mode to acquire an Image, one of the following messages displays.

Table 6-3: AEC Alarm Messages

Message	Reason and Corrective Action
Calculated Exposure Time Less Than Min	The AEC Sensor is over a thin or small area of breast tissue. Use the AEC Sensor in a position under the breast tissue. Alternately, use AutoTime and a lower kV.
Calculated Exposure Time Exceeds Max	Very dense breast tissue, an implant, pacemaker, or other anatomy changed the AEC sensor. Move the manual AEC sensor where there are no obstructions. Alternately, use AutoTime and a higher kV.

1. Select **OK** in the dialog box.
2. Select the **Reset Generator** button.

8.0 Menu Bar

The Menu Bar provides easy access to non-routine Acquisition Workstation options. These options are available while the Acquisition Workstation application runs unless a dialog box or Procedure is opened.

The table below shows the Menu Bar options when a patient procedure is *not* selected.

Table 6-4: Menu Bar Options and Functions

Menu	Options	Function
File	Exit	Exit from Acquisition Workstation
Edit	Standard Setup	Select the startup default values. (Section 8.1.1, page 60)
	View Order Editor	Change the view order. (Section 8.1.3, page 61)
	User Setup	Add, edit, delete a user profile. (Section 8.1.4, page 62)
	Outputs	Add, remove output devices. (Section 8.1.5, page 63)
Admin	Image Management	Retrieve the images from the image spool. (Section 8.2.1, page 65)
	Protect Patients	Protect or Unprotect patient records. (Section 8.2.4, page 67)
	Manage Queues	Find a specified job in the queue. (Section 8.2.5, page 68)
	Eject	Remove a disk from CD/DVD drive. (Section 8.2.6, page 69)
	Import	Copy images from a CD to the system. (Section 8.2.7, page 69)
	Retrieve Priors	View previously-acquired images. (Section 8.2.8, page 70)
	PPS Status*	Additional options (when MPPS Service is installed) for close of a procedure. (Section 8.2.9, page 70)
	Calibrate	Access calibration procedures. (Section 8.2.10, page 70)
	Test Patterns	Access test pattern procedures. (Section 8.2.11, page 71)
	DR Device Control	Access by Service personnel. (Section 8.2.12, page 71)
Available Disk Space	Displays status of hard disk drive space. (Section 8.2.13, page 72)	
Stereotactic**	En/Disable Stereo mode**	Engage or disengage the stereo mode. Available when the StereoLoc II is installed.
Info	About the Acquisition Workstation	(Section 8.3, page 72)
* This option appears when the MPPS Service is installed.		
** This option appears when the Stereotactic Service is installed.		

8.1 Edit Menu Function

8.1.1 Standard Setup Screens




Figure 6-13: The Setup Screen

1. Select **Edit>Standard Setup**.
2. Use the trackball or the keyboard Tab key to move between fields.
3. Use the trackball or arrow keys to display and select options.
4. Select the **Close** button to exit the dialog box.



Note...

When Manual Collimation is selected, the function that automatically changes the Collimation when a paddle is changed or moved turns off.

8.1.2 Set Technique Defaults

Select the Set Technique Defaults button to display the Technique Default Setup screen. You must set the default values for all modes.

1. Select **Manual** Exposure Mode and select the system startup defaults.
2. Change the mode to Auto-time and select the options.
3. Continue to change the mode and select the options for all remaining exposure modes including Magnification modes.

8.1.3 View Order Editor




Figure 6-14: View Order Editor

1. Select **Edit>View Order Editor**.
2. Select the procedure from the drop-down menus.
3. Select the Views for each Procedure:
4. Move the View icons to change the order.
5. Select **Close**. Your view order preferences load when you log in.

8.1.4 Add or Edit a User

A signed-in user can add or delete a user at the same or lower level. To access the Edit User options, select **Edit>User Setup**.




Figure 6-15: The Edit Users Dialog Box

Table 6-5: User Setup Options

Option	Steps
<ul style="list-style-type: none"> Add a new user 	<ol style="list-style-type: none"> Select New. Enter the requested information. The password must contain a minimum of 6 characters. Select the Accept button.
<ul style="list-style-type: none"> Edit a user 	<ol style="list-style-type: none"> Select your name from the list. Select Edit. Change the Initials, if necessary. If requested, enter the Password. The password must contain a minimum of 6 characters. Select the Accept button.
<ul style="list-style-type: none"> Delete a user 	<ol style="list-style-type: none"> Select the User to delete. Select the Delete button. Select the Yes button to the confirmation prompt.
<ul style="list-style-type: none"> Close the dialog box 	Select the OK button to close the Edit Users dialog box.

8.1.5 Outputs

To Add New Output Groups:

1. Select **Edit > Outputs**.
2. Select **New**.
3. Enter the name for this Output.
4. Select **OK**.
5. Select the new name in the **Edit Outputs** dialog box.
6. Select **Edit**.




Figure 6-16: Edit Output Device Dialog Box



Note...

Outputs are first created during installation.

7. Select the tab for the type of device.
 - CDRW = CD-RW
 - Laser Camera = Printer
 - Physician Display = Diagnostic Review Workstation, Technologist Review Workstation, CAD, etc.
 - Archive = PACS
8. Select the device name.
9. Select the **Add** button to move it to the right column.
10. To add another output device to the new group, select the tab for the device type.
11. Select the device name.
12. Select the **Add** button.
13. Select the **Accept** button when finished with the group.
14. Select the **OK** button.

To Edit Outputs:

1. Select **Edit > Outputs**.
2. Select an Output name. If the Output keeps its current name, go to step 7.
3. Select the **Copy** button.
4. Enter a new **Output name**.
5. Select the **OK** button.
6. Select the new name in the **Edit Outputs** dialog box.
7. Select the **Edit** button.
8. Select the tab for the type of device to add.
9. Select the device to add.
10. Select the **Add** button to move it to the right column.
11. Select the name of any output in the **Devices currently selected** column that is not needed in the group.
12. Select the **Remove Item** button.
13. Select the **Accept** button.
14. Select the **OK** button.

8.2 Admin Menu

8.2.1 Image Management

1. Select **Admin>Image Management**. The Choose a New Patient dialog box appears.

Name	ID	Sex	DOB
Patient, Two	00581382	F	05/06/1932

Figure 6-17: Choose a New Patient Dialog Box

2. Select the **Patient Name** tab or the **Patient ID** tab.
3. Enter the search information in the **Enter Patient Name** or **ID** box.
4. Select the **Search** button.
5. Select a patient from the list.
6. Select the thumbnail image of the image to resend or repreview.
 - You can select more than one to resend.
 - If there are more than four images, use the arrows on the side to scroll.
7. Select the Tab for **Repreview** or **Resend**.
8. Continue with directions in “The Resend Options,” page 66 or “The Repreview Option,” page 67.



Note...

If the commit function is enabled and archive has committed the image to storage, the line “Commit: Accepted” appears in the information about the image in the upper right pane. See Figure 6-18.

8.2.2 The Resend Options



Figure 6-18: Image Management Resend Options

1. Select the resend button.
 - To resend all images on the hard disk drive for this patient, select **Resend All Images**.
 - To resend selected images, select the thumbnails of the images to resend, then select the **Resend Selected** button.
2. Select the Output from the drop-down list.
3. If the selected Output has a printer, select the **Print Image Processing Options**.
4. Select the **Resend** button.
5. Select the **OK** button.



Note...

You cannot change any Patient or image object information from the resend function. This includes image, LUT modifications and DICOM information such as Patient name spelling, ID, and comments. To send new information you must make a new file with the Repreview function, save it, and Resend the new file.

8.2.3 The Repreview Option

1. Select the thumbnail image to review.
2. Select the **Repreview** tab, then select the **Repreview** button.
3. When you finish reviewing, select the **Cancel** button.
4. To create a new file to resend, make any required changes to the markers, comments, or use the **Edit View** function, then click **Save**. If you select a FOR PRESENTATION image, you cannot change the markers and ROI is not enabled.
5. Make any necessary changes in the **Contrast** or **Density**. These changes do not save when you save the image.
6. Select the **Save** button.



Note...

Make sure that you send the newly-saved image. The original is a possible selection.

8.2.4 Protect Patients

Use the Protect function to prevent automatic removal of patient records from the hard drive.

1. Select **Admin>Protect Patients**.
2. Search for the Patient to protect.
3. Select the Patient to protect from the search results.
4. Select the **Protect** button at the bottom of the screen.
5. When the confirmation dialog box appears, select the **Yes** button.
6. Select the **Close** button.



Note...

To create more space on the hard disk drive, unprotect patients who do not need protection.

To unprotect a patient record:

1. Select the patient.
2. Select the **Unprotect** button.

8.2.5 Manage Queues

Use Manage Queues to find problem jobs or a specified job in the queues.



Warning:

If you select Delete Job, you cannot undo the deletion. If it still needs to be sent, resend the job to a storage device later.




Figure 6-19: Manage Queues

Legend for Figure 6-19

1. View the list of available queues.
2. Enable or disable a job spooler.
3. Update the list of Queued Jobs.
4. View the list of the jobs in the queue.
5. Remove inactive jobs from the queue. Resend stopped jobs.

Table 6-6: How to Use Manage Queues

Option	Steps
<ul style="list-style-type: none"> • To find the status of a job: 	<ol style="list-style-type: none"> 1. In the Devices list, select the queue to view. 2. Select the OK button.
<ul style="list-style-type: none"> • To delete a job: 	<ol style="list-style-type: none"> 1. Select the queue in the Devices list. 2. Select the Suspend <name> Spooler button. 3. Select a job from the Queued Jobs list, or Ctrl+click to select many jobs. 4. Verify that the selected job is the job to delete. There is no “undo”. 5. Select the Delete From Queue button. 6. Select the OK button.

Table 6-6: How to Use Manage Queues

Option	Steps
<ul style="list-style-type: none"> To resend a stopped job: 	<ol style="list-style-type: none"> Select the queue in the Devices list. Select a job marked 'Stopped' from the Queued Jobs list, or Ctrl+click to select many jobs. Verify that the selected job is the job to resend. Select the Resend Stopped Job button. Select the OK button.

8.2.6 Eject

Select **Admin>Eject** to open the CD-RW drive drawer when there is a CD in the CD-RW drive.



Note...

When the CD-RW drive is empty, press the button on the CD-RW drive to open the drawer.

When there is a disk in the drive, the drive drawer opens only from the menu. To close the drive drawer, press the CD-RW drive button.

8.2.7 Import

To import the information from a CD:

- Open the CD-RW drive drawer.
 - Press the CD-RW drive button to open the tray if there is no CD in the drive.
 - Select **Admin>Eject** to open the drive if there is a CD in the drive.
- Put the CD that contains the images to import on the tray.
- Press the CD-RW drive button to close the tray.
- Select **Admin>Import**.
- In the dialog box, find the file (or files) to import.
 - Select a "+" to open a list of folders and/or files.
 - Select the empty box to select the file or folder.



Hint...

To deselect a file, click the checkmark. When you deselect a folder, the files below the folder are not deselected. You must clear each file.

- Select the **Import** button. The file is copied from the disk and is available in the **Select a Patient** screen.



Note...

*Do Not press the **Eject** button on the drive while the system reads from or writes to the disk.*

If you try to eject or read the CD while the import process finishes, you can cause the drive to stop.

8.2.8 Retrieve Priors

To find the Images in the Local Database, select the Local Exams tab and search for the Patient.

To Retrieve Priors from the PACS system:

1. Select **Admin > Retrieve Priors**.
2. Enter the Patient name in the **Name** field or ID in the **ID** field (wildcards allowed).
3. To limit the results to a date range, select the checkbox **Search by Date Range**, then enter the Date Range for the Prior Images. The date format must match the configuration setting, (mm/dd/yyyy or dd/mm/yyyy).
4. Select the **Query** button or press the **Enter** key. If the button is not enabled, you need more characters in either the **Name** field or **ID** field.
5. Select the Patient name (or names) to import.
 - You can select many Patients.
 - You can deselect Patients.
 - You can use the title bars to sort by field.
6. When the **Import** button activates, select **Retrieve and Import**. The images import into the Acquisition Workstation Local Database.
7. Select the **OK** button.
8. Select the **Close** button.

After the import completes, select the Local Exams tab to display the Patient Images. If there are more than 50 records in your query request results, adjust your search.

8.2.9 MPPS Status

When this optional feature is installed, a menu item is added to sites with the MPPS Service Class Provider. Workflow changes with three Close Procedure options instead of one option.

8.2.10 Calibrations

When Calibration is needed, the Select a Patient screen remains in Standby and the **Calibration Needed Alarm** message appears in the **System Message** box.

When you turn on the system, a message appears if the system needs any calibration. You cannot acquire the patient exposures until you perform the required calibration. All the review, resend, and administrative functions continue to work.



Note...

Clean the detector and acrylic block before you perform a Calibration.



Caution:

Make sure that you follow the directions on the screen when you calibrate a system.

If you are in the Calibration Mode screen by accident, select the **End Calibration** button. If you started the Calibration, but you did not accept images, you can reject the Image and exit the Calibration.



Note...

When you accept one image, do Not end the calibration before you acquire, view, and accept all the images.

8.2.11 Test Patterns

1. Select **Admin>Test Patterns**.
2. Select the Preview size (Figure 6-20, No. 1).
3. Select the **Test Pattern** from the **Pattern** field (No. 2).
4. From the **Output** area of the screen (No. 3):
 - Select the Output **Device**.
 - Select the Film **Size**.
 - Select the **Print True Size** checkbox, if needed.




Figure 6-20: Test Pattern Dialog Box

5. Select from the following options (No. 4):
 - To view the test pattern **Full Screen**, select the **Preview** button. To return to the dialog box, click any area on the Full Screen.
 - To send the Test Pattern to the selected output, select the **Send** button.
 - To remove any files from previous “sends” before you send the current images, select the **Clean Up** button.
 - To exit the Test Pattern option, select **Close**.

8.2.12 DR Device Control

This function is for use by service personnel.

8.2.13 Available Disk Space

This option displays the remaining number of images the hard disk drive and CD can hold.

The system automatically deletes some Images/Patients at specified levels to get back hard disk drive storage space for new Procedures. Automatic image removal does not include the following:

- The Images in a queue to an output device.
- The images for Protected Patients. (Make sure you unprotect Patients when you do not need their files.)
- The Images that are not successfully committed to an archive device. (Archive devices send commitments after Images go to long-term storage.)
- Rejected Images. The Managers need to remove rejected Images.



Caution:

The system will not acquire an image if the disk does not have enough available space. A message displays the capacity of the disk. Delete the images that are not needed to restore the disk space.

8.3 Info Menu

The Info menu displays the information about the Acquisition Workstation.

9.0 Special Modes Buttons

9.1 Non-imaging Mode Button

The Non-imaging mode is available to all users and allows the x-ray exposure without the acquisition of images. This mode is not available when you have a Patient selected. When you select the Non-imaging mode, the Image Receptor deactivates.

To access this mode, select the **Non-Imaging Mode** button.



Figure 6-21: Non-Imaging Mode Button



Caution:

Make sure that you protect the Image Receptor from excessive radiation by covering it with lead while using this mode.

To exit the Non-Imaging mode, select the **End Non-Imaging Mode** button.

9.2 The Review Button

This button is used during the Reject /Unreject process. See Chapter 8, Section 5.2, page 93.

9.3 The Simulate Capture Button

This function is for use by service personnel.

10.0 Generator Status

This area displays the current status of the system. See Figure 6-1, page 47.

11.0 System Messages

The System Messages area displays the system status. The displays in the message area are important when you restart the system after an Emergency Stop Switch is used.

12.0 Alarms



Figure 6-22: The Alarm Icon

Alarms are generally output problems that occur due to jobs that fail at the output device. The Alarm Icon shows the number of alarms. Select the icon to view the details. Erase the alarm to clear it.

Instructions for Use

Chapter 6 - The Select a Patient Screen
Alarms

Chapter 7 - The Patient Views Screen

When you select a Patient, the Views for the scheduled Procedure display on the Patient View screen in your preset order. The program automatically highlights the first View in the sequence. To choose a different View, select that View.




Figure 7-1: Patient View Screen

Table 7-1: Actions on the Patient Views screen

Topic	Information
Work with Views (1)	Section 1.0, page 76
Add/Edit the Patient Information	Chapter 6, Section 3.0, page 49
Add/Edit a Procedure (2)	Chapter 6, Section 5.0, page 51
Select the Image Output	Chapter 6, Section 6.0, page 53
Jobs in the Queue (3)	Displays the number of jobs in the outputs.
Change Exposure Techniques	Chapter 6, Section 7.0, page 54
Acquire the View	Section 2.0, page 77

1.0 Work with Views




Figure 7-2: Buttons Used on the View Screen

Table 7-2: View Options

Option	Actions
Add View (1)	<ol style="list-style-type: none"> 1. Select the Add View button. 2. Select an additional View from the list of Standard Mammography Views. To add many views, hold the Control key down and select the View. 3. Select the OK button. 4. Select the new view before you acquire the image.
Remove Unexposed (2)	<p>Unused views are automatically reclaimed, but you can delete the views with the Remove Unexposed Views button.</p> <ol style="list-style-type: none"> 1. Select the All button (or the Selected button then indicate which views to delete). 2. Select the Remove Unexposed button.
Implant Present Checkbox (3)	<p>Select the box if there is an implant present in the breast. The box should remain checked for all views in the procedure.</p>
Close Procedure (4)	<ol style="list-style-type: none"> 1. Select the Close procedure button. 2. Select the Yes button to close the procedure or the No button stay on the screen.

2.0 Acquire Images

The system is ready to acquire Images when the Generator Status indicates “Ready” (Figure 7-3, #1), and the X through the X-ray symbol (#2) in the Generator Status box is not displayed.




Figure 7-3: Ready for Exposure Indicators

Instructions for Use

Chapter 7 - The Patient Views Screen
Acquire Images

Chapter 8 - The Preview Image Screen

1.0 Introduction

The Preview Image screen appears when you acquire an image. The right side of the screen has image tools and buttons to accept or reject the image.



Figure 8-1: The Preview Image Screen

Legend for Figure 8-1

1. Acquired Image Display Area
2. Tools

2.0 The Preview Image

The Preview Image is a non-diagnostic, processed image for anatomy and position confirmation. The Patient Name, ID, and the View appear on a bar on the top of the Preview.

2.1 Progressive Preview

The first preview image appears for a quick check of the breast position. This image updates with an image of better quality to check for motion or image artifacts.

2.2 Film Label Area

The film label area which is the small rectangle in the upper, left corner is blank in the Preview.

A Digital Marker next to the label displays the Technologist initials and the marker for the selected View. You can move or change the marker.

2.3 Hanging Options

The Images in this manual are displayed in the Dorsal orientation (the Operator faces the patient).




Figure 8-2: Hanging Options

This option can affect how images display on your diagnostic workstation. Confirm the effect of the options with the vendor.

3.0 Tools for Image Enhancement and Annotation

- 1 Comments
- 2 Markers
- 3 Edit View
- 4 Implant Present
- 5 Zoom/Pan Thumbnail
- 6 Crosshair
- 7 Measurement
- 8 Window/Level
- 9 ROI size Drop-down
- 10 Set Center/Width
- 11 Display AEC Regions used
- 12 Exposure Index
- 13 Exposure Techniques
- 14 Reprocess View
- 15 Magnification
- 16 Full Zoom/Pan
- 17 Quick Zoom/Pan
- 18 ROI




Figure 8-3: The Preview Screen Tools

3.1 Comments

The comment text is inserted in the DICOM header of the Acquired Image. The diagnostic review workstation displays the comments.




Figure 8-4: Image Comments

1. Select the **Comments** button.
2. Select the **Additional Comment** checkbox.
3. Select the field below the checkbox and enter your comment.
4. Select the **OK** button to save the comment.



Note...

To remove all comments, select the **Clear Comments** button. To close the dialog box without a change to the comments, select the **Cancel** button.

3.2 The Edit View Button and Markers

Select the **Edit View** button to update the DICOM header and enable the image to hang correctly at the Diagnostic Review Workstation. The tabs and buttons allow you to change just the Marker.



Figure 8-5: The Markers in the Preview Pane

3.2.1 How to Use the Marker function

1. If the marker is wrong, select the **Clear All Markers** button
2. Select the Tab for the set of markers you need.
3. Select the marker you need. You can select a maximum of three markers from the available set for each image.
4. Move the marker to the correct position.

If you change a Marker, you change the label on the image, but not the DICOM header. You will not see the change on the Patient View page.

3.2.2 How to Use the Edit View button

Before You Accept the Image

1. Select the **Edit View** button.
2. Select the required view from the displayed list.
3. Select **OK**. The correct image displays on the preview screen and the marker and the direction the image hangs are correct.

After You Accept the Image

If the image is Accepted before you notice the wrong view, the incorrect image was sent to the selected output. Make the corrections, then resend the corrected image with the Image or Spool Management menu. See Chapter 6, Section 8.2.1, page 65.

1. Before You Make Any Changes:
 - If you printed the mislabeled image, find and discard that image.
 - If the image was sent to a diagnostic review workstation, delete the mislabeled image.
 - If the image was sent to the PACS, tell the PACS administrator to delete the mislabeled image.
2. Repreview the Image (Chapter 6, Section 8.2.3, page 67).
3. Select the Edit View button on the Preview screen.
4. Select the required view from the displayed list.
5. Select **Save**. A new thumbnail image appears. Resend this image.

3.3 Implant Present

Select the checkbox on the Patient View screen when an implant exists. When the checkbox is selected, the system marks future images with an Implant Present label. The checkbox on this screen enables you to change the selection before you accept the image.

3.4 Quick Zoom/Pan

Quick Zoom/Pan is active after the image is first displayed, before image processing completes.

1. Select the **Quick Zoom/Pan** button.
2. To pan through the image, drag with the right trackball button.
3. To turn off the function, select the **Quick Zoom/Pan** button.

3.5 Full Zoom/Pan




Figure 8-6: Full Zoom/Pan Image

The Full Zoom/Pan function becomes available after image processing is complete. The box in the thumbnail image on right side of the screen shows where the magnification view is in the complete image. You can move this box to any area in the image.

1. Select the **Full Zoom/Pan** button.
2. To move through the image, drag with the right trackball button.
3. To turn off the function, select the **Full Zoom/Pan** button again.

3.6 ROI




Figure 8-7: ROI Size Drop-down Menu

Method 1: Draw a custom box.

1. Select the **ROI** button.
2. From the drop-down menu, select **User Draw**.
3. Click in the image and drag to draw the size you need.

Method 2: Use a pre-defined box size.

1. Select the **ROI** button.
2. From the drop-down menu, select the size.
3. Click the area of interest in the image.

3.7

Crosshairs




Figure 8-8: Preview with Crosshairs

1. Install a localization paddle.
2. Acquire the image.
3. If necessary, select the **Zoom/Pan** button and move through the image.
4. Select the **Crosshair** button.
5. Select the area of interest in the image.
6. Drag the crosshairs to the needed position, or select another area.
7. If necessary, use the **Zoom/Pan** functions to follow the crosshair lines to the grid from the localization paddle. The crosshair lines remain when the **Zoom/Pan** function is enabled and disabled.

3.8 Measurement




Figure 8-9: A Measurement on the Preview Screen

1. Select the Measurement button.
2. Click the first point then drag to the second point. The distance displays when you release the button.
3. Select the “+” at either end of the line to clear. There lines are not saved with the image.

3.9

Magnification



Use the Magnification function to view part of the image at 2X the displayed resolution.

1. Select the **Magnification** button.
2. Move the box that appears at the area of interest.
3. To turn off the magnification function, select the **Magnification** button a second time.



Note...

You can use the magnification function on a Full Zoom image.



Figure 8-10: Preview with Magnification

3.10 Window/Level



1. Select the **Window/Level** button.
2. Select in the image and hold while you move the trackball.
 - Roll Up or Down to change the brightness (Level).
 - Roll Right or Left to change the contrast (Window).
3. To return to the original settings, double click the image.
4. Select the **Window/Level** button to turn off the feature.



Note...

This operation does not change the final, processed Image.




Figure 8-11: The Window/Level Settings

To enter a known **Window/Level** value:

1. Select the field (Center or Width) to change.
2. Enter a new number.
3. Select the **Set** button.

To return to the original setting, double click the image.

3.11 Display AEC Regions

To display the AEC area selected by the system, select the **Display AEC Regions** checkbox.

- Manual sensors show on the image as a white rectangle with rounded corners.
- The two Auto Sensor areas show on the image as 1 cm by 1 cm white squares.

To hide the AEC areas, deselect the **Display AEC Regions** checkbox.




Figure 8-12: AEC Regions

3.12 Exposure Index

Use the Exposure Index as a general guide to compare the selected Exposure Techniques to the quality of the acquired image.

- “Image Status: OK” indicates the image processing was applied.
- “Image Status: Raw” indicates that image processing was not applied.
- A number after the image status indicates a problem with the image processing. Reject the image. Contact Service Support.




Figure 8-13: Exposure Information

Legend for Figure 8-13

1. Exposure Index-pixel count
 2. Low exposure area—evaluate the image for excess noise. Re-acquire, if necessary.
 3. Acceptable, but below ideal exposure area
 4. Ideal exposure area
 5. High exposure area—the image is more than sufficient. You can reduce the radiographic technique for the next exposures.
 6. The Exposure Index line related to imaging of the ACR phantom
- For Manual and TEC modes, the Exposure Index calculates an area approximately one centimeter from the breast wall on the centerline of the image.
 - For all AEC modes, the Exposure Index calculates for the AEC areas used.

3.13 Exposure and Dose/Entrance Surface Exposure Information




Figure 8-14: Exposure Information

The exposure information is displayed in the Tools area. An estimate of the average glandular organ dose and the entrance surface exposure are displayed.

4.0 Accept or Reject

Select Accept or Reject to close the Preview screen. The sequence continues by highlighting the next view.

4.1 Accept

To accept the image, select the **Accept** button.

- The full resolution image with all related attributes is marked as accepted in the pool.
- The system transmits the image to the selected output devices.

4.2 Reject

To reject the image:

1. Select the **Reject** button.



Note...

Reject Analysis uses this information (required for Quality Control).




Figure 8-15: The Reject Reasons List

2. Select the reason for the reject.
3. Select the **Reject** button. The system moves the image to the Reject Bin. An asterisk displays next to the View label of a rejected image.



Figure 8-16: A Rejected Image

5.0 Additional information about Accepted/Rejected Images

After you reject an image, the Preview closes and you can:

- Re-acquire the rejected View.
- Select a different View.
- Close the Procedure.

5.1 The Image Repetition Information Dialog Box

When you accept a second image for the same View, the Image Repetition Information dialog box appears.

1. Select the reason for the repeat.
2. Close the dialog box.

5.2 How to Accept a Rejected Image

To accept a rejected image:

1. Reject the image that was taken as the replacement for the rejected image.
2. In the Patient View screen, select the **Review** button (below the Output box).
3. In the **Review** dialog box, select the tab for the View.
4. Select the **Rejected Image** thumbnail image to accept.
5. Select the **Comments** button to add or delete comments.
6. Select the **OK** button to close the **Comments** dialog box.




Figure 8-17: Set Accepted Button on Review Screen

7. Select the **Set Accepted** button.
8. Select the **OK** button to close the **Review** dialog box.




Figure 8-18: Accept Rejected Image

When you accept a rejected image, the Send Image to Output dialog box displays. Select Yes to send the images to the selected output group.



Note...

*If you select the **No** button, you can resend the images with the Image/Spool Management function.*

5.3 How to Reject an Accepted Image

To reject the Accepted Image:

1. In the **Patient View** screen, select the **Review** button.
2. In the dialog box, select the tab for the View to reject.
3. Select the Accepted Image.




Figure 8-19: Set Rejected Button

4. Select the **Set Rejected** button.
5. Select a reason for the reject. See Figure 8-15, page 92.
6. Select the **Reject** button.
7. Select the **OK** button. The displayed thumbnail image disappears from the Patient View screen.



Note...

*After you reject an Image with this method, the image is **not** removed from the output queues. To delete the image from the queues, see Table 6-6, page 68. The Rejected Image is added to the Reject Analysis Bin for tracking purposes.*

Chapter 9 - Maintenance and Cleaning

1.0 General Information

Call the Hologic Technical Support telephone number for the current list of recommended cleaning solutions.

Before each examination, clean and use a disinfectant on any part of the system which touches a patient. Give the attention to the paddles and the Image Receptor Device.



Caution:

Do not use any heat source (like a heating pad) on the image receptor.

Be careful with the plastic compression paddles. Inspect the paddles. Replace the paddle when you see damage.

1.1 For General Cleaning

Use a lint-free cloth or pad and apply a diluted dishwashing liquid.



Caution:

The fluid must not drip or run.

If more than soap and water is required, Hologic recommends any one of the following:

- 10% chlorine bleach and water with one part commercially available chlorine bleach (normally 5.25% chlorine and 94.75% water) and nine parts water
- Commercially available isopropyl alcohol solution (70% isopropyl alcohol by volume, not diluted)
- 3% maximum concentration of hydrogen peroxide solution

After any of the above solutions, use a pad and apply a diluted dishwashing liquid to clean any parts which touch the patient.



Caution:

Do not spray disinfectant on the system, because the moisture can enter into the system and damage the electronic components.



Warning:

If a paddle touches any possible infectious materials, call your Infection Control Representative for decontamination instructions.

1.2 To Prevent Possible Injury or Equipment Damage

Never use a corrosive solvent or abrasive detergents or polishes. Select a cleaning agent which will not damage plastics, aluminum, or carbon fiber.

- Do not use a strong detergent, abrasive cleaner, high alcohol concentration, or methanol at any concentration. If skin preparation contains a high alcohol concentration, allow time to dry before compression.
- Do not expose equipment parts to any steam or high temperature sterilization.

- Never allow any liquids to enter the internal parts of the equipment. Do not apply cleaning sprays or liquids to the equipment. Always use a clean cloth and apply the spray or liquid. If liquid enters the system, disconnect the electrical supply and schedule inspection by service personnel before the system is returned to use.



Caution:

Wrong cleaning methods can damage the equipment, decrease imaging performance, or increase the risk of electric shock.

- Always follow the instructions from the germicide manufacturer. The instructions include the methods of concentration, storage, application, contact time, wash requirements, protective clothing, shelf life, and for removal. The instructions describe the best and safest product use.

2.0 Care and Cleaning—Acquisition Workstation and Value Console

2.1 Display

Avoid touching the display screen.

Care should always be used when cleaning the outer surface of the LCD screen. Always use a clean soft, lint-free cloth. Microfiber cloths, available at most camera stores, are highly recommended.

- Never spray or pour a liquid directly onto the screen.
- Never apply excessive pressure to the screen.
- Never use detergents with fluorides, ammonia, alcohol, or abrasives.
- Never use bleach.
- Never use steel wool, or cloth woven with metal.
- Never use a sponge with abrasives.

There are many commercially available products that are specifically designed for cleaning LCD displays. Any of the products free of the ingredients mentioned above and used in accordance with the manufacturer's directions can be used.

2.2 Keyboard and Trackball or Mouse

Wipe clean the keyboard and trackball or mouse surfaces using a standard CRT wipe. If necessary, clean the keyboard with a vacuum. If liquids spill into the keyboard, call Service for a replacement.

3.0 Preventive Maintenance

Table 9-1: User Preventive Maintenance

Maintenance Task Description	Recommended Frequency						
	Each Use	Daily	Weekly	Biweekly	Monthly	Quarterly	Semi-annually
Clean & disinfect paddle	x						
Clean & disinfect breast platform	x						
Visually inspect all paddles for damage	x						
Recommended/ required start of operation verifications		x			x		
All daily recommended/ required calibration		x					x
Diagnostic Review Workstation Quality Control *			x				
Detector Flat Field Calibration *			x				
Artifact Evaluation *			x				
Phantom Image *			x				
Signal to Noise / Contrast to Noise Measurements *			x				
Geometry Calibration (Tomosynthesis Option) *							x
Compression Thickness Indicator *				x			
DICOM Printer Quality Control			x				
View boxes and Viewing Conditions			x		x		
Repeat/Reject Analysis*						x	
Compression *							x
Visual Checklist *					x		
Emergency Stops					x		
Clean array dust filter					x		

* Refer to Quality Control Manual

Instructions for Use

Chapter 9 - Maintenance and Cleaning
Preventive Maintenance

Appendix A—System Specifications

1.0 Dimensional Information

1.1 Tubestand (Gantry with C-arm)




Figure A-1: Tubestand Dimensions

Height	190.0 cm (74.8 in.) \pm 1.0 cm (0.4 in.); 213.0 cm (84.0 in.) \pm 1.0 cm (0.4 in.) at highest C-arm travel
Width	65.0 cm (25.6 in.) \pm 1.0 cm (0.4 in.) (C-arm at 0° position)
Depth	110.0 cm (43.3 in.) \pm 1.0 cm (0.4 in.)
Weight	274 kg (604 lb) \pm 20 kg (44.0 lb)

1.2 Acquisition Workstation




Figure A-2: Acquisition Workstation Dimensions

Overall Dimensions	189.2 cm (74.5 in.) (H) x 107 cm (42.1 in.) (W) x 76 cm (30.0 in.) (D) (maximum)
Weight	154 kg (340 lb) (maximum)

1.3 Acquisition Workstation with Dual Swivel Arms




Figure A-3: Acquisition Workstation with Dual Swivel Arms

Overall Dimensions 195.6 cm (77.0 in.) (H) x 107 cm (42.1 in.) (W) x 76 cm (30.0 in.) (D) (maximum)
Weight 218 kg (480 lb)

1.4 Selenia Value Console




Figure A-4: Value Console Dimensions

Overall Dimensions 189.4 cm (74.6 in.) (H) x 101.2 cm (39.8 in.) (W) x 60.5 cm (23.8 in.) (D) (maximum)
Weight 163.3 kg (360 lb) (maximum)

2.0 Operating Environment

2.1 General Operating Conditions

<i>Temperature Range</i>	20 °C (68° F) to 30 °C (86° F)
<i>Maximum rate of temperature change</i>	<10 °C/hr
<i>Relative Humidity Range</i>	10% to 80% non-condensing
<i>BTU output</i>	Typical Range 1700–2500 BTU/hr

3.0 Storage Environment

3.1 Tubestand

<i>Temperature Range</i>	-25 °C (-13° F) to +60 °C (140° F)
<i>Humidity</i>	Zero to 95% humidity—non-condensing (not packaged for outdoor storage)

3.2 Image Receptor

<i>Storage Temperature Range</i>	10 °C (50° F) to 40 °C (104° F)
<i>Maximum rate of temperature change</i>	<15 °C/hr
<i>Storage Humidity Range</i>	10% to 80% humidity— non-condensing (not packaged for outdoor storage)

4.0 Electrical Input

4.1 Tubestand

<i>Mains Voltage</i>	200/208/220/230/240 VAC nominal, true sinusoidal, (tap selectable at installation) ±10%, Single Phase, 50/60 Hz ± 5%, permanently wired
<i>Mains Impedance</i>	Maximum line impedance not to exceed 0.20 ohms for 220/230/240 VAC, 0.16 ohms for 200/208 VAC
<i>Maximum Power Consumption</i>	6.5 kVA for 5 second duration
<i>Standby Current</i>	2.0 A (maximum)
<i>Maximum Line Current</i>	35 A for 5 seconds (momentary rating per NEC 660.2)
<i>Recommended Branch Circuit Breaker Rating</i>	25 A minimum, refer to NEC 660.6

4.2 Acquisition Workstation

<i>Input Line Voltage</i>	100/120/200/208/220/230/240 VAC nominal, true sinusoidal, (tap selectable at installation), ± 10%
<i>Input Current</i>	8.0 A maximum @ 100/120 VAC 3.5 A maximum @ 200/208/220/230/240 VAC
<i>Frequency</i>	50/60 Hz ±5%
<i>Number of Phases</i>	Single

5.0 Acquisition Workstation Technical Information

<i>Operating System</i>	<i>Sun™ Solaris OS</i>
<i>Computer Memory</i>	<i>2 GB RAM minimum</i>
<i>Disk capacity</i>	<i>> 60 GB Image Storage</i>
<i>Storage Media</i>	<i>CD-RW Disks</i>
<i>Display Adapter Card</i>	<i>1600 x 1200 matrix minimum 8 bit gray scale display</i>
<i>Display</i>	<i>1600 x 1200 matrix minimum 450 cd/m2 nominal brightness</i>
<i>Network Interface</i>	<i>10/100 Base-T Ethernet</i>
<i>Remote Diagnostics</i>	<i>Internet</i>
<i>Graphical User Interface</i>	<i>X-ray exposure control Configurable mammographic Workflow Patient demographics Brightness and contrast control Magnification screen Pixel value readout QC test tools System Status Monitoring Error reporting Unattended archiving and printing</i>
<i>Radiation Shield Pb equivalence</i>	<i>Rated for a 0.5 mm Pb (lead) equivalence</i>

6.0 Tubestand Technical Information

6.1 C-arm

<i>Motorized C-arm Rotation Speed</i>	Variable speed (18° per second maximum). Rotation speed is service selectable: 50% to 100% in 5% increments. Motor Control provides soft start and dynamic braking.
<i>C-arm Rotation Range</i>	+195° +2° /-0° to 0° ±0.5° to -150° +0° /-2° (detent at 0°, ±90°). Rotation angle is displayed on both sides of Gantry.
<i>Vertical Travel</i>	74.5 cm (29.5 in.) total travel.
<i>Source-to-Image Distance (SID)</i>	66.0 cm (26.0 in.) ± 1.0 cm (0.4 in.)
<i>Source-to-Image Receptor Support Device Distance</i>	64.0 cm (25.2 in.) ± 1.0 cm (0.4 in.)
<i>Magnification Ratio</i>	1.8 x for objects 22.5 mm above the magnification platform breast support surface

6.2 Compression

<i>Manual Compression Force</i>	Limited to a maximum of 300 N +0/-89 N (67.4 lb +0/-20 lb) from 0° to +/-90° C-arm rotation. Not less than 169 N (38 lb) for a C-arm angle range greater than +150° and an angle less than -150°.
<i>Motorized Compression</i>	Functions in three operating modes, Pre-compression, Full-Range, Dual Compression, user selectable through software.
<i>Pre-Compression Force</i>	67 N +0/-22.3 N to 133.5 N ±22.3 N (15 lb +0/-5 lb to 30 lb ±5 lb)
<i>Full-Range Compression Force</i>	89 N ±22.3 N to 178 N ±22.3 N (20 lb ± 5.0 lb to 40.0 lb ± 5.0 lb)
<i>Dual Mode Compression</i>	Provides Pre-Compression force upon first activation of compression switch; then, if switch is activated within 2 seconds, the force is increased incrementally for each additional switch activation, up to the user selected FULL compression force.
<i>Compression Controls</i>	Up/Down controls on both sides of C-arm and on 2-position footswitch (Motorized). Handwheel on both sides of Compression Device (Manual).
<i>Compression Release</i>	Manual or Automatic. Motorized Release mode controlled by push-buttons on both sides of the C-arm. User selectable automatic release mode raises Compression Device upon exposure termination. All release functions are disabled if a Localization paddle is detected.
<i>Automatic Compression Release</i>	Moves the compression device upward a predetermined distance (10 cm). For Magnification Mode this may be less.

<i>Compression Down Motion Variable Speed</i>	<i>Service Selectable between approximately 10% through 100% of full speed.</i>
<i>Compression Force Display</i>	<i>Two LED Displays on the Compression Device show the compression force through the range of 10 lb to 67.4 lb (44.5 N to 300 N) in 1 lb (4.4 N) increments.</i>
<i>Compression Force Display Accuracy</i>	<i>±4.5 lb (±20 N)</i>
<i>Compression Thickness Display</i>	<i>Two LED Displays on Compression Device measures between 0 and 15 cm above image receptor in 0.1 cm increments. The display is visible from both sides of the patient.</i>
<i>Compression Thickness Accuracy</i>	<i>±0.5 cm for thicknesses between 0.5 cm and 15 cm</i>
<i>Compression Paddles</i>	<i>Compression paddles are transparent. The paddles are composed of polycarbonate. The paddles provide a parallel plane to the image receptor and do not deflect by more than 1 cm difference from any surface providing compression under 25 lb (111.1 N) of compression force (except for F.A.S.T. paddles.). The paddles are adjustable to provide the focal spot, compression device and image receptor alignment requirement.</i>

6.3 X-ray Tube: Molybdenum

<i>Focal Spot (NEMA / IEC)</i>	<i>Large (0.3 mm) Nominal Small (0.1 mm) Nominal</i>
<i>Tube Voltage</i>	<i>20 kV to 39 kV</i>
<i>Tube Current</i>	<i>Large Focal Spot = 100 mA between 25 and 32 kV Small Focal Spot = 30 mA between 25 and 32 kV</i>
<i>Anode Rotation</i>	<i>180 Hz (9600 RPM minimum)</i>
<i>Anode Angle</i>	<i>Bi-angular: Large focal spot at 16°, Small focal spot at 10°. X-ray tube angle at 6° to provide 22° (Large FS) and 16° (Small FS) anode to Image Receptor plane angle.</i>
<i>Anode Material</i>	<i>Molybdenum</i>
<i>X-ray Window</i>	<i>Beryllium 0.8 mm thickness (maximum)</i>

6.3.1 X-ray Tube Housing

<i>Over Temperature Protection Sensor</i>	<i>Internally connected in series with the stator common lead.</i>
<i>Maximum Temperature, Tube Housing Surface</i>	<i>55 °C (151° F)</i>
<i>Maximum Temperature, Tube Head Cover Surface</i>	<i>41 °C (105.8° F)</i>
<i>Safety Class</i>	<i>IEC 60601-1, Class I, IEC 60601-2-28</i>

6.3.2 X-ray Beam Filtration and Output

<i>Inherent Tube Filtration</i>	<i>0.0 mm Al equivalent</i>
<i>Added Filtration</i>	<i>Two-position filter changer mechanism to carry a 30-micron (0.03 mm) Molybdenum foil filter and a 30-micron (0.03 mm) Rhodium foil Filter. Type of filter is user selectable.</i>
<i>Beam Quality HVL for Mo/Mo Operation</i>	<i>At a given kilovolts, the measured HVL with the compression paddle in the x-ray beam is equal to or greater than the value of $kV/100+0.03$ (in units of mm of aluminum) but less than the value of $kV/100+0.12$ (in units of mm of aluminum).</i>
<i>Beam Quality HVL for Mo/Rh operation</i>	<i>At a given kilovolts, the measured HVL with the compression paddle in the x-ray beam is equal to or greater than the value of $kV/100+0.03$ mm Al (in units of mm of aluminum) but less than the value of $kV/100+0.19$ mm Al (in units of mm of aluminum).</i>
<i>Radiation Output</i>	<i>Equal to or greater than 800 mR/second for at least 3 seconds. Output is measured through the Compression Paddle, 4.5 cm above the breast support, 4.0 cm from the chest wall edge, using exposure techniques of Mo/Mo target/filter, large focal spot, 28 kV.</i>

6.4 X-ray Tube: Tungsten

<i>Focal Spot (NEMA / IEC)</i>	<i>Large (0.3 mm) Nominal Small (0.1 mm) Nominal</i>
<i>Tube Voltage</i>	<i>22 kV to 39 kV</i>
<i>Tube Current</i>	<i>Large Focal Spot = 100 mA between 25 and 32 kV Small Focal Spot = 30 mA between 25 and 32 kV</i>
<i>Anode Rotation</i>	<i>180 Hz (9500 RPM minimum)</i>
<i>Anode Angle</i>	<i>Bi-angular: Large focal spot at 16°, Small focal spot at 10°. X-ray tube angle at 6° to provide 22° (Large FS) and 16° (Small FS) anode to Image Receptor plane angle.</i>
<i>Anode Material</i>	<i>Tungsten</i>
<i>X-ray Window</i>	<i>Beryllium 0.8 mm thickness (maximum)</i>

6.4.1 X-ray Tube Housing

<i>Over Temperature Protection Sensor</i>	<i>Internally connected in series with the stator common lead.</i>
<i>Maximum Temperature, Tube Housing Surface</i>	<i>55 °C (151° F)</i>
<i>Maximum Temperature, Tube Head Cover Surface</i>	<i>41 °C (105.8° F)</i>
<i>Safety Class</i>	<i>IEC 60601-1, Class I, Type B, IEC 60601-2-28</i>

6.4.2 X-ray Beam Filtration and Output

<i>Inherent Tube Filtration</i>	<i>0.0 mm Al equivalent</i>
<i>Added Filtration</i>	<i>Two-position filter changer mechanism to carry a 60-micron (0.06 mm) Rhodium foil filter and a 60-micron (0.06 mm) Silver foil Filter. Type of filter is user selectable.</i>
<i>Beam Quality HVL for W/Rh Operation</i>	<i>At a given kilovolts, the measured HVL with the compression paddle in the x-ray beam is less than the value of kV/100 +0.03 (in units of mm of aluminum).</i>
<i>Beam Quality HVL for W/Ag operation</i>	<i>At a given kilovolts, the measured HVL with the compression paddle in the x-ray beam is less than the value of kV/100 +0.03 (in units of mm of aluminum).</i>
<i>Radiation Output</i>	<i>Equal to or greater than 230 mR/second for at least 3 seconds. Output is measured through the Compression Paddle, 4.5 cm above the breast support, 4.0 cm from the chest wall edge, using exposure techniques of W/Rh target/filter, large focal spot, 28 kV.</i>

6.5 X-ray Collimation

Available Collimated X-ray Fields:

Table A-1: Collimation Settings

Predefined Collimation Settings
24 x 29 cm
18 x 24 cm
15 x 15 cm
10 x 10 cm
7.0 x 8.5 cm

6.6 Light Field Indication

Light Field Lamp

Illuminates for 30 seconds, ±5 seconds, upon pressing a Light Field switch on either side of the x-ray tubehead or by pressing a Compression Down switch. Extinguishes automatically upon exposure initiation. A shatter shield is provided.

Light Field Illuminance

160 lux (minimum)—meets 21 CFR 1020.31 requirements. Lamp is adjustable to provide alignment of the light field to the x-ray field.

Light Field-to-X-ray Field Congruency

Within 2% of SID

6.7 X-ray Generator

Output Rating

3.2 kilowatt, maximum (isowatt), 100 mA @ 32 kV

Ripple

2% or less (typical), maximum 4%

Topology

Pulse width modulated High Frequency, active servo controlled

6.7.1 kV /mA Range

Table A-2: mA setting as a function of kV

Large Focal Spot	
kV	mA
20*	75 mA
21*	80 mA
22	85 mA
23	90 mA
24	95 mA
25-32	100 mA
33	85 mA
34-35	80 mA
36-37	75 mA
38-39	70 mA
*Molybdenum Tubes only	

Table A-3: mA setting as a function of kV

Small Focal Spot	
kV	mA
20	20 mA
21	22 mA
22	24 mA
23	26 mA
24	28 mA
25-32	30 mA
33-34	28 mA
35-37	26 mA
38-39	24 mA



Note...

All above values are for mAs values ≥ 40 mAs. For mAs values < 40 mAs, the mA value is adjusted by the mA Factor listed in the table below as a function of mAs.

Table A-4: mA Factor as a function of mAs

mAs	mA Factor
4-6	10%
8-10	20%
12-22	30%
24-38	60%
≥ 40	100%

Large Focal Spot Manual mAs Range:

TABLE 1 (default): 3 mAs through 400 mAs, 22 steps: 3, 4, 5, 6.4, 8, 10, 12.5, 16, 20, 25, 32, 40, 50, 64, 80, 100, 125, 160, 200, 250, 320, 400.



Note... Lower limit is 4 mAs when grid is in field.

TABLE 2 (user selected): 3 mAs through 400 mAs, 55 steps: 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32.5, 35, 37.5, 40, 42.5, 45, 47.5, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 110, 120, 130, 140, 150, 160, 170, 180, 190, 200, 220, 240, 260, 280, 300, 325, 350, 375, 400.



Note... Lower limit is 4 mAs when grid is in field.



Note... In AEC modes, in Large Focal spot, the mA setting adjusts downward, if necessary, to achieve the calculated mAs. mA can be adjusted as low as 10 mA.

Small Focal Spot Mag Manual mAs Range:

TABLE 1: 3 mAs through 100 mAs, 16 steps: 3, 4, 5, 6.4, 8, 10, 12.5, 16, 20, 25, 32, 40, 50, 64, 80, 100

TABLE 2: 3 mAs through 120 mAs, 38 steps: 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32.5, 35, 37.5, 40, 42.5, 45, 47.5, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 110, 120

6.7.2 Accuracy, Reproducibility, and Linearity

Reproducibility	< 0.05 coefficient of variation for 10 consecutive exposures (21 CFR).
Linearity	< 0.10 for adjacent mAs selections per the following: (X1-X2) is less than or equal to 0.10 (X1+X2) where X1 and X2 are average mR/mAs values for consecutive exposures (21 CFR).
mAs Accuracy	±5% or ±2 mAs, whichever is greater, from indicated, measured from the ground side of the tube circuit.
kV Accuracy	Within 1 kV of the indicated kV

7.0 Imaging System Technical Information

7.1 Image Receptor

<i>Fluid ingress</i>	<i>No fluid from incidental spillage on the top surface of the Image Receptor seeps inside</i>
<i>Deflection</i>	<i>Does not exceed 1.0 mm at maximum compression</i>
<i>Active Imaging Area</i>	<i>24 cm x 29 cm Nominal. The active image area is marked on the digital image receptor/breast platform cover.</i>
<i>Image Size, Screening and Diagnostic Exams</i>	<i>18 x 24 cm nominal; locations: center, left, right</i>
<i>Image Size, Diagnostic-Spot Compression</i>	<i>24 x 29 cm nominal; center location only</i>
<i>Image Size, Diagnostic-Magnification</i>	<i>18 x 24 cm nominal; locations: center, left, right</i>
<i>Image Pixel Size</i>	<i>18 x 24 cm nominal; center location only</i>
<i>Digital Image Receptor MTF - Nyquist frequency</i>	<i>70 μm</i>
<i>DQE at 0 c/mm</i>	<i>> 0.40</i>
<i>DQE at Nyquist frequency</i>	<i>50% or greater at 7.0 mR -0.0/+0.7 mR x-ray exposure</i>
<i>Saturation</i>	<i>15% or greater at 7 mR -0.0/+0.7 mR x-ray exposure</i>
<i>Dynamic Range</i>	<i>X-ray exposure level at which image pixels are saturated is not less than 1000 mR</i>
<i>Output Image</i>	<i>Linear response over at least 400:1 in x-ray exposure</i>
<i>Scatter Rejection</i>	<i>14 bit image data</i>
<i>Chest Wall Access</i>	<i>Lorad HTC™ high transmission cellular grid</i>
<i>Lateral Wall Access</i>	<i>The distance from the outside edge of the Image Receptor enclosure to the Active Image Area along the chest wall is less than 5 mm.</i>
<i>ACR Phantom Score at MGD = 2mGy</i>	<i>The distance from the outside edge of the detector enclosure to the active detector area along the edges perpendicular to the chest wall is less than 40 mm.</i>
<i>Image Preview Time</i>	<i>At least 5 fibers, 4 specks, 4 mass</i>
<i>Transmission Limit</i>	<i>The time between completion of an x-ray exposure and availability of the Preview image: less than 20 seconds.</i>
	<i>Within Federal Regulatory limit for screen-film mammography systems (21CFR 1020): 0.1 mR/h.</i>

8.0 Tissue Exposure Control (TEC) Mode—Enhanced Manual Mode

<i>The Breast Density Default Setting</i>	<i>Selectable via the Generator default screen on the Acquisition Workstation. Breast density defaults to the default setting at the beginning of each new study and when the exposure mode is first changed to TEC.</i>
<i>Breast Density Panel</i>	<i>Update rate once every 2 seconds maximum</i>
<i>Mammography Unit Status Polling</i>	<i>Between image capture cycles, Acquisition Workstation requests and displays current technique at a rate not slower than once every 2 seconds. In TEC mode, the rate increases to no slower than once per second, but not faster than once per 0.75 seconds. System response time for a single Generator Technique command is 1 second or less.</i>

9.0 Automatic Exposure Control (AEC)

<i>AEC Sensor</i>	<i>Positionable in seven locations, centered laterally in the image receptor support device. Position #1 is located 1 cm from the chest wall. Positions are spaced in 1.7 cm increments. The AEC region can be selected automatically by processing the pre-exposure image data to find the corresponding dense portion of the breast. When the AEC region auto-selection is enabled, the detector area that is scanned is dependent on the default X-ray collimation that is linked to the compression paddle ID. Change of collimation by the operator is ignored for AEC region auto-selection. LEDs for each AEC sensor position on the compression gantry indicate the selected position. The selected AEC region can be viewed after selecting the check box. The size and available position of the manual AEC regions are indicated at the x-ray input surface of the compression paddle.</i>
<i>AEC Auto-Filter Mode</i>	<i>The filter and kV are selected using the recommended table based exposure technique based on compressed breast thickness. The exposure is terminated at an mAs value as determined by the AEC system to yield a pixel count to which the unit has been calibrated.</i>
<i>AEC Auto-kV Mode</i>	<i>The Filter is not selectable. The filter defaults to Mo if the system is configured with Moly Tube. The filter defaults to Rh if the system is configured with a Tungsten Tube. Starting kV is determined by the Gantry. mAs is determined by the Digital Detector.</i>

<i>AEC Auto-Time Mode</i>	<i>The kV and filter are both user selectable. The exposure is terminated at a mAs value as determined by the AEC system to yield a pixel count to which the unit has been calibrated.</i>
<i>Exposure Termination</i>	<i>An exposure is aborted or prevented if: The predicted exposure time exceeds the x-ray tube limits. The predicted exposure time exceeds the safety backup timer setting. The predicted exposure time is less than 400 msec when the grid is in the x-ray field or less than 30 msec when the grid is out of the x-ray field. The user releases the exposure button prematurely. There is an error within the Selenia FFDM system.</i>
<i>AEC Reproducibility</i>	<i>AEC exposure reproducibility meets MQSA requirements as follows: For a 4 cm thick average breast phantom, coefficients of variation in mAs and digital values shall be less than 0.05 for four consecutive exposures of this phantom under each of the three AEC modes.</i>
<i>Exposure Adjustment</i>	<i>There is a “user-adjustable” scale factor that can be used to increase/decrease the final x-ray exposure (mAs) by 15% for each step from the default setting. A total of eight settings are available on the Acquisition Workstation, which allows the user to change the exposure from -55% of the default value to 160% of the default value.</i>

10.0 Hardcopy Film Printing Devices

This section lists the relevant technical characteristics for the hardcopy film printing device.

10.1 Image Engine

<i>Film sizes</i>	<i>20 cm x 25 cm (8 in. x 10 in.) 25 cm x 30 cm (10 in. x 12 in.), optional</i>
<i>Pixel size</i>	<i>≤70 μm</i>
<i>Grey level resolution</i>	<i>≥ 8 bits</i>
<i>Maximum Film Optical Density</i>	<i>≥ 3.5</i>

10.2 Interface

<i>Connectivity</i>	<i>DICOM 3.0 compatible</i>
<i>Data port</i>	<i>Ethernet</i>
<i>Transfer rate</i>	<i>Minimum: 10 Mbits/sec</i>
<i>Look up Tables (LUT)</i>	<i>Field programmable; Linear</i>

10.3 Printed Film

<i>Print to fit</i>	<i>If a 24 x 29 cm image is printed on a 20 x 25 cm (8 x 10 in.) film, the printer prints the image to fit onto a 20 x 25 cm (8 x 10 in.) film without cropping the image or without dropping rows and columns of data</i>
<i>Patient Information</i>	<i>Film meets ACR/MQSA requirements for patient demographics. Printed Information Is User Configurable.</i>

11.0 TechMate

<i>CPU</i>	<i>Dual Intel® Xeon® Processors</i>
<i>Operating System</i>	<i>Microsoft Windows®</i>
<i>Computer Memory</i>	<i>2 GB RAM minimum</i>
<i>Disk capacity</i>	<i>400 GB SATA minimum</i>
<i>Display</i>	<i>3 Mega Pixel, 21 in. LCD, 2048 x 1536 resolution</i>
<i>User Interface</i>	<i>Keyboard, Trackball</i>

Appendix B—The Mobile Selenia

This appendix provides information about Selenia systems installed in a mobile environment.

1.0 General Information

1.1 The Vertical Position Override Switch (VPOS)

The Mobile Selenia has a Vertical Position Override Switch (VPOS) on the back of the Gantry (see Figure B-1.) This switch can start the Selenia if the C-arm is lower than its normal operating position after rough-road travel.




Figure B-1: Location of Vertical Position Override Switch (VPOS)

1.1.1 When to use the VPOS

This condition shows on the first startup after transport. The Acquisition Workstation starts normally, but the Gantry does not startup. After Logon, an alert message appears that communication to the Generator is lost.

1.1.2 How to Use the Vertical Position Override Switch to Restart the Gantry

1. Press and hold the Vertical Position Override Switch—the Gantry starts. Continue to hold the switch pressed.
2. After approximately 5 seconds, lift the C-arm approximately 7.5 cm (3 inches) with the normal system controls. Release the VPOS after you lift the C-arm.
3. At the Acquisition Workstation, accept the prompt to restart the Generator. The Gantry turns off for a few seconds, then automatically restarts.
4. The system is ready for normal use.

2.0 Safety Conditions and Other Precautions



Warning:

The Selenia radiation shield is not approved for mobile use and is not provided. The coach manufacturer must provide radiation shielding that meets requirements.



Caution:

When shore power is not available, use mobile power sources that provide the equivalent performance (see “Appendix A”). When continuous true sinusoidal VAC power is supplied per the specifications for Selenia power input and loading, system function and performance is correct. The power source must provide a minimum of 35A at 208VAC for a minimum of 5 seconds and support this load one time every 30 seconds. Separate, dedicated circuits must supply power to the Acquisition Workstation and Gantry. Use an uninterruptible power supply with active line conditioner for each power circuit. Other circuits must supply all secondary mobile coach power. After the first installation and when the mobile coach is moved, verify that the electrical installation meets the specifications of system power input and IEC 60601-1 safety requirements.



Caution:

Control the temperature and humidity inside the vehicle at all times. When the unit is not in use, do not allow environmental parameters to be more than specified specifications.



Caution:

When in use, connect the Selenia to a PACS system or an acceptable hard copy printer for permanent archive.

3.0 Mobile Specifications

3.1 Shock and Vibration Limits

Vibration Limit	Not greater than 0.35 G (2 Hz to 200 Hz), as measured at mounting point of system to coach.
Shock Limit	Not greater than 1.0 G (½ sine pulse), as measured at mounting point of system to coach. An “air ride” coach suspension is recommended.

3.2 Coach Environment

3.2.1 Operating Environment

Temperature Range	20 °C (68° F) to 30 °C (86° F)
Maximum Rate of Temperature Change	<10 °C/hr
Relative Humidity Range	10% to 80% non-condensing
BTU output	Typical Range 1700–2500 BTU/hr

3.2.2 Non-operating/Transit Environment

Temperature Range	15 °C (59° F) through 35 °C (95° F)
Maximum Rate of Temperature Change	<15 °C/hr
Relative Humidity Range	10% to 80% non-condensing

3.3 Electrical Input

3.3.1 Tubestand

Mains Voltage	200/208/220/230/240 VAC nominal, true sinusoidal, (tap selectable at installation) ±10%, Single Phase, 50/60 Hz ± 5%, permanently wired
Mains Impedance	Maximum line impedance not to exceed 0.20 ohms for 220/230/240 VAC, 0.16 ohms for 200/208 VAC
Maximum Power Consumption	6.5 kVA for 5 second duration
Standby Current	2.0 A (maximum)
Maximum Line Current	35 A for 5 seconds (momentary rating per NEC 660.2)
Recommended Branch Circuit Breaker Rating	25 A minimum, refer to NEC 660.6

3.3.2 Acquisition Workstation

Input Line Voltage	100/120/200/208/220/230/240 VAC nominal, true sinusoidal, (tap selectable at installation), ± 10%
Input Current	8.0 A maximum @ 100/120 VAC 3.5 A maximum @ 200/208/220/230/240 VAC
Frequency	50/60 Hz ±5%
Number of Phases	Single

4.0 How to Prepare the System for Transport

Perform these steps before you transport the Mobile System:

1. Rotate the C-arm to 0 degrees (CC position).
2. Lower the C-arm to its lowest position.
3. Remove all power from the system.
4. Remove all system accessories (like the Face Shield, and Compression Paddles) before transport. Store all accessories correctly to make sure the components are not damaged.

5.0 How to Test the System Integrity After Transport

Instructions: Photocopy this checklist for use after the system is relocated.

Name _____ Date _____

System Serial Number _____

Table B-1: Integrity Checklist

Evaluated	Control Functions	Reference
Mechanical System Tests		
	Compression Up/Down	Chapter 3, "Compression Down," page 26
	Compression Release	Chapter 3, "Compression Release," page 26
	C-arm Rotation	Chapter 3, "C-arm and Gantry Rotation Switches," page 25
	C-arm Up/down	Chapter 3, "C-arm Up," page 27
	Collimator Override	Chapter 3, "Collimator Override," page 28
	Light Field Lamp	Chapter 3, "Light Field Lamp," page 27
	Smart Paddle System	Chapter 3, "Smart Paddle," page 28
	Emergency Off Switches	Chapter 2, "Emergency Off Switches," page 13
Quality Control Tests		
	Artifact Evaluation	Selenia QC Manual
	SNR/CNR Measurement	Selenia QC Manual
	Phantom Image Evaluation	Selenia QC Manual
	Compression Thickness	Selenia QC Manual

Index

A

accept

- image 91
- rejected image 93

Accept button 44–46, 91

accessories 1, 4

- available 33
- compression paddles 35–39
 - Smart Paddles 28, 36
- magnification platform 40

acquire images 46, 77

Acquisition Workstation

- circuit breaker location and function 12
- description 11

add

- additional view to procedure 45, 76
- new entry with function key 16
- new patient 49
- new procedure with new accession number 51
- user 62

AEC

- alarm messages 58
- Exposure Adjustment settings 57
- modes 57
- region, display on preview screen 90
- select Exposure Mode with keyboard 16
- sensor position 16, 20–21, 28
- sensor position settings 58
- using 57

alarms 73

- AEC messages 58

annotation tools list 81

archive, commit to storage function enabled 65

asterisk, rejected image 92

AUTO AEC sensor

- display region used 90
- use 58

average glandular organ dose display 91

B

bar code scanner 3, 14, 16

Barco display 17

breast density 16

- selecting in TEC 56

button

- compression release 15
- Eject, on CD-RW drive 69
- on and x-ray 15

button (software)

- Accept 49, 52, 91
- Add View 45, 76
- Brightness/Contrast (Window/Level) 89
- Cancel/Clear 49, 52
- Clean Up 71
- Clear All Markers 83
- Comments 82
- Crosshairs 86
- Edit View 83
- Full Zoom/Pan 84
- Implant Present 84
- Import Priors 70
- Magnification button 88
- Markers 83
- Measurement 87
- New Patient 49
- New Procedure 51
- Query Retrieve Priors 70
- Quick Zoom/Pan 84
- Reject 93
- Remove Unexposed Views 76
- Resend 66
- Resend Stopped Job 69
- Retrieve and Import 70
- Review 93
- ROI (Region of Interest) 85
- Set Accepted 93
- Set Rejected 94
- Simulate Capture 73
- Window/Level (Brightness/Contrast) 89

C

calibration messages 70

care and cleaning 96

C-arm

- components 11
- controls and displays 19–20
- Gantry rotation switches 19
- rotation switches 20, 25
- up and down 27

CD, CD-RW

- available space on 72
- ejecting disks 69
- importing images from 69

CD-RW Drive

- Eject button 69
- location 14

change Patient Information on Accepted Image 50

circuit breaker 4

- Acquisition Workstation and Value Console 14
- Gantry 12, 19

Clean Up button 71

cleaning, recommended 95

Clear All Markers button 83

Close

- Exam (Procedure) key 16
- Patient button 52
- Procedure 46, 76

collimator

- automatically set to paddle size 60
- override 28

Comments

- button 93
- entering or modifying 82

commit function, archive 65

complaints, product xvi

compliance 6–7

- label locations 8–9

compression

- automatic release 103
- controls and displays 21
- down/up 26
- force display accuracy 104
- functional tests 26
- release buttons 15
- thickness

- display 21

- display accuracy 104

compression paddles 35–39

- care 95
- cleaning 35
- FAST paddles 36
- installing and removing 35

contrast, change for printer 67

controls 12–21

- AEC sensor position 21
- automatic compression release 26, 103
- C-arm 20–21, 25, 27
- C-arm rotation 20, 25
- collimator override 20, 28
- compression 20–21, 26
- compression release buttons 15
- handwheels 21
- light field lamp 20, 27
- location and appearance 25
- manual compression 20

Crosshair device

- aligning to localization paddle 42
- installing and removing 41

Crosshairs 86

D

defaults

- exposure modes 54
- system 60
- system view order 61

definitions xvii

delete job from queue 68

delete user 62

determine job status 68

digital markers 80

dimensions 99

disk

- ejecting CD from drive 69
- space available 72

Dismiss button 24

display 14, 17

- AEC region used 90
- Barco description 17
- care and cleaning 96
- Totoku description 17

Dorsal/Ventral hanging 80

dose/entrance surface exposure information 91

E

edit

- patient data 50
 - on an accepted image 50
- patient data from MWL and PACS 50
- procedure accession number 52

edit user 62

Edit View button 83

Eject CD 69

electrical input 101

Emergency Off switches 14, 19

- location and function 13
- restart dialog box 29

enhancement tools list 81

enter new patient data 49

entrance surface exposure display 91

environment 101

equipment damage, prevention 95

exam room door indicators 5

exposure

- changing mode resets exposure techniques 55
- duration 5
- information 91
- ready indicators 77

exposure index 91

exposure modes

- AEC modes 57
- function and use 55
- manual 55
- TEC 56

exposure pane, location 54

exposure techniques

- AEC adjustments 57
- reset to system defaults 54
- select the mode 16
- select with function keys 16
- used 66

F

face shield 34

FAST paddles 36

files, select and unselect for import 69

film label area, location 80

film printer 11

filter, select with function key 16

flat screen display 17

focal spot, select 16

footswitches 3, 19, 21

- function 11

Full Zoom/Pan 84

function keys 15–16

functional tests 25–28

G

Gantry

- circuit breaker location and function 12
- components 11
- controls and displays 19
- description 11
- Emergency Off switches 19
- restart mobile with VPOS 115
- rotation switches for C-arm 19

Generator Status 77

grid in/out position toggle key 16

H

hanging options 80

Hologic technical support 95

I

image display area 79

image enhancement and annotation tools 81

Image Management

- Repreview 67
- Resend 66
- using 65–67

image processing display 91

image receptor

- aligning localization paddle with 42
- no heat on 3
- rotation switch 20
- storage temperature 101
- technical information 110

Image Repetition 93

image status 91

images

- comments 82
- film labeling area 80
- identifying 49
- import files from CD 69
- outputs 53
- Preview screen 79
- processed 91
- rejected 93
- remaining space for 72
- reprviewing 67
- resending 66
- Review dialog box thumbnails 93
- rules for automatic reclamation 72
- sizes 110
- tools 81
- viewing full resolution 88

Implant Present 76

button 84

import

image files 69

priors in background 70

indicator, exam door 5**infectious materials on equipment** 95**Info menu** 72**initials, where used** 62**initiating x-ray exposure** 15**integrity checklist, mobile** 118**intended uses** xv**interlocks** 6**international symbols** xix**J****job status, determine** 68**jobs in queue** 53, 75**K****keyboard** 15

care and cleaning 96

function keys, set exposure techniques 16

keypads 15**kV, select** 16**kV/mA range** 108**L****label locations** 8–9**laser film printer** 11**laser printer** 2**Launch dialog box** 24**light field lamp** 27**local database** 48**Local Exams, search** 48**localization paddle** 39

aligning crosshair device with 42

M**Magnification**

crosshair device 42

paddles 39

platform 40

Magnification button 88**maintenance, general** 95**Manage Queues** 68

delete (cancel) job status 68

determine job status 68

manual collimation, effects 60**Manual exposure mode** 55**Markers** 83

digital 80

mAs

focal spots manual range 109

select 16

Measurement button 87**Menu bar** 59**mobile information** 115–118**mobile system**

C-arm position switch 115

integrity checklist 118

safety considerations 116

specifications 117

temperature/humidity control requirements 116

tests, pre/post transport 118

Vertical Position Override Switch (VPOS) 115

modality worklist provider, query 48**monitor see display****monthly system tests** 29**mouse** 15

care and cleaning 96

N**New**

Patient 49

Procedure 51

New Entry function key 16**Non-Imaging Mode, use** 72**O****on and x-ray buttons** 15**operating conditions** 101, 117**Operating System (OS)** 24**organization, Instructions for Use** xx**outputs**

accepted images sent to 91

add new output group 63

change/select destinations 53

edit output group 64

send accepted/rejected images to 94

P**paddles**

identification 38

realigning front edge with image receptor 37

Panel Power dialog box 29

Patient

- add new 49
- edit data 50
- Face Shield 5
- protecting from reclamation 67
- unprotect 67

Patient Name, tab for local database 48

Patient View screen 43

pixel count 91

positioning, optimized for small breast with Smart Paddle 36

post start-up tests 25

power button, TechMate 18

power failure 1, 30

Power On

- button 15, 23
- indicator light 15

PPS status 70

- menu item 59

preventive maintenance 1, 97

Preview Image 43

- screen 79
- tools 79

printer 67

- copies of images 66
- print data on film 11
- test patterns 71

procedure

- add another
 - with different accession number 51
 - with same accession number 52
- close with function key 16

Progressive Preview 79

Protect Patients 67, 72

Q

Quality Control

- controls and functions tests 25
- Reject Analysis 92

query

- Local Exams 48
- Modality Worklist 44, 48
- retrieve priors 70
- retrieve priors, results 70

Query button 44

queues

- managing 68

Quick Zoom/Pan 84

R

radiation symbol, X-ray indicator 15

raw image status 91

Ready 77

recommended cleaning solutions 95

Region of Interest (ROI), button 85

region used to calculate exposure index 91

reject

- Accepted image 94
- image 92–93

Reject Analysis 92

Reject Bin 92

Reject button 46, 92, 94

rejected image, asterisk 92

Remove Unexposed Views 76

Repreview, using 67

reprint images 66

resend 66

Resend Stopped Job button 69

reset UPS 30

restart

- Gantry with mobile VPOS 115
- system after Emergency 29

restart system in Sleep mode 29

retract/engage grid 16

Retrieve Priors 70

Review button 93

Review dialog box

- thumbnails, accepted images 93
- view tab 93–94

review reject Accepted image 94

Rotation Angle displays 19

rotation switches

- C-arm and Gantry 25
- Gantry 19

S

safety

- isolation integrity 1
- mechanical 5
- mobile requirements 116
- radiation 5

screening exam sequence 44

screens

- Patient view 43
- Preview Image 43, 79
- Review dialog box 93
- Select a Patient 43
- View Order Editor 61

search

- hints for use 48
- Local Exams 44, 48
- retrieve priors 70

SecurView TechMate 11

select

- AEC position with keyboard 16
- AEC sensor position 58
- breast density in TEC mode 16
- exposure techniques 16
- grid position 16
- outputs 53

Select a Patient screen 43, 47

Selenia Display 14

send

- accepted rejected images to outputs 94
- images to an output-post procedure 66

sensor position, AEC 58

Set Rejected button 94

Set Technique Defaults screens 60

shield

- patient face 5
- radiation 5

shutdown, how to 32

Simulate Capture, button 73

Smart Paddle, functions 28

Smart paddles 36

SMPTE pattern 71

software version 72

specifications, mobile 117

spool

- accepted images 91
- full resolution image sent to 91

Standard Setup screens 60

status, image 91

stopped job in queue, resend 69

switch users 32

switches

- AWS and Value Console circuit breakers 12
- Gantry circuit breaker 19
- Vertical Position Override 115

symbols, international, defined xix

system

- controls and functional tests 25
- description 11
- restart in sleep mode 29
- tests after transport 118

system messages 73

system tests 29

T

TEC

- Exposure mode 56
- Override mode 57
- select Exposure Mode with keyboard 16

TechMate 11

- controls 18
- display 18
- power button 18, 23

technical information 102–110

- exposure modes 111
- printers 113
- TechMate 113

technologist review workstation 11

technologist, default view order 75

temperature (IR) 72

terms and definitions xvii

test patterns 71

tests

- after transport 118
- monthly system 29

thumbnails 93

tools, list of image enhancement and annotation 81

Totoku display 17

trackball 14–15

- care and cleaning 96

tubestand, components 11, 19

turning on system in Sleep mode 29

U

unprotect patients 67

UPS

- reset, TechMate 31
- resetting Selenia 30
- shut down Acquisition Workstation 30
- Value Console, location 14

UPS reset button, Acquisition Workstation 14

user

- add or edit 62
- delete 62

User Interface Screens 43

user setup function 62

User Training xvi

users, OS/Application names & passwords 24

V

Value Console

- circuit breaker location and function 12
- description 11

Vertical Position Override switch 115

view

- default order 75
- prior images on AWS 70
- rejected image 93
- selecting 75

View Order Editor screen 61

W

warnings, cautions and notes xix–6
defined xix

Window/Level 89

worklist 48, 50
tab 48

X

x-ray

- collimated fields 28, 107
- indicators 15

x-ray symbol 77

x-ray tube

- Molybdenum 104
- Tungsten 106

Instructions for Use

Index

X

Selenia®

Digital Mammography System