

# APEX™ Software



## Customer Release Notes

MAN-09980 Revision 001

**HOLOGIC®**

# APEX™ Software

Operating System for the Horizon® Series of Densitometers

## Customer Release Notes

For Software Version 5.6.1.3

Part Number MAN-09980

Revision 001

July 2023

### Product Support

USA: +1.877.371.4372

Asia: +852 37487700

Europe: +32 2 711 4690

Australia: +1 800 264 073

All Other: +1 781 999 7750

Email: [BreastHealth.Support@hologic.com](mailto:BreastHealth.Support@hologic.com)

© 2023 Hologic, Inc. Printed in the USA. This manual was originally written in English.

Hologic, InSight, Horizon, Discovery, and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners.

This product may be protected by one or more U.S. or foreign patents as identified at [www.Hologic.com/patent-information](http://www.Hologic.com/patent-information).



Hologic, Inc.  
600 Technology Drive  
Newark, DE 19702 USA  
1.800.447.1856  
[www.hologic.com](http://www.hologic.com)

Asia Pacific Hologic Hong Kong, Inc.  
7th Floor, Biotech Centre 2  
No. 11 Science Park West Avenue  
Hong Kong Science Park  
Shatin, New Territories  
Hong Kong

Australia /  
New Zealand Hologic (Australia) Pty Ltd.  
Suite 402, Level 3  
2 Lyon Park Road  
Macquarie Park NSW 2113  
Australia

## Table of Contents

<b>APEX Software Release Notes</b>	<b>4</b>
1.1 Version 5.6.1.3.....	4
1.2 Version 5.6.1.2.....	4
1.3 Version 5.6.1.1.....	5
1.4 Version 5.6.1.....	5
1.5 Version 5.6.0.7.....	6
1.6 Version 5.6.0.6.....	6
1.7 Version 5.6.0.5.....	6
1.8 Version 5.6.0.4.....	6
1.9 Version 5.6.0.3.....	6
1.10 Version 5.6.0.2.....	7
1.11 Version 5.6.0.1.....	7
1.12 Version 5.6.....	7
1.13 Version 5.5.3.1.....	8
1.14 Version 5.5.3.....	8
1.15 Version 5.5.2.2.....	8
1.16 Version 5.5.2.1.....	8
1.17 Version 5.5.2.....	8
1.18 Version 5.5.1.....	9
1.19 Version 5.5.....	9

# APEX™ Software Release Notes

APEX 5.X is the operating system for the Horizon® series of densitometers. For differences in results from previous versions of APEX, refer to the *Statement of Equivalency* (MAN-09016).



**Note**

Refer to the *Horizon User Guide* (MAN-08072-002) for operating instructions related to new features.

## 1.1 Version 5.6.1.3



**Notes**

APEX version 5.6.1.3 is for Windows® 10 operating system only.

Database sync (DBSYNC- APEX) is no longer compatible with Discovery systems.

**Enhancements/Changes**

- Add support for new Oriental motors for table and c-arm motion.

**Resolved Issues**

- Position of the 'Cancel' button on the 'Scan parameters' screen for Japanese language option.
- Duplicate SOP Instance UID in DICOM SR
- Failure restoring DxReport backup created with previous version
- Incorrect DxReport text for forearm scans w/ significant bone loss
- Intermittent issue of FRAX values less than 0 .1%
- Image Style dropdown not visible when ROI Coordinates selected
- CPLD version mismatch message displayed when CPLD is correct

## 1.2 Version 5.6.1.2



**Notes**

APEX version 5.6.1.2 is for Windows® 10 operating system only.

Database sync (DBSYNC- APEX) is no longer compatible with discovery systems.

**Enhancements/Changes**

- Add support for new HP® Z2 computer running Windows 10 LTSC 2019.

**Resolved Issues**

- Correct missing Radiation Dose Structured Reports (RDSR) fields.
- Correct intermittent issue where the C-arm does not move away from knees before lateral repositioning.

## 1.3 Version 5.6.1.1



### Notes

APEX version 5.6.1.1 is for Windows 10 operating system only.

Database sync (DBSYNC- APEX) is no longer compatible with discovery systems.

---

### Enhancements

- Resolve sporadic scan errors when running under Windows 10.
- Add MDD/MDR tier 1 language support.
- Add RDSR.

## 1.4 Version 5.6.1



### Note

APEX version 5.6.1 is for Windows 7 operating system only.

---

### Features

- Display of Horizon hip and spine scans in the discovery, legacy, and classic image styles.
- Export of patient database in a comma separated value (csv) file format.
- 3rd party developer access to scan images and DAP values.
- Custom buttons to launch user-defined applications.
- Whole body lean mass values *included in the* DxReport database.
- Pediatric reports to display height-adjusted Z-scores.
- Patient positioning guide for whole body scans that comply with NHANES recommendations.
- Improved support for localized languages.
- Change of default setting on body composition calibration to “Classic” calibration.

## 1.5 Version 5.6.0.7

### Features

- Added support for Windows 10-64.
- Resolved incorrect Ay motion: Table rising prior to C-arm being at the head position when performing IVA scans.

## 1.6 Version 5.6.0.6

- This version was not released.

## 1.7 Version 5.6.0.5

### Resolved Issues

- Incorrect Tx motion table continues to rise after aborting a radiographic uniformity test or whole body scan.
- Incorrect Tx motion table continues to rise after releasing the enable button during the lateral portion of an AP/lateral scan.
- Scans archived from a pre-5.6 version that had invalid data in the p-file could not be restored.
- Corrected reference curve database to use subtotal body lean mass in the pediatric lean mass vs BMC curve.
- DxReport does not display the 'Other' medications and medical conditions questionnaire answers.
- Patient names that contain Chinese characters do not display properly in patient lists.
- The database export filter 'new data only' exports all data.
- Incorrect error message "*Data in the following fields are either out of range or missing country code*" displayed when editing a patient on the patient confirmation screen.
- Updated application training scans (SE femur, whole body, and pediatric).

## 1.8 Version 5.6.0.4

- Enabled manufacturing flexibility.

## 1.9 Version 5.6.0.3

### Resolved Issues

- DICOM language option corrections.
- Corrected problem resulting in 'unable to find record' when analyzing legacy scan. Scans will need to be restored to be corrected.
- Correct abort when Auto-WB analysis fails.

## 1.10 Version 5.6.0.2

### Features

- Database synchronization added.

### Resolved Issues

- Extraneous database messages are no longer saved in the log.
- The call back list report now handles large databases.
- Detail worklist query not functioning.
- Local motion active message at startup on system running at 50 Hertz.
- DICOM Send rate of change reports generating message “*Attempt to scroll past end or beginning of data*”.
- Abort analyzing IVA scans on Japanese Windows 7.
- Reference curve source can be listed as 2012 BMDCS instead of 2012 BMDS/NHANES.
- Unable to import database with blank referring physician.
- Unable to import only QC data.
- Graph for whole body report for patient less than 20 correctly labeled as ‘Age vs. Sub Total’.

## 1.11 Version 5.6.0.1

### Resolved Issues

- DICOM extended dual hip rate of change report does not archive scan to PACS.
- If the DICOM queue is suspended when a scan is performed it does not restart until APEX is restarted.
- Reports not properly configured upon install.

## 1.12 Version 5.6

### Features

- Patient database converted from Access to SQL. Improves the speed of APEX when displaying lists of patients or scans.
- SE femur analysis and reporting enhanced to include the ability to annotate the image and add comments.
- Updated reference database for Caucasian and Black children with 2012 BMDCS data.
- New Pediatric reporting for hip, spine, and whole body exams. Allows for the configuration of these reports via the reporting tab under the System Configuration utility.
- Hologic Connect has been added to allow for remote service and support for those sites who have internet connection and sign up for the service.

### 1.13 Version 5.5.3.1

#### Resolved Issues

- Increased system stability and performance to mediate '*Device not ready for requested action*' or prolonged unresponsiveness of APEX.
- Corrected issues with DICOM filtered query/retrieve and disabling of MPPS.

### 1.14 Version 5.5.3

#### Features

- APEX reports T-score based on configurable ethnicity and gender in accordance with ISCD guidelines.
- APEX report configurable to report L2-L4 in accordance with JSBMR 2012 guidelines.
- BCA clinical report configurable to report weight in lbs.
- DICOM character set configurable.

### 1.15 Version 5.5.2.2

#### Features

- Adjustment of patient on/off position for Horizon-A.

### 1.16 Version 5.5.2.1

#### Features

- Support of lower pedestals in the Horizon-A.
- DICOM configuration modified to be able to generate any combination of hip rate of change and extended hip rate of change.
- Subscripts removed from DICOM SR.

### 1.17 Version 5.5.2

#### Features

- Support for the Horizon A.
- Updated infant whole body algorithm.
- HP whole body added to database export.
- Export of QC data enhanced.
- The adjusted LSC values displayed under the rate of change table.
- DICOM modified to not mark spine scan as archived with extended spine report.
- SOP instance UID for the DICOM SR modified to be different from DICOM send.
- Added fields to DICOM SR, including reference curve source, analysis method, height, weight, and ethnicity.
- Performed series sequence added to DICOM MPPS message.



- Checking added to ensure proper dual hip positioning.
- Computer power setting adjusted during scanning to prevent computer from going into sleep mode during scans.

### 1.18 Version 5.5.1

- This version was not released.

### 1.19 Version 5.5

#### Features

- Support for Horizon models Ci, C, Wi, W.
- SE Femur for the visualization of atypical femur fractures.
- Rulers available on the SE femur analysis for measurements of scan features.
- DxReport integration into APEX.
- Histogram of bone pixels to aid in spine analysis.
- Dual hip reports.
- Wide screen support.
- Legacy media access enhanced.
- HP whole body added to database export.
- Export of QC data enhanced.



#### Note

APEX 5.5.x can exchange data with APEX 3.5.x and 4.5.x. APEX 5.5.x can read data from APEX 4.0.x but APEX 5.5.x will update the data such that it will no longer be compatible with APEX 4.0.x.

---