
Brevera®

Breast Biopsy System

Система за биопсия на гърда Brevera
Sustav za biopsiju dojke Brevera
Système pro biopsii prsu Brevera
Brevera-brystbiopsisystem
Brevera-rintabiopsiajärjestelmä
Σύστημα βιοψίας μαστού Brevera
Brevera emlőbiopsziás rendszer
Brevera-brystbiopsisystem
System do biopsji piersi Brevera
Brevera Sistem de biopsie mamară
Систем за биопсију дојке Brevera

Système na biopsiu prsníka Brevera
Sistem za biopsijo dojki Brevera
Brevera bröstbiopsisystem
Brevera Borstbiopsiesystem
Brevera Sistema de biopsia de mama
Brevera Sistema de biopsia mamaria
Brevera Sistema di biopsia mammaria
Brevera Brustbiopsiesystem
Brevera Système de biopsie mammaire
Brevera 乳房生検システム

Hologic Biopsy Guidance System Setup Guide (GE Senographe® DS Stereotaxy and GE Senographe Essential Stereotaxy Systems) MAN-09401 Revision 002

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PLEASE READ PRIOR TO PERFORMING PROCEDURE:

This setup guide is designed to be used in conjunction with, not replace, the Brevera breast biopsy device instructions for use. Prior to performing the procedure, the physician must review and be familiar with the full operating instructions for the Brevera breast biopsy device and biopsy image guidance system. Consult the Brevera breast biopsy device instructions for use for full operational instructions, as well as any contraindications and safety information.

Please retain this setup guide. Note that while Hologic frequently updates the contents of this guide, medical information and technology can change rapidly. Therefore some of the information in this guide may be out of date. Use of these instructions with any product other than those specifically noted may result in injury to the patient and/or damage to the equipment.

Program the breast biopsy image guidance system for the Brevera breast biopsy device according to the needle setup instructions provided in the appropriate breast biopsy guidance system user manual. Specific parameters are listed below.

GE Senographe DS Stereotaxy and GE Senographe Essential Stereotaxy is Compatible with the following Brevera Components and Biopsy Accessories:

| Catalog Number | Description |
|------------------|---|
| BREVDRV | Brevera Driver |
| BREVADPTR | Brevera Adapter |
| BREVADPTRG | Brevera Adapter for GE Senographe® Stereotaxy |
| BREVSTYLBKKT | Stylet Bracket |
| BREV09 | Standard Stereotactic Biopsy Needle |
| BREVDISP09 | Standard Stereotactic Biopsy Needle |
| EVIVA_NG09R | Needle Guide |
| EVIVA_BUSHING_GE | Bushing for GE Senographe® Stereotaxy |

Setup of GE Senographe DS Stereotaxy and GE Senographe Essential Stereotaxy Systems **Vertical Approach** for Use with Brevera Biopsy

To Install the Brevera Adapter:

NOTE: The following is a general setup guide. Refer to the GE Senographe DS, GE Senographe Essential, GE Senographe DS Stereotaxy, or GE Senographe Essential Stereotaxy Operator Manuals as appropriate for more details.

1. Setup GE Senographe stereotaxy system for vertical approach using appropriate bracket(s). Refer to GE Senographe stereotaxy system operator manual for setup instructions.
2. Install the Brevera Adapter for GE, BREVADPTRG, onto the GE Senographe stereotaxy system. Check the Brevera Adapter components (reference Figure 1) to ensure they are secure. If the components are not secure call technical support.
3. Install the bushing, EVIVA_BUSHING_GE, into the lower needle guide holder of the GE Senographe stereotaxy system.
4. Install the disposable needle guide, EVIVA_NG09R, into the Bushing.
5. Install the Brevera Needle and Driver onto the Brevera Adapter for GE. Ensure the needle goes through the needle guide. Check to make sure the cable, tubing, and other parts of the Brevera biopsy device are clear of the x-ray tube during rotation.
6. Ensure the proper Brevera biopsy device is selected in the GE Senographe stereotaxy software prior to sending the target to the stereotactic system. The selected Brevera biopsy device should match the biopsy device that is being used in the procedure.
7. Perform the biopsy procedure per the Brevera breast biopsy device instructions for use.
8. Clean the Brevera Adapter for GE and Eviva Bushing for GE per the cleaning instructions found in this document.


Figure Legend

- 1. Adapter
- 2. Needle Guide
- 3. Bushing (hidden from view)
- 4. Needle Guide Holder

Figure 1: How to Install the Brevera Adapter

To Remove the Brevera Adapter:

- 1. Remove the Brevera Needle and Driver from the Brevera Adapter for GE.
- 2. Slide the Brevera Adapter for GE straight up and remove from mount.

Parameter Setup for use of Brevera Biopsy Device on the GE Senographe Stereotaxy System Vertical Approach


Program the GE Senographe stereotaxy system for Brevera biopsy device according to the needle setup instructions provided in the GE Senographe stereotaxy operator manual, in conjunction with the steps below.

When “Creating” needles, select “VAD with vertical holder” for type of needle. Use the dimensions provided in the table below. The diagram below the table shows a schematic of the dimensions.

Note: Do not select “CB” because this will cause the needle to move to an incorrect depth.

| Parameters | BREV09, BREVDISP09 | |
|------------------------------|------------------------|------------------------|
| | 20mm (Standard) | 12mm (Petite) |
| Recommended Name in Software | BREV0920, BREVDISP0920 | BREV0912, BREVDISP0912 |
| Needle Post-Fire A (mm) | 126.5 | 118.5 |
| Notch B (mm) | 19.8 | 13.6 |
| Tip C (mm) | 7.9 | 7.9 |
| Device D (mm) | 264 | 272 |
| Vertical Holder E (mm) | 3.4 | -4.6 |

All values entered above must be verified with the appropriate Brevera biopsy device and applicable tests to ensure the target is near the center of the biopsy device aperture.

Setup of GE Senographe DS Stereotaxy and GE Senographe Essential Stereotaxy Systems Lateral Approach for Use with Brevera Biopsy

Contact Technical Support to confirm compatibility of your biopsy device. Install the biopsy adapter as shown below, then install the biopsy device. (Refer to the user manual provided by the biopsy device manufacturer for installation and operating instructions.)

To Install the Brevera Adapter:

1. Setup GE Senographe stereotaxy system for lateral approach using appropriate bracket(s). Refer to GE Senographe stereotaxy operator manual for setup instructions.
2. Install the Brevera Adapter, BREVADPTR, onto the GE Senographe stereotaxy system. Check the Brevera Adapter components (reference Figure 2) to ensure they are secure. If the components are not secure call technical support.
3. Install the disposable needle guide, EVIVA_NG09R, into the GE Senographe lateral arm Needle Guide Holder.
4. Ensure the proper Brevera biopsy device is selected in the GE Senographe stereotaxy software prior to sending the target to the stereotactic system.
5. Advance the Brevera biopsy device to the Reference Position (edge of Biopsy Device Support nearest to breast) on the lateral adapter bar. See the table for this position.
6. Perform the biopsy procedure per the Brevera breast biopsy device instructions for use.
7. Clean the Brevera Adapter per the cleaning instructions found in this document.


Figure Legend

1. Adapter
2. Mount Screw
3. Needle Guide Holder
4. Needle Guide
5. Biopsy Device Support
6. Lateral Adapter Bar
7. Reference Position

Figure 2: How to Install the Brevera Adapter on the Stage

To Remove the Brevera Adapter:

1. Turn the Mount Screw counterclockwise to release the Brevera Adapter.

Parameter Setup for use of Brevera Biopsy Device on the GE Senographe Stereotaxy System Lateral Approach


Program the GE Senographe stereotaxy system for a lateral needle according to the instructions provided in the GE Senographe stereotaxy operator manual, in conjunction with steps below.

1. When “Creating” needles, select “LAT” for type of needle. The recommended name for the needle is shown in the table below.
2. For Length A reference table below to create the needle entry. The diagram below shows a schematic of the dimension.

NOTE: This length ensures that the needle will be aligned with the needle guide. This length is common for all Eviva, ATEC, and Brevera devices. No information about the length of the needle, notch (aperture) or end (dead space) is entered into the system; therefore, only one needle needs to be programmed for all current Eviva, ATEC, and Brevera devices. This distance does not control the insertion depth of the needle into the breast. The insertion depth will be controlled by advancing the biopsy device to the recommended position in the table below.

| Parameters | BREV09, BREVDISP09 | |
|---|------------------------|------------------------|
| | 20mm (Standard) | 12mm (Petite) |
| Recommended Needle Name | BREV0920, BREVDISP0920 | BREV0912, BREVDISP0912 |
| Type of Needle | LAT | LAT |
| Length A (mm) | 107.8 | 107.8 |
| Proper Device Reference Position on the Lateral Adapter Bar (mm) | 67 | 71 |

All values entered above must be verified with the appropriate Brevera biopsy device and applicable test to ensure the target is near the center of the biopsy device aperture. The Brevera biopsy device must be advanced to the proper device position on the lateral adapter bar. See above table.

Cleaning Instructions for Adapter Components

Remove adapter from biopsy guidance system before cleaning.

Use a lint-free cloth and apply a diluted, mild detergent to clean components and surfaces. To remove more difficult soils, use a soft bristle brush as needed.

Caution: Use the least possible amount of cleaning fluids. The fluids must not flow or run.

Disinfecting the Adapter

To disinfect the system's components and surfaces, Hologic recommends the following disinfecting solutions:

- 10% chlorine bleach and water with one part commercially available chlorine bleach (normally 5.25% chlorine and 94.75% water) and nine parts water. Mix this solution daily for best results
- Commercially available isopropyl alcohol solution (70% isopropyl alcohol by volume, not diluted)
- Commercially available wipes containing less than or equal to 2% quaternary ammonium compounds by weight in water

Apply disinfecting solutions to the system using a clean, lint-free cloth or wipe. Wipe the system with the wet cloth, observing appropriate contact (wet) times. Do not spray cleaning solutions directly on the console or other parts of the Brevera breast biopsy system.

Caution: To prevent damage to the electronic components, do not use disinfectant sprays on the system.

To Prevent Possible Injury or Equipment Damage

Do not use a corrosive solvent, abrasive detergent, or polish. Select a cleaning/disinfecting agent that does not damage the plastics, aluminum, or carbon fiber.

Do not use strong detergents, abrasive cleaners, high alcohol concentration, or methanol at any concentration.

Do not expose equipment parts to steam or high temperature sterilization.

Do not let liquids enter the internal parts of the equipment. Do not directly apply cleaning sprays or liquids to the equipment. Always use a clean, lint-free cloth and apply the spray or liquid to the cloth. If liquid enters the system, disconnect the electrical supply and examine the system before returning it to use.

Caution: Wrong cleaning methods can damage the equipment, decrease imaging performance, or increase the risk of electrical shock.

Always follow instructions from the manufacturer of the product you use for cleaning and disinfection. The instructions include the directions and precautions for the application and contact time, storage, wash requirements, protective clothing, shelf life, and disposal. Follow the instructions and use the product in the most safe and effective method.

Warranty

Except as otherwise expressly stated in the Agreement: i) Equipment manufactured by Hologic is warranted to the original Customer to perform substantially in accordance with published product specifications for one (1) year starting from the date of shipment, or if Installation is required, from the date of Installation ("Warranty Period"); ii) digital imaging mammography x-ray tubes are warranted for twenty-four (24) months, during which the x-ray tubes are fully warranted for the first twelve (12) months and are warranted on a straight-line prorated basis during months 13-24; iii) replacement parts and remanufactured items are warranted for the remainder of the Warranty Period or ninety (90) days from shipment, whichever is longer; iv) consumable Supplies are warranted to conform to published specifications for a period ending on the expiration date shown on their respective packages; v) licensed Software is warranted to operate in accordance with published specifications; vi) Services are warranted to be supplied in a workman-like manner; vii) non-Hologic Manufactured Equipment is warranted through its manufacturer and such manufacturer's warranties shall extend to Hologic's customers, to the extent permitted by the manufacturer of such non-Hologic Manufactured Equipment. Hologic does not warrant that use of Products will be uninterrupted or error-free, or that Products will operate with non-Hologic authorized third-party products. These warranties do not apply to any item that is: (a) repaired, moved, or altered other than by Hologic authorized service personnel; (b) subjected to physical (including thermal or electrical) abuse, stress, or misuse; (c) stored, maintained, or operated in any manner inconsistent with applicable Hologic specifications or instructions, including Customer's refusal to allow Hologic recommended Software upgrades; or (d) designated as supplied subject to a non-Hologic warranty or on a pre-release or "as-is" basis.

How Supplied

The Brevera adapter is supplied as a separately packaged accessory.

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.
250 Campus Drive
Marlborough, MA 01752 USA
Phone: 1-877-371-4372
www.hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

European Representative
Hologic BV
Da Vincilaan 5
1930 Zaventem
Belgium
Tel: +32 2 711 46 80

Symbols Used on Labeling

The following symbols may be found on the product labeling for the Brevera adapter:

| Symbol | Description | Standard |
|--|---|--|
| | Prescription use only | FDA 21 CFR 801.109 |
| | Authorized Representative in the European Community | ISO 15223-1, Reference 5.1.2 |
| | CE Mark European Conformity | MDR Regulation (EU) 2017/745 |
| | Translations in Box | Hologic |
| www.hologic.com/package-inserts | Consult instructions for use | ISO 15223-1, Reference 5.4.3 |
| | Follow instructions for use | IEC 60601-1, Reference No. Table D.2, Safety sign 10 (ISO 7010-M002) |
| | Medical device | ISO 15223-1, Reference 5.7.7 |
| | Catalog number | ISO 15223-1, Reference 5.1.6 |

| Symbol | Description | Standard |
|---|--|--|
|  | Contents | N/A |
| QTY | Quantity | Hologic |
|  | Manufacturer | ISO 15223-1, Reference 5.1.1 |
|  | Country of Manufacture | ISO 15223-1, Reference 5.1.11 |
| Patents | Patents | Hologic |
|  | Keep dry | ISO 15223-1, Reference 5.3.4 |
|  | Not safe for magnetic resonance imaging | ASTM F2503 Reference no. Table 2, Symbol 7.3.3;7.4.9.1; Fig. 9 |
|  | Non-sterile | ISO 15223-1, Reference 5.2.7 |
|  | Caution | ISO 15223-1, Reference 5.4.4 |
|  | Do not use if package is damaged | ISO 15223-1, Reference 5.2.8 |
| LOT | Batch code | ISO 15223-1, Reference 5.1.5 |
| SN | Serial number | ISO 15223-1, Reference 5.1.7 |
| YYYY-MM-DD | Date format: YYYY represents the year MM represents the month DD represents the day | Hologic |
|  | Country code for translation | ISO 3166 |

© 2023 Hologic, Inc.

Hologic, Affirm, Brevera, and Multicare are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners.

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