Artwork consists of:

- 302 8.26 inch x 11.69 inch (A4) sheets attached.

Artwork contains the following file(s):


File Name	Description
AW-16146-4320_004_02.pdf	View File

Artwork is black and white.

REV AUTHORED BY S.BHAMMAR	DATE 6/6/18	 		
REV DRAFTED BY N.FONSECA	DATE 12/16/20			
PROPRIETARY: This document contains proprietary data of Hologic, Inc. No disclosure, reproduction or use of any part thereof may be made except by written permission from Hologic.		TITLE TEXT, IFU, AQUILEX FLUID MANAGEMENT SYSTEM, 4320, NEW SCALE	DOCUMENT NUMBER AW-16146-4320	REV 004
REV. RELEASE DATE:	1/12/2021	SIZE A	SHEET 1 OF 1	

# Aquilex<sup>®</sup>

Fluid Control System

Instructions for Use **en**

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**HOLOGIC<sup>®</sup>**

**en**

Hologic and MyoSure are registered trademarks of Hologic, Inc. and its subsidiaries in the United States and other countries. Aquilex is a trademark of Hologic, Inc. and its subsidiaries in the United States and other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners.

These instructions for use contain information that is subject to copyright. All rights reserved. These instructions for use may not be photocopied, duplicated on microfilm, or otherwise copied or distributed, completely or in part, without the approval of W.O.M. WORLD OF MEDICINE GmbH.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in these instructions for use are associated with registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or issues concerning these instructions for use.

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### Federal Law (only for U.S. market)

**Caution: Federal law restricts this device to sale by or on the order of a physician.**

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Copyright © W.O.M. WORLD OF MEDICINE GmbH



---

Symbols

	Follow instructions for use (white image on a blue background)
	Consult instructions for use
	Type BF applied part
	Equipotentiality
IP 41	Degrees of protection provided by enclosures (IP-Code)
IP 21	Degrees of protection provided by enclosures (IP-Code)
	Alternating current
	Service
<b>REF</b>	Catalogue number
	Do not reuse
	Do not resterilize
<b>STERILE EO</b>	Sterilized using ethylene oxide
<b>LOT</b>	Batch code
<b>SN</b>	Serial number

	Manufacturer
	Date of manufacture (YYYY-MM-DD)
	Use by date (YYYY-MM-DD)
<b>QTY</b>	Quantity
	Not made with phthalates
	Made with phthalates
	Keep dry
	This way up
	Fragile
	Stacking limit by number
	Keep away from sunlight
	Protect from heat and radioactive sources
	Non-ionizing electromagnetic radiation
	Do not use if package is damaged

	<p>Authorized for Sale or use by Physician only</p>
	<p>Temperature limit</p>
	<p>Humidity limitation</p>
	<p>Atmospheric pressure limitation</p>
	<p>Waste management</p>
	<p>ON/OFF (push-push)</p>
	<p>Max. use 30 days</p>
	<p>Non sterile</p>
	<p>Scale connection</p>
	<p>Input/Output</p>
	<p>Caution</p>
	<p>Do not use this power output</p>
	<p>Certification Mark</p>

## Table of Contents

<b>1</b>	<b>Important Operator/User Notes</b>	<b>9</b>
<b>2</b>	<b>Safety Information</b>	<b>10</b>
<b>3</b>	<b>Purpose</b>	<b>11</b>
3.1	Warnings and Precautions	11
3.1.1	Warnings	11
3.1.2	Precautions	16
3.2	Description of the Aquilex® Fluid Control System	18
<b>4</b>	<b>Initial System Set-Up</b>	<b>19</b>
4.1	Preparing the System for Use	19
4.2	System Components	21
<b>5</b>	<b>System Operation</b>	<b>22</b>
5.1	Front of Irrigation Pump Unit	22
5.2	Rear of Irrigation Pump Unit	22
5.3	Fluid Monitoring Unit Set-up	23
5.3.1	Setting of the Container Scale	26
5.3.2	Connecting the Vacuum Tube	27
5.4	Turning On the Aquilex® System	28
5.5	Hanging the Fluid Bags	29
5.6	Using Tube Sets	30
5.7	Tube Overview	30
5.8	Connecting the Outflow Tube Set	31
5.8.1	Connecting Outflow Tube of Tissue Removal Handpiece (e.g. MyoSure®)	32
5.9	Inserting the Inflow Tube Set	33
5.10	Presetting the Intrauterine Pressure	34
5.11	Deficit Limit Setting	35
5.12	Using the Pump during Surgery	35
5.13	Changing Bags during Surgery	36
5.14	Changing Container during Surgery	37
5.15	Changing Instrument during Surgery	37
5.16	Total Inflow Volume Displayed	38
5.17	Turning System Off	38
<b>6</b>	<b>Function Test</b>	<b>39</b>
6.1	General Device and Setup Check	39
6.2	Scale Test	39
<b>7</b>	<b>Safety Functions</b>	<b>40</b>
<b>8</b>	<b>Care and Maintenance</b>	<b>41</b>
8.1	Cleaning the System	41
8.2	Maintenance Intervals	41
8.3	Maintenance by Authorized Service Technician	41
8.4	Replacing of the Fuse	42
<b>9</b>	<b>Annual Inspection</b>	<b>43</b>
9.1	Electrical Safety Test	43
9.2	Basic Function Tests	43
9.2.1	Scale Test	43
9.2.2	Flow Rate Test	44
9.2.3	Pressure Measuring Test	45
9.2.4	Fluid Deficit Measurement Test	46
9.2.5	Testing the Vacuum Pump	47
9.3	Determine the Software Version	47
<b>10</b>	<b>Error and Warning Messages</b>	<b>48</b>
<b>11</b>	<b>Technical Data</b>	<b>49</b>
<b>12</b>	<b>Electromagnetic Compatibility</b>	<b>51</b>
12.1	Electrical Connections	51
12.2	Guidance and Manufacturer's Declaration - Electromagnetic Emissions	51
12.3	Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity	52

<b>13</b>	<b>Accessory List.....</b>	<b>54</b>
<b>14</b>	<b>Warranty Information.....</b>	<b>55</b>
<b>15</b>	<b>Appendix.....</b>	<b>56</b>
15.1	Test log.....	56
	<b>Glossary.....</b>	<b>58</b>
	<b>Index.....</b>	<b>59</b>

## 1 Important Operator/User Notes

Read the Instructions for Use carefully and become familiar with the operation and function of the **Aquilex® Fluid Control System** (Aquilex® System) and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead to:

- life-threatening injuries to the patient,
- severe injuries of the surgical team, nursing staff or service personnel, or
- damage or malfunction of the system and/or accessories.

The Instructions for Use are only for the Aquilex® Fluid Control System, consisting of pump, container scale and bag scale.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of the product.

The words WARNING, CAUTION, and NOTE carry special meanings. Read these sections with special attention.

**Subject to technical changes**

**Please note**

---

### WARNING!

Warnings indicate risks to the safety of the patient or operator. Failure to follow warnings may result in injury to the patient or operator.

---

---

### CAUTION!

Warnings indicate risks to the equipment. Failure to follow cautions may result in damage to the system.

---

---

### NOTE!

These paragraphs contain information to clarify the instructions or provide additional useful information.

---

**Exclusion of liability**

**2 Safety Information**

Hologic is not liable for indirect, incidental and consequential damages, including, but not limited to, loss of profit. Any liability and applicable warranty become null and void if:

- the system and/or the accessories/ peripherals are improperly used, transported, stored, prepared, or maintained;
- the system and/or the accessories are improperly used, prepared, or maintained,
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories/peripherals;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedules are not adhered to.

Receipt of technical documentation from Hologic does not authorize individuals to perform repairs, adjustments, or alterations on or to the system or accessories.

**WARNING! Modification of the Aquilex® Fluid Control System is not permitted.**

**Authorized service technician**

Only an authorized service technician may perform repairs, adjustments, or alterations on the system or accessories/peripherals and use the service menu. Any violation will void any applicable warranty. Authorized service technicians are only trained and certified by the manufacturer.

**Normal Use**

The system may be used only as intended.

**Care and maintenance**

The service and maintenance of the device and its accessories/ peripherals has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use. Maintenance of the device may not be performed during surgery or operation of the device.

**NOTE! Service or maintenance work may not be carried out during surgery.**


**Contamination**

Before shipping, decontaminate the device and accessories/ peripherals in order to protect the service personnel. Follow the instructions listed in these instructions for use. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to refuse the repair of contaminated devices or accessories/peripherals.

**Waste management**

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. For disposal of the device and its accessories, please consult Hologic or an authorized disposal company, in compliance with legal or national regulations.

### 3 Purpose

The Aquilex® Fluid Control System is intended to provide fluid distension of the uterus during diagnostic and operative hysteroscopies and to monitor the volume differential between the irrigation fluid flowing into and out of the uterus.

The system may not be used to introduce fluids into the uterus when hysteroscopy is contraindicated. See the operator's manual of your hysteroscope for absolute and relative contraindications.

Relative contraindications to endometrial ablation:

Hysteroscopic endometrial ablation, whether by laser or electrosurgery, should not be undertaken before adequate training, preceptorship, and clinical experience. Additionally, tissue sampling is required prior to destruction of the endometrium. The following are clinical conditions that can significantly complicate hysteroscopic endometrial ablation:

- Adenomatous endometrial hyperplasia
- Uterine leiomyoma
- Severe adenomyosis
- Pelvic pain (subtle PID)
- Uterine anomalies
- Surgical skill (see above)
- Severe anemia
- Inability to circumnavigate the myoma (re: myoma size) - predominantly intramural myomas with small submucous components.

#### Intended Use

#### Contraindications

### 3.1 Warnings and Precautions

#### 3.1.1 Warnings

##### WARNING!

###### Distension media

When performing monopolar hysteroscopic electrosurgery, the distension medium must be electrically non-conductive. Examples include glycine, sorbitol and mannitol. Isotonic saline irrigation fluids may only be used when performing bipolar electrosurgical resective procedures.

##### WARNING!

###### Pressure

The pressure should be kept as low as possible to allow for a sufficient intrauterine distention and to reduce the forces that could allow fluid, ambient air, and/or gas to enter the circulatory system.

##### WARNING!

###### Intrauterine distention

Intrauterine distention is usually possible with pressure values between 35 to 70 mmHg. A pressure above 75 mmHg is required only in rare cases or if the patient has an excessively high blood pressure.

##### WARNING!

###### Fluid overload

There is a risk of irrigation fluid reaching the circulatory system of the patient's soft tissue by passing through the uterus. This can be affected by distention pressure, flow rate, perforation of the uterine cavity and duration of the hysteroscopic surgery. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Fluid deficit**

The fluid left in the patient must be monitored. The deficit is the total amount of fluid left in the patient or unaccounted for otherwise. Take notice of the measurement tolerance of the system (see Chapter Technical Data [▶ 49]). Estimating the fluid volume remaining in the patient is the physician's responsibility.

**WARNING!****Fluid intake and output surveillance**

Strict fluid intake and output surveillance should be maintained due to the risk of fluid overload. For healthy patients, the maximum fluid deficit of 1,000 ml is suggested when using a hypotonic solution (e.g. glycine, sorbitol and mannitol). If isotonic solutions (e.g. saline, Lactated Ringer's) are used, the fluid deficit should not exceed 2,500 ml.

**WARNING!****Serum sodium concentration**

It is also necessary to monitor the concentration of sodium in the blood of the patient to prevent electrolyte disturbances. Monitoring of the concentration of sodium in the blood must be performed by the physician and is not performed or supported by the system.

**WARNING!****Loss of deficit and inflow value**

The deficit display value is lost in case of a power loss or "brownout."

**WARNING!****Failure of scale connection**

If the message "Check Scale Connection" appears, the deficit must be calculated manually. The pump keeps displaying the last known deficit value determined prior to the failure of the scale connection.

**WARNING!****Container change during surgery**

A container change during surgery is only allowed, if the container holds at least 0,5 liters of fluid. Otherwise, the deficit value may be falsified. In this case, the manufacturer recommends manual deficit calculation.

**WARNING!****Hyponatremia**

Some distension fluids may lead to fluid overload and, consequently, hyponatremia with its attending sequelae. This can be affected by the distending pressure, flow rate, and duration of hysteroscopic procedure. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Pulmonary edema**

Hysteroscopic surgery is associated with a risk of developing pulmonary edema resulting from fluid overload with isotonic fluids. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Cerebral edema**

Hysteroscopic surgery is associated with a risk of developing cerebral edema resulting from fluid overload and electrolyte disturbances with hypoosmolar (non-ionic) fluids such as glycine 1.5 % and sorbitol 3.0 %. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Idiosyncratic reactions**

In rare cases, idiosyncratic reactions including

- intravascular coagulopathy
- allergic reaction including anaphylaxis

may occur while performing a hysteroscopy if a liquid distention medium is used.

**WARNING!****Hypothermia (monitoring body temperature)**

Continuous flow of distention fluids can lead to a lowering of the patient's body temperature during hysteroscopic surgery. Lower body temperatures can cause coronary and cardiovascular problems. Always monitor the patient's body temperature during the entire surgery. Make especially sure that the following, hypothermia promoting, operation conditions are avoided as best as possible:

- longer operating times
- use of cold irrigation fluid.

**WARNING!****Rupture of the fallopian tube secondary to tubal obstruction**

Distention of the uterus may lead to a tear of the fallopian tube should there be an obstruction or permanent occlusion. The rupture could lead to irrigation fluid flowing into the patient's peritoneal cavity, resulting in a fluid overload. It is critical to closely monitor the input and outflow of the distending liquid at all times.

**WARNING!****Air embolism**

An air embolism can be the result of air contained in the tube set or connected instrument reaching the patient. Ensure there is always fluid in the bag to prevent air from being pumped into the patient.

**WARNING!****Irrigation fluid bags**

The system is only intended for use with flexible fluid bags. Do not use glass containers as they might break. With rigid containers, fluid cannot flow quickly enough due to the vacuum being generated inside of the containers. Risk of implosion with rigid containers.

**WARNING!****Resetting the deficit display**

Filling the tubing with irrigation fluid and resetting the deficit display to zero are to be done at the physician's discretion.

**WARNING!****Function test**

The function test must be performed prior to each device use.

**WARNING!****Obvious defects**

Never use the system if it has suspected or confirmed defects, especially if these involve the power plugs or the mains power supply connection cables. In this case have the device repaired by authorized service personnel.

**WARNING!****System error**

Do not use the Aquilex® System if a defect is suspected or detected during the function test. This also applies to any obvious defects, especially defects on the power connector or plug and power cord.

**WARNING!****Disconnect the power cord**

Pressing the ON/OFF switch does not disconnect the system from the wall power outlet. This requires pulling the power cord located in the rear of the system.

**WARNING!****Technique and procedures**

Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.

**WARNING!****Check all factory settings**

Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.

**WARNING!****Original accessories**

For your own safety and that of your patient, use only Aquilex® accessories (see Chapter Accessory List [▶ 54]).

**WARNING!****Additional equipment**

Additional equipment connected to medical electrical devices must be demonstrated to be compliant with their respective IEC or ISO standards (IEC 60601-1, IEC 60950 or IEC 62368 for data processing equipment). Furthermore, all configurations must comply with the normative requirements for medical systems (see section 16 of the last valid edition of IEC 60601-1). Anyone who connects additional devices to medical electrical equipment is a system configurator and as such is responsible for the system's compliance with the normative requirements for systems. Please contact the technical service if you have additional questions.

---

**WARNING!****Not explosion-proof**

The system is not explosion-proof. Do not use in an area where flammable anesthetic gases are present.

---

---

**WARNING!****Risk of electrical shock**

To avoid risk of electrical shock, this system may only be connected to a supply mains with protective earth.

---

---

**WARNING!****Risk of electrical shock**

To prevent electrical shock, do not open this device. Never open this device yourself. Notify the authorized service technicians of any required repairs.

---

---

**WARNING!****Professional qualification**

The instructions for use do not include descriptions or instructions for surgical procedures/techniques. It is not suitable for training physicians in the use of surgical techniques. Medical peripherals and devices may be used only by physicians or medical assistants with the appropriate technical/medical qualifications working under the direction and supervision of a physician.

---

---

**WARNING!****Sterile media and accessories**

Always work exclusively with sterile substances and media, sterile fluids, and sterile accessories if so indicated.

---

---

**WARNING!****Replacement device and accessories**

In case the system or any of the accessories fail during a procedure, an alternative system and replacement accessories should be kept within easy reach to be able to finish the operation with the replacement components.

---

---

**WARNING!****Condensation / Water penetration**

Protect system from moisture. Do not use if moisture has penetrated the system.

---

---

**WARNING!****Replace the fuse**

Replace the fuse only with a fuse of the same type and rating (see Chapter Technical Data [ 49]).

---

**WARNING!****Electromagnetic emissions**

Electromagnetic emissions may increase and rise above the permissible limits if other equipment (e.g. MyoSure® Control Unit) is stacked onto or placed directly next to the Aquilex Fluid Control System. The user is responsible for monitoring the devices to make sure they function properly.

**WARNING!**

If the Aquilex® Fluid Control System is configured as part of a ME SYSTEM, the entire ME SYSTEM should be tested for compliance with IEC 60601-1-1, and any equipment used with the Aquilex® Fluid Control System should be Type BF.

**WARNING!**

If the leakage current of the configured ME SYSTEM exceeds the limits of IEC 60601-1-1, install an appropriately rated UL 2601-1/IEC 60601-1 approved isolation transformer and retest the system.

**WARNING!****Incorrect determination of fluid deficit**

Always use the hooks of the bag scale to hang the fluid bags to ensure an accurate determination of the fluid deficit. In addition, leave the empty fluid bags hanging on the bag scale until the end of surgery.

**WARNING!****Filter**

The vacuum tube with integrated filter is designed for max. 30 days. The vacuum tube may not be sterilized. Replace the vacuum tube sooner if it is obviously contaminated. The filter prevents body fluids from entering the interior of the device.

**3.1.2 Precautions****CAUTION!****Federal Law (only for U.S. market)**

Federal law restricts this device to sale by or on the order of a physician.

**CAUTION!****Indoor climate**

Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

**CAUTION!****Device setup**

Device should be positioned outside of the sterile area in such a way that

- it can be easily disconnected,
- it is easy to use and switch off and on,
- it allows an easy monitoring of the display values, device functions, and access to the control elements.

**CAUTION!****Modifying the system**

Modifying the system and/or its accessories is not permitted.

**CAUTION!****Combination of low set pressures and excessive vacuum pressures**

When using the Aquilex® Fluid Control System with tissue removal systems, e.g. MyoSure®, the combination of low set pressures and excessive vacuum pressures may result in a significant loss of intrauterine distension pressure which has the potential to affect the visibility of the surgical field. Conversely, when employing a high distension pressure, the deactivation of the tissue removal system can lead to pressure spikes that can exceed 150 mmHg.

**CAUTION!**

Do not use the covered power output at the rear of the irrigation pump unit.

**CAUTION!****Hysteroscope**

The system may only be connected with hysteroscopes designed for and featuring the technical specification permitting such a combined use. Any utilized hysteroscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600.

**CAUTION!****Electrical interference**

(see Chapter Electromagnetic Compatibility [▶ 51]): Electrical interference with other devices or instruments was considered when developing this system and none was detected during testing. However, if you still detect or suspect such interference, please follow these suggestions:

- Move the Aquilex® System, the other device, or both devices to a different location.
- Increase distance between devices used.
- Consult an electro-medical expert.

**CAUTION!****Wall outlet voltage**

Check to ensure the available wall outlet voltage matches the data listed on the label attached to the back of the pump. Incorrect voltage can cause errors and malfunctions and may destroy the system.

**CAUTION!****Transport**

The device is transportable. The roller wheels of the Fluid Monitoring Unit (cart/scale) are used for positioning at the place of use. To transport the device, remove all fluid bags from the hooks and make sure there are no containers or only completely emptied containers on the cart/scale. Inflow and outflow tubes must be completely removed. Make sure the power supply line does not touch the ground and there are no other objects located on the Aquilex® Fluid Control System. Always use the handle to move the system safely.

**CAUTION!****Accuracy of the deficit**

To avoid affecting the accuracy of the deficit calculation ensure that the first step of the canister change is to disconnect tubing from the full canisters. Reconnect tubing to the new canisters only if they are already inserted into the scale.

**CAUTION!****Combination of AQL-100PBS with AQL-100CBS**

The irrigation pump unit AQL-100PBS may only be used with the fluid monitoring unit AQL-100CBS as only this combination is approved to bear the MET NRTL mark.

**CAUTION!****Cleaning the system / Sterilization not allowed**

The pump and the cart/scale can be disinfected by wiping off the outer surfaces. Do not sterilize the pump and the cart/scale.

**Technical application scope of the system****3.2 Description of the Aquilex® Fluid Control System**

The intrauterine pressure can be adjusted on the front of the pump. It can be pre-set to a range between 40 and 150 mmHg. The maximum inflow rate is 800 ml/min and is reduced automatically by the pump once the pre-set intrauterine pressure setting has been reached.

The system has been designed to provide both fluid and vacuum systems that maximize the performance of tissue removal systems, e.g. MyoSure®.

**Suggested distension media**

The **Aquilex® Fluid Control System** can be used with hypotonic, electrolyte-free media (e.g., glycine 1.5 % and sorbitol 3.0 %) and isotonic, electrolyte containing media (e.g., saline 0.9 % and Lactated Ringer's).

**Pressure measuring and regulating**

The system operates with a completely non-contact pressure measurement of the irrigation medium. The contact-free pressure measurement is achieved by integrating the pressure chamber into the tubing system. The pressure chamber transmits the irrigation fluid pressure to the electronics of the device via a pressure sensor. The pressure control circuit continuously compares the desired preset intrauterine pressure with the actual intrauterine pressure. The function of this algorithm is to maintain the pre-set intrauterine pressure. Check for possible leaks if the pre-set intrauterine pressure cannot be achieved.

## 4 Initial System Set-Up

Always check all parts and accessories of the system when performing initial setup. If the system has obvious defects, contact Hologic Technical Support (Chapter Warranty Information [▶ 55]).

Place the system on a level surface and install in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in Chapter Technical Data [▶ 49].

### Initial system set-up

#### CAUTION!

##### Device setup

Device should be positioned outside of the sterile area in such a way that

- it can be easily disconnected,
- it is easy to use and switch off and on,
- it allows an easy monitoring of the display values, device functions, and access to the control elements.

#### CAUTION!

##### Indoor climate

Before switching on the device, sufficient time must have passed to adjust to the indoor climate.

### 4.1 Preparing the System for Use

#### CAUTION!

##### Possible malfunctions

The device Aquilex® Fluid Control System should not be used directly next to other devices as this could result in malfunctions. The device Aquilex® Fluid Control System was tested for compliance with IEC 60601-1-2 as a standalone system. Therefore, do not stack other devices (e.g. MyoSure® Control Unit) on the system or the Irrigation Pump Unit. In particular, do not place any other device than the AQL-100PBS on the trays of the AQL-100CBS. If usage in the manner described above is nevertheless required, this system and the other devices should be monitored to make sure they function properly.

#### CAUTION!

##### Combination of AQL-100PBS with AQL-100CBS

The irrigation pump unit AQL-100PBS may only be used with the fluid monitoring unit AQL-100CBS as only this combination is approved to bear the MET NRTL mark.

### Mains connection

#### CAUTION!

##### Mains connection

- Make sure the available mains voltage matches the data listed on the type label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.
- Make sure the connection data and technical specifications of the power supply comply with DIN VDE or national requirements. The mains connection cable may be plugged only into a properly installed, grounded safety wall socket (shockproof socket) (see DIN VDE 0100-710).
- Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

**Grounding contact**

The power connection must be equipped with a grounding contact. Use the original power cable to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.

The grounded, shockproof safety wall socket should be near the device and within easy reach. Disconnect the device from the mains power supply (pull cable out of the grounded safety wall socket) if the device is not being used for several days or an extended period. The device is ready for use as soon as all connections are established, and all cables have been plugged in.

**Only for U.S. operators**

Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 and IEC 60320-C13. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

**Potential equalization**


The equipotential bonding is used as a protective measure against the failure of the protective conductor according to requirements of IEC 60601-1 in the respectively valid version. The installation must be according to the relevant local safety regulations.

**Precautionary measures**

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

This system is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.

## 4.2 System Components

**Fig. 4-1 System Components**

- ① Bag scale
- ② Fluid bag
- ③ Irrigation Pump Unit
- ④ Trays
- ⑤ Container scale cable/connector
- ⑥ Container scale
- ⑦ Container holders
- ⑧ Container
- ⑨ Roller wheel base
- ⑩ Locking foot brake

The Aquilex® Fluid Control System is divided into two separate boxes for shipping:

**Box 1** contains:

- Irrigation Pump Unit
- Instructions for Use
- Power cord
- Aquilex® system vacuum tube set (low and high vacuum)
- 1000 g weight

**Box 2** contains:

- Fluid Monitoring Unit (cart with scale)
- Container rings


## 5 System Operation

Please make sure that the function test according to Chapter Function Test [▶ 39] has been performed prior to each device use.

### 5.1 Front of Irrigation Pump Unit

Fig. 5-1 Front of irrigation pump unit

- ① Pump display
- ② Intrauterine pressure display
- ③ Fluid deficit limit display
- ④ Deficit meter
- ⑤ Deficit display
- ⑥ Inflow tube holder
- ⑦ Roller wheel
- ⑧ Pressure sensor
- ⑨ Reset deficit button (Zero)
- ⑩ Decrease deficit limit
- ⑪ Increase deficit limit
- ⑫ Decrease intrauterine pressure setting
- ⑬ Increase intrauterine pressure setting
- ⑭ Intrauterine pressure setting display
- ⑮ ON/OFF switch
- ⑯ Pause/Resume button
- ⑰ Prime button


Please familiarize yourself with the layout of the individual elements on the front of the irrigation pump unit.

### 5.2 Rear of Irrigation Pump Unit

Fig. 5-2 Rear of irrigation pump unit

- ① Connection for low vacuum (white)
- ② Product Label
- ③ Performance data of the device
- ④ Adjustment controller for high vacuum
- ⑤ Connection for high vacuum (green)
- ⑥ Fuse holder
- ⑦ Covered power output
- ⑧ Power connection for pump
- ⑨ Potential equalization connection
- ⑩ Service interface
- ⑪ Connection for scale
- ⑫ Suction openings

Please familiarize yourself with the layout of the individual elements at the rear of the irrigation pump unit.

#### CAUTION!


**Do not use the covered power output at the rear of the irrigation pump unit.**

**WARNING!****Additional equipment**

Additional equipment connected to medical electrical devices must be demonstrated to be compliant with their respective IEC or ISO standards (IEC 60601-1, IEC 60950 or IEC 62368 for data processing equipment). Furthermore, all configurations must comply with the normative requirements for medical systems (see section 16 of the last valid edition of IEC 60601-1). Anyone who connects additional devices to medical electrical equipment is a system configurator and as such is responsible for the system's compliance with the normative requirements for systems. Please contact the technical service if you have additional questions.

**5.3 Fluid Monitoring Unit Set-up****CAUTION!****Combination of AQL-100PBS with AQL-100CBS**

The irrigation pump unit AQL-100PBS may only be used with the fluid monitoring unit AQL-100CBS as only this combination is approved to bear the MET NRTL mark.

**Fig. 5-3 Fluid monitoring Unit (cart with scale)**

- ① Bag scale
- ② Fluid bags
- ③ Bag deflector
- ④ Container scale cable/connector
- ⑤ Container scale
- ⑥ Holder for upper containers (Serres, Medela)
- ⑦ Holder for lower container (Abbott, Bemis®, Medi-Vac®, DeRoyal®)
- ⑧ Roller wheel base
- ⑨ Container

The fluid monitoring unit (cart/scale) consists of a weighing unit for fluid bags ①, a weighing unit for fluid containers ⑤, and a roller wheel base ⑧.


1. Remove the cart/scale from the cardboard shipping box.
2. Remove the pump and the power cords from the first cardboard box.
3. Loosen the handwheel ⑥ (Fig. Scale and pump connection [▶ 24]) and pull bag scale upwards to stop. The screw ⑦ (Fig. Scale and pump connection [▶ 24]) must be inserted into the provided opening. Secure the bag scale with the handwheel.
4. Depending on the type of container used, attach the container rings (included in the second box) to the upper ⑥ or lower ⑦ container holders (Fig. Fluid monitoring Unit (cart with scale) [▶ 23]).

5. Guide the power cord through the holes provided for this purpose and connect to pump ② (Fig. Scale and pump connection [▶ 24]) and plug into grounded, shockproof safety wall socket.
6. Attach the scale to the pump by connecting and screwing the connector of the bag scale together with the connector of the container scale ④ (Fig. Scale and pump connection [▶ 24]) and fix the connected cables below the lower pump tray by means of the provided wire clips. Connect and screw the second connector of the bag scale to the port on the back of the pump ⑤ (Fig. Scale and pump connection [▶ 24]).

Fig. 5-4 Scale and pump connection

- ① Bag deflector
- ② Power cord/pump connection
- ③ Covered power output
- ④ Plug connection, bag scale with container scale
- ⑤ Bag scale connector, connection to pump
- ⑥ Handwheel
- ⑦ Screw

**WARNING!****Scale error**

Ensure that nothing weighs down the scale during system start-up. Doing so may result in an inaccurate deficit value.

**WARNING!****Fluid deficit**

The fluid left in the patient must be monitored. The deficit is the total amount of fluid left in the patient or unaccounted for otherwise. Take notice of the measurement tolerance of the system (see Chapter Technical Data [▶ 49]). Estimating the fluid volume remaining in the patient is the physician's responsibility.

**WARNING!****Serum sodium concentration**

It is also necessary to monitor the concentration of sodium in the blood of the patient to prevent electrolyte disturbances. Monitoring of the concentration of sodium in the blood must be performed by the physician and is not performed or supported by the system.

**NOTE!****Automatic scale detection**

The pump's software automatically detects whether the container scale and bag scale are connected or only the container scale.

**NOTE!****Deficit accuracy**

The greater the consumption of irrigation fluid, the greater the deviation between the actual and the displayed deficit (see Chapter Technical Data [▶ 49]).

Try to collect all the fluid running out of the uterine cavity during the procedure to achieve the most exact deficit value possible.

**Precise balancing**

The container scale can be loaded with a weight of up to 25 kg (55 lbs). The maximum load of the bag scale is 12 kg (27 lbs). If the scales are loaded beyond these limits, the message **Scale Overload/Check Scale** is triggered and displayed. Three audible warning beeps are emitted as well (see Chapter Error and Warning Messages [▶ 48]).

**Scale capacity**

**WARNING!****Containers and fluid bags**

Make sure the containers and fluid bags hang freely, are not resting on something, and do not touch other objects except the bag deflectors. Failure to follow these instructions means the deficit cannot be calculated correctly.


**NOTE!****Scale recognition**

Connect the scale to the pump before turning the system on to ensure the system recognizes the scale.


### 5.3.1 Setting of the Container Scale


The container scale can be used with containers from different manufacturers.


Bemis® 3 liters


DeRoyal® Crystaline™ 2.1 l


Abbott 2 liters


Serres 2 & 3 liters

Medi-Vac® 3 liters

Medela 3 liters

Medi-Vac® Flex Advantage 3000 cc

**NOTE!****Container position**

Ensure containers are positioned properly in the respective holders.

**NOTE!****Containers with overflow protection**

Only use suction containers with overflow protection.

### 5.3.2 Connecting the Vacuum Tube

**CAUTION!****Combination of low set pressures and excessive vacuum pressures**


When using the Aquilex® Fluid Control System with tissue removal systems, e.g. MyoSure®, the combination of low set pressures and excessive vacuum pressures may result in a significant loss of intrauterine distension pressure which has the potential to affect the visibility of the surgical field. Conversely, when employing a high distension pressure, the deactivation of the tissue removal system can lead to pressure spikes that can exceed 150 mmHg.

Connect vacuum tube with hygiene filter to suction containers. The vacuum tube with hygiene filter must be replaced when dirty and after 30 days at the latest. The vacuum tube with hygiene filter should not be cleaned.

**Connection for low vacuum (white)**

- Connect vacuum tube with white connector to low vacuum port (white) ① Fig. Low vacuum tube [▶ 27]. This vacuum pump has a fixed vacuum pressure (~225 mmHg).
- Use the connecting tube (⑤ Fig. Low vacuum tube [▶ 27]) when two containers are serially connected to the same vacuum port.

**Fig. 5-5 Low vacuum tube**


- ① Connection for low vacuum (white)
- ② Hygiene filter
- ③ Vacuum tube
- ④ Container
- ⑤ Connecting tube

**Connection for high vacuum (green)**

- Connect vacuum tube set with the green connectors to the high vacuum port (green) ⑧ in Fig. High vacuum tube [▶ 28]. This vacuum can be adjusted to a maximum 500 mmHg using adjustment controller.
- Use the connecting tube (⑫ Fig. High vacuum tube [▶ 28]) when two containers are serially connected to the same vacuum port.


Fig. 5-6 High vacuum tube

- ⑥ Hygiene filter
- ⑦ Vacuum tube (green connectors)
- ⑧ Connection for high vacuum (green)
- ⑨ Container
- ⑩ Tissue trap
- ⑪ Adjustment controller
- ⑫ Connecting tube

#### 5.4 Turning On the Aquilex® System

Fig. 5-7 Front of Device

1. Press the ON/OFF switch. The displays and indicators light up and system turns on.
2. The system now performs a device self-test. The device self-test is used by the system to check whether the container and the bag scale are connected, among others. If the bag scale is not connected, the message **System OK. Bag Scale Not Connected** is displayed. In this case, please check the appropriate data connection of the scales with the pump according to Chapter Fluid Monitoring Unit Set-up [▶ 23].
3. If a tube set is in the inflow tube holder when the pump is switched on, the **pump display** (Fig. Front of irrigation pump unit [▶ 22] ①) shows the message **Remove Tube Set**. The device self-test resumes once the tube set is removed from the roller wheel.

If the device self-test is unsuccessful, the corresponding error messages are displayed (see Chapter Error and Warning Messages [▶ 48]).


The system has successfully completed the device self-test when a single audible beep is heard. The message **System OK** is displayed for 5 seconds followed by the message **Insert Tube Set**.

**WARNING!****Obvious defects**

Never use the system if it has suspected or confirmed defects, especially if these involve the power plugs or the mains power supply connection cables. In this case have the device repaired by authorized service personnel.

## 5.5 Hanging the Fluid Bags

**Fig. 5-8 Fluid bag suspension**

- ① Bag scale
- ② Fluid bag on bag hook
- ③ Bag deflector

**WARNING!****Distension media**

When performing monopolar hysteroscopic electrosurgery, the distension medium must be electrically non-conductive. Examples include glycine, sorbitol and mannitol. Isotonic saline irrigation fluids may only be used when performing bipolar electrosurgical resective procedures.

Hang one or two fluid bags filled with distension fluid appropriate for procedure.

**WARNING!****Irrigation fluid bags**

The system is only intended for use with flexible fluid bags. Do not use glass containers as they might break. With rigid containers, fluid cannot flow quickly enough due to the vacuum being generated inside of the containers. Risk of implosion with rigid containers.

### Tube set recognition

## 5.6 Using Tube Sets

The **Aquilex® Fluid Control System** is designed for use with sterile disposable inflow and outflow tube sets.

Each inflow tube set is equipped with tube set recognition technology. An RFID transponder detects the type of tube, whether it has been used, and its reliability automatically. The pump display indicates this information. This eliminates accidental reuse of tube sets on more than one patient (see Chapter Tube Overview [▶ 30]).

### WARNING!

#### Visual inspection of the tube set

Before the operation, perform a visual inspection of the tube set and its packaging. Damaged tube sets or tube sets from damaged packagings may not be used.

### WARNING!

#### Reprocessing of sterile disposable products

Reuse of inflow or outflow tube can cause an infection hazard for patients and/or users as well as impair of product functionality. Contamination and/or impaired functionality of the system can cause risk of injury, illness, or death. Do not re-process or reuse single-use inflow or outflow tube sets.

### NOTE!

#### Disposing of tubes and containers

Comply with national disposal and hygiene rules when disposing of tubes, collected fluid, and the containers.

## 5.7 Tube Overview

Three different tube sets are necessary to operate the system. The following table lists each type of tube set and its application.

Article number	Description
AQL-110	Tube set for irrigation, single-use, for <b>Aquilex® Fluid Control System</b>
AQL-111	Tube set for suction, single-use, for <b>Aquilex® Fluid Control System</b>
AQL-112	<b>Aquilex® Fluid Control System</b> complete tube set (inflow and outflow), disposable, sterilized using ethylene oxide
AQL-114	Tube set for vacuum incl. filter, 30-day use, for <b>Aquilex® Fluid Control System</b>

Table 1: Tube sets

## 5.8 Connecting the Outflow Tube Set

### CAUTION!


#### Combination of low set pressures and excessive vacuum pressures

When using the Aquilex® Fluid Control System with tissue removal systems, e.g. MyoSure®, the combination of low set pressures and excessive vacuum pressures may result in a significant loss of intrauterine distension pressure which has the potential to affect the visibility of the surgical field. Conversely, when employing a high distension pressure, the deactivation of the tissue removal system can lead to pressure spikes that can exceed 150 mmHg.

Fig. 5-9 Outflow tube set

- ① To low vacuum port (white)
- ② Container
- ③ Connecting tube
- ④ Patient port
- ⑤ Outflow tube set
- ⑥ Drape
- ⑦ Removable outflow channel or hysteroscope outflow sheath stopcock



Using low vacuum configuration of Fig. Outflow tube set [▶ 31], connect outflow tube set (Y-tube) to patient port (4) of second container. Yellow flexible connector attaches to drape (6). Yellow Luer fitting connects to stopcock (7) of removable outflow channel or hysteroscope outflow valve.

5.8.1 Connecting Outflow Tube of Tissue Removal Handpiece (e.g. MyoSure®)


Fig. 5-10 Port for tissue removal systems

- ① To high vacuum port (green)
- ② Container
- ③ Connecting tube
- ④ Specimen tissue port
- ⑤ Tissue trap
- ⑥ Vacuum tube of tissue removal handpiece (yellow)
- ⑦ Tissue removal handpiece

If intrauterine pathology is identified, the outflow tube of a tissue removal handpiece (6) is connected to the tissue trap (5) located in the second container.

## 5.9 Inserting the Inflow Tube Set


**Fig. 5-11 Tube set elements**

- ① Protective caps
- ② Bag spikes
- ③ Tube clamps
- ④ Y-connector
- ⑤ Inflow section
- ⑥ Roller wheel section
- ⑦ Pressure chamber with membrane and RFID transponder
- ⑧ Hysteroscope section
- ⑨ Luer lock connector (blue)
- ⑩ Roller wheel connector

(See Fig. Tube set elements [▶ 33]) The inflow tube set consists of three tube sections, a Y-connector (4) and two bag spikes (2). The three tube sections are: Roller wheel section (6), inflow section (5), and hysteroscope section (8). The bag spikes (2) are used to connect the tube sections to the fluid bags.

The Luer lock connector (9) connects the hysteroscope tube with the hysteroscope.


**Fig. 5-12 Inserting the tube set**

- ① Fluid bags
- ② Bag spikes
- ③ Bag clamps
- ④ Y-connector
- ⑤ Inflow tube
- ⑥ Roller wheel tube
- ⑦ Pressure chamber with membrane and RFID transponder
- ⑧ Hysteroscope tube
- ⑨ Luer lock connector (blue)

**Open outer packaging****Connect to hysteroscope****Inserting the tube set****Connect the fluid bags****Fig. 5-13 Attach roller wheel tube**

- ⑤ Inflow tube
- ⑥ Roller wheel tube
- ⑦ Pressure chamber
- ⑧ Hysteroscope tube
- ⑪ Roller wheel
- ⑫ Holder for inflow tube
- ⑬ Alignment arrows

- Inflow tube set - To be carried out by non-sterile technician:
  - Open outer packaging of the inflow tube set.
  - A sterile nurse then removes the inner tube set package and opens it.
- To be carried out by sterile technician:
  - Keep the blue Luer lock connector (9) in the sterile area and hand the tube end with the bag spikes (2) to the non-sterile technician.
  - Connect the blue Luer lock connector (9) with the hysteroscope inflow stopcock. Open stopcock.
- To be carried out by non-sterile technician:
  - Ensure system is turned on.
  - Close the clamps (3) on the inflow tube below the bag spikes (2).
  - Insert the inflow tube set into the inflow tube holder. Insertion of the roller wheel tube is depicted in Fig. Attach roller wheel tube [▶ 34].
  - Carefully insert the pressure chamber (7) into the lower notch of the inflow tube holder (12) until you feel resistance. Align pressure chamber and inflow tube holder using arrows (see Fig. Attach roller wheel tube [▶ 34]).
  - When inserting the roller wheel tube, make sure not to damage the membranes of the pressure chamber. Insert the pressure chamber (7) only if chamber is not pressurized.
  - Place the roller wheel tube (6) around the roller wheel (11).
- When connecting or removing the tube to or from the irrigation fluid bags, always grasp the bag spike at the provided handle. Observe aseptic technique when inserting the spike(s) into the fluid bag(s). The surgeon must select a distension fluid suitable for the type of procedure.

**5.10 Presetting the Intrauterine Pressure****Intrauterine pressure setting**

The intrauterine pressure setting can be adjusted while the system is in operation. Use the ▲ and ▼ buttons (Fig. Front of irrigation pump unit [▶ 22]). The pressure setting can be adjusted to between 40 to 150 mmHg in steps of 5 mmHg.

The intrauterine pressure is shown on the intrauterine pressure display (2).

**Safety threshold**

If when scrolling with the ▲ button (Fig. Front of irrigation pump unit [▶ 22]) the safety threshold of 100 mmHg is reached, an audible warning signal is emitted. Release the ▲ button for one second and scroll again to set higher values up to 150 mmHg.

**CAUTION!****Risk of intravasation**

If the intrauterine pressure does not react to an increase in the pressure setting during the procedure, a perforation of the uterine cavity might be the cause. This results in an increased risk of intravasation. Examine the uterine cavity for injuries.

**5.11 Deficit Limit Setting**

The deficit limit can be adjusted while the system is in operation. Use the ▲ and ▼ buttons (see Fig. Front of irrigation pump unit [▶ 22]). The deficit limit can be adjusted to between 600 and 2500 ml in increments of 100 ml. The deficit limit is shown on the deficit limit display (3). The deficit meter is designed to help the operator track the deficit volume. The color of the deficit meter changes as the deficit limit is approached. The deficit limit set by the operator is marked with a red LED on the top of the deficit meter. If during surgery the actual deficit rises, the LEDs will light up sequentially representing the actual deficit volume until the deficit limit is reached (see section **Deficit Limit** in Chapter Safety Functions [▶ 40]).

**Setting the deficit limit****5.12 Using the Pump during Surgery**

- Open bag clamps on the fluid bags (3 Fig. Inserting the tube set [▶ 33]).
- Fully open hysteroscope inflow stopcock.
- If drainage stopcock is available: Fully close drainage stopcock.
- Keep the hysteroscope at the height of the patient and above the drape to ensure the fluid can be collected. Do not insert the hysteroscope into the uterus at this time.
- Press the button **Prime** (17 Fig. Front of irrigation pump unit [▶ 22]).
- Pump will run for approximately 20 seconds to purge air from tubes and run the automatic lumen calibration function.
- Pump will display **Calibration Running**.

**Setting the deficit limit**

The pump is equipped with an automatic lumen calibration function. The system determines the flow resistance of the hysteroscope. This resistance is used to calculate the pump pressure necessary to maintain the pre-set intrauterine pressure. To overcome this resistance, the pump allows a pressure of up to 80 mmHg during calibration. This will be shown in the display of the actual intrauterine pressure. If calibration fails due to high resistance, calibration is repeated with a permissible pressure of up to 150 mmHg. If calibration is then still not completed, the pump displays **Prime Fail - Open Stopcock Clamps**.

**Automatic lumen calibration**

The automatic lumen calibration starts once the **Prime** button is pressed.

- Three audible warning signals are heard once the automatic lumen calibration is finished. The pump display will show **Prime Successful Close Stopcock** for 5 seconds followed by **System Operating**.
- Close hysteroscope inflow stopcock to stop fluid flow. Once all fluid has been removed from the drape, zero the deficit display.
- Check to see if fluid has leaked in the area of the pressure chamber. If you find leaked irrigation fluid in the area of the pump, change the tube set and retry the automatic lumen calibration.

**NOTE!****Correct lumen calibration and deficit calculation**

The calibration must be performed outside of the patient to ensure a correct lumen calibration and deficit calculation.

**NOTE!****To stop the pump**

The pump continues to operate after automatic lumen calibration is complete. The pump should be stopped by closing hysteroscope inflow stopcock.

**CAUTION!****New lumen calibration**

A new lumen calibration must be carried out with each change of the hysteroscope during surgery (see Section Changing Instrument during Surgery [▶ 37]).

**System operation****Completing system operation**

- Open stopcock and guide the hysteroscope with fluid flowing into the uterus.
- Adjust intrauterine pressure setting as necessary to obtain adequate distension and visualization.
- After the medical procedure, if system operation is complete close the hysteroscope inflow stopcock.
- Wait until the entire fluid volume from the drape and the tube set has been collected in both containers.
- Press the **Pause/Resume** button.
- Note the deficit volume indicated on the deficit display. This is the total fluid volume that was absorbed by the patient.

**WARNING!****System error**

Do not use the Aquilex® System if a defect is suspected or detected during the function test. This also applies to any obvious defects, especially defects on the power connector or plug and power cord.

**WARNING!****Failure of scale connection**

If the message “Check Scale Connection” appears, the deficit must be calculated manually. The pump keeps displaying the last known deficit value determined prior to the failure of the scale connection.

**NOTE!****Change the containers and bags during surgery**

It is possible to change the containers and bags during surgery without losing the previously measured deficit.

**Bag change during surgery****5.13 Changing Bags during Surgery**

The device system detects automatically the replacement of a bag. Brief fluctuations in the deficit calculation (< 10 s) may occur when a bag is replaced. The replacement of a bag is indicated with the message **Bag Change, Please Proceed**.

- Close the clamp of the empty bag.
- Hang a new fluid bag to a hook of the bag scale.
- Reconnect the new fluid bag with the inflow tube set.

**WARNING!****Touching the bags and the hooks**

Touching the bags and the hooks as well as vibrations of the balancing system should be avoided during surgery to prevent false detection of the bag change and not negatively affect the accuracy of the deficit calculation.

**WARNING!****Empty fluid bags**

The empty fluid bags should remain on the hooks so that the measurement accuracy is not diminished.

**WARNING!****Changing fluid bags**

Fluid bags should be changed quickly to avoid affecting the accuracy of the deficit calculation.

**5.14 Changing Container during Surgery**

The device system detects automatically the replacement of a container. The pump will stop immediately, and the deficit display is locked to insure an accurate deficit count is maintained. Brief fluctuations in the deficit calculation (< 10 s) may occur when a container is replaced. The replacement of a container is indicated with the message **Container Change, Press Resume**.

**Container change during surgery**

- Disconnect tubing from full containers.
- Immediately remove full containers from scale.
- Install new containers.
- Reconnect tubing to new containers.
- Press **Pause/Resume** button to resume procedure.

**WARNING!****Container change during surgery**

A container change during surgery is only allowed, if the container holds at least 0,5 liters of fluid. Otherwise, the deficit value may be falsified. In this case, the manufacturer recommends manual deficit calculation.

**WARNING!****Touching containers and their holders**

Touching the containers and their holders as well as vibrations of the balancing system should be avoided during surgery to prevent false detection of the container change and not negatively affect the accuracy of the deficit calculation.

**WARNING!****Change of the containers**

Containers should be changed quickly to avoid affecting the accuracy of the deficit calculation.

**CAUTION!****Accuracy of the deficit**

To avoid affecting the accuracy of the deficit calculation ensure that the first step of the container change is to disconnect tubing from the full containers. Remove full containers from the scale immediately after that.

**5.15 Changing Instrument during Surgery**

- Pause the pump by pressing **Pause/Resume** button.
- Press the **Prime** button for 2 seconds.
- Change the instrument.

**Instrument change during surgery**

- Fully open hysteroscope inflow stopcock.
- Keep the hysteroscope at the height of the patient and above the drape to ensure the fluid can be collected. Do not insert the hysteroscope into the uterus at this time.
- Press the **Prime** (17) button (Fig. Front of irrigation pump unit [▶ 22]).
- Pump is running to carry out the automatic lumen calibration. The pump display depicts **Calibration Running**.
- Three beeps are heard once the automatic lumen calibration is finished.
- The pump display depicts **Prime Successful Close Stopcock** for 5 seconds, followed by **System Operating**.
- Close the stopcock for the hysteroscope inflow to stop the inflow.

### Total inflow volume displayed

#### 5.16 Total Inflow Volume Displayed

If a manual check of the fluid deficit is desired, the total fluid volume supplied from the fluid bags can be obtained by simultaneously pressing and holding both the up and down arrows ((10) and (11) Fig. Front of irrigation pump unit [▶ 22]). The number in the deficit display is the total inflow fluid volume in ml. Once one or both of these buttons is released, the deficit display will return to the fluid deficit value.

### Turning off

#### 5.17 Turning System Off

Press the **ON/OFF** switch to turn pump off. The displays and indicators are no longer illuminated.

#### **WARNING!**

**Disconnect the power cord**

**Pressing the ON/OFF switch does not disconnect the system from the wall power outlet. This requires pulling the power cord located in the rear of the system.**

## 6 Function Test

en

### WARNING!

#### Function test

The function test must be performed prior to each device use.

### 6.1 General Device and Setup Check

1. Perform a visual check of the devices. Do not use the system in case of obvious damage.
2. Check the rollers of the roller wheel to make sure they move easily and smoothly.
3. Turn on the device, check whether power switch and indicators and displays light up accordingly.
4. The device self-test must run successfully; there are no displayed error messages (see Section Turning On the Aquilex® System [▶ 28]).
5. The bags with the irrigation fluid must hang freely and may not touch the scale.
6. Check that all tube connections (vacuum/inflow/outflow) were made correctly and are intact.
7. Check to make sure all tube connections are free of mechanical stresses and are routed without snagging. The tube connections may not touch the scale. Non-observance can lead to a distortion of the deficit calculation.
8. The automatic lumen calibration has been carried out successfully; there are no displayed error messages (see Section Using the Pump during Surgery [▶ 35]).
9. Check for leaking irrigation fluid in the area of the pressure chamber.

### 6.2 Scale Test

You need the following for this test: 1000 g weight (included with each pump).

1. Switch system on.
2. Once the display depicts **Insert Tube Set**, press **Pause/Resume** and **Zero** at the same time.
3. The pump display depicts the message **Scale Test**.  
The first option is the test of the container scale.
4. Place the 1000 g weight on the container scale.
5. The displayed value should be 1000 g. The tolerance range is  $\pm 20$  g.
6. If the displayed value is outside the tolerance range, the scale must be calibrated by a service technician.
7. Remove weight from container scale.
8. Press the ▲ button (see Fig. Front of Irrigation Pump Unit [▶ 22], (11) Increase the deficit limit) to start the test of the bag scale.
9. The test of the bag scale must be carried out on both sides of the bag scale (using left and right bag hook).
10. Place the 1000 g weight on the right bag hook of the bag scale.
11. The display of the fluid deficit limit depicts the weight. The displayed value should be 1000 g. The tolerance range is  $\pm 20$  g.
12. If the displayed value is outside the tolerance range, the scale must be calibrated by a service technician.
13. Repeat the steps 10. to 12. with using the left bag hook of the bag scale.
14. Remove weight from bag scale.

(Switch to the container scale test by pressing the ▼ button (see Fig. Front of Irrigation Pump Unit [▶ 22], (10) Decrease deficit limit).)

15. Press the **Pause/Resume** button to end the test.

Enter the results into the test log in section Test log [▶ 56]. The test is successfully completed when the results are within the permissible tolerance limit.

## 7 Safety Functions

The electronic components continuously monitor the proper function of the system. System malfunctions are indicated with audible warning signals, error messages, and/or the blocking of system functions. A table listing a summary of possible error and warning messages is provided in Error and Warning Messages [▶ 48].

<b>Intrauterine pressure 10 mmHg above preset intrauterine pressure setting</b>	If the intrauterine pressure exceeds the intrauterine pressure setting by 10 mmHg for longer than 5 seconds, the pressure reduction function is activated. The roller wheel will move forward or backward a few times during the pressure reduction process. If the pressure cannot be reduced, the message <b>Overpressure/ Open Stopcock</b> is displayed, and three audible warning signals are emitted.
<b>Intrauterine pressure &gt; 150 mmHg</b>	The message <b>Maximum Pressure</b> is displayed, and the pump will emit 3 audible warning signals once the intrauterine pressure exceeds 150 mmHg. The maximum permissible pressure has now been reached.
<b>Intrauterine pressure &gt; 200 mmHg</b>	If the intrauterine pressure exceeds 200 mmHg for longer than 5 seconds, the roller wheel stops, and the message <b>Overpressure/Check Stopcock</b> is displayed. Three short continuous audible warning signals are emitted until the pressure is reduced. Once the intrauterine pressure falls below 200 mmHg, the audible warning signals stop and the pump wheel resumes turning automatically.
<b>Check Tube Set Installation</b>	If the inflow tube set is not inserted properly into the roller wheel, pressing the <b>Prime</b> button results in a short audible warning signals and <b>Check Tube Set Installation</b> is displayed. The roller wheel does not start to turn.
<b>Pressure measuring system malfunctions</b>	If a malfunction is detected in the pressure measurement electronics, <b>Sensor Error</b> is displayed, and five short audible warning signals are emitted. The roller wheel stops turning.
<b>Scale overload</b>	If the maximum permissible weight of the scale is exceeded (bag scale or container scale, a continuous audible warning signals is emitted, and <b>Scale Overloaded Check Scale</b> is displayed. The warning stops once the excess weight is removed from the scale.
<b>Loading/unloading scale while in operation</b>	The following messages are displayed when a container is removed from the scale or a new fluid bag is added during operation of the system: <b>Container Change, Press Resume or Bag Change, Please Proceed.</b>
<b>Pressure setting at restart</b>	If the last pre-set intrauterine pressure setting is greater than 80 mmHg, this value is reset to the default value of 80 mmHg.
<b>Deficit limit</b>	Each additional deficit increase by 100 ml above the selected deficit threshold triggers 3 audible warning signals that are repeated while the pump continues to operate. The message <b>Deficit Limit Exceeded</b> is displayed.
<b>Deficit rate &gt;300 ml/min</b>	If the deficit rate exceeds 300 ml/min, 3 audible warning signals are emitted, and the message <b>High Fluid Loss Check Leakage</b> is displayed. If no obvious source of high fluid loss can be identified, an assessment of potential cervical or uterine perforation should be made.
<b>Serious system defect</b>	Five short audible warning signals are emitted, and the message <b>Motor Error</b> is displayed. It is possible that errors can occur in the startup sequence prior to the enabling of the pump display. In this situation, the pump display will remain blank.

## 8 Care and Maintenance

### NOTE!

**Service and maintenance work may not be carried out during surgery.**

The service and maintenance of the system and its accessories has to be carried out as per instructions to ensure the safe operation of the system. For the protection of the patient and the operating team, check that the system is complete and functional before each use.

Special care is necessary when servicing, maintaining, and storing the system and its accessories to maintain the functionality of the equipment and any attached devices.

### 8.1 Cleaning the System

1. Use the **ON/OFF** switch to turn off the system.
2. Remove the power cord.
3. Wipe the surface of the system with a soft cloth moistened with a disinfectant (for example, Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the system.

### CAUTION!

**Cleaning the system / Sterilization not allowed**

**The pump and the cart/scale can be disinfected by wiping off the outer surfaces. Do not sterilize the pump and the cart/scale.**

### 8.2 Maintenance Intervals

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. This inspection has to be carried out once a year. The tests are described in Chapter Annual Inspection [▶ 43].

Regular inspections will assist in early detection of possible malfunctions. This helps to preserve the device and increases its safety and service life.

### 8.3 Maintenance by Authorized Service Technician

An authorized service technician has to inspect and service the device at appropriate intervals to ensure its safety and functionality. The minimum service interval is two years, depending on frequency and duration of use. If the service interval is not maintained, the manufacturer does not assume any liability for the functional safety of the device.

A sticker located on the rear panel of the device states the latest date for the next service or maintenance check.

Authorized service technicians are only trained and certified by the manufacturer.

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians.

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

Care and maintenance

Manufacturer's specifications

Two-year maintenance interval

Authorized trained personnel

Unauthorized personnel

Liability

Technical documents

## Certification

Ask the service technician for a certificate after he or she has inspected the unit or performed any service tasks. This certificate lists the type and scope of the service as well as the date and name of the servicing company together with the signature of the service technician.

**WARNING!****Risk of electrical shock**

To prevent electrical shock, do not open this device. Never open this device yourself. Notify the authorized service technicians of any required repairs.

**WARNING!****Modification of the device**

This device may not be modified without the permission of the manufacturer.

**WARNING!****Modified device**

If the device is modified, suitable examinations and tests must be carried out to ensure further safe use of the device.

**CAUTION!****Replacing the fuse**

Before replacing the fuse, check the values of the fuse to be inserted according to Technical Data [▶ 49].

The fuse may be defective and needs to be replaced if:

- one or more of the pump displays does not light up,
- the system does not function.

Check whether

- the power cord is properly connected to the power cord connection (Figure Opening the fuse holder [▶ 42]) and to a grounded safety wall outlet,
- the wall outlet has power.

**WARNING!****Unplug the power cord from the system before checking the fuse.**

The system does not have to be opened to replace the fuse.

1. Turn system off.
2. Disconnect system from wall power outlet.
3. The fuse holder is located on the back of the pump, next to the male connection.
4. Remove both fuse holders as depicted in Fig. Opening the fuse holder [▶ 42], using small flathead screwdriver.
5. Pull out the fuse holders.
6. Check the fuses.
7. Insert new fuses. Use only the specified type of fuse (see Technical Data [▶ 49]).
8. Insert the fuse holders.


Fig. 8–1 Opening the fuse holder

## 9 Annual Inspection

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. These inspections have to be carried out on an annual basis. Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

The following tests are designed specifically for trained personnel or a hospital technician. The operation of the device as well as its functionality and serviceability are easily checked. Each test conducted has to be documented with date and signature in the test log.

### Manufacturer's specification

### Inspection tests

#### WARNING!

##### Parameters and tolerances exceeded

**If the specified parameters and tolerances are exceeded, the system must be returned to Hologic for evaluation.**

### 9.1 Electrical Safety Test

1. Perform a visual inspection. Make sure that
  - the fuse corresponds with the specifications indicated by the manufacturer,
  - labels and stickers on device are legible,
  - the mechanical condition of the device allows for its safe use,
  - the device is clean to ensure proper and safe functionality.
2. Carry out the measurements for the ground leakage current, short-circuit current/housing leakage current, and the protective conductor resistance as per IEC 62353 in the current version or according to the applicable national standard.

Measuring the protective conductor resistance should be carried out only according to IEC 60601-1 / EN 60601-1 using a current of 25 A AC.

### 9.2 Basic Function Tests

The basic function tests check displays, buttons, and overall system performance.

You need the following for this test:

- Aquilex® inflow tube set
- Fluid bag containing at least 1.5 l of fluid, to prevent the bag from running empty during the test procedure
- Measuring cup with marked scale (1 Liter)
- Stopwatch
- Precision weights (e.g., Ohaus 1 kg 49016-11 or 41000-00 or equivalent)
- 1000 g weight (included in the scope of delivery of the fluid monitoring unit)

#### NOTE!

##### Function test fails

**If the device does not work as described and the test fails the device must be sent to the service.**


#### 9.2.1 Scale Test

The scale test must be carried out as part of the annual inspection. Enter the results into the test log in section Test log [▶ 56]. The test is successfully completed when the results are within the permissible tolerance limit.

See Chapter Scale Test [▶ 39] for details.

Fig. 9–1 Flow rate test

## 9.2.2 Flow Rate Test

## Testing the flow rate

The test setup is depicted in Fig. Flow rate test [▶ 44].

1. Switch system on. (See Chapter Turning On the Aquilex® System [▶ 28])
2. Insert tube set into pump and close clamps on the bags.
3. Hang the fluid bags onto the fluid bag hooks.
4. Insert tap spikes into the fluid bags and open clamps on the fluid bags.
5. Insert hysteroscope tube into the measuring cup.
6. Set intrauterine nominal pressure to 150 mmHg.
7. Press the **Prime** button
8. The roller wheel starts to rotate in order to remove the air from the tubes and to carry out the automatic lumen calibration.
9. When the automatic lumen calibration is completed, press **Pause/Resume** button (about 20 seconds).
10. Empty the measuring cup.
11. Place hysteroscope tube back into the measuring cup.
12. Press the **Pause/Resume** button.
13. Press **Pause/Resume** button after one minute. The measuring cup should contain approx. 800 ml +/- 60 ml of fluid.


Enter the results into the Test log [▶ 56] in Section Appendix [▶ 56]. The test is successfully completed when the results are within the permissible tolerance limit.

### 9.2.3 Pressure Measuring Test

The test setup is depicted in Fig. Set up of pressure measuring test.

Fig. 9–2 Set-up of pressure measuring test

h Height of the water line

The pressure test checks the pressure chamber, the pressure sensor, and the accuracy of the pressure measurement, to ensure that all elements function correctly. For this test, an inflow tube set and a water-filled container is required. The height of the water column (hydrostatic pressure) is used to test the pressure transducer.

1. Hang the inflow end of the tube with the tap spikes for the fluid bags into a container filled with water.
2. Fill the tube end completely with water by running the pump with the Prime button. Let the pump run until the calibration is complete. Press the **Pause/ Resume** button to stop the roller wheel. The intrauterine actual pressure display depicts 0 mmHg.
3. Seal the end of the hysteroscope tube (with finger on the tip of the luer connector).
4. Hold the water level at the end of the hysteroscope tube (h) 30 cm above the pressure chamber. The water column exerts a hydrostatic pressure on the pressure transducer.
5. Remove your finger from the end of the hysteroscope tube.
6. The intrauterine actual pressure display should depict 20 mmHg ( $\pm 5$  mmHg).
7. Change the water column height by changing the height of the water filled end of the tube set. The value depicted on the intrauterine actual pressure display should change accordingly.


Enter the results into the test log in section Test log [▶ 56]. The test is successfully completed when the results are within the permissible tolerance limit.

### 9.2.4 Fluid Deficit Measurement Test

The test setup is depicted in Fig. Test setup for fluid deficit measurement [▶ 46]. It is critical that the collection container is positioned **on** the scale as shown in Fig. Test setup for fluid deficit measurement [▶ 46].

1. Carry out the scale test (see Chapter Scale Test [▶ 39]).
  2. Hang the fluid bag (holding at least 1.5 l of fluid) onto one of the hooks of the bag scale.
  3. Position one collection container properly in the respective holders of the container scale. Use one of the places at the back of the scale.
  4. Turn the system on.
  5. Once the display depicts **Insert Tube Set**, insert a tube set into the pump.
  6. Connect the tube set with the fluid bag and open the clamps.
  7. Fix the other end of the tube set in the recesses at the back of the pump tray (see Fig. Test setup for fluid deficit measurement [▶ 46]) so that the tube end with the luer lock hangs freely in the container.
  8. Place the middle tube section on the pump so that it does neither touch container scale nor the bag scale (see Fig. Test setup for fluid deficit measurement [▶ 46]).
  9. Press the **Prime** button.
  10. Once the display depicts **Prime Successful Close Stopcock**, press the **Pause/Resume** button to stop fluid flow.
  11. Zero the deficit display.
  12. Press the **Pause/Resume** button to start fluid flow.
  13. Allow the system to run for 1 minute (use a stop watch).
  14. Press the **Pause/Resume** button to stop fluid flow.
  15. The displayed fluid deficit should be 0 ml. The permissible tolerance is  $\pm 50$  ml.
- The test is successfully completed when the results are within the permissible tolerance limit.

Fig. 9–3 Test setup for fluid deficit measurement

### 9.2.5 Testing the Vacuum Pump


This test is not designed as a performance test for the measurement of negative pressure. The test only indicates whether or not the vacuum pump is operational.

1. Check to be sure at least one port is opened into the container.
2. Press the **Pause/Resume** button.
3. Place a finger lightly over the suction openings (Fig. Vacuum pump exhaust ports [▶ 47], ①) and check whether you can feel a sucking air flow.

Enter the results into the test log in section Test log [▶ 56]. The test is successfully completed when a stream of air can be demonstrated.

**Fig. 9–4 Vacuum pump exhaust ports**

① Suction openings

### 9.3 Determine the Software Version

The software version of the pump may be required for further analysis.

Follow these steps to determine the software version:

1. Turn device on and wait for the device self-test.
2. Press the Pause/Resume button for at least 2 seconds.
3. The software version of the pump is depicted in the "Deficit Limit" display and consists of a 5-digit numerical sequence.

Enter the determined software version into the test log in Section Test log [▶ 56].

## 10 Error and Warning Messages

The messages appear in the pump display and audible warning signals are heard. Audible warning signals indicating a warning or operating message are sounded several times.

Message in the pump display	audible warning signals	Method
Check Tube Set Installation	1 audible warning signal	Remove and re-insert tube set. If the message reappears, insert a new tube set.
Tube Set Over Usage Limits	1 audible warning signal	Tube set detection indicates that the tube set has been used already. Insert new tube set.
Check Flow Path, Stopcock, Clamps	3 audible warning signals	Flow path is blocked. Check that bag clamps and hysteroscope stopcock are open. Make sure the tube set is not blocked.
Incorrect Tube Set	1 audible warning signal	Replace tube set. The tube set does not match the type approved for the <b>Aquilex® Fluid Control System</b> .
Pump Paused, Press Resume	1 audible warning signal	Pause/Resume button has been activated. Press <b>Pause/Resume</b> button again to resume surgery.
Overpressure Open Stopcock	3 audible warning signals	Most commonly triggered when hysteroscope stopcock is closed while pump is operating at peak flow rate. Open hysteroscope stopcock or remove other closure to relieve pressure.
Overpressure Check Stopcock	5 audible warning signals, constantly repeating until pressure is reduced	Pressure has exceeded the 200 mmHg safety limit and must be reduced. The most common cause is a closed hysteroscope stopcock while the pump operates at the highest flow rate. Open inflow stopcock on hysteroscope or remove blockage clamping off the inflow tube set.
Maximum Pressure	No audible signal	The message "Maximum Pressure" is depicted when the intrauterine pressure is greater than 150 mmHg.
Deficit Limit Reached	3 audible warning signals	Actions are at the discretion of the physician.
Deficit Limit Exceeded	3 audible warning signals	Actions are at the discretion of the physician. If necessary, perform manual deficit control.
Pressure Threshold	1 audible warning signal	For a hysteroscopy, pressures above 100 mmHg are usually not necessary. Careful monitoring of the fluid deficit is recommended.
Check Scale Connection	3 audible warning signals	Check scale connection. Reconnect scale, restart device. If the message reappears, contact Hologic.
Remove Tube Set for System Check	1 audible warning signal	Make sure all tube sets are removed from the roller wheel during the system test. Remove tube set and wait until you hear an audible warning signal and the message "Insert Tube Set" appears.
Prime Fail-Open Stopcock, Clamps	3 audible warning signals	Check bag clamp(s) and hysteroscope inflow stopcock are open. Press the "Prime" button to restart.
High Fluid Loss Check Leakage	3 audible warning signals	Actions are at the discretion of the physician. If necessary, perform manual deficit control.
Scale Overloaded/Check Scale	3 audible warning signals	Weight on the <b>container scale</b> exceeds 25 kg (55 lbs) or weight on the <b>bag scale</b> exceeds 12 kg (26.5 lbs). The weight must be reduced. System function will resume when the excess weight has been removed.
Communication Error	5 audible warning signals	Contact Hologic Technical Support.
Calibration Error	5 audible warning signals	Contact Hologic Technical Support.
Sensor Error	5 audible warning signals	Contact Hologic Technical Support.
Motor Error	5 audible warning signals	Contact Hologic Technical Support.
Low Vac Failed Use Alternative	3 audible warning signals	An alternative low vacuum source must be used in order to continue the procedure. Contact Hologic Technical Support.
High Vac Failed Use Alternative	3 audible warning signals	An alternative high vacuum source must be used in order to continue the procedure. Contact Hologic Technical Support.
Vac Systems Out Use Alternative	3 audible warning signals	An alternative low vacuum source must be used in order to continue the procedure. Contact Hologic Technical Support.

## 11 Technical Data

Type designation	Aquilex Fluid Control System (REF: AQL-100S) consisting of: Irrigation Pump Unit (REF: AQL-100PBS) Fluid Monitoring Unit (REF: AQL-100CBS)	
Manufacturer Information	W.O.M. WORLD OF MEDICINE GmbH Salzufer 8, 10587 Berlin	
Software version	See user manual/operator instructions to determine software version (9.3, Determine the Software Version)	
Mains voltage range [V]	100 - 240 V~	
Supply frequency range [Hz]	50/60 Hz	
Fuse designation	2 x T 3,15 AH, 250 V, UL-recognized	
Internal voltage supply	No	
Power consumption:	Current [A]	Power consumption [VA]
Normal operation 100 V/60 Hz	0.5 A	50 VA
Normal operation 240 V/50 Hz	0.3 A	72 VA
Max. allowed load of the additional socket / multiple sockets [A or VA]	1.6 A	
Protection class (I, II, III)	I	
Application part type (B, BF, CF):	BF (with Aquilex Fluid Control System inflow tube set) <sup>1</sup>	
Defibrillator protection (yes/no)	No	
Protection class (IP code)	IP41 (AQL-100PBS), IP21 (AQL-100CBS)	
Classification (I, IIa, IIb, III) in accordance with the Medical Device Directive 93/42/EEC, Appendix IX:	IIb	
Conformity with the following standards: (in the currently valid version)	IEC 60601-1 / EN 60601-1 / CAN/CSA-C22.2 NO. 60601-1 / AAMI ANSI ES 60601-1 IEC 60601-1-2 / EN 60601-1-2 / AAMI ANSI IEC 60601-1-2	
Operating conditions [°C] [°F], [%], [kPa]	10 to 40 °C/50 to 104 °F 30 to 70 % rel. humidity 70 to 106 kPa air pressure 3000 m max. altitude above sea level for device use	
Possible use with explosive anesthetic gases	This device is not designed for use with flammable anesthetic gases (Class AP) or flammable anesthetic gases with oxygen (Class APG).	
Storage conditions [°C] [°F], [%], [kPa]	5 to 40 °C/41 to 140 °F 5 to 85 % rel. humidity 70 to 106 kPa air pressure	
Transport conditions [°C] [°F], [%], [kPa]	-20 to 70 °C/-4 to 158 °F 5 to 90 % rel. humidity at 30 °C/86 °F 70 to 106 kPa air pressure	
Max. sound level [dB]	< 80 dB(A) (with acoustic signals)	
Max. load of container scale	25 kg/55,1 lbs, 4 containers	
Max. load of bag scale	12 kg/ 26.5 lbs, 2 bags 5 l each	
Maximum flow rate [l/min]	0.8 l/min ± 10 %	
Maximum negative suction pressure	-67 kPa	
Adjustable values		
	Pressure range [mmHg]	40 – 150 mmHg, 5.3 – 20 kPa
	Deficit limit [l]	0.6 - 2.5 l
	Suction pressure	0, LOW: -30 kPa (fixed value), HIGH: -40 to -67 kPa
Measurement range		
	Pressure [mmHg]	0 – 500 mmHg

	Deficit [ml]	-995/+9995 ml
Accuracy		
	Pressure [mmHg]	± 10 mmHg
	Deficit [ml]	< 1 l: ± 60 ml, > 1 l: ± 6 % deficit accuracy related to provided inflow volume, but maximum 300 ml
	Suction pressure	± 20 %
	Total inflow volume [ml]	± 10 %
Dimensions	Width x Height x Depth [mm], [in]	300 mm x 140 mm x 300 mm / 11.8 in x 5.5 in x 11.8 in (AQL-100PBS), 670 mm x 1390 mm x 670 mm / 26.4 in x 54.8 in x 26.4 in (AQL-100CBS)
Weight [kg], [lbs]		5.8 kg, 12.8 lbs (AQL-100PBS) 14.2 kg, 31.3 lbs (AQL-100CBS)
Interfaces:		
	ON/OFF signal for components	Irrigation Pump Unit: 1 x scale connection (flanged socket/round connection socket with 5 pins) 1 x service port (RS232 socket DSUB9/D-SUB9) Fluid Monitoring Unit/Container Scale: 1 x data connection (flanged female / round socket with 5 pins) for connection to bag scale. Fluid Monitoring Unit/Bag Scale: 2 x data ports (flanged socket/round connection socket with 5 pins): Connector to connect to pump, socket to connect with container scale
	Mains power socket	IEC 60320-1 C14
	RFID Transponder Technology	Working frequency: 13.5609 MHz, Transmitting power: -7.51 dBμA/m at 10 m
Essential Performance		Pressure build-up in the body cavity, control and measurement, limit: max. 150 mmHg (Normal Condition), 200 mmHg for max. 5 seconds (Single Fault Condition). Suction Pressure: up to -67 kPa (Normal Condition), 0 kPa by no function or -84 kPa by limit due to design constraints (Single Fault Condition) Deficit measurement: < 1 l: ± 60 ml, > 1 l: ± 6 % deficit accuracy related to provided inflow volume, but maximum 300 ml Single Fault Condition 1: no function of deficit measurement Single Fault Condition 2: < 1 l: greater deviation than ± 60 ml, > 1 l: greater deviation than ± 6 % deficit accuracy related to provided inflow volume or more than 300 ml

<sup>1</sup> A tube set is not an application part in terms of the standard. However, it meets all the technical requirements for an application part.

## 12 Electromagnetic Compatibility

### CAUTION!

#### Accessories

To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device **Aquilex® Fluid Control System** must be used only with the accessories listed in Chapter Accessory List [▶ 54].

#### Precautionary measures

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

The device is to be used only for the purposes described in the Instructions for Use and is intended for use in environments in Professional Healthcare Facility Environment. This applies even if individual requirements meet the conditions for deviating electromagnetic environments. During installation and commissioning as well as during operation of the device, the compliance with the notes and instructions for EMC must be strictly observed.

To ensure the basic safety and essential functionality in relation to electromagnetic interference over the life of the device, the device must be restarted after 24 hours so that a diagnostic self-test can be performed. The maintenance intervals indicated in Chapter Maintenance Intervals [▶ 41] must also be observed.

This device complies with the electromagnetic compatibility (EMC) requirements for medical electrical devices as defined by IEC 60601-1-2. The limits used in testing provide a basic level of safety against typical electromagnetic interference likely to occur in professional health care facilities. Nevertheless, it can happen that individual performance features are no longer available or only to a limited extent due to the presence of EM interference.

### 12.1 Electrical Connections

The following are ESD precautionary measures:

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

#### ESD (Electrostatic Discharge) precautionary measures

Hospital employees should be informed about and trained in ESD precautionary measures.

### 12.2 Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The device **Aquilex® Fluid Control System** is intended for use in an environment as described below. The user/operator of the device **Aquilex® Fluid Control System** should make sure the device is operated within such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device <b>Aquilex® Fluid Control System</b> uses RF energy solely for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device <b>Aquilex® Fluid Control System</b> is suitable for use in all establishments, other than public establishments and those directly connected to the public mains network that supplies buildings used for public purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

### 12.3 Guidelines and Manufacturer's Statement - Electromagnetic Interference Immunity

The device **Aquilex® Fluid Control System** is intended for use in an electromagnetic environment as described below. The user/operator of the device **Aquilex® Fluid Control System** should make sure the device is operated within such an environment.

Electromagnetic interference immunity tests	IEC 60601 test level	Compliance levels	Electromagnetic environment guidelines
Electrostatic discharges (ESD) IEC 61000-4-2	± 8 kV contact ± 2, ± 4, ± 8, ± 15 kV air	± 8 kV contact ± 2, ± 4, ± 8, ± 15 kV air	Floors should be made from wood or concrete or covered with ceramic tiles. If the floor covering consists of synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/bursts according to IEC 61000-4-4	± 2 kV for power lines Repetition frequency 100 kHz	± 2 kV for power lines Repetition frequency 100 kHz	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Surges according to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	The quality of the supply voltage should be the same as the voltage of a typical business or hospital environment.
Blackouts, brownouts, and fluctuations of the power supply according to IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	
	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	
	0 % UT; 250/300 cycles	0 % UT; 250/300 cycles	
Supply frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields of the mains power frequency should comply with the typical values of business and hospital environments.
Conducted RF IEC 61000-4-6	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0.15 MHz to 80 MHz 6 Vrms in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Transients RF IEC 61000-4-3	3 V/m 0.15 MHz to 80 MHz 80 % AM by 1 kHz	3 V/m 0.15 MHz to 80 MHz 80 % AM by 1 kHz	

**IMMUNITY to proximity fields from RF wireless communications equipment****Test method IEC 61000-4-3**

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM $\pm 5$ kHz deviation 1 kHz sinus	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100-5800	WLAN WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

**WARNING!****Portable HF communication equipment**

Portable HF communication equipment can affect the performance characteristics of the device Aquilex® Fluid Control System. Such equipment must therefore comply with a minimum distance of 30 cm (regardless of all calculations) from the device Aquilex® Fluid Control System, its components and cables.

### 13 Accessory List

The following accessories and peripherals are available:

Article	Order No.
Tube set for irrigation, single-use, for <b>Aquilex® Fluid Control System</b>	AQL-110
Tube set for suction, single-use, for <b>Aquilex® Fluid Control System</b>	AQL-111
<b>Aquilex® Fluid Control System</b> complete tube set (inflow and outflow), disposable, sterilized using ethylene oxide	AQL-112
Tube set for vacuum incl. filter, 30-day use, for <b>Aquilex® Fluid Control System</b>	AQL-114
<b>Aquilex® Fluid Control System</b> container rings	AQL-200
<b>Aquilex® Fluid Control System</b> power cord (US)	AQL-215
<b>Aquilex® Fluid Control System</b> power cord (UK)	AQL-216
<b>Aquilex® Fluid Control System</b> power cord (EU)	AQL-217
<b>Aquilex® Fluid Control System</b> spare calibration weight	AQL-218

## 14 Warranty Information

Hologic warrants to the original purchaser of the **Aquilex® Fluid Control System** that it shall be free of defects in material and workmanship when used as intended under normal surgical conditions and in conformance with its instructions for use and maintenance instructions. The obligation of Hologic under this warranty shall be limited to the repair or replacement, each at no charge, at the option of Hologic within one year from the date of purchase. Alternatively, Hologic may elect to repay or credit the original purchaser an amount equal to the purchase price of the defective equipment.

THIS WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES EXPRESSED OR IMPLIED INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR USE AND ALL OTHER OBLIGATIONS AND LIABILITIES ON THE PART OF HOLOGIC. HOLOGIC'S ENTIRE WARRANTY RESPONSIBILITY IS EXPRESSLY LIMITED TO REPAIR OR REPLACEMENT (AT HOLOGIC'S OPTION AND IN THE FORM ORIGINALLY SHIPPED) OF PRODUCT OR CORRECTION OF SERVICE SUBJECT TO ANY CLAIM, OR, AT HOLOGIC'S ELECTION, REPAYMENT OF, OR CREDITING CUSTOMER WITH, AN AMOUNT EQUAL TO THE HOLOGIC PRICE, FEE OR CHARGE THEREFOR. SUCH LIMITED WARRANTY IS GIVEN SOLELY TO THE ORIGINAL PURCHASER AND IS NOT GIVEN TO, NOR MAY IT BE RELIED UPON BY, ANY THIRD PARTY, INCLUDING, WITHOUT LIMITATION, CUSTOMERS OF PURCHASER. THIS WARRANTY IS VOID UPON TRANSFER OF PRODUCT BY PURCHASER TO ANY ENTITY WHO HAS LESS THAN FIFTY (50) PERCENT OWNERSHIP IN THE PRODUCT. THIS WARRANTY SHALL NOT APPLY TO AN AQUILEX® SYSTEM OR TO THE AQUILEX® FLUID CONTROL SYSTEM WHICH HAS BEEN SUBJECT TO ACCIDENT, NEGLIGENCE, ALTERATION, ABUSE, OR MISUSE, OR THAT HAS BEEN REPAIRED, MOVED, OR ALTERED BY ANYONE OTHER THAN AN AUTHORIZED HOLOGIC SERVICE PERSON. HOLOGIC MAKES NO WARRANTY WHATSOEVER WITH REGARD TO ACCESSORIES OR PARTS USED IN CONJUNCTION WITH THE AQUILEX® FLUID CONTROL SYSTEM NOT SUPPLIED AND/OR MANUFACTURED BY HOLOGIC. THE TERM "ORIGINAL PURCHASER", AS USED IN THE WARRANTY, SHALL BE DEEMED TO MEAN THAT PERSON OR ORGANIZATION AND ITS EMPLOYEES, IF APPLICABLE, TO WHOM THE AQUILEX® SYSTEM WAS SOLD BY HOLOGIC.

### Technical Support and Product Return Information

Contact Hologic Technical Support if the **Aquilex® Fluid Control System** fails to operate as intended. If product is to be returned to Hologic for any reason, Technical Support will issue a Returned Materials Authorization (RMA) number. Return **Aquilex® Fluid Control System** according to the instructions provided by Technical Support. Be sure to clean the **Aquilex® Fluid Control System** with a clean damp cloth and germicide or isopropyl alcohol before returning it and include all accessories in the box with the returned unit.

Hologic and its distributors and customers in the European Community are required to comply with the Waste Electrical and Electronic Equipment (WEEE) Directive (2012/19/EC). Hologic is dedicated to meeting country specific requirements related to the environmentally sound treatment of its products. Hologic's objective is to reduce the waste resulting from the disposal of its electrical and electronic equipment. Hologic realizes the benefits of subjecting such WEEE to potential reuse, treatment, recycling or recovery to minimize the amount of hazardous substances entering the environment. Hologic customers in the European Community are responsible for ensuring that medical devices marked with the following symbol, indicating that the WEEE Directive applies, are not placed into a municipal waste system unless authorized to do so by local authorities.

Contact Hologic Technical Support to arrange for proper disposal of the **Aquilex® Fluid Control System** in accordance with the WEEE Directive.

### *Hologic Technical Support*


#### **United States and Canada:**

Phone: 1 800 442 9892 (toll-free) or 1 508 263 2900

Fax: 1 508 229 2795

#### **Authorized European Representative:**

Phone: +32 2 255 17 74

### 15 Appendix

#### 15.1 Test log

Date	Result	Comment	Signature

### 15.2 Return form

Please fill out this form when returning the device:

**Name of owner:**

**Sales partner:**

**Address of person returning unit:**

**Street:**

**House number:**

<input type="text"/>	<input type="text"/>
----------------------	----------------------

**ZIP/Postal code:**

**City:**

<input type="text"/>	<input type="text"/>
----------------------	----------------------

**Country:**

**IMPORTANT!**

**Serial number (see identification plate):**

**Device type:**

**Description of defect:**

<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>
<input type="text"/>

\_\_\_\_\_  
**Contact**

\_\_\_\_\_  
**Signature**

\_\_\_\_\_  
**Date**

## Glossary

### **Contamination Soiling**

Pollution of rooms, water, foods, objects, or persons due to microorganisms or radioactive materials, biological poisons or chemical agents

### **Contraindication**

Circumstances (e.g., age, pregnancy, certain illness, medication) prohibiting the use of an otherwise indicated measure (contrary to an indication)

### **Embolism**

Sudden capillary blockage due to embolus

### **Flow rate**

Quantity (in ml) of irrigation fluid flowing through tube set per minute

### **Hypervolemia**

An increased volume of circulating blood

### **Hyponatremia**

A low concentration (< 130 mmol/l) of sodium in the patient's bloodstream

### **Hysteroscope**

Endoscope to look inside the uterus

### **Intrauterine pressure**

Pressure in uterine cavity

### **Intravasation**

Entry of foreign matter into a blood vessel

### **Saline**

Isotonic saline solution, i.e., one liter (l) contains 9.0 grams of sodium chloride.

---

**Index**

<b>A</b>		<b>U</b>	
Authorized service technician	10	Unauthorized personnel	41
Authorized trained personnel	41		
<b>C</b>		<b>W</b>	
Care and maintenance	10	Waste management	10
Certification	42		
Check Tube Set Installation	40		
Contamination	10		
<b>D</b>			
Deficit limit	40		
Deficit rate >300 ml/min	40		
<b>E</b>			
Exclusion of liability	10		
<b>G</b>			
Grounding contact	20		
<b>I</b>			
Inspection tests	43		
Intrauterine pressure > 150 mmHg	40		
Intrauterine pressure > 200 mmHg	40		
Intrauterine pressure 10 mmHg above preset intrauterine pressure setting	40		
<b>L</b>			
Liability	41		
Loading/unloading scale while in operation	40		
<b>M</b>			
Mains connection	19		
Manufacturer's specification	43		
Manufacturer's specifications	41		
<b>O</b>			
Only for U.S. operators	20		
<b>P</b>			
Potential equalization	20		
Pressure measuring system malfunctions	40		
Pressure setting at restart	40		
<b>S</b>			
Scale overload	40		
Serious system defect	40		
<b>T</b>			
Technical documents	41		
Two-year maintenance interval	41		

**List of Figures**


Fig. 4-1	System Components.....	21
Fig. 5-1	Front of irrigation pump unit.....	22
Fig. 5-2	Rear of irrigation pump unit .....	22
Fig. 5-3	Fluid monitoring Unit (cart with scale).....	23
Fig. 5-4	Scale and pump connection .....	24
Fig. 5-5	Low vacuum tube.....	27
Fig. 5-6	High vacuum tube.....	28
Fig. 5-7	Front of Device .....	28
Fig. 5-8	Fluid bag suspension .....	29
Fig. 5-9	Outflow tube set.....	31
Fig. 5-10	Port for tissue removal systems.....	32
Fig. 5-11	Tube set elements .....	33
Fig. 5-12	Inserting the tube set.....	33
Fig. 5-13	Attach roller wheel tube.....	34
Fig. 8-1	Opening the fuse holder.....	42
Fig. 9-1	Flow rate test.....	44
Fig. 9-2	Set-up of pressure measuring test.....	45
Fig. 9-3	Test setup for fluid deficit measurement.....	46
Fig. 9-4	Vacuum pump exhaust ports .....	47

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Estas instrucciones de uso contienen informaciones protegidas por el derecho de propiedad (copyright), que forma parte de los derechos de autor. Todos los derechos están protegidos. Sin autorización de W.O.M. WORLD OF MEDICINE GmbH, estas instrucciones de uso no podrá ser ni total ni parcialmente reproducidas ni divulgadas por medio de fotocopia, microfilm u otros medios y procedimientos.
Deze gebruiksaanwijzing bevat auteursrechtelijk beschermde informatie waar copyright op bestaat. Alle rechten voorbehouden. Het is verboden om deze gebruiksaanwijzing zonder toestemming van W.O.M. WORLD OF MEDICINE GmbH geheel of gedeeltelijk door middel van fotokopieën, microfilm of met andere middelen te vermenigvuldigen of te verspreiden.