
Brevera®

Breast Biopsy System

Brevera Σύστημα βιοψίας μαστού
Brevera 乳房生検システム
Brevera Borstbiopsiesysteem
Brevera-bröstbiopsisystem
Brevera Rintabiopsiajärjestelmä
Brevera Brystbiopsisystem
Brevera Brystbiopsisystem
Brevera Sustav za biopsiju dojke
Système pro biopsii prsu Brevera
Brevera Emlőbiopsziás rendszer

Brevera Sistema de biopsia de mama
Brevera Sistem de biopsie mamară
Brevera Système na biopsiu prsníka
Brevera Sistem za biopsiju dojke
Brevera System do biopsiji piersi
Brevera Система за биопсия на гърда
Brevera Sistem za biopsiju dojke
Brevera Sistema de biópsia mamária
Brevera Sistema di biopsia mammaria
Brevera Brustbiopsiesystem
Brevera Système de biopsie mammaire

Hologic Biopsy Guidance System Setup Guide

MAN-09400 Revision 002

Σύστημα καθοδήγησης βιοψίας Hologic Οδηγός ρύθμισης
Hologic 生検ガイダンスシステム セットアップガイド
Hologic-biopsiegeleidingssysteem installatiegids
Konfigurationsguide till Hologic-vägledningssystem för biopsi
Hologic-biopsiaohjausjärjestelmä Asetusopas
Hologic veiledningssystem for biopsi Konfigurasjonsveiledning
Hologic-styresystem til biopsi Installationsvejledning
Sustav za navodenje biopsije Hologic Vodič za postavljanje
Průvodce nastavením systému Hologic pro řízenou biopsii
Hologic biopsziás irányítórendszer Beállítási útmutató
Sistema de guía de biopsia Hologic Guia de configuración
Sistem de ghidare pentru biopsie Hologic Ghid de configurare
Navádzací systém na biopsiu Hologic Sprievodca nastavením
Sistem za vodenje biopsije Hologic Priročnik za nastavitve
System kierowanej biopsji Hologic Instrukcja konfiguracji
Система за насочвана биопсия на Hologic Ръководство за настройка
Систем за биопсијско навођење компаније Hologic Приручник за подешавање
Sistema de orientação da biópsia Hologic Guia de Configuração
Sistema di guida per biopsia Hologic Guida all'impostazione
Hologic Biopsieführungssystem Konfigurationshandbuch
Système de guidage de biopsie Hologic Guide de configuration

HOLOGIC®

Brevera® System

Setup Guide for Hologic Affirm® and Multicare® Platinum Products

PLEASE READ PRIOR TO PERFORMING PROCEDURE:

This setup guide is designed to be used in conjunction with, not replace, the Brevera breast biopsy device instructions for use. Prior to performing the procedure, the physician must review and be familiar with the full operating instructions for the Brevera breast biopsy device and biopsy image guidance system. Consult the Brevera breast biopsy device instructions for use for full operational instructions, as well as any contraindications and safety information.

Please retain this setup guide. Note that while Hologic frequently updates the contents of this guide, medical information and technology can change rapidly. Therefore some of the information in this guide may be out of date. Use of these instructions with any product other than those specifically noted may result in injury to the patient and/or damage to the equipment.

Program the breast biopsy image guidance system for the Brevera breast biopsy device according to the needle setup instructions provided in the appropriate breast biopsy guidance system user manual. Specific parameters are listed below.

Hologic Affirm and Multicare Platinum Products are Compatible with the following Brevera Components and Biopsy Accessories:

Catalog Number	Description
BREVDRV	Brevera Driver
BREVDPTR	Brevera Adapter
BREV09	Standard Stereotactic Biopsy Needle
BREVDISP09	Standard Stereotactic Biopsy Needle
EVIVA_NG09L	Needle Guide

Setup of Hologic Affirm Biopsy Guidance System with or without the Lateral Arm Upright Biopsy Accessory and Affirm Prone Biopsy System for Use with Brevera Biopsy

To Install the Brevera Adapter:

- Align the outer holes in the adapter with the guide pins on the mount.
- Align the center hole with the mount screw.
Check the Brevera Adapter components (reference Figure 1) to ensure they are secure. If the components are not secure call technical support.
- Turn the knob clockwise on the mount to attach the adapter.
- Place needle guide in needle guide holder.
Check to make sure the cable, tubing, and other parts of the Brevera biopsy device are clear of the x-ray tube during rotation.


Figure Legend

1. Adapter
2. Mount Screw
3. Mount Knob
4. Needle Guide Holder

Figure 1: How to Install the Brevera Adapter

To Remove the Brevera Adapter:

1. Turn the knob counterclockwise on the mount to release the adapter.
2. Pull the adapter away from the mount.

Parameter Setup for use of Brevera Biopsy Device on the Affirm Upright Biopsy Guidance System

Parameters	Safety Margin Parameters (mm)	
	BREV09, BREVDISP09	
	20mm (Standard)	12mm (Petite)
Aperture	19.81	13.56
Dead Space	7.87	7.87
Diameter	4.17	4.17
Length	202.74	202.74
Offset	37.34	37.34
Prefire Pullback	-5.23	-0.35
Stroke	23.00	15.00
Front Edge	86.36	86.36
Body Size	16.00	16.00

Parameter Setup for use of Brevera Biopsy Device on the Affirm Prone Biopsy System

Parameters	Safety Margin Parameters (mm)	
	BREV09, BREVDISP09	
	20mm (Standard)	12mm (Petite)
Aperture	19.81	13.56
Dead Space	7.87	7.87
Diameter	4.17	4.17
Length	202.74	202.74
Offset	37.34	37.34
Prefire Pullback	-5.23	-0.35
Stroke	23.00	15.00
Needle Front	86.36	86.36
Needle Width	55.88	55.88

Setup of Hologic Multicare Platinum Prone Biopsy System for Use with Brevera Biopsy

Contact Technical Support to confirm compatibility of your biopsy device. Install the biopsy adapter as shown below, then install the biopsy device. (Refer to the user manual provided by the biopsy device manufacturer for installation and operating instructions.)

To Install the Brevera Adapter:

1. Align the adapter with the pins on the Stage platform, and put the adapter on the platform.
2. Turn the knob clockwise to tighten the adapter.
Check the Brevera Adapter components (reference Figure 2) to ensure they are secure. If the components are not secure call technical support.
3. Place needle guide in needle guide holder.


Figure Legend

- 1. Adapter
- 2. Mount Knob
- 3. Needle Guide Holder

Figure 2: How to Install the Brevera Adapter on the Stage

To Remove the Brevera Adapter:

1. Turn the knob counterclockwise on the mount to release the adapter.
2. Pull the adapter away from the mount.

Parameter Setup for use of Brevera Biopsy Device on the Multicare Platinum Biopsy System (PRONE)

Parameters	Safety Margin Parameters (mm)	
	BREV09, BREVDISPO9	
	20mm (Standard)	12mm (Petite)
Mode	Needle Core Biopsy	Needle Core Biopsy
Stroke	23	15
Differential	-5.2	-0.4
Safety Equation	$Z + 6 \leq \text{compression}$	$Z + 3 \leq \text{compression}$
Minimum Compression	28	22

Cleaning Instructions for Adapter Components

Remove adapter from biopsy guidance system before cleaning.

Use a lint-free cloth and apply a diluted, mild detergent to clean components and surfaces. To remove more difficult soils, use a soft bristle brush as needed.

Caution: Use the least possible amount of cleaning fluids. The fluids must not flow or run.

Disinfecting the Adapter

To disinfect the system’s components and surfaces, Hologic recommends the following disinfecting solutions:

- 10% chlorine bleach and water with one part commercially available chlorine bleach (normally 5.25% chlorine and 94.75% water) and nine parts water. Mix this solution daily for best results
- Commercially available isopropyl alcohol solution (70% isopropyl alcohol by volume, not diluted)
- Commercially available wipes containing less than or equal to 2% quaternary ammonium compounds by weight in water

Apply disinfecting solutions to the system using a clean, lint-free cloth or wipe. Wipe the system with the wet cloth, observing appropriate contact (wet) times. Do not spray cleaning solutions directly on the console or other parts of the Brevera breast biopsy system.

Caution: To prevent damage to the electronic components, do not use disinfectant sprays on the system.

To Prevent Possible Injury or Equipment Damage

Do not use a corrosive solvent, abrasive detergent, or polish. Select a cleaning/disinfecting agent that does not damage the plastics, aluminum, or carbon fiber.

Do not use strong detergents, abrasive cleaners, high alcohol concentration, or methanol at any concentration.

Do not expose equipment parts to steam or high temperature sterilization.

Do not let liquids enter the internal parts of the equipment. Do not directly apply cleaning sprays or liquids to the equipment. Always use a clean, lint-free cloth and apply the spray or liquid to the cloth. If liquid enters the system, disconnect the electrical supply and examine the system before returning it to use.

Caution: Wrong cleaning methods can damage the equipment, decrease imaging performance, or increase the risk of electrical shock.

Always follow instructions from the manufacturer of the product you use for cleaning and disinfection. The instructions include the directions and precautions for the application and contact time, storage, wash requirements, protective clothing, shelf life, and disposal. Follow the instructions and use the product in the most safe and effective method.

Warranty

Except as otherwise expressly stated in the Agreement: i) Equipment manufactured by Hologic is warranted to the original Customer to perform substantially in accordance with published product specifications for one (1) year starting from the date of shipment, or if Installation is required, from the date of installation ("Warranty Period"); ii) digital imaging mammography x-ray tubes are warranted for twenty-four (24) months, during which the x-ray tubes are fully warranted for the first twelve (12) months and are warranted on a straight-line prorated basis during months 13-24; iii) replacement parts and remanufactured items are warranted for the remainder of the Warranty Period or ninety (90) days from shipment, whichever is longer; iv) consumable Supplies are warranted to conform to published specifications for a period ending on the expiration date shown on their respective packages; v) licensed Software is warranted to operate in accordance with published specifications; vi) Services are warranted to be supplied in a workman-like manner; vii) non-Hologic Manufactured Equipment is warranted through its manufacturer and such manufacturer's warranties shall extend to Hologic's customers, to the extent permitted by the manufacturer of such non-Hologic Manufactured Equipment. Hologic does not warrant that use of Products will be uninterrupted or error-free, or that Products will operate with non-Hologic authorized third-party products. These warranties do not apply to any item that is: (a) repaired, moved, or altered other than by Hologic authorized service personnel; (b) subjected to physical (including thermal or electrical) abuse, stress, or misuse; (c) stored, maintained, or operated in any manner inconsistent with applicable Hologic specifications or instructions, including Customer's refusal to allow Hologic recommended Software upgrades; or (d) designated as supplied subject to a non-Hologic warranty or on a pre-release or "as-is" basis.

How Supplied

The Brevera adapter is supplied as a separately packaged accessory.

For More Information

For technical support or reorder information in the United States, please contact:

Hologic, Inc.
250 Campus Drive
Marlborough, MA 01752 USA
Phone: 1-877-371-4372
www.hologic.com

International customers, contact your distributor or local Hologic Sales Representative:

European Representative
Hologic BV
Da Vincilaan 5
1930 Zaventem
Belgium
Tel: +32 2 711 46 80

Symbols Used on Labeling

The following symbols may be found on the product labeling for the Brevera adapter:

Symbol	Description	Standard
	Prescription use only	FDA 21 CFR 801.109
	Authorized Representative in the European Community	ISO 15223-1, Reference 5.1.2
	CE Mark European Conformity	MDR Regulation (EU) 2017/745
	Translations in Box	Hologic
 www.hologic.com/package-inserts	Consult instructions for use	ISO 15223-1, Reference 5.4.3

Symbol	Description	Standard
	Follow instructions for use	IEC 60601-1, Reference No. Table D.2, Safety sign 10 (ISO 7010-M002)
	Medical device	ISO 15223-1, Reference 5.7.7
	Catalog number	ISO 15223-1, Reference 5.1.6
	Contents	N/A
	Quantity	Hologic
	Manufacturer	ISO 15223-1, Reference 5.1.1
	Country of Manufacture	ISO 15223-1, Reference 5.1.11
	Patents	Hologic
	Keep dry	ISO 15223-1, Reference 5.3.4
	Not safe for magnetic resonance imaging	ASTM F2503 Reference no. Table 2, Symbol 7.3.3;7.4.9.1; Fig. 9
	Non-sterile	ISO 15223-1, Reference 5.2.7
	Caution	ISO 15223-1, Reference 5.4.4
	Serial number	ISO 15223-1, Reference 5.1.7
YYYY-MM-DD	Date format: YYYY represents the year MM represents the month DD represents the day	Hologic
	Country code for translation	ISO 3166

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© 2023 Hologic, Inc.