Siemens MAMMOMAT 3000 Nova with Opdima® and Siemens MAMMOMAT Novation DR with Opdima® Setup Guide MAN-09293 Revision 001

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PLEASE READ PRIOR TO PERFORMING PROCEDURE:

This setup guide is designed to be used in conjunction with, not replace, the Eviva® breast biopsy device instructions for use and/or the Siemens MAMMOMAT 3000 Nova with Opdima® operator manual and/or the Siemens MAMMOMAT Novation DR with Opdima® operator manual. Prior to performing the procedure, the physician must review and be familiar with the full operating instructions for the Eviva breast biopsy device and Siemens system(s). Consult the Eviva breast biopsy device instructions for use and the Siemens system(s) operator manual for full operational instructions, as well as any contraindications and safety information. Verification testing was conducted to ensure accuracy of the Eviva biopsy device on the Siemens system(s). However, the setup numbers provided herein are for guidance purposes only. The setup numbers must be verified onsite with the corresponding Eviva biopsy device prior to performing any procedure.

Please retain this setup guide. Note that while Hologic frequently updates the contents of this guide, medical information and technology can change rapidly. Therefore some of the information in this guide may be out of date. Use of these instructions with any product other than those specifically noted may result in injury to the patient and/or damage to the equipment.

Setup of Siemens MAMMOMAT 3000 Nova with Opdima® and Siemens MAMMOMAT Novation DR with Opdima® for Eviva® Biopsy Device


Program the Siemens system(s) for Eviva biopsy device according to the needle setup instructions provided in the Supplement to the Instructions for Use Mammomat 3000 Opdima manual, in conjunction with steps below.

1. Press the Needle Setup button in the service dialog
2. Press Add
3. Select needle type “Core” from the option menu
4. Enter needle length (L1) and aperture length (L2) in the text fields and press OK

Device	Eviva 0910-20 Eviva 1210-20	Eviva 0910-12	Eviva 0910-12T
Recommended Name in Software	EvivaXX1020v	EvivaXX1012	Eviva XX1012T
L1 (mm)	134	130	134
L2 (mm)	18	10	14

Performing a Stereotactic Biopsy Using the Eviva Biopsy Device on the Siemens MAMMOMAT 3000 Nova with Opdima and Siemens MAMMOMAT Novation DR with Opdima

1. Attach the Locating Bar to the Eviva Adapter
2. Install the Eviva Adapter with Locating Bar onto the Siemens system(s).

3. Install the appropriate disposable needle guide, Eviva_NG09R or Eviva_NG12R, onto the Siemens system(s).
4. Ensure the proper Eviva biopsy device is selected in the Siemens system(s) software prior to sending the target to the stereotactic system.

NOTE: The selected Eviva biopsy device should match the biopsy device that is being used in the procedure.

5. Perform the biopsy procedure per the Eviva breast biopsy device instructions for use and per the Supplement to the Instructions for Use Mammomat 3000 Opdima Manual.
6. Perform the biopsy procedure per the Eviva breast biopsy device instructions for use.
7. Clean the *Eviva Adapter* per the cleaning instructions found on page 3 of this document.

Cleaning Instructions for Eviva Adapter

1. Upon completion of the biopsy procedure, immediately rinse the *Eviva Adapter* with warm water.
2. Spray all components with a cleaner

The following cleaning products are **NOT RECOMMENDED** as they may damage the exposed components:

- Bleach based cleaning agents
 - Hydrogen Peroxide
 - Any cleaners/disinfectants with a pH of less than 4 or greater than 10
3. Allow the necessary time for the chosen cleaner to disinfect the components.
 4. Rinse with warm water.
 5. Wipe dry with cloth or paper towel.

Symbols:

The following symbols may be found on the product labeling for the Eviva adapter:

Symbol	Description	Standard
	Authorized Representative in the European Community	ISO 15223-1, Reference 5.1.2
	Prescription use only	21 CFR 801.109
	Caution	ISO 15223-1, Reference 5.4.4
 www.hologic.com/package-inserts	Consult instructions for use	ISO 15223-1, Reference 5.4.3 Hologic
	Follow instructions for use	IEC 60601-1, Reference No. Table D.2, Safety sign 10 (ISO 7010-M002)
	CE Mark with notified body reference number	MDR Regulation (EU) 2017/745
	CE Mark European Conformity	MDR Regulation (EU) 2017/745
	Medical device	ISO 15223-1, Reference 5.7.7
	Translations in Box	Hologic
	Do not use if package is damaged	ISO 15223-1, Reference 5.2.8
	Catalog number	ISO 15223-1, Reference 5.1.6
	Quantity	Hologic

Symbol	Description	Standard
	Manufacturer	ISO 15223-1, Reference 5.1.1
	Country of Manufacture with Date of Manufacture	ISO 15223-1, Reference 5.1.11
YYYY-MM-DD	Date format: YYYY represents the year MM represents the month DD represents the day	Hologic
	Keep dry	ISO 15223-1, Reference 5.3.4
	Serial number	ISO 15223-1, Reference 5.1.7
	Patents	Hologic
	Non-sterile	ISO 15223-1, Reference 5.2.7
	Not safe for magnetic resonance imaging	ASTM F2503 Reference no. Table 2, Symbol 7.3.3;7.4.9.1; Fig. 9
	Country code for translation	ISO 3166

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© 2024 Hologic, Inc. All Rights Reserved. HologicおよびEvivalは、米国および/またはその他の国においてHologic, Inc.および/またはその子会社の商標および/または登録商標です。本文書で使用されているその他すべての商標、登録商標、および製品名は、そのそれぞれの所有者の財産です。
© 2024 Hologic, Inc. All Rights Reserved. Hologic과 Eviva는 미국 및/또는 기타 국가에서 Hologic, Inc. 및/또는 그 회사의 등록 상표입니다. 본 문서에서 사용된 기타 모든 상표, 등록 상표 및 제품명은 각 소유자의 재산입니다.
Hologic, Inc.