



InSight[®] FD Mini C-arm Technical Specifications Manual
MAN-10058-001 Revision 001

InSight[®] FD

Mini C-arm Imaging System

Technical Specifications Manual

Part Number: MAN-10058-001

Revision 001

July 2023

HOLOGIC[®]

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Customer Support

Toll Free in USA: +1.800.321.4659

Email: SkeletalHealth.Support@hologic.com

In Europe, South America, or Asia, contact your local dealer or distributor.

© 2023 Hologic, Inc. Printed in the USA. This manual was originally written in English.

Hologic, Fluoroscanner, InSight, MegaView, and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the U.S. and/or other countries. Any other product and company names mentioned herein are the trademarks of their respective owners.

This product may be protected by one or more U.S. or foreign patents as identified at www.Hologic.com/patent-information.

Table of Contents

1 Product Specifications	1
2 InSight FD Entrance Exposure Rate Data	5
3 Iso-Exposure Rates	8
4 Symbols	17
5 InSight FD System Labels	19
5.1 Manufacturer Serial Number Label	19
5.2 X-ray System Warning Label	20
5.3 Electric Shock Caution/Danger Label.....	20
5.4 Class 1 Laser Product Label	21
5.5 Caution on Incline Label.....	21
5.6 X-ray Controller Compliance Label.....	21
5.7 X-ray Source Compliance Label	22
5.8 Collimator Compliance Label.....	22
5.9 C-arm Compliance Label.....	23
5.10 AC Input Power Fuse Label.....	23
5.11 Printer Fuse Label.....	23
5.12 Extremities Only Label	24
5.13 No Pushing Warning Label.....	24
5.14 ISO 7010-M002 Label.....	25
5.15 Field of View	26
5.16 Computer On/Standby Label.....	26
5.17 Ground Labels.....	26
5.18 Back Panel Labels	27

List of Figures

Figure 1: Phantom Positioning	8
Figure 2: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	9
Figure 3: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	10
Figure 4: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	11
Figure 5: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	12
Figure 6: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	13
Figure 7: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	14
Figure 8: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	15
Figure 9: Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h	16
Figure 10: System Main Label	19
Figure 11: X-ray System Warning Label	20
Figure 12: Electric Shock Warning Label.	20
Figure 13: Class 1 Laser Product Label	21
Figure 14: Caution on Incline Label	21
Figure 15: X-ray Controller Compliance Label	21
Figure 16: X-ray Source Compliance Label	22
Figure 17: Collimator Compliance Label	22
Figure 18: C-arm Compliance Label	23
Figure 19: Input Power Fuse Label	23
Figure 20: Printer Fuse Label	23
Figure 21: Extremities Only Label	24
Figure 22: No Pushing Warning Label*	24
Figure 23: ISO 7010-M002 Label	25
Figure 24: Field of View	26
Figure 25: Computer On/Standby Label	26
Figure 26: Ground Labels	26
Figure 27: Back Panel Labels	27

List of Tables

Table 1: Power Requirements	1
Table 2: Environmental Requirements	1
Table 3: Weights/Dimensions	1
Table 4: Space Requirements	1
Table 5: X-ray Specifications	2
Table 6: Performance	2
Table 7: Imaging Specifications	3
Table 8: Laser Performance Specifications	4
Table 9: Equipment Classification	4
Table 10: Mode: Auto, Standard Dose	5
Table 11: Mode: Auto, High Resolution	5
Table 12: Mode: Auto IQ, Standard Dose	6
Table 13: Mode: Auto, Low Dose	6
Table 14: Mode: Cine Record, Auto	6
Table 15: Mode: Cine Record, Auto IQ	7
Table 16: Mode: Cine Record, Low Dose	7
Table 15: Symbols	17

1 Product Specifications

Table 1 Power Requirements

Parameter	Value
AC Input Voltage	Nominal 100/120/220/230/240 VAC selectable, Single-Phase
Wattage	750 watts maximum
AC Input Frequency	Nominal 50/60 Hz
Power Failure Protection	Images software protected
Power Cord Length	20 ft (6 m)

Table 2 Environmental Requirements

Parameter	Value
Temperature Range (Operating)	+15.5 °C to +32 °C (+60 °F to +90° F)
Temperature Range (Storage)	-20 °C to +50 °C (-4 °F to +122 °F)
Relative Humidity (Operating)	20 – 80% non-condensing
Relative Humidity (Storage and Transport)	20 – 85% non-condensing
Atmospheric Pressure	500 -1060 hPa (375 - 795 mm Hg)

Table 3 Weights/Dimensions

Height	66 in. (168 cm)
Width	28 in. (71 cm)
Depth	32 in. (81 cm)
Weight	475 lb (215 kg)

Table 4 Space Requirements

Required Floor Space:	5 ft x 6.6 ft (1.5 m x 2.0 m)
Door Size	32 in. (81.2 cm) minimum
Ground Clearance:	3.875 in. (9.8 cm) minimum
Floor Capacity	151.1 lb/ft ² (737.2 kg/m ²)
Arm Length	51.8 in. (132 cm) fully extended from the pivot to mid point of detector

Table 5 X-ray Specifications

Source - Image receptor distance	44 cm (17.5 in.)
X-ray source	Grounded tungsten anode X-ray tube with custom designed high voltage generator
X-ray tube window	0.005 in. (0.127 mm) Beryllium
Beam filtration	Aluminum with minimum thickness 3.0 mm
Focal spot	0.045 mm (0.0018 in.) @ 7.5 watts
Field of View (FOV)	Operator selectable collimation of Limited Field or Full Field. Limited Field uses a 4.4 in. x 4.4 in. collimation mode, which magnifies a small area. Full Field uses a 5.7 in. x 4.5 in. collimation mode, which displays large area.
Rated peak tube potential	75 kVp
Tube kVp range	40 to 75 kVp
Tube current range	0.020 to 0.100 mA
Accuracy of displayed values for kV/mA	kVp \pm 5% of displayed value mA \pm 8% \leq 0.035 mA of displayed value mA \pm 5% $>$ 0.035 mA of displayed value
Maximum Duty cycle at 75 kVp@0.100 mA	50%
Maximum Output	0.100 mA @75 kVp

Table 6 Performance

Parameter	Value
Acquisition Frame Rate	Standard Resolution image / Cine 30 FPS High Resolution image* and Low Dose image / Cine are approximately 50% of the Standard Resolution image rate. *Cine recording not available in High resolution mode
Fluoroscopic Time Alarm Range	0 seconds – 5 minutes, (default is 5 minutes)
Image Storage Capacity	10000 Images (maximum)* *Fewer images can be stored at high resolution mode
Cine Loop Capacity	40 minutes maximum
Laser Alignment Guide Timeout	60 seconds (press and release) No timeout (press and hold for 2 seconds)

Table 7 Imaging Specifications

Parameter	Value
CMOS detector	75 micron pixel array 2k x 1.5k resolution Full +/-90 degree rotation Auto-tracking collimation
C-arm dimensions	Source to detector distance: 17.5 in. (44 cm) Free space: 13.5 in. (34.3 cm) Depth of arc: 20 in. (50.6 cm) Orbital rotation: 120 degrees Horizontal travel: 23 degrees Vertical travel: 23 in. (58 cm) Pivot rotation: 380 degrees Panning motion: 320 degrees
Image processing	Digital, with up to 32 frames of video averaging and last frame freeze 16-bit image processing
Video system	24 in. (61 cm) LED touch screen monitor with a native resolution of 1920 x 1200 and an aspect ratio of 16:10. Solid state video camera HDMI output
Resolution (using Nuclear Associates test tool 07-539, located at the minimum SSD in normal (Full Field) mode)	7.0 lp/mm, minimum

Table 8 Laser Performance Specifications

Parameter	Condition	Value
Laser Power	Device on	IEC Class 1
Laser Wavelength	Device on	635 nm
Projected Line Width	At CMOS detector surface	0.05 ± 0.03 in.
Projected Line Length	At center of CMOS detector surface	7.00 ± 0.50 in.
Line Intersection Location Accuracy	At center of CMOS detector surface	± 0.125 in.
Auto-off delay Auto-off disabled	Press and release Press and hold for 2 seconds	Off in 60.0 ± 2.0 seconds Laser remains on until turned off (press again)
Input Power Voltage	N/A	5.00 ± 0.25 volts
Input Power Current	Device off	0.010 mA, max.
Input Power Current	Device on	0.100 mA, max.
Maximum Output	Device on	0.20 mW
Pulse Duration	Device on	CW (Continuous Wave)

The InSight[®] FD Mini C-arm Imaging System is classified* as follows:

Table 9 Equipment Classification

Type of protection against electric shock	CLASS 1
Degree of protection against electric shock	TYPE B
Laser per IEC-60825-1	CLASS 1
Degree of protection against harmful ingress of water	IPX0
Mode of operation	Continuous Operation

The InSight FD Mini C-arm Imaging System is also classified* as:
 Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

* Medical Electrical Equipment - UL 60601-1.

2 InSight FD Entrance Exposure Rate Data

Typical entrance exposure rates, measured at a point in center of the CMOS detector, 5 cm above surface of the detector.

Technique Controls:

- Patient Entry Air Kerma Rate (AKR) and Exposure Rate
- kVp and mA
- FS to AK reference location: 40 cm
- Measurements made at 5cm above IA

Table 10 Mode: Auto, Standard Dose

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.056	0.493	0.056
Limited			50	0.058	0.502	0.057
Full	19 mm Al	Average	63	0.084	1.338	0.152
Limited			64	0.086	1.429	0.162
Full	2 mm Pb	Maximum	74	0.099	2.067	0.235
Limited			74	0.099	2.186	0.248

Table 11 Mode: Auto, High Resolution

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.099	0.879	0.100
Limited			50	0.099	0.915	0.104
Full	19 mm Al	Average	67	0.099	1.832	0.208
Limited			68	0.099	1.882	0.214
Full	2 mm Pb	Maximum	74	0.097	2.210	0.251
Limited			74	0.099	2.210	0.259

Table 12 Mode: Auto IQ, Standard Dose

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.094	0.839	0.095
Limited			49	0.094	0.823	0.094
Full	19 mm Al	Average	67	0.099	1.833	0.208
Limited			67	0.099	1.823	0.207
Full	2 mm Pb	Maximum	74	0.099	2.218	0.252
Limited			74	0.099	2.196	0.250

Table 13 Mode: Auto, Low Dose

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.022	0.193	0.022
Limited			50	0.023	0.207	0.024
Full	19 mm Al	Average	63	0.033	0.524	0.060
Limited			64	0.034	0.565	0.064
Full	2 mm Pb	Maximum	74	0.040	0.868	0.099
Limited			74	0.040	0.865	0.098

Table 14 Mode: Cine Record, Auto

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.056	0.491	0.056
Limited			49	0.056	0.488	0.056
Full	19 mm Al	Average	63	0.084	1.329	0.151
Limited			64	0.086	1.430	0.163
Full	2 mm Pb	Maximum	74	0.099	2.196	0.250
Limited			74	0.099	2.209	0.251

Table 15 Mode: Cine Record, Auto IQ

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.094	0.835	0.095
Limited			49	0.094	0.830	0.094
Full	19 mm Al	Average	67	0.099	1.790	0.203
Limited			67	0.099	1.822	0.207
Full	2 mm Pb	Maximum	74	0.099	2.185	0.248
Limited			74	0.099	2.173	0.247

Table 16 Mode: Cine Record, Low Dose

Field Size	Phantom	Patient Size	kVp	mA	Measured Values	
					Converted to mGy/min at AK Ref Pt.	Exposure R/min
Full	3.3 mm Al	Small	49	0.022	0.196	0.022
Limited			49	0.023	0.212	0.024
Full	19 mm Al	Average	63	0.033	0.514	0.058
Limited			64	0.034	0.564	0.064
Full	2 mm Pb	Maximum	74	0.040	0.855	0.097
Limited			74	0.040	0.858	0.097

3 Iso-Exposure Rates


Scatter radiation measurements around the Fluoroscanner[®] InSight[®] FD Mini C-arm were made with a hand phantom and a foot phantom placed on the imaging assembly as scatter sources.

Figure 1 Phantom Positioning




Phantom	Hand
Field	Full
Mode	Auto, Hi Res
Technique	kV 52 mA 0.099

Figure 2 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h




Phantom	Hand
Field	Limited
Mode	Auto, Hi Res
Technique	kV 51 mA 0.99

Figure 3 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h




Phantom	Hand
Field	Full
Mode	Low Dose, Std
Technique	kV 51 mA 0.024

Figure 4 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h




Phantom	Hand
Field	Limited
Mode	Low Dose, Std
Technique	kV 51 mA 0.024

Figure 5 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h




Phantom	Foot
Field	Full
Mode	Auto, Hi Res
Technique	kV 53 mA 0.099

Figure 6 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h




Phantom	Foot
Field	Limited
Mode	Auto, Hi Res
Technique	kV 54 mA 0.099

Figure 7 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h




Phantom	Foot
Field	Full
Mode	Low Dose, Std
Technique	kV 54 mA 0.026

Figure 8 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h



Phantom	Foot
Field	Limited
Mode	Low Dose, Std
Technique	kV 53 mA 0.025

Figure 9 Contours: 0.1, 0.25, 0.5, 1, 2, 4 and 8 mR/h



4 Symbols

Table 15 Symbols


























	Warning: Electricity		Caution		
	X-ray source assembly		Radiation Filter		
	Type B applied part		Alternating current		
	Earth (ground)		Protective earth (ground)		
	Fuse		Equipotentiality		
	USB connector		Electrostatic sensitive device		
	Temperature limit		Humidity limitation		
	Computer On switch		Computer Standby switch		
	Date of Manufacture		Manufacturer		
	Equipment to be disposed in compliance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment.	<table border="1" data-bbox="813 1562 951 1612"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	Authorized representatives in the European Community
EC	REP				

Table 15 Symbols


	Catalog number		Serial number
	Follow instructions for use		No pushing
	Combined weight of the equipment and its safe working load		CE Mark

5 InSight FD System Labels

5.1 Manufacturer Serial Number Label

The Manufacturer Serial Number Label is on the back of the base cabinet.


Figure 10 System Main Label



5.2 X-ray System Warning Label

The X-ray System Warning Label is on the front of the base cabinet below the Computer On/Standby switch.

Figure 11 X-ray System Warning Label



5.3 Electric Shock Caution/Danger Label

The Electric Shock Caution/Danger Label is located near the fasteners for panels covering high voltage components.


Figure 12 Electric Shock Warning Label.



5.4 Class 1 Laser Product Label

Class 1 Laser Product complies with FDA performance standards for laser products except for deviations pursuant to Laser Notice number 50 dated July 26, 2001.

Figure 13 Class 1 Laser Product Label



5.5 Caution on Incline Label

The Caution on Incline Label is on the back of the unit, above the handle.


Figure 14 Caution on Incline Label



5.6 X-ray Controller Compliance Label

The X-ray Controller Compliance Label is on the X-ray Controller assembly and on the back of the base cabinet.


Figure 15 X-ray Controller Compliance Label



5.7 X-ray Source Compliance Label

The X-ray Source Compliance Label is on the X-ray source module and on the back of the base cabinet.

Figure 16 X-ray Source Compliance Label



5.8 Collimator Compliance Label

The Collimator Compliance Label is on the X-ray source module (near the collimator assembly).


Figure 17 Collimator Compliance Label



5.9 C-arm Compliance Label

The C-arm Compliance Label is on the bottom of the High Voltage Power Supply cover.


Figure 18 C-arm Compliance Label



5.10 AC Input Power Fuse Label

The AC Input Power Fuse Label is on the back of the base cabinet near the AC power connector.


Figure 19 Input Power Fuse Label



5.11 Printer Fuse Label

The Printer Fuse Label is located at the back of the printer compartment near the printer outlet.


Figure 20 Printer Fuse Label



5.12 Extremities Only Label

The Extremities Only Label is located on the back of the base cabinet.

Figure 21 Extremities Only Label



5.13 No Pushing Warning Label

The No Push Warning Label is located on the back of the base cabinet.

Figure 22 No Pushing Warning Label*



*Do not push the system with the casters locked or obstructed and the C-arm extended.

5.14 ISO 7010-M002 Label

Figure 23 ISO 7010-M002 Label



The ISO 7010-M002 Label, refer to instruction manual/booklet, is located on the main label.



Note *On ME Equipment* “Follow instructions for use.”*

**ME EQUIPMENT is equipment*

1. *provided with not more than one connection to a particular SUPPLY MAINS; and*
2. *intended by its MANUFACTURER to be used:*
 - a. *in the diagnosis, treatment or monitoring of a PATIENT; or*
 - b. *for compensation or alleviation of disease, injury or disability.*

having an APPLIED PART or transferring energy to or from the PATIENT or detecting such energy transfer to or from the PATIENT...”.

Note: *The above definition aligns well with the definition for medical electrical devices in the European Medical Device Directive.*

5.15 Field of View

Field of View is represented by raised rectangles on the underside of the High Voltage Power Supply. When the Field of View selection lever is positioned at the large rectangle full field of view is selected. When the Field of View selection lever is positioned at the small rectangle reduced field of view is selected.

Figure 24 Field of View



5.16 Computer On/Standby Label

The Computer On/Standby Label is on the front of the base cabinet above the Computer On/Standby switch

Figure 25 Computer On/Standby Label



5.17 Ground Labels

Each system grounding point, depending upon type of ground, has one of the Ground Labels shown below.

Figure 26 Ground Labels

Earth (ground)



Protective earth (ground)




Equipotentiality



5.18 Back Panel Labels

The Ethernet and DVI labels are located in the middle of the system back panel below their respective connectors.

Figure 27 Back Panel Labels



HOLOGIC®



Hologic, Inc.
600 Technology Drive
Newark, DE 19702 USA
1.800.447.1856

Australia

Hologic (Australia) Pty Ltd.
Suite 402, Level 3
2 Lyon Park Road
Macquarie Park NSW 2113
Australia
1.800.264.073



Hologic, BV
Da Vincilaan 5
1930 Zaventem
Belgium
Tel: +32.2.711.46.80
Fax: +32.2.725.20.87

CE
2797