

CTB-01235 Rev 001

Date: November 8, 2023

Author: Service Engineering

Product: Selenia Dimensions/
3Dimensions

Subsystem: Mobile Systems

Subject: Transportation of Mobile Selenia Dimensions /
3Dimensions Systems

Purpose

The purpose of this document is to provide guidance on the transportation of the Selenia Dimensions and 3Dimensions Mobile systems.

Scope

This bulletin is applicable to all mobile Selenia Dimensions and 3Dimensions systems when preparing the system for travel per *Appendix C: Use of Mobile System* of the Selenia Dimensions and 3Dimensions User Manuals.

Technical Bulletin (cont.)

Reason

The current versions of *Appendix C: Use of Mobile Systems* of the Selenia Dimensions and 3Dimensions User Manuals instructs the user to place the C-Arm of the Gantry at 0° (CC Position) prior to transportation by the mobile vehicle.

The Selenia Dimensions and 3Dimensions User Manuals will be updated to instruct the user to rotate the C-Arm of the Gantry to 180° and move it to its lowest position prior to transportation.

Figure 1 shows an example as to how a C-Arm should be positioned at 180° and lowered prior to transportation.



Figure 1 – Example of C-Arm Rotated to 180° and Moved to Lowest Position

Technical Bulletin (cont.)

By placing the C-Arm at 180° and lowering it closest to the ground, this puts the center of gravity for the C-Arm at its lowest position so that it is not top-heavy, and that the weight distribution is more balanced and stable during transportation.

Please follow this procedure prior to transporting. These instructions will be incorporated into the next edition of the Selenia Dimensions and 3Dimensions User Manuals.