

## ThinPrep® Stain Protocol using Rinse II and Bluing II Solutions on the SAKURA TISSUE-TEK DRS™ 2000 Autostainer

Use this protocol ONLY when ThinPrep Rinse II and Bluing II Solutions (P/Ns ASY-04875, ASY-04876) are components of the ThinPrep Stain.

Use this protocol for staining slides fixed with 95% alcohol.

At the Start Process screen, be sure “Mode” (F3) is set to “Continuous”.

Remove all wash nozzles and insert plugs. Use only Solution Reservoirs in all stations.

**Maintain bath heights to completely cover the slides at full immersion.** (Note: Bath fill volume = 650 ml)

**Caution:** Change Distilled / Instrument Feed Water bath at Station 9 after every **eight racks** of slides.

**NOTE:** Stations 21 and 22 must first be assigned as Wash stations prior to assigning other reagent stations.

Required Setup Parameters (System Setup Screen):					
Mix Amplitude = 30mm, Mix Frequency = 10, Mix Speed = 1, Basket Lifting Speed = 1					
Step	Station	Solution	Duration	Time	Mix
1	S (27)	Start (95% reagent or ethyl alcohol)	—	—' —"	
	S (26)	Start (95% reagent or ethyl alcohol)	—	—' —"	
2	25	70% reagent or ethyl alcohol	Infinity	01' 00"	+
3	24	50% reagent or ethyl alcohol	Infinity	01' 00"	+
4	23	Distilled or Instrument Feed Water (IFW) <sup>1</sup>	Infinity	01' 00"	+
5	8	<b>ThinPrep Nuclear Stain</b>	Exact	07' 00"	+
6	9	Distilled or Instrument Feed Water <sup>1</sup> (Change after 8 racks.)	Exact	00' 10"	+
7	10	<b>ThinPrep Rinse II Solution</b>	Exact	01' 00"	+
8	11	Distilled or Instrument Feed Water <sup>1</sup>	Infinity	00' 30"	+
9	12	<b>ThinPrep Bluing II Solution</b>	Exact	00' 30"	+
10	13	Distilled or Instrument Feed Water <sup>1</sup>	Infinity	00' 30"	+
11	1	Distilled or Instrument Feed Water <sup>1</sup>	Infinity	00' 30"	+
12	2	95% reagent or ethyl alcohol	Infinity	00' 30"	+
13	3	<b>ThinPrep Orange G Solution</b>	Exact	02' 00"	+
14	4	95% reagent or ethyl alcohol	Exact	00' 15"	+
15	5	95% reagent or ethyl alcohol	Exact	00' 15"	+
16	6	<b>ThinPrep EA Solution</b>	Exact	04' 00"	+
17	7	95% reagent or ethyl alcohol	Exact	01' 00"	+
18	20	95% reagent or ethyl alcohol	Exact	01' 00"	+
19	19	100% reagent or ethyl alcohol	Infinity	00' 30"	+
20	18	100% reagent or ethyl alcohol	Infinity	00' 30"	+
21	17	100% reagent or ethyl alcohol	Infinity	00' 30"	+
22	16	Xylene or other Hologic approved clearing agent <sup>2</sup>	Infinity	01' 00"	+
23	15	Xylene or other Hologic approved clearing agent <sup>2</sup>	Infinity	03' 00"	+
24	E (14)	End (Xylene)	Infinity	—' —"	+

Remove slides to a separate clearing bath then coverslip with the appropriate Hologic approved mountant.<sup>2</sup>

<sup>1</sup> Hologic specification for IFW:  $\geq 1.0$  megohm-cm Resistivity **OR**  $\leq 1.0$   $\mu$ Siemens/cm Conductivity (Refer to Clinical and Laboratory Standards Institute (CLSI) document C3-A4, 2006)

<sup>2</sup> See ThinPrep Stain User's Manual, Section 1, "STAINING" or contact Hologic for current list of Hologic approved clearing agents and mounting media.



Hologic, Inc. • 250 Campus Drive • Marlborough, MA 01752 USA • 1-800-442-9892 • [www.hologic.com](http://www.hologic.com)



Hologic BV • Da Vincilaan 5 • 1930 Zaventem • Belgium

## Solution Stations

### Sakura Tissue-Tek DRS™ 2000 Autostainer

Using ThinPrep Rinse II and Bluing II Solutions

For slides fixed with 95% alcohol

<i>Bath Station Number</i> →	25
Solution →	70% Alcohol
Step Number →	<b>2</b>

**Upper Level**

1	2	3	4	5	6	7
Distilled Water	95% Alcohol	<b>OG Solution</b>	95% Alcohol	95% Alcohol	<b>EA Solution</b>	95% Alcohol
<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>

8	9	10	11	12	13
<b>Nuclear Stain</b>	Distilled Water	<b>Rinse II Solution</b>	Distilled Water	<b>Bluing II Solution</b>	Distilled Water
<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>

14	15	16	17	18	19	20	21	22	23	24	25	26	27	DRY
Xylene	Xylene	Xylene	100% Alcohol	100% Alcohol	100% Alcohol	95% Alcohol			Distilled Water	50% Alcohol	70% Alcohol	95% Alcohol	95% Alcohol	
<b>END 24</b>	<b>23</b>	<b>22</b>	<b>21</b>	<b>20</b>	<b>19</b>	<b>18</b>			<b>4</b>	<b>3</b>	<b>2</b>	<b>START 1</b>	<b>START 1</b>	

**Lower Level**