

# ThinPrep® Expansion Stain Protocol for the Sakura Tissue-Tek® Prisma® and Prisma® Plus Autostainers using 95% alcohol fixed slides

Use for staining slides with Small Solution Reservoirs (Sakura P/N 6145).

Use this protocol only for staining slides fixed with 95% alcohol.

**Caution:** Change the Distilled / Instrument Feed water bath at Step 5 (Station 48) after every **NINE** racks of slides. Maintain bath heights to completely cover the slides at full immersion. (Note: Bath fill volume is 280 ml, except for Stations 47 through 52 which are standard, 820 ml baths.)

<b>Required Setup Parameters:</b> Mix (Amplitude = 30 mm, Frequency = 10, Speed = 1), Basket Lifting Speed = 1				
Step	Solution Name	Time	Delay	Mix
1	Start Station (Reagent or Ethyl alcohol: 70%)	--:--:--	-----	-----
2	Reagent or Ethyl alcohol: 50%	0:01:00	**	on
3	Distilled or Instrument Feed Water (IFW) <sup>1</sup>	0:01:00	**	on
4	<b>ThinPrep Nuclear Stain</b>	0:07:00	==	on
5	Distilled or Instrument Feed Water <sup>1</sup> (Replace after every 9 racks.)	0:00:10	==	on
6	<b>ThinPrep Rinse Solution</b>	0:01:00	==	on
7	Distilled or Instrument Feed Water <sup>1</sup>	0:00:30	==	on
8	<b>ThinPrep Bluing Solution</b>	0:00:30	==	on
9	Distilled or Instrument Feed Water <sup>1</sup>	0:00:30	**	on
10	Reagent or Ethyl alcohol: 50%	0:00:30	**	on
11	Reagent or Ethyl alcohol: 95%	0:00:30	**	on
12	<b>ThinPrep Orange G Solution</b>	0:02:00	==	on
13	Reagent or Ethyl alcohol: 95%	0:00:15	==	on
14	Reagent or Ethyl alcohol: 95%	0:00:15	==	on
15	<b>ThinPrep EA Solution</b>	0:04:00	==	on
16	Reagent or Ethyl alcohol: 95%	0:01:00	==	on
17	Reagent or Ethyl alcohol: 95%	0:01:00	==	on
18	Reagent or Ethyl alcohol: 100%	0:00:30	**	on
19	Reagent or Ethyl alcohol: 100%	0:00:30	**	on
20	Reagent or Ethyl alcohol: 100%	0:00:30	**	on
21	Xylene or other Hologic approved clearing agent <sup>2</sup>	0:01:00	**	on
22	Xylene or other Hologic approved clearing agent <sup>2</sup>	0:03:00	**	on
23	End Station (Xylene or other Hologic approved clearing agent <sup>2</sup> )	--:--:--	-----	-----

Remove slides to a separate clearing bath then coverslip with the appropriate Hologic-approved mountant. <sup>2</sup>

\*\* Unlimited, == Exact (no delay)

<sup>1</sup> Hologic specification for IFW:  $\geq 1.0$  megohm-cm Resistivity OR  $\leq 1.0$   $\mu$ Siemens/cm Conductivity (Refer to Clinical and Laboratory Standards Institute (CLSI) document C3-A4, 2006)

<sup>2</sup> See ThinPrep Stain User's Manual, Section 1, "STAINING" or contact Hologic for current list of Hologic approved clearing agents and mounting medias.



Hologic, Inc. • 250 Campus Drive • Marlborough, MA 01752 USA • 1-800-442-9892 • www.hologic.com




Hologic BV • Da Vincilaan 5 • 1930 Zaventem • Belgium

# Expansion Solution Stations

## Sakura Tissue-Tek Prisma and Prisma Plus Autostainers

For staining slides fixed with 95% alcohol



1 <b>THINPREP ORANGE G</b> 12	3 95% ALCOHOL 11	5 50% ALCOHOL 10	7 (W4)	8 (W3)	9 (W2)	10 (W1)	13 (D1a)
2 95% ALCOHOL 13	4 95% ALCOHOL 14	6 DISTILLED WATER 9	11	12	54	55	14 (D1b)
15 95% ALCOHOL 16	17 <b>THINPREP EA</b> 15	19 <b>THINPREP BLUING</b> 8	21	23	25	27	29 (D2a)
16 95% ALCOHOL 17	18 100% ALCOHOL 18	20 DISTILLED WATER 7	22	24	26	28	30 (D2b)
31 100% ALCOHOL 20	33 100% ALCOHOL 19	35 <b>THINPREP RINSE</b> 6	37	39	41	43	45*
32 XYLENE 21	34 XYLENE 22	36 <b>THINPREP NUCLEAR STAIN</b> 4	38	40	42	44	46*
53  LINK	47 (E1) XYLENE 23 / End	48 (E2) POST- NUCLEAR DISTILLED WATER 5	49 (E3) DISTILLED WATER 3	50 (S3) 50% ALCOHOL 2	51 (S2) 70% ALCOHOL 1 / Start	52 (S1) 70% ALCOHOL 1 / Start	

E# = Exit Station

W# = Running Water

S# = Start Station

D# = Drying Station

\* Configuration shown is for Prisma Model 6130 or Prisma Plus Model 6170. In Models 6131 and 6171, Stations 45 & 46 are replaced with heat stations.