

With Aptima[®] mRNA-based HPV testing,
the result comes straight
from the messenger.

The Aptima[®] HPV Assay:

Identifying the presence and threat of a high-risk HPV infection.

The Aptima[®] HPV 16 18/45 Genotype Assay:

Enhanced detection for improved management of HPV positive patients.

Worth it.

Aptima[®] HPV
Assay

Aptima[®] HPV 16 18/45 Genotype
Assay



The Aptima® HPV assay targets E6/E7 mRNA.

Don't just sense a presence. Sense a threat.

Nearly all sexually active men and women will have an HPV infection at some point in their lives. Very few will go on to develop cancer.²


The Aptima® HPV assay targets high-risk HPV mRNA.³

Studies have shown mRNA identifies the presence and activity of a high-risk HPV infection.^{3,4}

E6 E7

E6/E7 mRNA expression is indicative of the HPV infections most likely to lead to disease.^{3,4}

mRNA and Cervical Disease Progression⁵




Because HPV DNA levels may decrease as infections progress toward cancer, some HPV DNA tests may provide false-negative results in more than 10% of the most severe cervical disease cases.⁶

“The optimal screening strategy should identify those cervical cancer precursors likely to progress to invasive cancers (maximizing the benefits of screening) and avoid the detection and unnecessary treatment of transient HPV infection and its associated benign lesions that are not destined to become cancerous (minimizing the potential harms of screening).”¹

Test Design

Aptima® HPV targets oncogenic activity of the HPV virus by detecting E6/E7 mRNA.⁷ Tests that target only the L1 gene are detecting an area that is not needed for disease progression and that can be deleted during integration. There is a risk of up to 10-15% of the most severe disease cases being missed with a L1-based DNA HPV test.⁸⁻¹⁰

Why E6/E7




Deletions may occur during integration
Opening of circular DNA during integration

HPV genome – genotype 16 example

HPV Integration

- HPV DNA must linearize to integrate into human DNA.
- L1 region can be deleted.
- HPV assays that only target the L1 region are at risk for false negative results.¹¹

HPV Detection Strategies




Maximizing the Benefits

With intervals between recommended screenings for cervical cancer extended, identifying those patients at risk becomes increasingly important. Excellent sensitivity means minimizing false-negative test results. The Aptima HPV assay, which targets mRNA, has shown the same excellent sensitivity as DNA-based tests:

The Aptima® HPV assay provides the same excellent sensitivity you have come to expect from DNA-based tests.

Sensitivity†

Clinical Sensitivity for ≥ CIN3 ^{7,16-30}




*Clinical Sensitivity for ≥ CIN2

Specificity†

The Aptima HPV assay has been shown to deliver fewer false-positive test results compared with DNA-based tests.

Clinical Specificity for CIN2 + ^{2-4,6-7,16-26}



†This chart is a representation of clinical data from multiple published sources. The clinical studies represented within these sources were conducted using different study designs with various assays.

Aptima® HPV assay showed

24%


fewer false-positive test results compared to a DNA-based test.⁷

- Minimizing difficult patient conversations
- Minimizing the potential for over-treatment

While Minimizing Potential Harms

Minimizing false-positives helps clinicians target the right patients for colposcopy. In the NILM arm of the CLEAR trial, the Aptima HPV assay showed 24%† fewer false-positive test results compared to a DNA-based test.⁷

Fewer False-Positive Results^{7†}



†The graph above represents data adapted from the Aptima HPV Assay Package Insert Table 22.⁷



Aptima® HPV 16 18/45 Genotype Assay


Enhanced detection for improved management of HPV positive patients.

The Next-Generation Genotype Test

Result for type 16 with separate combined result for HPV types 18 and 45.

Reflex positive Aptima HPV assay results to genotyping for types 16, 18 and 45. Identification of these types as part of reflex testing may identify up to 94% of all cervical adenocarcinomas.²

HPV Genotypes in Invasive Cervical Cancer¹



HPV types 16, 18, and 45 are associated with up to

94%


of HPV-related cervical adenocarcinomas.³¹

HPV type 45:

- Is uncommon and only prevalent in 0.4% of women with normal cytology.³¹
- Is the third most common HPV type in invasive cervical cancer.^{31,32}
- Types 16, 18 and 45 show higher carcinogenic potential relative to all other high-risk HPV types.^{30,33}
- The addition of HPV type 45 in the Aptima HPV 16 18/45 genotype assay identifies more women at risk for adenocarcinoma, with minimal impact to colposcopy rates.³¹

Worth it.





References: 1. Saslow D, et al. American Cancer Society, American Society for Colposcopy and Cervical Pathology, and American Society for Clinical Pathology Screening Guidelines for the Prevention and early Detection of Cervical Cancer. *Am J Clin Pathol* 2012;137:516-542. 2. Genital HPV Infection - CDC Fact Sheet. <http://www.cdc.gov/std/hpv/hpv-factsheet-march-2014-press.pdf>. Published 2014. Accessed August 24, 2015. 3. Tinelli A, et al. HPV viral activity by mRNA HPV molecular analysis to screen the transforming infections in precancer cervical lesions. *Curr Pharm Biotechnol*. 2009;10(8):767-771. 4. Cuschieri K, et al. Human Papillomavirus Type Specific DNA and RNA Persistence—Implications for Cervical Disease Progression and Monitoring. *J Med Virol*. 2004;73(1):65-70. 5. Doorbar J. Molecular biology of human papillomavirus infection and cervical cancer. *Clin. Sci. (Lond)*. 2006 May;110(5):525-41. 6. Wright C, et al. The ATHENA human papillomavirus study: design, methods & baseline results. *Am J Obstet Gynecol*. 2012;206(1):46.e1-46.e11. (Study included cobas® HPV, Hybrid Capture® 2 assay) 7. Aptima HPV Assay [package insert, AW-11141-001 Rev 003 (EN)]. San Diego, CA; Hologic, Inc., 2015. Table#22. 8. De Sanjose. *Lancet Oncol*. 2010;11:1048-56. 9. Wheeler CM, et al. *J Natl Cancer Institute* 2009;101:475-487. 10. Coutlee F, et al. *J. Med. Virol*. 2011; 83:1034-1041. 11. Morris, *Clin Chem Lab Med* 2005; 43(11):1171-1177. 12. Cervista HPV HR CE Package Insert, 15-3053 Rev. 105. 13. cobas c4800 package insert #05641268001-0N. 14. Abbott RealTime High Risk HPV package insert #2N09 Rev 49-2028/R3 15. Hybrid Capture® 2 High-Risk HPV DNA Test package insert #L00665 Rev. 2 2007. 16. Szarewski A, et al. Comparison of predictors for high-grade cervical intraepithelial neoplasia in women with abnormal smears. *Cancer Epidemiol Biomarkers Prev*. 2008;17(11):3033-3042. 17. Reuschenbach M, et al. Performance of p16INK4a-cytology, HPV mRNA, and HPV DNA testing to identify high grade cervical dysplasia in women with abnormal screening results. *Gynecol Oncol*. 2010;119(1):98-105. 18. Clad A, et al. Performance of the Aptima high-risk human papillomavirus mRNA assay in a referral population in comparison with Hybrid Capture 2 and cytology. *J Clin Microbiol*. 2011;49(3):1071-1076. 19. Ratnam S, et al. Aptima HPV E6/E7 mRNA test is as sensitive as Hybrid Capture 2 assay but more specific at detecting cervical precancer and cancer. *J Clin Microbiol*. 2011;49(2):557-564. 20. Ovestad IT, et al. Comparison of different commercial methods for HPV detection in follow-up cytology after ASCUS/LSIL, prediction of CIN2-3 in follow up biopsies and spontaneous regression of CIN2-3. *Gynecol Oncol*. 2011;123(2):278-283. 21. Szarewski A, et al. Comparison of seven tests for high-grade cervical intraepithelial neoplasia in women with abnormal smears: the Predictors 2 study. *J Clin Microbiol*. 2012; 50(6):1867-1873. 22. Eaton, et al. Comparison of the Aptima HPV assay and the cobas HPV test in an ASC-US population [abstract]. Paper presented at: 28th International Papillomavirus Conference. 23. Cuschieri K, et al. Clinical performance of RNA and DNA based HPV testing in a colposcopy setting: Influence of assay target, cut off and age. *J. Clin. Virol*. 2014;59(2):104-108. 24. Cubie HA, et al. Evaluation of commercial HPV assays in the context of post-treatment follow-up: Scottish Test of Cure Study (STOCS-H). *J Clin Pathol*. 2014;67(6):458-463. 25. Wu R, et al. Human papillomavirus messenger RNA assay for cervical cancer screening: the Shenzhen Cervical Cancer Screening Trial I. *Int J Gynecol Cancer*. 2010;20(8):1411-1414. 26. Monsonego J, et al. Evaluation of oncogenic human papillomavirus RNA and DNA tests with liquid-based cytology in primary cervical cancer screening: the FASE study. *Int J Cancer*. 2011;129(3):691-701. 27. Iftner, et al. Comparison of Aptima and HC2 in a routine screening trial in Germany with follow up [abstract]. Paper presented at: 28th International Papillomavirus Conference; Nov 30-Dec 6, 2012; San Juan, Puerto Rico. 28. Cuzick J, et al. Comparing the performance of six human papillomavirus tests in a screening population. *Br J Cancer*. 2013;108:908-913. 29. Nieves L, et al. Primary cervical cancer screening and triage using an mRNA human papillomavirus assay and visual inspection. *Int J Gynecol Cancer*. 2013;23(3):513-518. 30. Guan P, et al. Human papillomavirus types in 115,789 HPV-positive women: a meta-analysis from cervical infection to cancer. *Int. J. Cancer*. 2012;131(10):2349-2359. 31. de Sanjose S, et al. Human papillomavirus genotype attribution in invasive cervical cancer: a retrospective cross-sectional worldwide study. *Lancet Oncol*. 2010;11(11):1048-1056. 32. Hopenhayn C, et al. Prevalence of Human Papillomavirus Types in Invasive Cervical Cancers From 7 US Cancer Registries Before Vaccine Introduction. *J of Low Genit Tract Dis*. 2014;18(3):182-9. 33. Tjalma WA, et al. Differences in human papillomavirus type distribution in high-grade cervical intraepithelial neoplasia and invasive cervical cancer in Europe. *Int. J. Cancer*. 2013;132(4):854-867.

Aptima® HPV
Assay

Aptima® HPV 16/18/45 Genotype
Assay

Diagnostic Solutions | Hologic.com | euinfo@hologic.com

PB-00221-EUR-EN Rev. 001. © 2015 Hologic, Inc. Hologic, Science of Sure, Aptima Combo 2, Aptima, ThinPrep and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to euinfo@hologic.com.

