

HOLOGIC®

SPECIMEN COLLECTION
AND TRANSFER GUIDE

Expand your infectious disease
testing on the Panther® and
Panther Fusion® systems.



GROW ON
PANTHER

Consolidated testing with proven performance		PRD-03546 Aptima® Multitest Swab Specimen Collection Kit						301040 Aptima Urine Specimen Collection Kit		301041 Aptima Unisex Swab Specimen Collection Kit		302657 Aptima Cervical Specimen Collection and Transport Kit	PRD-05110 / 301154c Aptima Specimen Transfer Kit		
		Vaginal Clinician-Collected	Vaginal Self-Collected	Rectal	Throat	Penile Meatal	Skin Lesions Anogenital or Oral	Urine Male	Urine Female	Endocervical	Male Urethral	Endocervical	Cervical ThinPrep®	Cervical Surepath	Skin Lesions Anogenital or Oral in Virus Transport Medium
PRD-05576 Aptima Combo 2® 100 Test	PRD-05571 Aptima Combo 2 250 Test	Y	Y	Y	Y			Y	Y	Y	Y	Y			
302925 Aptima CT 100 Test	302927 Aptima NG 100 Test	Y	Y					Y	Y	Y	Y	Y			
		Store at 2-30°C for up to 60 days.						Store at 2-30°C for up to 30 days.		Store at 2-30°C for up to 60 days.		Store at 2-30°C for up to 30 days before transfer. After transfer store at 2-8°C for up to 30 days or 15-30°C for up to 14 days.			
PRD-03374 Aptima Mgen 100 Test (no calibrators)	PRD-03919 Aptima MG 100 Test (with calibrators)	Y	Y			Y		Y	Y	Y	Y	Y			
		Store at 2-30°C for up to 60 days.						Store at 2-30°C for up to 30 days.		Store at 2-30°C for up to 60 days.		Store at 2-30°C for up to 30 days before transfer. After transfer store at 2-30°C for up to 60 days.			
303163 Aptima TV 250 Test	303209 Aptima TV 100 Test	Y							Y	Y		Y			
		Store at 2-30°C for up to 60 days.						Store at 2-30°C for up to 30 days.		Store at 2-30°C for up to 60 days.		Store at 2-30°C for up to 30 days before transfer. After transfer store at 2-8°C for up to 30 days or 15-30°C for up to 14 days.			
PRD-03568 Aptima HSV 1 and 2 100 Test						Y								Y	
		Store at 2-30°C for up to 60 days.										Store at 2-8°C for up to 3 days before transfer. After transfer store at 2-30°C for up to 30 days.			
303093 Aptima HPV 250 Test	302929 Aptima HPV 100 Test										Y	Y	Y		
303236 Aptima HPV 16 18/45 Genotype 100 Test															
										Store specimens at 2-30°C for up to 60 days.		Transfer specimen within 105 days of collection. Prior to transfer specimens should be stored at 2-30°C with no more than 30 days at temperatures above 8°C. Transferred specimens may be stored at 2-30°C for up to 60 days.		Transfer specimen within 7 days of collection. Prior to transfer specimens should be stored at 2-25°C. Transferred specimens may be stored at 2-25°C for up to 7 days.	
PRD-05186 Aptima BV 100 Test	PRD-05189 Aptima CV/TV 100 Test	Y	Y												
		Store at 2-8°C for up to 30 days and 15-30°C for up to 30 days													

Features and benefits:

- Flexible, can be used to collect samples from a variety of specimen collection sites
- Eliminate manual uncapping, recapping to reduce human error and repetitive motion injuries
- Penetrable cap technology
- Guanidine free-non-toxic media
- Non-Hazardous lysis buffer inactivates common respiratory viruses including SARS-CoV-2, reducing risk in the lab

Simplified workflow with the direct load capture cap collection kit	PRD-03546 Aptima® Multitest Swab Specimen Collection Kit		PRD-06951 Direct Load Capture Cap Collection Kit	
	Nasal Swab	Oropharyngeal	Nasal Swab	Oropharyngeal Swab
PRD-06419 Aptima SARS-CoV-2 250 Test	Y	Y	Y	Y
PRD-06815 Aptima SARS-CoV-2/Flu 250 Test	Store at 2-30°C for up to 6 days.		Store at 2-30°C for up to 6 days.	

Sample transfer	PRD-04339 Panther Fusion® Specimen Lysis Tubes PRD-06554 Hologic Specimen Lysis Tube							
	Naso-pharyngeal VTM/UTM	Nasal VTM/UTM	Oropharyngeal VTM/UTM	Mid-turbinate VTM/UTM	Naso-pharyngeal Wash / Aspirate VTM/UTM	Nasal Aspirates	Saliva samples in Minimum Essential Media	
PRD-06419 Aptima SARS-CoV-2 250 Test	Y	Y	Y	Y	Y	Y	Y	
	Store at 2-8°C for up to 96 hours before transfer. After transfer store at 2-30°C for up to 6 days.						Store at 15-30°C for up to 12 hours before transfer to MEM. Samples in MEM can be stored at 15-30°C for up to 2 hours before transfer. After transfer store at 2-30°C for up to 6 days.	
PRD-06815 Aptima SARS-CoV-2/Flu 250 Test	Y	Y	Y	Y	Y	Y		
PRD-04328 Panther Fusion Flu A/B/RSV assay 96 Test	Y							
PRD-04329 Panther Fusion Paraflu assay 96 Test	Y							
PRD-04330 Panther Fusion Adv/hMPV/RV assay 96 Test	Y							
PRD-04868 Panther Fusion Bordetella assay 96 Test	Y							
	Store at 2-8°C for up to 96 hours before transfer. After transfer store at 15-30°C for up to 6 days or 2-8°C for up to 3 months.							



Multitest



Urine




Unisex



Cervical



Transfer



Hologic BV, Da Vincilaan 5,
1930 Zaventem, Belgium

Notified Body number wherever applicable
EC Representative information wherever applicable

Diagnostic Solutions | Hologic.com | euinfo@hologic.com

SS-01170-EUR-EN Rev 002 © 2021 Hologic, Inc. All rights reserved. Hologic, Aptima, Panther, Panther Fusion, ThinPrep and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshow, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to euinfo@hologic.com.

This document does not replace the Instructions for Use. For a listing of indications, warnings and precautions, specimen collection, storage and handling requirements please refer to the package insert for collection devices and Aptima assays at www.hologic.com

Aptima®