

ShearWave™ Elastography: a reliable and outperforming diagnostic tool for liver fibrosis assessment in chronic hepatitis. A literature review.

Jae Young Lee, MD Associate
Professor, Department of Radiology
Seoul National University Hospital

Address for correspondence:
Jae Young Lee, M.D.
Seoul National University Hospital
101 Daehak-ro, Jongno-gu
Seoul, 110-744, Korea
E-mail: leejy4u@snu.ac.kr

Introduction

Ultrasound (US) imaging plays a major role in the diagnosis, the regular follow-up, and the therapeutic management of chronic liver disease. Its use covers a wide spectrum of clinical applications, such as:

- Analyzing liver parenchyma echotexture and assessing risk of chronic liver disease (such as changes in the size of individual segments or liver dysmorphism and signs of portal hypertension),
- Detecting and characterizing nodules in the cirrhotic liver (and in particular identifying any suspicious lesion such as hepatocellular carcinoma (HCC)),
- Guiding while performing the percutaneous focal treatment (such as radiofrequency-ablation, cryotherapy, etc...) of lesions such as HCC,
- Evaluating therapeutic response.

Conventional US imaging is limited by the subjective nature and the variability in assessing the hepatic parenchyma echotexture alteration and liver dysmorphism, and thus is unable to accurately differentiate hepatic fibrosis stages. However, quantification of hepatic fibrosis is of critical importance in chronic hepatitis not only for diagnosis, but also for antiviral treatment decision-making. Two end-points are clinically relevant: detection of significant fibrosis, which is an indication for antiviral treatment, and detection of cirrhosis, which is an indication for specific monitoring of complications related to portal hypertension and of an increased risk of developing HCC¹.

ShearWave™ Elastography (SWE™) is an ultrasound-based elastography technique that has the ability to map and measure liver stiffness². It has been implemented on a complete ultrasound imaging system, the Aixplorer®, and therefore might address the limitations of conventional US imaging to characterize liver fibrosis. This modality could also become part of the routine examination of liver nodules (e.g. HCC) in cirrhotic contexts.

SWE has three advantages over other methods that perform liver stiffness measurements. Because it is integrated into a diagnostic ultrasound system, the use of grayscale images to guide SWE acquisitions (for example, to avoid large arteries) might increase the repeatability of stiffness measurements³. Also, it should benefit from improved separation of stiffness levels, i.e. fibrosis stages, thanks to the use of shear waves with greater bandwidths⁴. Finally, it provides a real-time, two-dimensional, quantitative, color-coded map of liver tissue stiffness. The spatial heterogeneity of liver stiffness can be visualized and the Q-Box™ (region of interest) size used for a measurement can be selectively placed and/or adjusted to target a homogenous part of the liver parenchyma. As a result, physiological variations of liver fibrosis can be averaged out. Its real time aspect also ensures that excessive liver motion is avoided.

We have reviewed the clinical results that have been reported in the literature up to September 2013 and we are providing an interpretation of these reports, taking into account our own experience of SWE in the assessment of liver fibrosis.

ShearWave™ Elastography has a low technical failure rate

Hudson et al recently investigated the reproducibility of SWE in healthy volunteers and demonstrated that 98% of SWE images were quantifiable in liver segments 6 and 8, whereas this percentage decreased to 83% in segment 2/3⁵.

When performed on the right liver lobe through the intercostal space on liver segments 6 and 8, SWE demonstrated a failure rate ranging from 2% to 3%. This low rate can be positively compared to the failure rate of Transient Elastography (TE), usually reported as ranging between 2.4-9.4%⁶. The difference between the 2 techniques may be due to the fact that SWE measurements are not impacted by the presence of ascites. Shared limiting factors for both techniques include narrow intercostal spaces and obesity^{3,7-8}. Ferraioli et al excluded patients with ascites from their study population and therefore could observe that the technical failure rate (2.5%) for both techniques was due to the same patients conditions, i.e. narrow intercostal spaces in 2 patients and a BMI>32 kg/m² in 1 patient³. However, SWE may be

less impacted by obesity as extra pressure on the probe reduces the thickness of the fatty layer between the probe and the rib cage, and the depth of SWE measurements can be adapted to go down to 10-12 cm.

The experience of Leung showed different conclusions: SWE was successful in 449/454 (98.9%) subjects including patients and healthy volunteers, while TE was successful in 407/454 (89.6%). Similarly to Ferraioli's experience, common reasons for technical failures between SWE and TE were obesity and narrow intercostal spaces. In addition, the inability of patients to perform an optimal breath suspension was also a factor for technical failure⁹.

Reproducibility of ShearWave™ Elastography

Intra-observer Reproducibility

In a prospective study to investigate the reproducibility of SWE™ measurements in normal livers, Ferraioli et al demonstrated that the intraclass correlation coefficients (ICC) for the intraobserver agreement were close to perfect for measurements performed on the same day by both an expert and a novice operator¹⁰. Similarly, Hudson reported that the reproducibility of SWE measurements was almost perfect on liver segments 6 and 8, with the maximum reproducibility obtained for measurements performed on the same scanning session by the most experienced operator⁵. Intra-operator reproducibility of SWE measurements on liver segment 2/3 was also very good (>0.60), although lower than that on segments 6 and 8.

The fact that the lowest agreement (although showing a good ICC > 0.60) was obtained for measurements performed by operators with less experience in ultrasound imaging and between 2

different scanning sessions, suggests that SWE operators must ensure good imaging technique, in order to reproduce a given scanning imaging plane over time.

Using the Bland-Altman statistical analysis, Ferraioli demonstrated that the mean differences between measurements within a scanning session or between scanning sessions were 0.01 kPa for the expert and -0.01 kPa for the novice, and 0.06 kPa for the expert and 0.26 kPa for the novice, respectively¹⁰.

Leung et al also reported intra-observer reproducibility data with ICCs ranging from 0.86 to 0.98 for 3 different operators⁹.

	Expert		Novice	
	Same day	Different days	Same day	Different days
Ferraioli ¹⁰	0.95	0.84	0.93	0.65
Hudson ⁵	0.91	0.63	0.92	0.84

Table 1. Within-session and between-session intra-operator ICCs of SWE measurements performed on right liver lobe segments.

Intra-observer Reproducibility

Ferraioli reported an inter-observer ICC of 0.88, indicating almost perfect agreement between 2 operators on liver segments 6 and 8, with a mean difference between measurements reported to be -0.12 kPa by the Bland-Altman analysis¹⁰. Similarly, Leung's experience showed an ICC of 0.85 (95% CI: 0.70-0.94) for the inter-observer reproducibility.

In Hudson's study, the inter-operator ICCs were 0.78 and 0.76 for segments 6 and 8, respectively. As in previous experience, the inter-operator agreement was poorer for

measurements performed on liver segment 2/3 (ICC=0.65)⁵. Liver measurements performed by 2 operators were found to be not statistically different in segments 6 and 8 ($p=0.16$ and $p=0.20$, respectively), whereas they led to significantly different measurements in segment 2/3 ($p=0.02$). However, the analysis showed that only 1% to 8% of the variance was due to the operator.




Figure 1. A 62 year-old man with chronic hepatitis. Mean elasticity was 8.0 kPa with ShearWave Elastography. Biopsy confirmed a METAVIR F2 liver fibrosis.




Figure 2. A 68 year-old woman with liver cirrhosis. Mean elasticity was 31.6 kPa with ShearWave Elastography. Surface nodularity was seen, representative of cirrhosis. Biopsy confirmed liver cirrhosis.

Diagnostic Performance of ShearWave™ Elastography in Chronic Hepatitis Patients

In the first study reporting the diagnostic performance of the supersonic shear imaging (SSI) technique (on which SWE™ is based) to evaluate liver fibrosis in patients with Hepatitis C, good correlation was found ($r=0.8296$) between the elasticity measured with SWE and TE, although a mean offset of 2.40 kPa was observed between the 2 techniques⁷. In addition to the measurement of liver stiffness, SWE also provides information on the heterogeneity of liver stiffness, which cannot be assessed with TE. As shown in Table 2, the Receiver Operating Characteristic (ROC) analysis of liver stiffness measurements performed with SWE showed areas under the ROC curve (AUROC) all greater than 0.95 for the diagnosis of significant fibrosis (METAVIR F \geq 2), severe fibrosis (METAVIR F \geq 3), and liver cirrhosis (METAVIR F4). The authors also assessed a better accuracy of SWE over TE, on the basis of several criteria such as misclassification rates, Youden's

index, specificity at 95% of sensitivity, and sensitivity at 95% of specificity.

As was further demonstrated by Ferraioli et al in patients with hepatitis C (Table 2), the use of different cut-off values for SWE and TE favored SWE in the assessment of early fibrosis stages³. The use of different cut-off values for both techniques is supported by the fact that the Young's modulus (corresponding to the liver stiffness) is derived from the shear group velocity, which is measured from the broadband mechanical excitation (60 Hz–600 Hz) generated using the acoustic radiation force for SWE, whereas TE elasticity values are assessed using an external vibrator acting at a single frequency of 50 Hz²⁻⁴. Therefore, the stiffness assessed by SWE is based on higher frequency vibrations and integrates both elasticity and viscosity properties⁷.

AUROC ¹	CLD ²	F \geq 2		F \geq 3		F=4	
		SWE	TE	SWE	TE	SWE	TE
Bavu ⁷	Hep C	0.95	0.85 P=0.005	0.96	0.86 P=0.001	0.97	0.94 P=0.15
Ferraioli ³	Hep C	0.92	0.84 P=0.002	0.98	0.96 P=0.14	0.98	0.96 P=0.48
Leung ⁹	Hep B	0.88	0.78 P=0.01	0.933	0.83 P=0.01	0.98	0.92 P=0.04

Table 2. Summarized performances of SWE and TE in published studies.1 Area under the ROC curve; 2 Chronic liver disease.

In results published by the Liver Fibrosis Study Group in Pavia, Italy, the optimal cutoff values for the diagnosis of METAVIR F \geq 2, F \geq 3 and F4 were 7.1 kPa (1.5 m/s), 8.7 kPa (1.7 m/s), and 10.4 kPa (1.9 m/s), respectively, with sensitivities of 87.5 to 97.3 % and specificities of 87.5 to 96.8 %³. Leung et al reported their experience on 226 chronic hepatitis B carriers, and demonstrated a better correlation to fibrosis METAVIR scores for stiffness measurements performed with SWE ($r=0.81$) as compared to TE ($r=0.58$)⁹. The values for AUROCs on this population of chronic hepatitis B infected patients are reported in Table 2. The calculated cut-off values were consistent with those found in other studies on chronic hepatitis C patients: 6.5 kPa (METAVIR F \geq 1), 7.1 kPa (METAVIR F \geq 2), 7.9 kPa (METAVIR F \geq 3), and 10.1 kPa (METAVIR F4)⁹. Interestingly, this group also studied the diagnostic value of spleen stiffness in staging liver fibrosis. AUROCs and the optimal cut-offs elasticity values (in kPa) for spleen stiffness for METAVIR F \geq 1, F \geq 2, F \geq 3, and F4 were found to be 0.81 and 19.4, 0.82 and 19.8, 0.83 and 20.6, and

0.84 and 22.0, respectively. The authors could not demonstrate any improvement of the AUROCs by combining the stiffness information from the liver and the spleen. However, the authors concluded that spleen stiffness measured with SWE may serve as an ancillary parameter to detect advanced fibrosis.

Another recent paper reported the results of a study performed without using biopsy as the gold standard¹¹. This study used the cut-off values defined for TE for both the TE and the SWE techniques. Therefore, it showed decreased performances of SWE™ as compared to TE. However, when analyzed separately, these results showed that the AUROC of SWE in differentiating F1 from F2, and F2 from F3, were higher than those of TE: 0.590 versus 0.574, and 0.600 versus 0.509, respectively.

SWE demonstrated a good correlation with other elastography techniques such as acoustic radiation force impulse imaging (ARFI) in patients with Hepatitis C, although diagnostic performances were not assessed due to the limited number of patients⁸.

Conclusion

This review of the results of ShearWave™ Elastography in assessing liver fibrosis and cirrhosis shows that SWE, only available on the Aixplorer®, has better performances in identifying early stages of liver fibrosis (especially METAVIR F_{≥2}) and similar performances in assessing liver cirrhosis, as compared to other elastography techniques currently available. However, specific stiffness cut-off values should be used, due to the inherent technical differences. SWE provides a 2D quantitative map of liver stiffness, thus the spatial heterogeneity of liver stiffness can be visualized in real time and easily averaged to better analyze the overall fibrosis state.

In our experience, as well as in the literature review, this map has proven to be very useful to avoid artifacts arising from pulsating vessels, reverberation, or motion. As a consequence, SWE has demonstrated almost perfect intra-observer reproducibility and a very good inter-observer reproducibility. As it is available on a premium ultrasound imaging system, which encompasses other imaging modes such as gray scale imaging, Doppler modes and contrast-enhanced ultrasound imaging, Aixplorer and SWE offer a complete diagnostic tool to assess chronic liver diseases.




Figure 3. A 52-year-old man with chronic hepatitis. a. On B-mode ultrasound, coarse echotexture, an ultrasound feature representing chronic hepatitis. b. Mean elasticity was 15.4 kPa on ShearWave Elastography. Biopsy confirmed a METAVIR F3 liver fibrosis.




Figure 4. A 61-year-old man with liver cirrhosis. a. Shrunken liver with surface nodularity and a large amount of ascites was clearly seen on B-mode ultrasound. b. On ShearWave Elastography, mean elasticity was 42.9 kPa.

No correlation to biopsy was available for this patient.

Appendix. Staging Methods for Liver Fibrosis

Chronic liver diseases are known to be diffuse, heterogeneous^{12,13}, and usually combine hepatocyte and/or cholangiocyte necrosis or apoptosis with inflammation and interstitial fibrosis. The extension of the latter may result in alterations of the hepatic architecture and the appearance of regeneration nodules, which define cirrhosis.

Liver Biopsy

The outcome of liver biopsy has traditionally been considered as the standard of reference for assessing liver fibrosis severity in patients with chronic liver diseases, and especially those with chronic hepatitis. Liver biopsy can be performed percutaneously, or by a transvenous route in case of hemostasis disorder. One of the main advantages of biopsy is that it provides additional information about the inflammatory reaction, the level of steatosis. Nevertheless, it has several drawbacks:

- It is an invasive technique, which is associated with morbidity (3%, including 0.6% severe complications) and mortality (approximately 1%).
- It is expensive, requiring a day of hospitalization¹⁴,
- It can lead to false outcomes: the biopsy core sample is not very large (<25 mm in length and 1 mm in diameter) and may not be representative of the liver fibrosis heterogeneity. Therefore, the diagnosis of fibrosis seems to be underestimated in 10 to 30% of cases¹⁵
 - Although the histo-pathological outcome of liver biopsy has been standardized by the use of scoring systems such as the METAVIR or the Ishak scores, these semi-quantitative methods show an inter-observer variability¹⁶. Indeed, the percentage of fibrotic areas that can be measured for successive intermediate METAVIR scores are very similar to each other: 2.0±0.1% for F0, 3.4±0.3% for F1, 5.8±0.7% for F2, 14.7±0.8 for F3, 25.1±1.4% for F4¹⁷.

• As a consequence, for biopsy cores under 20 mm in length, there is an increased risk (>30%) of misclassification of intermediate METAVIR stages¹⁷. Also the assessment of liver fibrosis by several pathologists can show a very high discordance rate (>60%) for intermediate METAVIR stages¹⁸.

- Liver biopsy is not ideal for repeated assessment of disease progression¹⁸.

Both the progression and the regression of hepatic fibrosis over time could be of clinical significance. Recent research has demonstrated a reduction in liver fibrosis with treatment, even in advanced stages^{19,20}.

Therefore new, non-invasive techniques to assess hepatic fibrosis have been an important focus of research in hepatology for the last 10 years. Currently available methods rely on two different approaches: a "biological" approach based on the dosage of serum biomarkers of fibrosis^{21,23}, and a "physical" approach based on the measurement of liver stiffness^{24,25}.

Non-invasive Staging

Although the large number of publications over the past decade confirms the growing interest regarding these new non-invasive methods, specific limitations must be considered. As an example, in most studies the fibrosis level was derived from the liver biopsy METAVIR score, which suffers from its own limitations described above for intermediate stages. Other examples preventing an accurate assessment of intermediate stages are the variability of the measurements and the limited sampling used for the measurement²⁶.

Serum Biomarkers

Serum markers are used to calculate a fibrosis score from the measurements of biological parameters. Several tests are available to the clinicians depending on the etiology of the underlying chronic liver disease. The FibroMeter^{®22}, the Hepascore^{®23} and the FibroTest[®] are amongst the most used blood tests. The latter combines the dosage of 5 serum markers (α_2 -macroglobulin, haptoglobin, apolipoprotein A1, total bilirubin, γ -glutamyltranspeptidase) with an adjustment for sex, age and body mass index (BMI). It has been extensively studied and has demonstrated a diagnostic accuracy ranging from 70% to 85%²¹. However, it has limitations in cases of hyperbilirubinemia, hemolysis, inflammation or concomitant illness. All serum markers and blood tests share similar strengths and limits. They are not routinely available in most hospital settings, therefore limiting their clinical use.

Elastography Techniques

Conventional imaging techniques provide anatomical, hemodynamic and perfusion information, which are valuable in the context of focal diseases, but are of limited benefit in diffuse chronic liver diseases. The elasticity (or, equivalently, stiffness) of body tissues varies greatly and is a parameter that can be coded to differentiate tissues and also lesions in surrounding tissues²⁷. Many disease processes produce changes in tissue elasticity. Tumors (especially malignant) are generally harder than surrounding normal tissue. Interstitial fibrosis, which appears in some diffuse diseases (liver cirrhosis, renal failure...), also causes a change of elasticity^{28,29}. As a result, additional information on the viscoelastic properties of the organs or tumors is of great interest to the clinicians.

Elasticity imaging of the human body is a fairly new modality currently being evaluated. It proposes replacement of subjective palpation with imaging the elastic properties of tissue. Static elastography is currently available on many ultrasound diagnostic imaging devices. However, it does not provide quantitative values of the elastic properties of tissues. Elastography imaging is also being developed in MRI (Magnetic Resonance Elastography, MRE, or elasto-MR)^{28,30}.

Three other techniques, based on the properties of shear waves, have been developed in the last decade to quantitatively measure the elastic properties of tissues. Indeed, the speed of a shear wave propagating in a medium is directly related to the longitudinal modulus of elasticity of the biological tissue; the tissue elasticity modulus can then be derived from this measurement. Accordingly, the shear wave speed in stiff or "hard" tissue will be greater than in a softer region.

One-Dimensional Transient Elastography

The first technique, called Transient Elastography (TE) is a one-dimensional non-invasive, non-imaging, bedside method to evaluate liver fibrosis by measuring liver stiffness²⁴. This technique is dedicated to liver fibrosis assessment and allows/permits the diagnosis of cirrhosis and significant fibrosis.

The shear wave is generated by an external low frequency vibrator (50 Hz), which strikes the patient's skin. This external pitch is sufficient to produce a shear wave whose propagation is measured by a one-dimensional ultrasound system (approximately 5 MHz) and provides an average elasticity measurement. This technique is currently commercially available (FibroScan[®], Echosens[™], Paris, France).

It has been widely studied and validated in clinical practice to measure the elasticity of the liver parenchyma in a cylindrical volume sample²⁴. The measurement is typically performed intercostally on the right liver lobe and covers a small (30 - 40 mm) region of interest (from a given depth). The outcome is a value that corresponds to the average elasticity in the single explored cylinder. The measurement is typically repeated 10 times and the median is considered to be the representative elasticity value.

When hepatic elasticity (liver stiffness) measured with TE produces values greater than 12.5-14.5 kPa, cirrhosis could be diagnosed with a high positive predictive value^{32,35}. Significant fibrosis could be suggested by TE when elasticity values would be greater than 7.1-8.7 kPa^{32,34}. Among all the non-invasive approaches that have been developed and evaluated to stage liver fibrosis in the last decade, TE is the only tool that has successfully entered clinical practice, particularly in Europe, and is now reimbursed in some countries.

However, there is considerable variation in the performances reported for TE to predict significant fibrosis in the literature (AUROCs of 0.75 to 0.91), most probably due to the known limitations of the technique³⁶. In fact, the majority of failed TE exams (between 2.4% and 20%) were reported to originate from variability within the acquisitions^{36,37}. These limitations are:

- Low volume of parenchyma explored,
- Absence of ultrasound imaging to guide the measurement,
- Spatial distribution of liver elasticity is not provided.
- Measurement/technical difficulties in obese patients and those with ascites,
- Lack of specificity for the distinction of significant fibrosis level,
- Learning curve required to acquire correctly, without imaging guidance.

ARFI-Based Techniques

The second technique, Acoustic Radiation Force Impulse (ARFI) quantification is also a one-dimensional technique but has been integrated onto a conventional ultrasound imaging system²⁵. Unlike TE, it relies on the mechanical excitation of tissue by providing localized, bursting, acoustic radiation force. This results in the propagation of a shear wave away from the region of excitation. Using conventional beamforming architecture, beams are continuously transmitted until the passing shear wave front is detected. Like TE, ARFI-based systems are commercially available to measure tissue stiffness. Limitations include:

- No elasticity map of tissue,
- The elasticity measurement is not real time,
- The elasticity measurement cannot be performed retrospectively,
- Only one acquisition can be acquired at a time,
- The evaluated area of parenchyma is a small pre-determined size and cannot be modified,
- Only the average elasticity in the ROI is calculated, without any information on standard deviation,
- The depth of the ROI is restricted due to transducer limitations; limiting the frequency and magnitude of push pulses prevents excessive heating.

ShearWave[™] Elastography

The third quantitative imaging technique is ShearWave[™] Elastography (SWE[™]) and has been implemented on the Aixplorer[®] ultrasound imaging system². SWE allows two-dimensional, real time, quantitative imaging of tissue elasticity in combination with conventional ultrasound grayscale imaging. This technique has been validated for the characterization of breast lesions³⁸⁻⁴⁰ and thyroid nodules⁴¹⁻⁴⁵, for the staging of liver fibrosis^{3,7,9}, for the diagnosis of liver nodules⁴⁶ and for the detection and characterization of prostate cancer^{47,48}.

SWE relies on the measurement of the shear wave propagation speed in soft tissue. Like ARFI-based techniques, it does not require an external vibrator to generate the shear wave and it is based on the generation of a radiation force in the tissue to create the shear wave. However, in SWE, several focal points are generated almost simultaneously, in a line perpendicular to the surface of the patient's skin. This creates a conical shear wave front around the focal point, which sweeps the image plane on both sides⁴⁹.

The progression of the shear wave is then captured by UltraFast[™] imaging: the very rapid acquisition of ultrasound images (up to 20,000 images per second). The process takes only a few milliseconds. The high-speed acquisition is necessary to capture the shear wave as its propagation velocity ranges from 1 to 10 m/s. A comparison of two consecutive ultrasound images allows the measurement of displacements induced by the shear wave and creates a "movie" showing the propagation of the shear wave whose local speed is linked to local tissue elasticity. The propagation speed of the shear wave is then estimated from the movie that is created and a real-time, two-dimensional color map is displayed. The color codes either the elasticity of the medium in kilopascals (kPa) or the shear wave speed in meters per second (m/s). This color map is displayed on top of the anatomic grayscale (or B-mode) image⁴⁹.

Using a region-of-interest quantification tool (ROI) called the Q-Box[™], local tissue elasticity or shear wave velocity can be measured retrospectively over an area of interest ranging from 1 to 700 mm². Since each pixel in the color-coded map corresponds to a tissue elasticity measurement, the stiffness of the tissue is locally assessed. Additionally, the automatic standard deviation calculation provides relevant information on the stiffness value distribution within the region of interest.

References

- 1- Diagnosis, Management and Treatment of hepatitis C: an update. Ghany MG, Strader DB, Thomas DL, Seeff LB. *Hepatology* 2009; 49: 1335–1374.
- 2- Quantitative viscoelasticity mapping of human liver using supersonic shear imaging: preliminary in vivo feasibility study. Muller M, Gennisson JL, Deffieux T, Tanter M, Fink M. *Ultrasound Med Biol*. 2009 Feb;35(2):219-29.
- 3- Accuracy of real-time shear wave elastography for assessing liver fibrosis in chronic hepatitis C: A pilot study. Ferraioli G, Tinelli C, Dal Bello B, Zicchetti M, Filice G, Filice C. *Hepatology*. 2012 Dec;56(6):2125-33
- 4- Shear wave spectroscopy for in vivo quantification of human soft tissues visco-elasticity. Deffieux T, Montaldo G, Tanter M, Fink M. *IEEE Trans Med Imaging*. 2009 Mar;28(3):313-22.
- 5- Inter- and Intra-Operator Reliability and Repeatability of Shear Wave Elastography in the Liver: A Study in Healthy Volunteers. Hudson JM, Milot L, Parry C, Williams R, Burns PN. *Ultrasound Med Biol*. 2013 Jun;39(6):950-5.
- 6- Non-invasive evaluation of liver fibrosis using transient elastography. Castera L, Forns X, Alberti A. *J Hepatol* 48:835-847
- 7- Noninvasive in vivo liver fibrosis evaluation using supersonic shear imaging: a clinical study on 113 hepatitis C virus patients. Bavu E, Gennisson JL, Couade M, et al. *Ultrasound Med Biol* 37:1361-1373
- 8- Use of Aixplorer for detection of liver fibrosis or cirrhosis in patients with hepatitis C. Karlas T, Tröltzsch M, Wiegand J, Keim V. Oral Communication WFUMB 2011
- 9- Quantitative Elastography of Liver Fibrosis and Spleen Stiffness in Chronic Hepatitis B Carriers: Comparison of Shear-Wave Elastography and Transient Elastography with Liver Biopsy Correlation. Leung VY, Shen J, Wong VW, Abrigo J, Wong GL, Chim AM, Chu SH, Chan AW, Choi PC, Ahuja AT, Chan HL, Chu WC. *Radiology*. 2013 Aug 2.
- 10- Reproducibility of real-time shear wave elastography in the evaluation of liver elasticity. Ferraioli G, Tinelli C, Zicchetti M, Above E, Poma G, Di Gregorio M, Filice C. *Eur J Radiol*. 2012 Nov;81(11):3102-6
- 11- Liver fibrosis evaluation using real-time shear wave elastography: Applicability and diagnostic performance using methods without a gold standard. Poynard T, Munteanu M, Luckina E, Perazzo H, Ngo Y, Royer L, Fedchuk L, Sattoune F, Pais R, Lebray P, Rudler M, Thabut D, Ratziu V. *J Hepatol*. 2013 Jan 12. doi:10.1016/j.jhep.2012.12.021.
- 12- Estimation of collagen content of liver specimens. Variation among animals and among hepatic lobes in cirrhotic rats. Gascon-Barre M, Huet PM, Belgiojorno J, Plourde V, Coulombe PA. *J Histochem Cytochem* 1989; 37:3:377-381.
- 13- Observer error and sampling variability tested in evaluation of hepatitis and cirrhosis by liver biopsy. Soloway RD, Baggenstoss AH, Schoenfield LJ, Summerskill WH. *Am J Dig Dis* 1971;16:1082-1086.
- 14- Practices of liver biopsy in France: results of a prospective nationwide survey. Cadranel JF, Rufat P, Degos F, for the Group of Epidemiology of the French Association for the Study of the Liver (AFEFL). *Hepatology*. 2000 Sep;32(3):477-81.
- 15- The role of laparoscopy in the diagnosis of cirrhosis. Poniachik J, Bernstein DE, Reddy KR, Jeffers LJ, Coelho-Little ME, Civantos F, Schiff ER. *Gastrointest Endosc*. 1996 Jun;43(6):568-71.
- 16- Checkmate to liver biopsy in chronic hepatitis C? Trifan A, Stanciu C. *World J Gastroenterol*. 2012 Oct 21;18(39):5514-20.
- 17- Sampling variability of liver fibrosis in chronic hepatitis C. Bedossa P, Dargère D, Paradis V. *Hepatology*. 2003 Dec;38(6):1449-57.
- 18- Sources of variability in histological scoring of chronic viral hepatitis. Rousselet MC, Michalak S, Dupré F, Croué A, Bedossa P, Saint-André JP, Calès P; *Hepatitis Network* 49. *Hepatology*. 2005 Feb;41(2):257-64.
- 19- Reversibility of cirrhosis in HIV/ HBV coinfection. Mallet VO, Dhalluin-Venier V, Verkarre V, Correas JM, Chaix ML, Viard JP, Pol S. *Antivir Ther*. 2007;12(2):279-83.
- 20- Reversibility of liver fibrosis. Ramachandran P, Iredale JP. *Ann Hepatol*. 2009 Oct-Dec;8(4):283-91.
- 21- Biochemical markers of liver fibrosis in patients with hepatitis C virus infection: a prospective study. Imbert-Bismut F, Ratziu V, Pieroni L, Charlotte F, Benhamou Y, Poynard T. *Lancet* 2001;357:1069-1075.
- 22- FibroMeters: a family of blood tests for liver fibrosis. Calès P, Boursier J, Oberfi F, Hubert I, Gallois Y, Rousselet MC, Dib N, Moal V, Macchi L, Chevailler A, Michalak S, Hunault G, Chaigneau J, Sawadogo A, Lunel F. *Gastroenterol Clin Biol*. 2008 Sep;32(6 Suppl 1):40-51.
- 23- Noninvasive diagnosis of hepatic fibrosis in chronic hepatitis C. Stauber RE, Lacker C. *World J Gastroenterol* 2007;13 (32):4287-94.
- 24- Transient elastography: a new noninvasive method for assessment of hepatic fibrosis. Sandrin L, Fourquet B, Hasquenoph JM, Yon S, Fournier C, Mal F, Christidis C, Ziol M, Poulet B, Kazemi F, Beaugrand M, Palau R. *Ultrasound Med Biol*. 2003 Dec; 29(12): 1705-13
- 25- Shear-Wave Generation using acoustic radiation for in vivo and ex vivo results. Nightingale K, Mcleavy S, Trahey G. *Ultrasound Med Biol* 2003;12:1715-1723
- 26- Prospective comparison of transient elastography, Fibrotest, APRI, and liver biopsy for the assessment of fibrosis in chronic hepatitis C. Castera L, Vergniol J, Foucher J, Le Bail B, Chanteloup E, Haaser M, Darriet M, Couzigou P, De Ledinghen V. *Gastroenterology* 2005;128:343-350.
- 27- Biophysical bases of elasticity imaging. Sarvazyan AP et al. In: Jones JP, ed. *Acoustical Imaging 21*. New York: Plenum Press, 223–240, 1995.
- 28- Advanced MRI methods for assessment of chronic liver disease. Taouli B, Ehman RL, Reeder SB. *AJR Am J Roentgenol*. 2009 Jul; 193(1):14-27.
- 29- Viscoelastic shear properties of in vivo breast lesions measured by MR elastography. Sinkus R, Tanter M, Xydeas T, Catheline S, Bercoff J, Fink M. *MagnReson Imaging*. 2005 Feb; 23(2): 159-65.
- 30- Assessment of hepatic fibrosis with magnetic resonance elastography. Yin M, Talwalkar JA, Glaser KJ, Manduca A, Grimm RC, Rossman PJ, Fidler JL, Ehman RL. *Clin Gastroenterol Hepatol*. 2007 Oct; 5(10): 1207-1213.
- 31- Quantitative assessment of breast lesion viscoelasticity: Initial clinical results using supersonic shear imaging. Tanter M, Bercoff J, Athanasiou A, Deffieux T, Gennisson JL, Montaldo G, Muller M, Tardivon A, Fink M. *Ultrasound Med Biol*. 2008 Sep;34(9):1373-86.
- 32- Non Invasive assessment of liver Fibrosis in chronic Hepatitis C. Castera L. *Hepatol Int* (2011) 5:625–634.
- 33- Performance of transient elastography for the staging of liver fibrosis: a meta-analysis. Friedrich-Rust M, Ong MF, Martens S, Sarrazin C, Bojunga J, Zeuzem S, et al. *Gastroenterology* 2008;134(4):960–74.
- 34- Elastography for the diagnosis of severity of fibrosis in chronic liver disease: a meta-analysis of diagnostic accuracy. Tsochatzis EA, Gurusamy KS, Ntaoula S, Cholongitas E, Davidson BR, Burroughs AK. *J Hepatol*. 2011 Apr;54(4):650-9.
- 35- Noninvasive assessment of liver fibrosis by measurement of stiffness in patients with chronic hepatitis C. Ziol M, Handra-Luca A, Kettaneh A, Christidis C, Mal F, Kazemi F, de Ledinghen V, Marcellin P, Dhumeaux D, Trinchet JC, Beaugrand M. *Hepatology* 2005;41:48-54.
- 36- Pitfalls of liver stiffness measurement: a 5-year prospective study of 13,369 examinations. Castéra L, Foucher J, Bernard PH, Carvalho F, Allaix D, Merrouche W, Couzigou P, de Ledinghen V. *Hepatology*. 2010 Mar;51(3):828-35.
- 37- Ultrasound-based Hepatic Elastography Origins, Limitations, and Applications. Cohen EB, Afdhal NH. *J Clin Gastroenterol* 2010;44:637–645.
- 38- Breast lesions: quantitative elastography with supersonic shear imaging- preliminary results. Athanasiou A, Tardivon A, Tanter M, Sigal-Zafrani B, Bercoff J, Deffieux T, Gennisson JL, Fink M, Neuschwander S. *Radiology*. 2010 Jul; 256(1): 297-303
- 39- Shear wave elastography for breast masses is highly reproducible. Cosgrove DO, Berg WA, Doré CJ, Skyba DM, Henry JP, Gay J, Cohen-Bacrie C, the BE1 Study Group. *Eur Radiol*. 2012 May;22(5):1023-32.
- 40- Shear-wave Elastography Improves the Specificity of Breast US: The BE1 Multinational Study of 939 Masses. Berg WA, Cosgrove DO, Doré CJ, Schäfer FKW, Svensson WE, Hooley RJ, Ohlinger R, Mendelson EB, Balu-Maestro C, Locatelli M, Tourasse C, Cavanaugh BC, Juhav V, Stavros AT, Tardivon A, Gay J, Henry JP, Cohen-Bacrie C, and the BE1 Investigators. *Radiology*. 2012 Feb;262(2):435-49.
- 41- Quantitative assessment of shear-wave ultrasound elastography in thyroid nodules: diagnostic performance for predicting malignancy. Kim H, Kim JA, Son EJ, Youk JH. *Eur Radiol*. 2013 Sep;23(9):2532-7.
- 42- Quantitative shear wave elastography as a prognostic implication of papillary thyroid carcinoma (PTC): elasticity index can predict extrathyroidal extension (ETE). Park YJ, Kim JA, Son EJ, Youk JH, Park CS. *Ann Surg Oncol*. 2013 Aug;20(8):2765-71.
- 43- A threshold value in Shear Wave elastography to rule out malignant thyroid nodules: a reality? Veyrieres JB, Albarel F, Lombard JV, Berbis J, Sebag F, Oliver C, Petit P. *Eur J Radiol*. 2012 Dec;81(12):3965-72.
- 44- Shear wave elastography may add a new dimension to ultrasound evaluation of thyroid nodules: case series with comparative evaluation. Slapa RZ, Piwowonski A, Jakubowski WS, Bierca J, Szopinski KT, Slowinska-Szrednicka J, Migda B, Mlosek RK. *J Thyroid Res*. 2012;2012:657147.
- 45- Shear wave elastography: a new ultrasound imaging mode for the differential diagnosis of benign and malignant thyroid nodules. Sebag F, Vaillant-Lombard J, Berbis J, Griset V, Henry JF, Petit P, Oliver C. *J Clin Endocrinol Metab*. 2010 Dec;95(12):5281-8.
- 46- Evaluation of shearwave elastography for the characterisation of focal liver lesions on ultrasound. Guibal A, Boularan C, Bruce M, Vallin M, Pilleul F, Walter T, Scoazec JY, Boublay N, Dumortier J, Lefort T. *Eur Radiol*. 2013 Apr;23(4):1138-49.
- 47- Ultrasound elastography of the prostate: state of the art. Correas JM, Tissier AM, Khairoune A, Khoury G, Eiss D, Hélénon O. *Diagn Interv Imaging*. 2013 May;94(5):551-60.
- 48- Shear wave ultrasound elastography of the prostate: initial results. Barr RG, Memo R, Schaub CR. *Ultrasound Q*. 2012 Mar;28(1):13-20.
- 49- Supersonic shear imaging: a new technique for soft tissue elasticity mapping. Bercoff J, Tanter M, Fink M. *IEEE Trans Ultrason Ferroelectr Freq Control*. 2004 Apr; 51(4): 396-409.

United States / Latin America

250 Campus Drive
Marlborough, MA 01752 USA
Tel: +1.877.371.4372
Sales: +1.781.999.7453
Fax: +1.781.280.0668

HQ / Europe, Middle East and Africa:

+33 (0)4 42 99 24 24
North America: +1 (954) 660 3528
China: +86 10 85861023/2951/2971
contacts@supersonicimagine.com
www.hologic.com

WP-00232 Rev 001 (9/18) US/International © 2018 Hologic, Inc. All rights reserved. Printed in the USA. Specifications subject to change without notice. Hologic, Aixplorer, ShearWave, and associated logos or trademarks are registered trademarks of Hologic, Inc. and/or subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners. Views and opinions expressed herein by third parties are theirs alone and do not necessarily reflect those of Hologic.

This is a general information tool for medical professionals and is not a complete representation of the product(s)' Instruction for Use (IFU) or Package Insert, and it is the medical professionals' responsibility to read and follow the IFU or Package Insert. The information provided may suggest a particular technique or protocol however it is the sole responsibility of the medical professional to determine which technique or protocol is appropriate. At all times, clinicians remain responsible for utilizing sound patient evaluation and selection practices, and for complying with applicable local, state, and federal rules and regulations regarding accreditation, anesthesia, reimbursement, and all other aspects of in-office procedures. In no event shall Hologic be liable for damages of any kind resulting from your use of the information presented.

Indications for Use: The SuperSonic Imagine Aixplorer®, Aixplorer® Ultimate, Aixplorer MACH 30® ultrasound diagnostic systems and transducers are intended for general purpose pulse echo ultrasound imaging, Doppler fluid flow analysis of the human body, and soft tissue elasticity imaging. The Aixplorer®, Aixplorer® Ultimate, Aixplorer MACH 30® ultrasound diagnostic systems are indicated for use in the following applications, for imaging and measurement of anatomical structures: Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Vascular, Peripheral Vascular, OB-GYN, Pelvic, Pediatric, Trans-rectal, Trans-vaginal, Urology, Neonatal/Adult Cephalic and Non-invasive Cardiac. In addition, the SuperSonic Imagine Aixplorer®, Aixplorer® Ultimate, Aixplorer MACH 30® ultrasound diagnostic systems and associated transducers are intended for: measurements of abdominal anatomical structures; measurements of broadband shear wave speed, and tissue stiffness in internal structures of the liver and the spleen; measurements of brightness ratio between liver and kidney; visualization of abdominal vascularization, microvascularization and perfusion; quantification of abdominal vascularization and perfusion. The shear wave speed and stiffness measurements, the brightness ratio, the visualization of vascularization, microvascularization and perfusion, the quantification of vascularization and perfusion may be used as an aid to clinical management of adult and pediatric patients with liver disease. They are intended for use by a licensed personnel qualified to direct the use of the medical ultrasound devices. CE certificate no. 26415, FDA cleared - K180572.