
Instructions For Use

ATEC[®]

Breast Biopsy and Excision System

CE 0086

Hologic UK Ltd.
Unit 2, Link 10, Napier Way
Crawley, West Sussex, RH10 9RA, UK

ATEC® Breast Biopsy Device

Instructions for Use (IFU)

Please read all information carefully. Failure to properly follow the instructions may lead to unintended consequences.

Important: This package insert is designed to provide instructions for clinical use (IFU) for the ATEC® breast biopsy device to be used with the ATEC® breast biopsy and excision system. It is not a reference to surgical techniques.

Indications

The ATEC® Breast Biopsy and Excision System is indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities. The ATEC Breast Biopsy and Excision System is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality. The extent of histologic abnormality cannot be reliably determined from its mammographic appearance. Therefore, the extent of removal of the imaged evidence of an abnormality does not predict the extent of removal of histologic abnormality, e.g. malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedure.

Contraindications

The ATEC breast biopsy device is for diagnostic use only and is NOT indicated for therapeutic use. The ATEC breast biopsy device is contraindicated for those patients who, based on physician's judgment, may be at increased risk or develop complications associated with core removal or biopsy. Patients receiving anticoagulant therapy or who may have bleeding disorders may be at increased risk of procedural complications.

Device Description

The ATEC Breast Biopsy and Excision System is comprised of an ATEC breast biopsy device, appropriate console, and associated accessories. The ATEC breast biopsy device is a sterile, latex free, single use device intended for use with visual imaging guidance, such as ultrasound, stereotactic and magnetic resonance imaging (MRI). The ATEC breast biopsy device is a vacuum-assisted biopsy device which is used to remove breast tissue in a minimally invasive manner designed for biopsy and diagnosis only. The device is NOT intended to be used outside the standards, specifications and limitations as outlined in this IFU and other related materials. The device is single use and the entire handpiece is fully disposable. The device primarily consists of a hollow needle with a side aperture and a sharpened inner cannula that rotates at high speeds and extends across the aperture to acquire targeted tissue. Throughout this process, vacuum is created inside the device to assist in pulling tissue into the aperture. Saline is supplied through the device to lavage the cavity and deliver tissue to the tissue filter. A compatible biopsy site marker may be used in conjunction with the ATEC breast biopsy device following tissue removal or biopsy to mark the location of the biopsy site for future procedures.

Instructions for Use

Refer to the “ATEC Breast Biopsy and Excision System Operator’s Manual” for console and handpiece instructions for use.

Warnings and Precautions

- **Rx Only** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- The biopsy procedure should be performed only by persons having adequate training and familiarity with this procedure. Consult medical literature relative to techniques, complications and hazards prior to performance of any minimally invasive procedure.
- The ATEC breast biopsy device should be used only by physicians trained in percutaneous biopsy procedures.
- Sound professional judgment should be used when the ATEC breast biopsy device is used on patients with breast implants.
- A complete and comprehensive preoperative medical history and physical examination are suggested. Radiographic evaluation and laboratory tests may be included.
- The MRI ATEC breast biopsy device is intended for use only with an ATEC MRI compatible console and footswitch.
- Prior to use of the ATEC breast biopsy device, inspect the protective packaging and device to verify that neither has been damaged during shipment. If it appears that the packaging has been compromised, do not use the device.
- As with any medical procedure, please ensure that users wear appropriate personal protective equipment to guard against potential contact with bodily fluids.
- Avoid operator or instrument contact with the sheathed needle portion of the ATEC breast biopsy device.
- Minimally invasive instruments and accessories manufactured or distributed by companies not authorized by Hologic, Inc. may not be compatible with the ATEC breast biopsy device. Use of such products may lead to unanticipated results and possible injury to the user or patient.
- Instruments or devices which come into contact with bodily fluids may require special disposal handling to prevent biological contamination.
- Dispose of all opened instruments whether used or unused.
- Do not resterilize and/or reuse the ATEC breast biopsy device. Resterilization and/or reuse may compromise the integrity of the instrument. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned and sterilized devices.
- The tubing set connecting the biopsy needle and console of this device contains di-(2-ethylhexyl) phthalate; DEHP. In accordance with European Commission Directive 67/548/EEC, it is noted here that

DEHP may impair fertility; it also may cause harm to the unborn child. This device is not contraindicated for use in nursing or pregnant women or children, however residual risks may exist from use of DEHP in these populations. Sound medical judgment should be used.

How Supplied

The ATEC breast biopsy device is supplied sterile for single patient use. Discard into an appropriate container after use. The biopsy device is sterilized by gamma radiation to sterility assurance level (SAL) 10^{-6} .

For More Information

For more information on the ATEC breast biopsy device, U.S. and Canadian customers can contact the Hologic Customer Support Department at: 1-877-887-8767 or CSupport@hologic.com. Customer complaints can be directed to the Complaint Hotline at: 1-888-355-7876 or quality.reporting@hologic.com.

International customers, please contact your local distributor.

© 2013 Hologic Inc. All rights reserved. Hologic, ATEC, and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries, in the United States and/or other countries.

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