



Date: 5/23/2019

Author: Software Engineering

System: Apex

Subsystem: Microsoft

Windows OS 7

Subject: Microsoft Windows 7 OS Patch release May 2019

Purpose:

To release the list of validated Microsoft critical patches on the Apex systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0 " which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.5 – v5.6.0.5
Discovery	APEX v4.0 – v4.6.0.3
Workstation	APEX v4.0 – v5.6.0.5

List of patches/KB's validated for May 2019:

KB4499175

KB4499164

KB4498206

KB4495606

KB4495588

List of all patches/KB's validated on the Apex System:

KB976932	KB4499175	KB4499164	KB4498206
KB4495606	KB4495588	KB4493448	KB4490628
KB4490128	KB4489885	KB4486564	KB4486459
KB4483187	KB4480960	KB4474419	KB4471328
KB4467106	KB4462915	KB4457145	KB4457044
KB4457035	KB4343899	KB4338823	KB4284867
KB4103712	KB4096418	KB4093108	KB4074587
KB4056897	KB4054521	KB4048960	KB4041678
KB4019990	KB3177467	KB3172605	KB3161949
KB3159398	KB3156016	KB3150220	KB3139914
KB3139398	KB3138910	KB3138612	KB3127220
KB3126587	KB3124275	KB3122648	KB3110329
KB3109560	KB3109103	KB3108664	KB3108371
KB3101722	KB3097989	KB3093513	KB3086255
KB3078601	KB3075226	KB3074543	KB3072305
KB3071756	KB3067903	KB3061518	KB3060716
KB3059317	KB3055642	KB3046269	KB3046017
KB3045685	KB3042058	KB3037574	KB3035126
KB3031432	KB3030377	KB3023215	KB3022777
KB3021674	KB3020388	KB3011780	KB3010788
KB3004375	KB2992611	KB2991963	KB2984972
KB2978742	KB2977292	KB2973201	KB2973112
KB2972100	KB2968294	KB2943357	KB2937610
KB2931356	KB2923545	KB2912390	KB2911501
KB2894854	KB2894844	KB2893294	KB2884256
KB2871997	KB2868038	KB2864202	KB2862335
KB2862330	KB2862152	KB2861698	KB2857650
KB2852386	KB2847927	KB2840631	KB2834140
KB2830477	KB2813430	KB2813347	KB2807986
KB2786081	KB2770660	KB2742599	KB2736422
KB2729094	KB2727528	KB2705219	KB2698365

KB2690533	KB2670838	KB2667402	KB2654428
KB2639308	KB2631813	KB2621440	KB2620704
KB2604115	KB2579686	KB2574819	KB2564958
KB2560656	KB2552343	KB2533552	KB2532531
KB2506212	KB2491683	KB2479943	