

# Marker Compatibility for MRI Breast Biopsy



Biopsy Device	Catalog Number (SKU)	Gauge	Length	Shape name or Shadow?	Metal	Bio-absorbable Carrier	Deployment Method	# per box
ATEC MRI	Smark-M-ss1	12G	13cm	Buckle	SS	Yes	Introducer	10
	Smark-M-ss2			Infinity	SS	Yes	Introducer	10
	Smark-M-ss3			Stoplight	SS	Yes	Introducer	10
	TriMark TD 13MR	9G		Cork	Ti	No	Introducer	10
	TriMarkTD-2S-13-MR			Hourglass	Ti	No	Introducer	10

<p><b>SecurMark® Marker</b></p> <p>The bio-absorbable suture-like net that carries a bio-compatible permanent marker, is designed to <b>minimize movement</b> within the biopsy cavity and <b>enhance ultrasound visualization</b> for three to four weeks post-biopsy.<sup>1</sup></p>	<p><b>TriMark® Marker</b></p> <p><b>Reliable</b> and <b>simple</b> bio-compatible <b>Titanium metal markers</b>, optimal for thin-breasted patients and superficial lesions.</p>

<sup>1</sup>. Based on Clinical Evaluation Report DHM-08585. Visibility is dependent on surrounding tissue, experience may vary.

[www.hologic.com](http://www.hologic.com) | [BSHSalesSupportUS@Hologic.com](mailto:BSHSalesSupportUS@Hologic.com) | + 1.800.442.9892

**Hologic Headquarters**

**United States / Latin America**

250 Campus Drive

Marlborough, MA 01752

**USA Tel:** +1.508.263.2900

**Sales:** +1.781.999.7453

**Fax:** +1.781.280.0668

**Email:** [info@hologic.com](mailto:info@hologic.com)

**Hologic Europe**

Everest (Cross Point)

Leuvensesteenweg

250A 1800

Vilvoorde, Belgium

**Tel:** +32.2.711.4680

**Fax:** +32.2.725.2087

**Hologic Asia Pacific Limited**

Unit 01-03A, 13/F

No. 909 Cheung Sha Wan

Road Cheung Sha Wan

Kowloon, Hong Kong

**Tel:** +852 3748 7700

**Fax:** +852 3526 0723

**Hologic Australia**

Suite 402, Level 4

2 Lyon Park Road

Macquarie Park NSW 2113

Australia

**Tel:** +61.2.9888.8000

**Fax:** +61.2.9870.7555