

| | | |
|--|---|---|
| IHE Integration Statement | | Date |
| | | 21 April 2011 |
| Vendor | Product Name | Version |
| Hologic, Inc. | Physician's Viewer | 7.0 |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: | | |
| Integration Profiles Implemented | Actors Implemented | Options Implemented |
| Scheduled Workflow | Image Display | Query / Retrieve Images |
| | | Image Storage |
| Internet address for the vendor's IHE information: http://www.hologic.com/cc/iheinteg.htm | | |
| Links to Standards Conformance Statements for the Implementation | | |
| HL7 | http://www.hologic.com/cc/dicom.htm | |
| DICOM | http://www.hologic.com/cc/dicom.htm | |
| Links to general information on IHE | | |
| In North America: www.rsna.org/ihe | In Europe: www.ihe-europe.org | In Japan: www.jira-net.or.jp/ihe-j |