

<b>1 IHE Integration Statement</b>		<b>1.1 Date</b>
		10 April 2012
<b>1.2 Vendor</b>	<b>1.3 Product Name</b>	<b>1.4 Version</b>
Hologic, Inc.	Fluoroscan InSight	5.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>1.5 Options Implemented</b>
Scheduled Workflow	Acquisition Modality	Modality Image Storage
		Storage Commitment
		Modality Worklist Query
		Modality Performed Procedure Step
<b>Internet address for the vendor's IHE information:</b> <a href="http://www.hologic.com/en/product-support/mini-carm/insight/">http://www.hologic.com/en/product-support/mini-carm/insight/</a>		
<b>1.6 Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="http://www.hologic.com/en/product-support/mini-carm/insight/">http://www.hologic.com/en/product-support/mini-carm/insight/</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://rsna.org/ihe.aspx">http://rsna.org/ihe.aspx</a>	In Europe: <a href="http://www.ihe.net/">http://www.ihe.net/</a>	In Japan: <a href="http://www.ihe-j.org/">http://www.ihe-j.org/</a>