

IHE Integration Statement for SecurView® DX

IHE Integration Statement			
Vendor	Product	Version	Date
Hologic, Inc.	SecurView DX Diagnostic Workstation	6-x through 10.x	October 2016
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display	None	
Access to Radiology Information	Image Display	None	
Mammography Image	Image Display	None	
Mammography Image	Print Composer	None	
Consistent Presentation of Images	Evidence Creator	None	
Portable Data for Imaging	Portable Media Creator	None	
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	http://www.hologic.com/support/securview-diagnostic-workstations		
Links for general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	