

Introduction

This document modifies the DICOM Printer Quality Control test in the Selenia® digital mammography system Quality Control Manual. The addendum clarifies that the DICOM Printer Quality Control test is not required if the facility does not use a printer.

The following change is applied to MAN-01476 *Selenia Quality Control Manual*, Revision 002. Please place this document with your product manual for future reference.

Addendum

Chapter 3, section 1.0 *DICOM Printer Quality Control, Frequency* is amended to include the following statement:



Note

If the facility does not own or use a printer, the DICOM Printer Quality Control test is not required. If the facility is using the printer sporadically, the facility must perform this test prior to using the printer to print clinical images.
