

1.1 Introduction

This addendum modifies your Selenia® Dimensions® User and Service Manuals for an added length Radiation Shield. The other measurements are not changed. To see the changed length, refer to item A in the following figures.

Place this addendum with your manuals for future reference.

1.1.1 Premium and Standard Workstations

Premium Acquisition Workstation




Figure 1: Premium Acquisition Workstation Dimensions system

A.	Height	222.1 cm (87.45 inches)
B.	Width	92.7 cm (36.5 inches)
C.	Depth	58.5 cm (23 inches)
	Weight	154 kg (340 pounds)

Standard Acquisition Workstation




Figure 2: Standard Acquisition Workstation Dimensions system

A.	Height	217.2 cm (85.51 inches)
B.	Width	107 cm (42 inches)
C.	Depth	76 cm (30 inches)
	Weight	219 kg (462 pounds)