

IMPORTANT: PLEASE READ

Recommended DVD R/W Media

It is important to use only good quality media in the DVD R/W drive supplied with your Hologic system. Only brands recommended by the manufacturer of the DVD R/W drive should be used.

Proper Handling of DVD R/W Media

Always read and follow the handling precautions shown on the packaging for your new DVD R/W media.

To protect the recording surface of the disc from scratches, fingerprints, dust, and smears, never touch the recording surface or place it face down on a hard surface. Hold the disc by the outer edges or by placing a finger into the hole.

Handle the disc carefully so that fingerprints, scratches, or dust will not get on the recording surface.

DO NOT stick labels or protective sheets onto the disc, or coat the disc surface. If tape or a label inadvertently gets put on the disc, do NOT try to remove it, as this may damage the disc. Use your backup media to make another copy of the information. If this is the only copy of the data, with the label/tape still attached, attempt to copy the information from the disc to your hard drive and then to another DVD R/W disc.

Use a felt-tip pen to write on the disc label area. Do not use hard-tipped writing instruments such as ballpoint pens.

Store the disc in a location that is not exposed to direct heat or sunlight, and has an ambient temperature of between -5°C (23°F) and 55°C (131°F) and a humidity level of 10% to 90% RH, with no condensation.