



IHE Integration Statement



SuperSonic[®] MACH 40 Ultrasound System V3.x software

Revision: A
Date: November 02, 2020

Doc Ref: RD.DD.530 (RevA)
Software Version: V3.x

Introduction

Integrating the Healthcare Enterprise® (IHE) is an initiative by healthcare professionals and industry to improve the way computer systems in healthcare share information. Systems that support IHE Integration Profiles work together better, are easier to implement, and help care providers use information more effectively. The goal is efficient delivery of optimal patient care.

IHE Integration Statements are documents prepared and published by vendors to describe the intended conformance of their products with the IHE Technical Framework. Integration Statements tell consumers the specific IHE capabilities a given product is designed to support in terms of the key concepts of IHE: Actors and Integration Profiles.

IHE Integration Statement

IHE Integration Statement			
Vendor	Product Name	Date	November 02, 2020
Hologic SuperSonic Imagine	SuperSonic MACH 40	Version	V3.x
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles	Actors	Options	
Radiology domain			
Scheduled Workflow Scheduled Workflow.B	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
Patient Information Reconciliation	Acquisition Modality		
Evidence Documents	Acquisition Modality		
Portable Data for Imaging	Portable Media Creator	Web Content	
Cardiology domain			
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
Evidence Documents	Acquisition Modality	Echo Evidence	
Image-Enabled Office	Acquisition Modality	Broad Worklist Query Patient Based Worklist Query PPS Exception Management	
IT infrastructure domain			
Consistent Time	Time Client		
Internet address for vendor's IHE information: http://www.supersonicimagine.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	Not Applicable		
DICOM	http://www.supersonicimagine.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org	