

Marker Compatibility for Stereotactic-Guided Breast Biopsy

Biopsy Device	Marker Name	Shape	Metal	Bio-absorbable Material (PLA/PGA)	Deployment Method
ATEC® Device	SecurMark® for ATEC® Biopsy Site Marker				
	SMark-ATEC-36-09	"mini cork"	Titanium	Yes	Handpiece (long shaft)
	SMark-ATEC-2s-36-09	"top hat"	Titanium	Yes	Handpiece (long shaft)
	SMark-A-36-09-ss1	"buckle"	Stainless Steel	Yes	Handpiece (long shaft)
	SMark-ATEC-36-12	"mini cork"	Titanium	Yes	Handpiece (long shaft)
	SMark-ATEC-13-09	"mini cork"	Titanium	Yes	Needle Hub (short shaft)
	SMark-ATEC-2s-13-09	"top hat"	Titanium	Yes	Needle Hub (short shaft)
	SMark-ATEC-13-12	"mini cork"	Titanium	Yes	Needle Hub (short shaft)
	SMark-ATEC-2s-13-12	"top hat"	Titanium	Yes	Needle Hub (short shaft)
	SMark-A-13-12-ss2	"infinity"	Stainless Steel	Yes	Needle Hub (short shaft)
	TriMark® Biopsy Site Marker				
	TriMarkTD-3609	"cork"	Titanium	No	Handpiece (long shaft)
	TriMarkTD-2s-3609	"hourglass"	Titanium	No	Handpiece (long shaft)
	TriMarkTD-3612	"cork"	Titanium	No	Handpiece (long shaft)
	TriMarkTD-2s-3612	"hourglass"	Titanium	No	Handpiece (long shaft)
	TriMarkTD-1309	"cork"	Titanium	No	Needle Hub (short shaft)
	TriMarkTD-2s-1309	"hourglass"	Titanium	No	Needle Hub (short shaft)
	TriMarkTD-1312	"cork"	Titanium	No	Needle Hub (short shaft)
	TriMarkTD-2s-1312	"hourglass"	Titanium	No	Needle Hub (short shaft)
	Tumark® Flex Marker for ATEC Device				
351250 (12G)	"U"	Nitinol	No	Handpiece (long shaft)	
351252 (9G)	"U"	Nitinol	No	Handpiece (long shaft)	
Brevera® and Eviva® Devices	SecurMark® for Eviva® Biopsy Site Marker				
	SMark-Eviva-13	"mini cork"	Titanium	Yes	Introducer
	SMark-Eviva-2s-13	"top hat"	Titanium	Yes	Introducer
	SMark-E-13-ss1	"buckle"	Stainless Steel	Yes	Introducer
	SMark-E-13-ss2	"infinity"	Stainless Steel	Yes	Introducer
	SMark-E-13-ss3	"stoplight"	Stainless Steel	Yes	Introducer
	SMark-Eviva-10	"mini cork"	Titanium	Yes	Introducer
	SMark-Eviva-2s-10	"top hat"	Titanium	Yes	Introducer
	TriMark® for Eviva® Biopsy Site Marker				
	TriMark-Eviva-10	"cork"	Titanium	No	Introducer
	TriMark-Eviva-2s-10	"hourglass"	Titanium	No	Introducer
	TriMark-Eviva-13	"cork"	Titanium	No	Introducer
TriMark-Eviva-2s-13	"hourglass"	Titanium	No	Introducer	

Note: Not all markers are available worldwide.

Hologic.com | info@hologic.com | +1.866.636.8837 | +1.508.263.2795 | +1.781.999.7300

Breast and Skeletal Health

MISC-01695 Rev. 008 (8/17) © 2017 Hologic Inc. Printed in USA. All rights reserved. Specifications are subject to change without prior notice. Hologic, ATEC, Brevera, Eviva, SecurMark, The Science of Sure, TriMark and associated logos are trademarks or registered trademarks of Hologic and/or Hologic subsidiaries in the United States and/or other countries. Tumark and associated logos are trademarks and/or registered trademarks of Somatex Medical Technologies GmbH. Hologic is an authorized distributor of Tumark Professional X, Q, Vision and Flex for ATEC biopsy site markers. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to womenshealth@hologic.com.