



Date: 10/5/2018

Author: Software Engineering

System: Apex

Subsystem: Microsoft

Windows OS

Subject: Microsoft Windows OS Patch release September 2018

Purpose:

To release the list of validated Microsoft critical patches on the Apex systems.

Note: If you are using the following versions of APEX, you must ensure your system(s) have been patched with the "Sentinel Run-Time 7.80.0.0 " which can be located in the Windows Control Panel "Program and Features". If needed you can install the patch by ordering "APEX Service Dongle Patch Release" (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.5 – v5.6.0.5
Discovery	APEX v4.0 – v4.6.0.3
Workstation	APEX v4.0 – v5.6.0.5

List of patches/KB's validated for September 2018:

KB4457426

KB4457145

KB4457144

KB4457044

KB4457035

KB4343899

KB4339284

List of all patches/KB's validated on the Apex System:

KB976932	KB4457426	KB4457145	KB4457144
KB4457044	KB4457035	KB4343899	KB4339284
KB4338823	KB4284867	KB4130978	KB4103712
KB4096418	KB4093108	KB4074587	KB4056897
KB4054521	KB4049068	KB4048960	KB4041678
KB4019990	KB3177467	KB3172605	KB3161949
KB3159398	KB3156016	KB3150220	KB3139914
KB3139398	KB3138910	KB3138612	KB3127220
KB3126587	KB3124275	KB3122648	KB3110329
KB3109560	KB3109103	KB3108664	KB3108371
KB3101722	KB3097989	KB3093513	KB3086255
KB3078601	KB3075226	KB3074543	KB3072305
KB3071756	KB3067903	KB3061518	KB3060716
KB3059317	KB3055642	KB3046269	KB3046017
KB3045685	KB3042058	KB3037574	KB3035126
KB3031432	KB3030377	KB3023215	KB3022777
KB3021674	KB3020388	KB3011780	KB3010788
KB3004375	KB3003743	KB2992611	KB2991963
KB2984972	KB2978742	KB2977292	KB2973201
KB2973112	KB2972100	KB2968294	KB2943357
KB2937610	KB2931356	KB2923545	KB2912390
KB2911501	KB2894854	KB2894844	KB2893294
KB2892074	KB2884256	KB2871997	KB2868038
KB2864202	KB2862335	KB2862330	KB2862152
KB2861698	KB2857650	KB2852386	KB2847927
KB2840631	KB2834140	KB2830477	KB2813430
KB2813347	KB2807986	KB2786081	KB2770660
KB2742599	KB2736422	KB2729094	KB2727528

KB2705219

KB2698365

KB2690533

KB2670838

KB2667402

KB2654428

KB2639308

KB2631813

KB2621440

KB2620704

KB2604115

KB2579686

KB2574819

KB2564958

KB2560656

KB2552343

KB2533552

KB2532531

KB2506212

KB2491683

KB2479943