

Introducing the NEW 3 in 1 Omni™ Hysteroscope

Taking SEE & TREAT to the next level!

3 HYSTEROSCOPES IN 1

Convenience to SEE & TREAT pathology with a single hysteroscope

ROD LENS


Quality visualization throughout procedures

LONG 200 mm WORKING LENGTH

Facilitates easier access and treatment for larger patients

SMALL DIAMETER

Easier insertion and greater patient comfort



Diagnostic through to operative - all in ONE kit


Base Scope

3.7 mm Diagnostic Sheath

5.5 mm Sheath and Outflow Channel

6 mm Sheath and Outflow Channel

Light Guide Adapters



MyoSure
Tissue Removal Device

THE OMNI HYSTEROSCOPE MYOSURE COMPATIBILITY	MYOSURE MANUAL DEVICE	MYOSURE LITE DEVICE	MYOSURE REACH DEVICE	MYOSURE XL DEVICE
5.5 mm Sheath	✓	✓	✓	
6 mm Sheath	✓	✓	✓	✓

ORDERING

Description	Order Number
The Omni™ Hysteroscope Standard Kit	60-250-1
The Omni™ Hysteroscope Instrument Tray	60-903
Rod Lens Scope Seals – 10 per pack	40-902

Contact a Hologic sales representative to learn more.

For detailed benefit and risk information, please refer to the IFU.

IMPORTANT SAFETY INFORMATION

The Omni™ hysteroscope is intended to provide viewing of the cervical canal and uterine cavity for the purpose of performing diagnostic and operative procedures. The Omni hysteroscope is not appropriate for patients who are exhibiting acute pelvic inflammatory disease or exhibit the following conditions: inability to distend uterus, cervical stenosis, cervical/vaginal infection, uterine bleeding or menses, known pregnancy, invasive carcinoma of the cervix, recent uterine perforation, medical contraindications, or intolerance to anesthesia.

SSL-00352-CAN-EN Rev.001 ©2019 Hologic, Inc. Hologic, MyoSure, Omni, The Science of Sure, and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries.

Omni
Hysteroscope