



Radiopaque tab to easily position the bag using a grasper


Made of ripstop nylon cells used in parachutes and hot air balloons, boosting strength and stopping runaway tears

High contrast stitching and color markers for an easy-to-see and easy-to-orientate opening

Long monofilament drawstring for easy opening and keeping open, with extended length for external control of bag closure

Internally stitched seams can withstand great force

Polyurethane coating fully contains biological fluids



EcoSac

Tissue Retrieval System

Comprehensive laparoscopic solutions from node biopsy to an enlarged spleen

A unique design that eliminates the need for a plastic introducer and applicator, reducing shelf space requirements and eliminating plastic waste

A material difference in retrieval bags: EcoSac Tissue Retrieval System

USER FRIENDLY

Seamlessly reopen and reuse throughout a procedure, with a smooth drawstring closure and nylon that becomes transparent when wet for easy viewing of the specimen.

ECO-AND-ECONOMICALLY FRIENDLY

Reduce packaging, storage space, waste, and costs without the need for a plastic introducer and applicator.



Ordering Information

Product, Volume (mL)	Product number	QTY	Diameter	Length	Procedure
5 MM TROCAR					
EcoSac, 60 mL	EMP60ECO-R-10	10/pack	70 mm	100 mm	Node biopsy, appendix, ectopic, small gallbladder; also recommended for pediatric use
EcoSac, 80 mL	EMP80ECO-R-10	10/pack	70 mm	125 mm	<i>Suitable for robotic procedures</i>
EcoSac, 200 mL	EMP200ECO-R-10	10/pack	70 mm	230 mm	Ectopic, oophorectomy, endometriosis, salpingectomy, biopsy, small gallbladder, adrenalectomy, nephrectomy, small fibroids <i>Suitable for robotic procedures</i>
10 MM TROCAR					
EcoSac, 100 mL	EMP70ECO-US-10	10/pack	57 mm	100 mm	Node biopsy, appendix, ectopic, small gallbladder, ovary, fibroids, adrenal glands
EcoSac, 200 mL	EMP100ECO-US-10	10/pack	70 mm	104 mm	Gallbladder, ovary, fibroids, prostate, adrenal glands
12 MM TROCAR					
EcoSac, 550 mL	EMP140ECO-US-10	10/pack	115 mm	185 mm	Large fibroid, uterus, kidney, spleen, lung, colon, thoracic
EcoSac, 1,000 mL	EMP160ECO-US-5	5/pack	120 mm	240 mm	Malignant kidney, enlarged spleen
EcoSac, 1,200 mL	EMP170ECO-US-5	5/pack	170 mm	180 mm	Large tumors, kidney, colon, lung, liver
EcoSac, 2,000 mL	EMP180ECO-US-5	5/pack	170 mm	240 mm	Colon tumors, large spleen, lung, liver, large fibroid, uterus
EcoSac, 3,100 mL	EMP230ECO-US-5	5/pack	165 mm	340 mm	Gynecological procedures, thoracic and general surgery procedures (e.g., large resections and tumors)
15 MM TROCAR					
EcoSac, 6,000 mL	EMP600ECO-5	5/pack	320 mm	410 mm	Large/very large fibroids, myomectomy, large spleen

Contact your Hologic sales representative for more information.

All Rights Reserved. Hologic and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. EcoSac is manufactured by Espiner Medical and distributed exclusively in the United States by Hologic.

Hologic, Inc. | 250 Campus Dr., Marlborough, MA 01752 | 508 263 2900

MISC-09899 Rev. 001 | 07/2024 | © 2024 Hologic, Inc.



Hologic.com