

# More comfortable mammograms. More satisfied patients.

Hologic addresses one of the main reasons why women might avoid having a potentially life-saving mammogram: discomfort.<sup>1</sup> Our SmartCurve™ breast stabilization system delivers an exceptional patient experience without compromising image quality, workflow or dose. The system includes three paddles to accommodate the majority of women and breast sizes.



The SmartCurve™ system is clinically proven to deliver a more comfortable mammogram without increasing dose.<sup>2</sup>



Curved compression surface mirrors the shape of a woman's breast.




95% of women surveyed would recommend facilities that use the SmartCurve™ system.<sup>2</sup>



# The SmartCurve™ system will change the way women look at breast screening.

The fear of pain stops many women from making regular breast screening a healthcare priority, potentially delaying life-saving cancer detection.

The SmartCurve™ system was designed to address this challenge, transforming the breast screening experience for you and your patients. Now, with the addition of the SmartCurve™ Mini paddle it will be easier to accommodate small-breast patients, including men, and patients with displaced implants.



The 3 pack comes with a large (24x29 cm), small (18x24 cm), and mini (18x8 cm) paddle\*\*.

- Compression surface conforms to the shape of the breast and is curved at the chest wall and sides to reduce pinching.
- Shown to **improve comfort in 93% of patients** who reported moderate to severe discomfort with standard compression paddles.<sup>2</sup>
- **Processing software** takes paddle geometry into account, so there's no impact on image quality.<sup>2</sup>
- There's **no change in positioning technique or workflow** for technologists.
- Can be used with a **MammoPad® breast cushion** for additional cushioning and increased tissue capture.<sup>3</sup>

**Proof. Not Promises.**

## Product information

The SmartCurve™ system is intended for screening and diagnostic mammography and should not be used for biopsy procedures or for contrast enhanced (I-View™) procedures. It is not recommended for cleavage views, rolled views, or mosaic views of very large breasts. The standard screening paddles are always included with the system and should be used for cases where the SmartCurve™ system is not appropriate. SmartCurve™ is compatible with breast implants.

The SmartCurve™ system comes standard with all 3Dimensions™ purchases, and is available to existing Selenia® Dimensions® systems as an upgrade.

\*\*Requires minimum Dimension 1.11/2.2 software.

## Ordering details

Part Number	Description
DIM-ASY-SCS-3PACK	SmartCurve™ Breast Stabilization system 3-pack upgrade option for Selenia® Dimensions® systems.
DIM-ASY-MINI	SmartCurve™ Mini small breast paddle only.

## References

1. Padoan M, Ferrante D, Pretti G, Magnani C. Study of socio-economic characteristics, diagnosis and outcome of women participating or not participating in mammogram screening. *Ann Ig.* 2014;26(6):518–526. 2. Data On File: DHM-06039 Rev 002 Bartoshevich J, Orefice T, Mansoor S, et al Internal Study comparing Hologic's flat paddle to the SmartCurve paddle (18x24cm) (2017). 3. Coryell T. Increasing Mammography Tissue Acquisition through Positioning Training and Use of a Foam Breast Cushion. National Consortium of Breast Centers, 16th Annual National Interdisciplinary Breast Conference, Las Vegas, NV (2006)

SS-00523-EUR-EN Rev.005 (03/21) Hologic Inc. ©2021 All rights reserved. Hologic, Dimensions, MammoPad, Selenia, SmartCurve, The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

[www.3dimensionsmammography.eu](http://www.3dimensionsmammography.eu) | [euinfo@hologic.com](mailto:euinfo@hologic.com) | +32.2.711.4680

**No Compromise. No Comparison.**