



KNOW YOUR IMAGE

FLUOROSCAN[®]
InSight-FD




FULLY AUTOMATIC

AUTO ULTRA
.1505 cGy.cm²



PULSE HIGH AUTO KF 16
.14785 cGy.cm²




SIMILAR POWER SETTINGS

ULTRA AUTO
kV 60 mA 52
.1505 cGy.cm²



MANUAL
kV60 mA 52
.12013 cGy.cm²



PEDAL HELD UNTIL IMAGE STABILIZES

ULTRA AUTO
kV 52 mA 60
.1505 cGy.cm²



PEDAL HELD DOWN KF 16
kV 60 mA 74
.65014 cGy.cm²



FD LOW ULTRA VS AUTO

LOW DOSE AUTO ULTRA
52kV .024mA
.1083 cGy.cm²



LOW DOSE AUTO
52kV .024mA
.0658 cGy.cm²



LOW DOSE VS PULSE

LOW DOSE ULTRA
.1083 cGy.cm²



K16 HIGH PULSE AUTO
.14785 cGy.cm²




K16 7.5pps .08844 cGy.cm²



FIELD OF VIEW COMPARISON

- ORTHOSCAN FD pulse²
15cm x 15cm
- FLUOROSCAN InSight FD¹
12cm x 15cm
- GE OEC Elite MiniView³
12.7 cm circle



Hologic is committed to addressing the continuum of skeletal health care, from the OR to in-office imaging. Just as we are leading the way in women's health imaging, we are raising the bar in skeletal health imaging.

www.hologic.com | info@hologic.com | +1.800.442.9892

1. MAN-05309-001 Rev 005 InSight FD technical specifications manual **2.** Orthoscan FD pulse specifications retrieved from <https://www.orthoscan.com/products/surgical/fd-pulse-mini-c-arm/> on Mar 28 2019 **3.** GE OEC Elite MiniView Specifications retrieved from 510 (k) K160131 on Mar 29 2018 <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm?id=K160131> **4.** Images obtained from Hologic InSight FD serial number 30-0118-12M, manufactured 01-2018 and OrthoScan FD Pulse serial number 5J0291 manufactured 10-2014.

PB-00742 Rev 001 (03/19) © 2019 Hologic Inc. All rights reserved. Hologic, FluorSCAN, MegaView, The Science of Sure, and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.

FLUOROSCAN[®]
InSight-FD