

TruNode[®] single-use gamma probe

With the TruNode[®] system, the single-use, wireless gamma probe from Hologic, you know with certainty.

The TruNode[®] system is the first high performance, sterile, single-use, wireless gamma probe that has been optimized for the most common radio-guided surgical procedures. The TruNode[®] system is effective, in applications such as sentinel lymph node mapping and radioactive seed localization.



TruNode[®] system features:

Industry-Leading Resolution

The TruNode[®] system features a 10% energy resolution, enabling the user to detect faint nodes closer to the injection site and facilitating precise hotspot localization.

Fingertip Control

Toggle between three modes of detection (high-sensitivity Scan, high-resolution Point and I-125 Seed) with a simple to use, fingertip controls from the sterile field.

Three Detection Modes

With its three detection modes, the TruNode[®] system optimizes sentinel node, parathyroid and tumor localization procedures.

Minimal Incision Size

With a 10mm operative diameter, the TruNode[®] system is designed to be minimally invasive and easy to maneuver.

Sterile

The TruNode[®] probe is fully sterile, eliminating drapes and reprocessing in the OR.

Voice Reporting

The TruNode[®] system's compact User Feedback Unit provides both audible and visual user feedback, reducing the need to monitor a display.

Simple Placement

The TruNode[®] system has a minimal footprint as the User Feedback Unit clips to any IV pole in the OR.

Single-Patient Use

The TruNode[®] system was designed to reduce the risk of cross-contamination and surgical site infections.

TruNode[®]
Single-use Gamma Probe

A comprehensive family of solutions

The TruNode[®] single-use gamma probe fits right into Hologic's full range of breast conserving solutions. Hologic has a complete set of products for breast conserving surgery to help you improve efficiency, optimize the patient experience and deliver the best outcomes possible. Talk to your representative about the systems that can best suit your and your patients' needs, ensuring your patients have the best possible care when they need it most.

BioZorb[®] 3D bioabsorbable marker



LOCalizer[®] wire-free guidance system



Be sure to talk with your Hologic representative about which breast conserving solution is right for you and your patients.

www.hologic.com | BSHInsideSales@hologic.com | + 1.800.442.9892

**Hologic Headquarters
United States / Latin America**
250 Campus Drive
Marlborough, MA 01752
USA Tel: +1.508.263.2900
Sales: +1.781.999.7453
Fax: +1.781.280.0668
Email: info@hologic.com

Hologic Europe
Everest (Cross Point)
Leuvensesteenweg
250A 1800
Vilvoorde, Belgium
Tel: +32.2.711.4680
Fax: +32.2.725.2087

Hologic Asia Pacific Limited
Unit 01-03A, 13/F
No. 909 Cheung Sha Wan Road
Cheung Sha Wan
Kowloon, Hong Kong
Tel: +852 3748 7700
Fax: +852 3526 0723

Hologic Australia
Suite 402, Level 4
2 Lyon Park Road
Macquarie Park NSW 2113
Australia
Tel: +61.2.9888.8000
Fax: +61.2.9870.7555

MISC-05510 Rev.001 (1/19) ©2019 Hologic Inc. All rights reserved. Hologic, BioZorb, TruNode, The Science of Sure and associated logos are trademarks and/or registered trademarks of Hologic, Inc., and/or its subsidiaries in the United States and/or other countries. Hologic is an exclusive distributor and licensee of the LOCalizer product and trademark. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative.