

**Date:** May 21, 2024

**Author:** Software Engineering

**System:** TRIDENT HD

**Subsystem:** Windows 10 & Software Security Patches

**Subject:** Windows 10 & Software Security Patch release April 2024

---

**Purpose:**

To release the list of validated Microsoft critical patches on the TRIDENT HD System.

**List of patches/KB's validated for April 2024:**

KB5037016                      KB5036899

**List of all patches/KB's validated on the TRIDENT HD System:**

KB890830	KB5027219	KB5016058
KB5037016	KB5026363	KB5015808
KB5036899	KB5025228	KB5014702
KB5035962	KB5023788	KB5014553
KB5035855	KB5023697	KB5014026
KB5034862	KB5022838	KB5013952
KB5034767	KB5022289	KB5012596
KB5034119	KB5021235	KB5012170
KB5033373	KB5021126	KB5011570
KB5032391	KB5021045	KB5011495
KB5032197	KB5019964	KB5010359
KB5031362	KB5018411	KB5009546
KB5030504	KB5017396	KB5008207
KB5030213	KB5017305	KB5007192
KB5029376	KB5017095	KB5006669
KB5029242	KB5016884	KB5005698
KB5028169	KB5016622	KB5005573

KB5005043	KB4576750	KB4512574
KB5004238	KB4571694	KB4512517
KB5003638	KB4565912	KB4509091
KB5003197	KB4565511	KB4507460
KB5002464	KB4562561	KB4505225
KB5002406	KB4561616	KB4498947
KB5001402	KB4557397	KB4494440
KB5001347	KB4556813	KB4494175
KB5001078	KB4550994	KB4493470
KB5000803	KB4550929	KB4489882
KB5000685	KB4540723	KB4487026
KB4601392	KB4540670	KB4485447
KB4601318	KB4537764	KB4346087
KB4598243	KB4535680	KB4132216
KB4593226	KB4534271	KB4091664
KB4589210	KB4530689	KB4022619
KB4586830	KB4525236	KB4018073
KB4583465	KB4521858	KB2565063
KB4580346	KB4520724	KB2467173
KB4577586	KB4519998	
KB4577015	KB4516044	