



Date: *September 15, 2021*

Author: Software Engineering

System: BREVERA / CORLUMINA

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release August 2021

Purpose:

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

List of patches/KB's validated for August 2021:

KB5005112 KB5005030

List of all patches/KB's validated on the BREVERA / CORLUMINA System:

KB890830	KB5000822	KB4577668
KB5005112	KB5000685	KB4577667
KB5005030	KB4601393	KB4577586
KB5004244	KB4601345	KB4570333
KB5003711	KB4601060	KB4570332
KB5003646	KB4598480	KB4561600
KB5003243	KB4598230	KB4557397
KB5003171	KB4592440	KB4552924
KB5001407	KB4589208	KB4539571
KB5001404	KB4587735	KB4535680
KB5001342	KB4586793	KB4494174
KB5000859	KB4580325	