



Date: September 12, 2024
Author: Software Engineering
System: Apex
Subsystem: Microsoft Win Pro Embedded 10 LTSP 2016
Subject: Microsoft Win Pro Embedded 10 LTSP 2016 Patch release August 2024

Purpose:

To release the list of validated Microsoft critical patches on the Apex Win Pro Embedded 10 LTSP 2016 systems.

Note: If you are using the following versions of APEX, if you are required to update to the latest Hologic validated Windows updates, you must ensure your system(s) have been patched with the “Sentinel Run-Time 7.80.0.0” which can be located in the Windows Control Panel “Program and Features”. If needed you can install the patch by ordering “APEX Service Dongle Patch Release” (ASY-12610). Please contact Hologic Support if you have any questions.

System Type	Software Versions
Horizon	APEX v5.6.0.7
Workstation	APEX v5.6.0.7

List of patches/KB’s validated for August 2024:

KB5041773 KB5041576

List of all patches/KB’s validated on the Apex System:

KB890830	KB5035962	KB5032197
KB5041773	KB5035855	KB5032196
KB5041576	KB5034862	KB5031984
KB5040562	KB5034767	KB5031362
KB5040434	KB5034127	KB5031361
KB5039334	KB5034119	KB5030504
KB5039214	KB5033904	KB5030214
KB5037763	KB5033373	KB5030213
KB5037016	KB5033371	KB5029931
KB5036899	KB5032391	KB5029247

KB5029242	KB5017305	KB5006672
KB5029185	KB5017095	KB5006669
KB5028169	KB5016884	KB5005698
KB5028168	KB5016623	KB5005573
KB5027222	KB5016622	KB5005568
KB5027219	KB5016058	KB5005112
KB5027131	KB5015811	KB5005043
KB5026363	KB5015808	KB5005030
KB5026362	KB5014702	KB5004244
KB5025229	KB5014692	KB5004238
KB5025228	KB5014553	KB5003711
KB5023788	KB5014026	KB5003646
KB5023702	KB5013952	KB5003638
KB5023697	KB5013941	KB5003243
KB5022840	KB5013641	KB5003197
KB5022838	KB5012647	KB5003171
KB5022511	KB5012596	KB5002542
KB5022289	KB5012170	KB5002520
KB5022286	KB5012128	KB5002497
KB5021237	KB5011570	KB5002464
KB5021235	KB5011503	KB5002406
KB5021126	KB5011495	KB5001407
KB5021045	KB5010359	KB5001404
KB5020866	KB5010351	KB5001402
KB5020627	KB5009557	KB5001347
KB5019966	KB5009546	KB5001342
KB5019964	KB5008873	KB5001078
KB5018419	KB5008218	KB5000859
KB5018411	KB5008207	KB5000822
KB5017396	KB5007206	KB5000803
KB5017315	KB5007192	KB5000685

KB4601393	KB4580346	KB4561616
KB4601392	KB4580325	KB4561608
KB4601345	KB4577668	KB4561600
KB4601318	KB4577667	KB4558998
KB4601060	KB4577586	KB4558997
KB4598480	KB4577015	KB4557397
KB4598243	KB4576750	KB4556813
KB4598230	KB4571694	KB4552924
KB4593226	KB4570333	KB4551853
KB4592440	KB4570332	KB4550994
KB4589210	KB4566424	KB4550929
KB4589208	KB4565912	KB4549949
KB4587735	KB4565511	KB4549947
KB4586830	KB4565349	KB4540723
KB4586793	KB4562562	KB4540670
KB4583465	KB4562561	KB4539571
KB4538461	KB4519998	KB4503308
KB4537764	KB4519338	KB4498947
KB4537759	KB4516115	KB4494440
KB4535680	KB4516044	KB4494175
KB4534273	KB4514366	KB4494174
KB4534271	KB4512937	KB4493470
KB4532691	KB4512578	KB4489882
KB4530715	KB4512577	KB4487026
KB4530689	KB4512574	KB4485447
KB4525236	KB4512517	KB4346087
KB4523205	KB4511553	KB4132216
KB4523204	KB4509091	KB4091664
KB4521862	KB4507460	KB4022619
KB4521858	KB4506998	KB4018073
KB4520724	KB4505225	KB2565063

KB2467173