



Date: February 11, 2025

Author: Software Engineering

System: Path Vision

Subsystem: Windows 10 LTSC 2019 (v1809) & Software Security Patches

Subject: Windows 10 LTSC 2019 (v1809) & Software Security Patch release January 2025

Purpose:

To release the list of validated Microsoft critical patches on the Path Vision System.

Note: Core Vision and BioVision are on similar product line as Path Vision and both use the same executable so patches validated on Path Vision system can be installed on Core Visions system and BioVision as well.

List of patches/KB's validated for January 2025:

KB5050008

KB5049608

List of all patches/KB's validated on the Path Vision System:

KB890830

KB5039879

KB5032196

KB5050008

KB5039217

KB5031984

KB5049608

KB5037765

KB5031361

KB5048661

KB5036896

KB5030214

KB5046615

KB5036604

KB5029931

KB5044277

KB5035849

KB5029247

KB5044016

KB5034768

KB5028168

KB5043050

KB5034127

KB5027222

KB5041578

KB5033904

KB5027131

KB5040430

KB5033371

KB5026362

KB5025229

KB5019966

KB5013641

KB5023702

KB5018419

KB5012170

KB5022840

KB5017315

KB5005112

KB5022511

KB5016623

KB4589208

KB5021237

KB5015811

KB4577586

KB5020866

KB5014692

KB4535680

KB5020627

KB5013941

KB4494174