



**Date:** August 06, 2021

**Author:** Software Engineering

**System:** BREVERA / CORLUMINA

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release July 2021

---

**Purpose:**

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Note: Below listed patches for this month are cumulative and include “Print Nightmare” fixes. Additional testing to test “Print Nightmare” fixes has been performed for this particular patch validation.

**List of patches/KB’s validated for July 2021:**

KB5004244

**List of all patches/KB’s validated on the BREVERA / CORLUMINA System:**

KB890830	KB5000822	KB4586793
KB5004244	KB5000685	KB4580325
KB5003711	KB4601393	KB4577668
KB5003646	KB4601345	KB4577667
KB5003243	KB4601060	KB4577586
KB5003171	KB4598480	KB4570333
KB5001407	KB4598230	KB4570332
KB5001404	KB4592440	KB4561600
KB5001342	KB4589208	KB4557397
KB5000859	KB4587735	KB4552924

KB4539571

KB4535680

KB4494174