



**Date:** *May 28, 2021*

**Author:** Software Engineering

**System:** BREVERA / CORLUMINA

**Subsystem:** Windows 10 LTSC 2019 (v1809) & Software Security Patches

**Subject:** Windows 10 LTSC 2019 (v1809) & Software Security Patch release May 2021

**Purpose:**

To release the list of validated Microsoft critical patches on the BREVERA / CORLUMINA System.

Note: Install SFW-03452 twice to ensure all patches are applied. The first execution of the patch disk installs the newly required Service Stack and a second run of the patch disk is required to install the latest cumulative patch for May

**List of patches/KB's validated for May 2021:**

KB5003243	KB5003171	KB5001407
-----------	-----------	-----------

**List of all patches/KB's validated on the BREVERA / CORLUMINA System:**

KB890830	KB4601393	KB4580325
KB5003243	KB4601345	KB4577668
KB5003171	KB4601060	KB4577667
KB5001407	KB4598480	KB4577586
KB5001404	KB4598230	KB4570333
KB5001342	KB4592440	KB4570332
KB5000859	KB4589208	KB4561600
KB5000822	KB4587735	KB4557397
KB5000685	KB4586793	KB4552924

KB4539571

KB4535680

KB4494174