

IHE Integration Statement for Selenia Dimensions

IHE Integration Statement			
Vendor	Product	Version	Date
Hologic, Inc.	Selenia Dimensions Acquisition Workstation (AWS)	1.0.x	December 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Assisted Acquisition Protocol Setting	
		PPS Exception Management Option	
		Broad Worklist Query	
		Patient Based Worklist Query	
Patient Information Reconciliation	Acquisition Modality	None	
Mammography Image Profile	Acquisition Modality	None	
Mammography Image Profile	Print Composer	None	
Access to Radiology Information	Image Display	None	
Portable Data for Imaging	Portable Media Creator	None	
Links to Standards Conformance Statements for the Implementation			
HL7	Not applicable		
DICOM	www.hologic.com/product-support/digital-mammography/dimensions/		
Links for general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	